

(200)
R290

UNITED STATES DEPARTMENT OF THE INTERIOR
GEOLOGICAL SURVEY

A Program for Mass Spectrometer Control and
Data Processing Analyses in Isotope
Geology--Written in BASIC for
an 8K Nova 1210 Computer

By

J. S. Stacey
U.S. Geological Survey

and

J. Hope
Bernard Price Institute of Geophysics
University of Witwatersrand
Johannesburg, South Africa

Open-file report 75-127
1975

This report is preliminary and has not
been edited or reviewed for conformity
with U.S. Geological Survey standards
and nomenclature.



CONTENTS

	Page
Abstract-----	1
Introduction-----	2
Mass spectrometer control system-----	3
Functions of the program-----	4
Method of Computation-----	6
Loading and running the program from paper tape-----	9
Modifications to demonstrate the operations of the program-----	11
BASIC software description-----	13
Program structure-----	13
Parameter symbols used in INPUT subroutine-----	16
Computations during data acquisition-----	17
Baseline and rubidium corrections to peak heights-----	17
Computation of ratios-----	17
Assembly language driver program-----	19
References-----	20
Appendixes-----	21
A. Basic program listing.-----	21
B. Input - output listing for a "STRONTIUM" analysis using synthetic data.-----	32
C. Assembly language driver listing-----	36
D. Interface hardware diagrams and connections using the NOVA general interface board.-----	40

Illustrations

Figure 1. Diagram for ratio computations in a 5-cycle block using 3 peaks in the main switching sequence-----	18a
--	-----

Tables

	Page
Table 1. Connections for input register-A general purpose interface board-----	46
2. Connections for input register-B general purpose interface board-----	47
3. Connections for output register general purpose interface board-----	48

A PROGRAM FOR MASS SPECTROMETER CONTROL AND DATA PROCESSING
ANALYSES IN ISOTOPE GEOLOGY--WRITTEN IN BASIC
FOR AN 8K NOVA 1210 COMPUTER

By

J. S. Stacey
U.S. Geological Survey

and

J. Hope
Bernard Price Institute of Geophysics
University of Witwatersrand
Johannesburg, South Africa

ABSTRACT

A system is described which uses a minicomputer to control a surface ionization mass spectrometer in the peak switching mode, with the object of computing isotopic abundance ratios of elements of geologic interest. The program uses the BASIC language and is sufficiently flexible to be used for multiblock analyses of any spectrum containing from two to five peaks. In the case of strontium analyses, ratios are corrected for rubidium content and normalized for mass spectrometer fractionation. Although almost any minicomputer would be suitable, the model used was the Data General Nova 1210 with 8K memory. Assembly language driver program and interface hardware descriptions for the Nova 1210 are included.

INTRODUCTION

This program was written while J. S. Stacey was visiting the Bernard Price Institute of Geophysics at the University of Witwatersrand in Johannesburg, South Africa. Although there was not sufficient operational time available to incorporate all possible levels of sophistication, the program does offer a very flexible and comprehensive package which could be used by almost any minicomputer. This article may therefore be of interest to people in other isotope geology laboratories. The main concepts were developed over a period of several years at the U.S. Geological Survey's laboratory in Denver, Colorado, using a variety of instruments and computer systems and we particularly wish to acknowledge the work of R. Terrezas and E. E. Wilson in this regard. The program described here is, however, the first we have written in BASIC language--which is well suited for laboratory use. Changes in BASIC are very easy to make and the format is easy to understand. The version of BASIC used "BASIC with CALL package," was written for the Data General Nova 1210, by Perseus Computing and Automation (Pty) Ltd., P. O. Box 3987, Pretoria, South Africa. It facilitates communication between BASIC and Nova's assembly language and considerably simplified the writing of the driver program to control the mass spectrometer. The program will handle any spectrum containing up to five peaks in a switching mode. Included therefore are analyses of uranium, thorium, rubidium, lead, and strontium. In the case of strontium, corrections are made for rubidium and the ratios are normalized for fractionation.

Features which could easily be incorporated are "peak top and baseline filtering" and the facility to reject data during an analysis (Stacey and others, 1971; 1972). Sensitivity range switching for the vibrating reed electrometer could easily be added if desired.

Actually, since this report was initiated, the program described here has been further developed and adapted to run on the Hewlett Packard HP9830 calculators now used in Denver for mass spectrometer control. Some simplification has resulted by writing the strontium option as a separate program from that for the other elements. The new listings are available from our Denver laboratory.

These new programs have retained the structure of the original program, and rather than outmoding the material in this report they serve to emphasize its fundamental simplicity and flexibility.

MASS SPECTROMETER CONTROL SYSTEM

The mass spectrometer at the Bernard Price Institute of Geophysics is provided with magnetic field control by a Bell Inc., type 640, incremental gaussmeter. 5 independent peak channels are offered by 5 potentiometers. Each potentiometer can be selected separately by its own relay operated from the computer. Two other potentiometers are provided so that magnetic field offsets can be adjusted for all peaks simultaneously. These offset positions are selected by two relays and correspond to "baseline below peak" (offset reduces field), "on peak" (zero offset), "baseline above peak" (offset increases field) respectively. The circuitry arrangement for this is shown in Appendix D. Altogether, a total of 7 relays are driven from the computer--5 for the peaks and 2 for the offsets.

Analog signals from the Carey Model 31 vibrating reed electrometer (1 volt for full scale), are measured by a voltage-to-frequency converter and counter combination (Hewlett Packard HP221 B and HP5326A respectively). Normally signals are integrated over a period of one second. Therefore almost once per second the computer is required to read a six-character BCD measurement from the counter. The A-register on the Nova general purpose interface board is a 16-bit device and is arranged to read the 4 least significant digits. The 2 most significant digits are read by 8 of the 12 bits of the B-register.

The C-register forms the "interblock print-time counter" which counts the number of measurements made by the counter(that is, the number of seconds taken) during data printout between blocks in a mass spectrometer analysis. Both the B and C registers were assembled on the Nova general purpose interface board.

FUNCTIONS OF THE PROGRAM

Any type of analysis involving up to ^{five} peaks may be made. Therefore analyses of U, Th, Rb, Pb, and Sr are included. In analyses where three or more peaks are measured, #3 peak is always the reference peak--that is, all ratios are measured with respect to this peak. Thus, in a lead analysis for instance, Pb^{206} may conveniently be assigned to #3 channel and the others as below:

^{204}Pb to #1 or #2 channels

^{206}Pb to #3 channel

^{207}Pb to #4 channel

^{208}Pb to #5 channel

For strontium analyses it is important to assign the isotopes as indicated below because of the ^{85}Rb correction and the normalizing computations.

^{84}Sr to #1 or #2 channels

^{85}Rb to #2 or #1 channels

^{86}Sr to #3 channel

^{87}Sr to #4 channel

^{88}Sr to #5 channel

In strontium analyses the ^{85}Rb correction peak is measured only at the beginning and end of each block--as are the baselines for each of the peaks in the spectrum.

A block comprises baseline measurements for each peak taken below only, above only, or above and below each peak as requested by the operator. For strontium analyses the ^{85}Rb peak is then measured with its 'below' baseline. This is followed by measurement of the main peaks switched through in sequence as many times as has been requested by the operator at the beginning of the analysis. In strontium analyses, after the main peak switching, the ^{85}Rb peak and its baseline are measured and the block is ended with a final set of baselines for each of the main peaks in the same sequence as at the beginning of the block.

If more than one block is to be measured, then while the teletype is printing ratios at the end of a block, the C-register built into the interface keeps track of time. Time is clocked by the number of measurements made by the main signal measuring counter (at the rate of one per second in our case). Time for the printout between blocks is read from the C-register and is added to the total time elapsed at the start of the next block of data. This enables the baselines at the end of the n th block to be used as the first baselines in the $(n + 1)$ th block and so on throughout the analysis.

The time spent on each peak (that is, the number of measurements made at one per second) can be different for each peak, as also can be the times (skip times) allowed for switching transients at the beginning of each peak. Each baseline measurement is given 8 seconds longer than its own peak measurement and 2 extra seconds are skipped to allow switching transients to fully decay. Similarly, extra periods are given for the ^{85}Rb peak over those ascribed to the ^{87}Sr peak.

As far as the computer is concerned there is no limit to the time which can be requested for a peak--but blocks may not consist of more than 10 switching cycles through the main peaks and an analysis may not comprise more than 10 blocks.

Between blocks, individual ratios and their times are printed, as also are their means and standard deviations. The data are filtered by rejecting ratios more than 2 standard deviations from the mean. The rejected ratios and times are printed, followed by the new mean and standard deviation. The process is repeated until no ratios remain outside the 2 sigma limit.

At the end of an analysis a table is printed of all the mean ratios for the blocks and the table ends with the final means and standard deviations for the whole analysis.

For strontium analyses, rubidium corrections are made for each block. Because of memory storage limitations with an 8 K memory--only one normalized $^{87}\text{Sr}/^{86}\text{Sr}$ ratio is computed for each block from the mean $^{87}\text{Sr}/^{86}\text{Sr}$ and $^{88}\text{Sr}/^{86}\text{Sr}$ ratios ($^{86}\text{Sr}/^{88}\text{Sr}$ is set to 0.1194). No spiking calculations are included in the program when the sample contains ^{84}Sr .

METHOD OF COMPUTATION

Baselines for each peak are stored as straight lines (slope and intercept, m_b and c_b respectively). These lines join the mean of the measurements at the beginning and end of each block. They have the form

$$y_b = m_b t + c_b.$$

Peak height information is also stored as straight lines (slope and intercept, m_p and c_p respectively). These lines join consecutive measurements of the same peak and are represented by

$$y_p = m_p t + c_p.$$

Peaks are corrected for baselines simply by subtracting the slope and intercept parameters in the following manner,

$$y'_p = m'_p t + c'_p$$

where $m'_p = m_p - m_b$
and $c'_p = c_p - c_b$.

In a strontium analysis, a similar procedure is used to correct the ^{87}Sr peak for the ^{87}Rb , based on the line parameters computed from the ^{85}Rb measurements,

$$y'_{87} = m'_{87} t + c'_{87}$$

where $m'_{87} = m_{87} - \frac{m_{85}}{2.59}$
and $c'_{87} = c_{87} - \frac{c_{85}}{2.59}$.

In each switching cycle a ratio is computed for each peak relative to the reference peak. Thus:

$$r(t_{\text{ref}}) = \frac{y_p(t_{\text{ref}})}{y_{\text{ref}}(t_{\text{ref}})}$$

$$= \frac{m'_p(t_{\text{ref}}) + c'_p}{m'_{\text{ref}}(t_{\text{ref}}) + c'_{\text{ref}}}.$$

An additional ratio for each peak relative to the reference peak is computed at the occurrence of each peak.

$$r(t_p) = \frac{y_p(t_p)}{y_{\text{ref}}(t_p)}$$

$$= \frac{m'_p(t_p) + c'_p}{m'_{\text{ref}}(t_p) + c'_{\text{ref}}}.$$

If there are N switching cycles in a block, a total of $2(N-1)$ ratios is obtained for each peak per block. This method of computation helps to eliminate the effects of curvature in the decay curve of the ion beam signal, without the redundancy of computing all ratios at the occurrence of each peak.

LOADING AND RUNNING THE PROGRAM FROM PAPER TAPE

After loading the CALL BASIC interpreter, press the STOP switch and check that the data switches are set to X7777. Load the assembled driver tape called "PACKAGE.LDA" in the teletype reader and press RESET and START. The driver tape will be read into core. Restart BASIC via address 00002 on the data switches and load the main control program "NEW" from the teletype reader.

The following conversation should then take place between the operator and the computer. Underlined parameters are typed by the operator. Assume we wish to run a strontium sample spiked with ^{84}Sr .

RUN

PKS? 4 4 main peaks in the switching sequence

? 1 ? 3 ? 4 ? 5 Assign ^{84}Sr to #1 channel
 ^{86}Sr to #3 channel
 ^{87}Sr to #4 channel
 ^{88}Sr to #5 channel

TIMES

? 20 ? 20 ? 20 ? 15 Total times per peak - including skip time
 Assign 20 seconds to peaks #1, #3, and #4
 Assign 15 seconds to peak #5

SKIP

? 4 ? 4 ? 4 ? 3 Ignore first 4 seconds for peaks #1, #3, and #4
 Ignore first 3 seconds for peak #5
 Total useful time per peak (TIME-SKIP).

B/L S? 0

Type 0 for baselines to be measured both sides of peak

Type 1 for baselines to be measured below peaks

Type 2 for baselines to be measured above peaks

RB#? 2

Assign ^{85}Rb to channel #2

If your reply is 0 the computer assumes this is not a strontium analysis being run

C/S? 5

Any number of CYCLES from 2 to 10 may be requested here

BLKS? 5

Any number of BLOCKS from 1 to 10 may be asked for here

After typing RETURN the computer will begin the analysis.

MODIFICATIONS TO DEMONSTRATE THE OPERATION OF PROGRAM

following

The/list of modifications below will enable the program to be run with the ordinary BASIC interpreter and without the driver program. A set of synthetic data will be generated having the following properties:

- (1) All main peak heights will be measured as 11
- (2) All baselines will be measured as 1
- (3) The #2 peak will be measured as 3.59

When baselines are subtracted all peaks except #2 will have a height of 10 and that of #2 will be 2.59.

If a strontium analysis is called with ^{85}Rb assigned to channel #2 and the other peaks assigned as described earlier, then the ^{87}Sr peak will be corrected by the amount

$$\left(\frac{1}{2.59}\right) \times ^{85}\text{Rb} = 1$$

$$\therefore ^{87}\text{Sr corrected} = 9$$

$$\text{and } ^{87}\text{Sr}/^{86}\text{Sr} = 0.9, \quad ^{88}\text{Sr}/^{87}\text{Sr} = 1.0$$

$$\text{and } ^{84}\text{Sr}/^{86}\text{Sr} = 1.0 \text{ if a spiked run is asked for.}$$

The $^{87}\text{Sr}/^{86}\text{Sr}$ ratio will be normalized by setting $^{86}\text{Sr}/^{88}\text{Sr} = 0.1194$

The statements to modify the program are:

```
710 LET R8 = 0
```

```
825 LET X = 11
```

```
910 LET R8 = 0
```

```
8095 and 8029 LET Y2 = 0
```

```
876 LET L(I,J1) = 1
```

```
891 IF A(I)<>R GOTO 893
```

```
892 LET H(I,J) = 3.95
```

If you wish to observe the complete switching sequence for an analysis, the following PRINT statements will print peak numbers and mean measurements as the analysis proceeds:

```
881 PRINT "881"; A(I); L(I,J1); M(I,J1)
```

```
896 PRINT "896"; A(I); H(I,J); Q(I,J)
```

Appendix B shows an input-output listing for a "strontium" analysis using the synthetic data described here.

BASIC SOFTWARE DESCRIPTION

PROGRAM STRUCTURE

	Line Numbers
<u>Initializing</u>	20 - 499
Dimensions, Set Counters, zeros arrays	
<u>Data Acquisition</u>	
"INPUT"	500 - 698
Acquires sequence control data for the analysis from the operator (see p. 16 for parameter symbols)	
"PEAK SWITCH"	700 - 740
Subroutine - does just that.	
"READ"	800 - 885
Routine for reading counter	
"RUBIDIUM"	900 - 950
Switches and reads rubidium peaks	
"BASELINE"	1000 - 1080
Switches and reads baselines	
<u>Data Processing</u>	
"RB LINE"	1140 - 1190
Slopes and intercepts computed for rubidium peak	
"BASELINE LINES"	1200 - 1298
Slopes and intercepts computed for baselines	

"PEAK LINES" 1300 - 1450

Slopes and intercepts computed and
stored for peaks

"RATIOS-REFERENCE" 1600 - 1720

Ratios for each peak with respect to
the reference peak computed at each
time of occurrence of the reference
peak

"RATIOS-PEAKS" 1800 - 1870

A ratio for each peak with respect to
the reference peak is computed at each
of the times of its own occurrence.

Headings printed for ratios; arrays
loaded and printed 3120 - 3340

"FPRINT" 3500 - 3960

Means and standard deviations printed.
Individual ratios then filtered by
rejecting those with deviations more
than 2σ from the mean. Rejected ratios
and their times are printed. New means
and σ s printed. The process is repeated
until no ratios lie outside 2σ from the
recomputed mean.

"NORMALIZE" .4000 - 4040

Normalizes the mean $^{87}\text{Sr}/^{86}\text{Sr}$ for the
block, assuming $^{88}\text{Sr}/^{86}\text{Sr} = 0.1194$.

"APRINT"

5000 - 5310

Headings and table printed for mean ratios from each block at end of analysis. Analysis may comprise up to 10 blocks.

"MAIN"

8000 - 8500

Routine for controlling overall sequence of program

On the Nova computer BASIC software, all arrays contain a "zeroth" position, and these positions have been utilized in our program to save memory. Thus, most counting parameters such as N, I, and B2 start at "zero" and are one less than their physical counterparts.

N is (N1-1)

T[I] Total time to be spent on Ith peak

S[I] Total SKIP time on Ith peak

B1	Baselines to be measured	"Below each peak"	B1 = 1
		"Both sides of each peak"	B1 = 0
		"Above each peak"	B1 = 2

R Channel number for rubidium-85

If $R = 0$, the analysis is NOT strontium.

C1 Number of switching cycles per block

B2 Number of blocks (less one)

COMPUTATIONS DURING DATA ACQUISITION

During data acquisition, all baseline and peak height measurements are accumulated and stored as mean values of the appropriate groups of incoming data--together with their mean times. They are stored as $L(I,J)$ and $M(I,J)$ respectively for baselines and $H(I,J)$, $Q(I,J)$ respectively for peak heights and times (lines 875-895 in the program). Note that accumulation does not start until after points to be "skipped" at the beginning of each group have been ignored.

BASELINE AND RUBIDIUM CORRECTIONS TO PEAK HEIGHTS

After the acquisition of data for a block is completed, parameters (slope and intercept) for lines joining mean baseline measurements for each peak at end of the block are calculated ($M1$ and $L1$, lines 1280, 1285). Slope and intercept for each line joining consecutive measurements of the same peak are computed as $Q1$ and $P1$ - lines 1320, 1331. These last two parameters are corrected for baseline by subtracting $M1$ and $L1$ (lines 1380, 1390) and stored as $G(I,K)$ and $H(I,K)$. I is the peak number (less one) and K is the line number (starting at $K = 0$). There are $C1-1$ lines for each peak where $C1$ is the number of switching cycles within a block.

Rubidium corrections to ^{87}Sr peaks are made in the same manner as baseline corrections - at lines 1360, 1370.

COMPUTATION OF RATIOS

Ratios are calculated for all peaks with respect to the reference peak at each occurrence of the reference peak. (Usually/^{the} reference peak channel number is 3-- / ^{that is,} $R2 = 2$). These ratio calculations are made in subroutine RATIOS-REF (1600 - 1620).

Within each switching cycle, another ratio for each peak with respect to the reference peak is calculated--at the occurrence of that peak (RATIOS-PEAKS, lines 1800-1870).

In each of the above instances two cases arise, depending on whether the peak number is greater than or less than that of the reference peak (that is, whether $I > R2$ or $I < R2$). Consider calculation of ratios at the reference peaks in subroutine RATIOS-REF. Let $I < R2$ so that $R2=2$ and $I=1$ as in figure 1. Starting at the FIRST occurrence of the reference peak ($C1-1$) ratios will be calculated at times $Q(R2,J)$ for $J=0$ to $J=(C1-2)$. The peak line parameters $H(I,J)$, $G(I,J)$ for the $I=1$ peak and $H(2,J)$, $G(2,J)$ for the reference peak are used in each calculation. For peaks where $I > R2$ ($I=3$ as in figure 1) the calculations are similar except that we start at the SECOND occurrence of the reference peak. Thus, ($C1-1$) ratios are computed at times $Q(2,J1)$ for $J1=(C1-1)$

Subroutine RATIOS-PEAK calculates individual ratios at the occurrence of the peaks themselves in a similar manner. ($C1-1$) ratios are computed for each peak, the first ratio being calculated at the first or second occurrence of that peak--depending on whether $I > R2$ or $I < R2$ respectively.

From the two subroutines, a total of $2*(C1-1)$ ratios are calculated for each peak in the main switching sequence, with respect to the reference peak.

● MEAN PEAK HEIGHT

○ INTERPOLATED VALUE

USED IN COMPUTING RATIOS.




FIG. 1 Diagram for ratio Computations in a 5-cycle block using 3 peaks in the main switching sequence.

ASSEMBLY LANGUAGE DRIVER PROGRAM

three
The complete driver comprises / subroutines labelled "START", "SIG",
and "BCNT".

START operates the peak SWITCH and baseline offset relays in the mass spectrometer, under control of the main BASIC program. It selects the relays to be closed from the statement: "710 CALL 0, D, A[I]-1" in which the '0' parameter calls SIG itself. The D parameter (value 8, 9 or 10) determines the selection of baseline below, baseline above or peak top. A[I] is the peak relay number. These last / ^{two} parameters are read by the assembly language statement GET. START then sets up the necessary bit pattern to operate the relays via the relay driver circuits in the interface unit.

The BASIC statement "825 CALL I, X" calls subroutine SIG and when the main counter is DONE, 6 BCD digits are read from the interface (4 from the ^{the} A register and 2 from/B register). This string is transferred into the BASIC program as a decimal number X by the statement 'PUTD'. The counter is reset by SIG to permit the next reading to accumulate in the counter.

When the final subroutine BCNT is called by "8095 CALL 2, 1" the interblock print-time register BCTR is initiated. It keeps track of time by counting the number of measurements made by the main counter. At the end of printout between blocks, statement ^{"8233} / CALL 2, -1, Y2" stops BCTR. The primary count from BCTR is returned to the BASIC program as parameter Y2 by the assembly language statement PUT.

REFERENCES

- Stacey, J. S., Wilson, E. E., Peterman, Z. E., and Terrazas, R., 1971,
Digital recording of mass spectra in geologic studies--I: Canadian Jour.
Earth Sci., v. 8, no. 3, p. 371-377.
- Stacey, J. S., Wilson, E. E., and Terrazas, R., 1972, Digital recording
of mass spectra in geologic studies--II: Canadian Jour. Earth Sci.,
v. 9, no. 7, p. 824-834.

APPENDIX A

BASIC PROGRAM LISTING

INITIALIZE .

```

20 DIM A(5),T(5),S(5),P(5,10),Q(5,10),L(5,5),M(5,5),G(5,10)
21 DIM H(5,10),W(5,10),E(18,2),B(5),C(5)
25 LET B3=-1
122 FOR I=0 TO 4
123     IF B3>= 0 GOTO 140
125 LET A(I)=0
126 LET T(I)=0
127 LET S(I)=0
129 LET B(I)=0
130 LET C(I)=0
140 FOR J=0 TO 9
142 LET G(I,J)=0
143 LET W(I,J)=0
146 LET Q(I,J)=0
147 LET H(I,J)=0
148     IF B3<> 0 GOTO 160
149     LET P(I,J)= 0
150     IF J>=4 GOTO 160
156     LET L(I,J)= 0
157     LET M(I,J)= 0
160 NEXT J
165 NEXT I
170 FOR I=0 TO 17
175 FOR K=0 TO 1
180 LET E(I,K)=0
182 NEXT K
183 NEXT I
185 IF B3<>-1 GOTO 8029
190 LET Z2= 0
191 LET J1= 0
192 LET J= 0
196 LET F=-1
197 LET I2= 0
499 GOTO 8000

```

```

500 REM      I/P - - - - - INPUT
505 PRINT " PKS";
510 INPUT N
515 LET N1=N
516 LET N=N-1
520 PRINT
530 FOR I=0 TO N
540 INPUT A(I)
550 NEXT I
555 PRINT
560 PRINT " TIMES"
570 FOR I=0 TO N
580 INPUT T(I)
590 NEXT I
600 PRINT
610 PRINT " SKIP"
620 FOR I=0 TO N
630 INPUT S(I)
640 NEXT I
650 PRINT
651 PRINT " B/LS";
652 INPUT B1
653 PRINT
655 PRINT " RB#";
656 INPUT R
657 IF R= 0 GOTO 664
658 LET A(N1)=R
660 LET S(N1)=S(N-1)+1
662 LET T(N1)=T(N-1)+8
664 PRINT
665 PRINT " C/S";
670 INPUT C1
675 PRINT
680 PRINT " BLKS";
690 INPUT B2
692 LET B2=B2-1
695 PRINT
698 RETURN

```

```

700 REM SW- - - - -PEAK SWITCH
705 FOR I=0 TO N
710 CALL 0,D,A(I)-1
720 GOSUB 800
730 NEXT I
740 RETURN
800 REM RD- - - - -READ
805 LET Y= 0
806 LET Z1=Z2
810 LET K= 0
812 IF D=10 GOTO 820
814 LET U=2+S(I)
815 LET V=8+T(I)
818 GOTO 825
820 LET U=S(I)
822 LET V=T(I)
825 CALL 1,X
828 LET Z2=Z2+1
830 LET K=K+1
840 IF K<=U GOTO 825
850 LET Y=Y+X
860 IF K<V GOTO 825
870 IF D=10 GOTO 890
875 LET L(I,J1)=Y/(K-U)
880 LET H(I,J1)=Z1+(K+U)/2
882 IF F<>-1 GOTO 935
885 RETURN
890 LET H(I,J)=Y/(K-U)
893 IF H(I,J)>=0 GOTO 895
894 LET H(I,J)=0
895 LET Q(I,J)=Z1+(K+U)/2
898 IF F<>-1 GOTO 946
899 RETURN

```

```

900  REM          RB- - - - - - - - - - -RUBIDIUM
901  LET I=N1
902  IF J<> 0 GOTO 952
903  LET F= 0
904  LET D=8
906  LET S1=J1
908  LET J1=F
910  CALL 0,D,R-1
920  GOTO 805
935  LET D=10
940  LET S2=J
942  LET J=F
944  GOTO 910
946  LET J1=S1
948  LET J=S2
949  LET F=-1
950  RETURN
952  LET F=1
954  GOTO 904
1000  REM          B/LS - - - - - - - - - -BASELINES
1003  IF J= 0 GOTO 1010
1006  LET J1=J9
1010  IF B1=2 GOTO 1020
1015  LET D=8
1017  GOTO 1030
1020  LET D=9
1030  GOSUB 700
1035  LET J1=J1+1
1040  IF B1<> 0 GOTO 1070
1050  LET D=D+1
1060  GOSUB 700
1062  LET J1=J1+1
1070  IF J<> 0 GOTO 1080
1075  LET J9=J1
1080  RETURN

```

```

1140 REM          RB/L - - - - - - - - - -RB LINE
1145 IF R= 0 GOTO 1200
1150 LET I=N1
1155 LET K= 0
1160 GOTO 1268
1170 LET S1=M1
1175 LET S2=L1
1178 IF B3= 0 GOTO 1320
1180 GOTO 1316
1185 LET S1=Q1-S1
1190 LET S2=P1-S2
1200 REM          B/L LS- - - - - - - - - -BASELINE LINES
1201 IF R= 0 GOTO 1210
1203 IF B1<> 0 GOTO 1210
1205 LET L(N, 0)=L(N,1)
1207 LET L(N,2)=L(N,3)
1210 FOR I=0 TO N
1230     IF B1<> 0 GOTO 1268
1235     LET K= 0
1240     FOR K1=0 TO 2 STEP 2
1242         IF B3= 0 GOTO 1250
1247         LET K1=2
1248         LET K=1
1250         LET L(I,K)=(L(I,K1)+L(I,K1+1))/2
1255         LET M(I,K)=(M(I,K1)+M(I,K1+1))/2
1260         LET K=1
1265     NEXT K1
1268     IF B3= 0 GOTO 1280
1270     LET L(I, 0)=B(I)
1275     LET M(I, 0)=C(I)
1280     LET M1=(L(I, 0)-L(I,1))/(M(I, 0)-M(I,1))
1285     LET L1=L(I, 0)-M1*M(I, 0)
1290     LET B(I)=L(I,1)
1295     LET C(I)=M(I,1)
1298     IF I=N1 GOTO 1170
1300 REM          PK/LS - - - - - - - - - -PEAK LINES
1310 FOR K=0 TO C1-2
1312     IF R= 0 GOTO 1320
1314     IF B3= 0 GOTO 1320
1316     LET H(N1, 0)=H8
1318     LET Q(N1, 0)=Q8
1320     LET Q1=(H(I,K)-H(I,K+1))/(Q(I,K)-Q(I,K+1))
1330     LET P1=H(I,K)-Q1*Q(I,K)
1331     IF I<>N1 GOTO 1340
1332     LET H8=H(N1,1)
1333     LET Q8=Q(N1,1)
1335     GOTO 1185
1340     IF R= 0 GOTO 1380
1350     IF I<>N-1 GOTO 1380
1360     LET Q1=Q1-S1/2.59
1370     LET P1=P1-S2/2.59
1380     LET G(I,K)=Q1-M1
1390     LET H(I,K)=P1-L1
1395 NEXT K
1399 NEXT I

```

```

1410 IF R= 0 GOTO 1450
1420 IF N<>3 GOTO 1450
1430 LET R2=I
1440 GOTO 1600
1450 LET R2= 0
1600 REM R/RF - - - - -RATIOS REFERENCE
1620 FOR I=0 TO N
1630 IF I=R2 GOTO 1870
1650 FOR J= 0 TO C1-2
1660 IF I<R2 GOTO 1700
1670 LET J1=J+1
1690 GOTO 1715
1720 LET J1=J
1715 LET H1=G(R2,J)*Q(R2,J1)+H(R2,J)
1720 LET H2=(G(I,J)*Q(R2,J1)+H(I,J))/H1
1800 REM R/PKS
1815 IF I>R2 GOTO 1840
1820 LET J2=J+1
1830 GOTO 1848
1840 LET J2=J
1848 LET H1=G(I,J)*Q(I,J2)+H(I,J)
1850 LET W(I,J)=H1/(G(R2,J)*Q(I,J2)+H(R2,J))
1855 LET H(I,J)=H2
1860 NEXT J
1870 NEXT I
3120 FOR I= 0 TO N
3140 IF I=R2 GOTO 3950
3150 PRINT
3155 PRINT
3160 PRINT A(I);"/";A(R2);"RATIOS-";B3+1
3177 PRINT
3178 LET K= 0
3180 FOR J= 0 TO C1-2
3190 IF I>R2 GOTO 3280
3210 LET E(K, 0)=H(I,J)
3220 LET E(K,1)=Q(R2,J)
3230 LET E(K+1, 0)=W(I,J)
3240 LET E(K+1,1)=Q(I,J+1)
3250 PRINT E(K, 0);E(K,1)
3260 PRINT E(K+1, 0);E(K+1,1)
3270 GOTO 3330
3280 LET E(K, 0)=W(I,J)
3290 LET E(K,1)=Q(I,J)
3300 LET E(K+1, 0)=H(I,J)
3310 LET E(K+1,1)=Q(R2,J+1)
3320 GOTO 3250
3330 LET K=K+2
3335 NEXT J
3340 PRINT

```



```

3500 REM          FPT- - - - -FPRINT
3502 LET K1=2*(C1-1)-1
3505 LET J=-1
3511 LET S5= 0
3512 LET T5= 0
3513 LET S7= 0
3520 FOR K= 0 TO K1
3530 LET S5=S5+E(K, 0)
3540 LET T5=T5+E(K,1)
3550 NEXT K
3555 LET K=K1+1
3560 LET S5=S5/K
3570 LET T5=T5/K
3580 FOR K= 0 TO K1
3590 LET S6=E(K, 0)-S5
3610 LET S7=S7+S6*S6
3640 NEXT K
3650 LET S6= SQR (S7/K1)
3655 IF B3>B2 GOTO 5080
3657 PRINT S5;S6;T5
3660 PRINT
3662 PRINT
3665 IF S6=0 GOTO 3925
3670 IF J=K1 GOTO 3925
3680 IF K1<5 GOTO 3925
3790 FOR K= 0 TO K1
3800 LET D5=E(K, 0)-S5
3810 IF ABS (D5)<2*S6 GOTO 3860
3830 PRINT E(K, 0);E(K,1);D5;S6
3840 PRINT
3850 GOTO 3900
3860 LET J=J+1
3870 IF K=J GOTO 3900
3880 LET E(J, 0)=E(K, 0)
3890 LET E(J,1)=E(K,1)
3900 NEXT K
3910 IF K1=J GOTO 3925
3915 LET K1=J
3920 GOTO 3511
3925 LET P(I,B3)=S5
3940 LET Q(I)=S6
3945 LET W(I)=T5
3950 NEXT I
3960 RETURN
4000 REM          NORM- - - - -NORMALIZE
4010 LET P(N1,B3)=P(1,B3)*2/(1+.1194*P(2,B3))
4015 PRINT "87/86-N IS";P(N1,B3)
4020 PRINT
4025 PRINT
4030 PRINT
4040 RETURN

```

```

5000 REM APT - - - - -APRINT
5001 PRINT
5002 PRINT
5003 PRINT,
5004 LET N2=N
5005 IF R= 0 GOTO 5010
5006 IF N=3 GOTO 5010
5007 LET N2=N1
5010 FOR I= 0 TO N2
5020 IF I=R2 GOTO 5110
5030 FOR K2= 0 TO B2
5040 LET E(K2, 0)=P(I,K2)
5050 NEXT K2
5060 LET K1=B2
5070 GOTO 3505
5080 LET Q(I)=S6
5090 LET H(I)=S5
5095 IF R= 0 GOTO 5105
5096 IF I2=2 GOTO 5110
5098 IF N=2 GOTO 5102
5100 PRINT"84/86","87/86","88/86",
5101 GOTO 5109
5102 PRINT"87/86","88/86","87/86-N",
5103 GOTO 5109
5105 PRINT A(I);"/";A(R2),
5109 LET I2=2
5110 NEXT I
5115 IF I>=5 GOTO 5119
5117 PRINT
5119 PRINT
5135 FOR K= 0 TO B2
5137 PRINT,
5140 FOR I= 0 TO N2
5153 IF I=R2 GOTO 5170
5162 PRINT P(I,K),
5170 NEXT I
5175 IF I>=5 GOTO 5190
5180 PRINT
5190 NEXT K
5201 PRINT
5202 PRINT " MEANS",
5210 FOR I= 0 TO N2
5225 IF I=R2 GOTO 5240
5230 PRINT H(I),
5240 NEXT I
5245 IF I>=5 GOTO 5255
5250 PRINT
5255 PRINT " ST.DEVS.",
5260 FOR I= 0 TO N2
5275 IF I=R2 GOTO 5290
5280 PRINT Q(I),
5290 NEXT I
5300 PRINT
5310 RETURN

```

```

8000  REM  MAIN - - - - - - - - - - -MAIN
8002  IF B3<>-1 GOTO 8010
8004  LET B3= 0
8010  GOSUB 500
8020  GOSUB 1000
8025  FOR B3=0 TO B2
8026  IF B3=0 GOTO 8032
8027  GOTO 122
8029  CALL 2,-1,Y2
8030  LET Z2=Z2+Y2
8032  IF R=0 GOTO 8045
8035  GOSUB 900
8045  LET D=10
8055  FOR J=0 TO C1-1
8060  GOSUB 700
8070  NEXT J
8080  IF R= 0 GOTO 8090
8085  GOSUB 900
8090  GOSUB 1000
8095  CALL 2,1
8100  GOSUB 1140
8105  IF R=0 GOTO 8130
8110  IF N=3 GOTO 8130
8115  GOSUB 4000
8130  NEXT B3
8132  LET B3=B3+1
8135  IF B2=0 GOTO 8145
8140  GOSUB 5000
8145  FOR J=1 TO 10
8150  PRINT
8155  NEXT J
8160  GOTO 25
8165  END

```

APPENDIX B
INPUT - OUTPUT LISTING
FOR A
"STRONTIUM" ANALYSIS
WITH "SYNTHETIC DATA".

RUN

PKS?3 - - - - - 3 peaks in main switching sequence

?3 86Sr → #3
?4 -Peak channel Nos. 87Sr → #4
?5 88Sr → #5

TIMES - - - - - for each peak

?15

?18

?10

SKIP - - - - - for each peak

?4

?4

?3

B/L5?0 - - - - - baselines - both sides of each peak

RB#?2 - - - - - ⁸⁵Rb → #2 peak channel

C/S?3 - - - - - 3 switching cycles

BLKS?2 - - - - - 2 blocks

(1) (2) (3) (4) - - - - - Column No.

881	3	1	14.5		
881	4	1	39		
881	5	1	60.5	Baselines	(1) Print statement #
881	3	1	81.5	Main Peaks	
881	4	1	106		(2) Peak channel #
881	5	1	127.5		
881	2	1	154.5	Peak and Baseline	(3) Mean measurements
896	2	3.59	183.5	Rb	1 assigned to Baselines
896	3	11	203.5		
896	4	11	220		11 assigned to main peaks
896	5	11	233.5	Synthetic data	
896	3	11	246.5		3.59 assigned to ⁸⁵ Rb
896	4	11	263	Main Peaks	
896	5	11	276.5		
896	3	11	289.5		(4) Mean times for peak
896	4	11	306		
896	5	11	319.5		measurements.
881	2	1	343.5	Peak and Baseline	
896	2	3.59	372.5	Rb	
881	3	1	397.5		
881	4	1	422		
881	5	1	443.5	Baselines	
881	3	1	464.5	Main Peaks	
881	4	1	489		
881	5	1	510.5		

4 / 3 RATIOS- 1

RATIO TIME

.9 220
.9 246.5
.9 263
.9 289.5

.9	.5376993E-9	254.75
Mean	Standard	Mean
Ratio	Deviation	Time

5 / 3 RATIOS- 1

1 233.5
1 246.5
1 276.5
1 289.5

1 0 261.5

87786-N IS 1.608004

2nd BLOCK STARTS

881 2	1	537.5	Peak and Baseline Rb
896 2	3.59	566.5	
896 3	11	586.5	
896 4	11	603	
896 5	11	616.5	
896 3	11	629.5	Main Peaks
896 4	11	646	
896 5	11	659.5	
896 3	11	672.5	
896 4	11	689	
896 5	11	702.5	Peak and Baseline Rb
881 2	1	726.5	
896 2	3.59	755.5	
881 3	1	780.5	
881 4	1	805	
881 5	1	826.5	Baselines Main Peaks
881 3	1	847.5	
881 4	1	872	
881 5	1	893.5	

4 / 3 RATIOS- 2

.9 603

.9 629.5

.9 646

.9 672.5

.9 .5376993E-9 637.75

5 / 3 RATIOS- 2

1 616.5

1 629.5

1 659.5

1 672.5

1 0 644.5

87/86-N IS 1.608004

FINAL TABLE OF RATIOS

	87/86	88/86	87/86-N
	.9	1	1.608004 ----- 1st Block
	.9	1	1.608004 ----- 2nd Block
MEANS	.9	1	1.608004
ST. DEVS.	0	0	0

PKS? - - - - - START NEXT ANALYSIS

APPENDIX C
ASSEMBLY LANGUAGE DRIVER
LISTING

;PERSEUS COMPUTING AND AUTOMATION (PTY) LTD
;P. O. BOX 3987, PRETORIA, SOUTH AFRICA
;COPYRIGHT (O) 1971.

;SYMBOL DEFINITIONS FOR 'CALL' PACKAGE
;(NO HIGH SPEED READER OR PUNCH)
;USE WITH SINGLE-USER BASIC KEY-B3 AND WHEN ASSEMBLING
DRIVER PROGRAM.

000300	DUSR	BUS=300	
006156	DUSR	WATT=JSR @156	
007620	DUSR	MASK=7770	
007771	DUSR	BEGIN=10127	
003170	DUSR	INIT=3170	
007000	DUSR	GET=JSR @0,2	
007001	DUSR	PUT=JSR @1,2	
003002	DUSR	RETURN=JMP @2,2	
007003	DUSR	PUTD=JSR @3,2	
007004	DUSR	GETD=JSR @4,2	
007005	DUSR	GETF=JSR @5,2	
007006	DUSR	PUTF=JSR @6,2	
	EOT		
000265		ERRE=265	;ERROR TERMINATE
000267		ERRS=267	;ERROR CONTINUE

DRIVER PROGRAM LISTING.

```

---      000000      BCTR = 60
          000060      CLR = 60
          000060      DRIVER = 60
          000300      .LUC .BUS
000300 010136      FREE
          007771      .LUC BEGIN
07771 000003      3
07772 007775      START
07773 010033      SIG
07774 010124      BCNT
07775 037050      START: GET
07776 054430      STA 3,STORE
07777 102520      SUBEL 0,0
10000 124000      COM 1,1
10001 125405      LOOP: INC 1,1,SNR
10002 000403      JMP DONE
10003 101120      MOVEL 0,0
10004 000775      JMP LOOP
10005 024422      DONE: LDA 1,CURR
10006 034422      LDA 3,LOW
10007 117415      AND# 0,3,SNR
10010 000412      JMP BLINE
10011 034420      LDA 3,HIGH
10012 167400      BACK: AND 3,1
10013 107000      ADD 0,1
10014 044413      STA 1,CURR
10015 065060      DOA 1,DRIVER
10016 034410      LDA 3,STORE
10017 175405      INC 3,3,SNR
10020 000404      JMP FIN
10021 000754      JMP START
10022 034406      BLINE: LDA 3,LOW
10023 000767      JMP BACK
10024 024406      FIN: LDA 1,BGN
10025 003002      RETURN
10026 000000      STORE: 0
10027 032000      CURR: 2000
10030 000377      LOW: 377
10031 177400      HIGH: 177400
10032 010032      BGN: .+1
10033 054773      SIG: STA 3,STORE
10034 034447      LDA 3,MASK
10035 056456      STA 3,@DIGIT
10036 060260      NIUC CLR
10037 063660      SKPDN CLR
10040 000777      JMP .-1
10041 060560      DIAS 0,CLR
10042 000402      JMP .+2
10043 061560      BREG: DIAS 0,CLR
10044 024437      NEXT: LDA 1,MASK
10045 107400      AND 0,1
10046 010445      ISZ DIGIT
10047 046444      STA 1,@DIGIT
10050 101200      MOVN 0,0
10051 010434      ISZ CN11
10052 000776      JMP .-2


```

10053	034431		LDA 3,RS11
10054	054431		STA 3,CN11
10055	010432		ISE CN12
10056	000766		JMP NEXT
10057	034427		LDA 3,RS12
10060	054427		STA 3,CN12
10061	010427		ISE CN13
10062	000761		JMP BREG
10063	034421		LDA 3,RS11
10064	054423		STA 3,CN12
10065	034421		LDA 3,RS12
10066	054422		STA 3,CN13
10067	034737		LDA 3,STONE
10070	010076	LEAD:	FEED
10071	020777		LDA 0,LEAD
10072	007003		PUTD
10073	024417		LDA 1,RS1
10074	044417		STA 1,DIGI1
10075	003002		RETURN
10076	022415	FEED:	LDA 0,@DIGI1
10077	024412		LDA 1,ANSCI
10100	123000		ADD 1,0
10101	014412		DSE DIGI1
10102	001400		JMP 0,3
10103	000017	MASK:	17
10104	177774	RS11:	177774
10105	177774	CN11:	177774
10106	177776	RS12:	177776
10107	177774	CN12:	177774
10110	177776	CN13:	177776
10111	000060	ANSCI:	60
10112	010114	RS1:	+.2
10113	010113	DIGI1:	STACK
	000010	STACK:	.BLK 10
10124	007000	BCNT:	GET
10125	125415		INC #1,1 SER
10126	000400		JMP ON
10127	066460		DIC 1,BCTR
10130	102400		SUB 0,0
10131	007000		PUT
10132	060360		NIOF BCTR
10133	000402		JMP OUT
10134	060360	ON:	DIC 0,BCTR
10135	003002	OUT:	RETURN
	101036		FREE=.
	003170		.END .INIT

APPENDIX D

INTERFACE HARDWARE DIAGRAMS AND CONNECTIONS USING NOVA GENERAL INTERFACE BOARD

PEAK SWITCH RELAY CIRCUIT





RL6, RL7 MADE	BELOW PEAK	BIT 9
RL6 MADE	ON PEAK	BIT 11
RL6, RL7 OPEN	ABOVE PEAK	BIT 10 (not used)

RV1 - RV5 1KΩ 10-turn

RV6, RV7 500Ω 10-turn

INCREASE
RESISTANCE AND FIELD,
CLOCKWISE





B-REGISTER

C-REGISTER

TABLE 1--CONNECTIONS FOR INPUT REGISTER-A GENERAL PURPOSE INTERFACE BOARD.

SIGNAL NAME	WRAP PINS		SLOT 9 PIN #	MOTHER BOARD EDGE CONNECTOR PIN #
M.S.B. Digit 4 Ctr.	12A	41	A67	11
	12	49	A73	8
	11	49A	A75	7
L.S.B. Digit 4 Ctr.	10A	57	A79	19
M.S.B. Digit 3 Ctr.	22	67	A81	20
	21	68A	A83	22
	20A	68	A85	24
L.S.B. Digit 3 Ctr.	20	69	A87	26
M.S.B. Digit 2 Ctr.	28A	70A	A89	27
	28	70	A91	3
	27	71	B6	29
L.S.B. Digit 2 Ctr.	26A	72A	B11	30
M.S.B. Digit 1 Ctr.	36	72	B13	31
	36A	76	B15	32
	35	83	B19	33
L.S.B. Digit 1 Ctr.	34	89	B25	35

TABLE 3--CONNECTIONS FOR OUTPUT REGISTER GENERAL PURPOSE INTERFACE BOARD.

SIGNAL NAME	WRAP PIN FROM TO	SLOT 9 PIN #	MOTHER BOARD EDGE CONNECTOR PIN #
RLY 1	24 131	B34	38
RLY 2	24A 132A	B36	39
RLY 3	25 132	B38	40
RLY 4	26 133	B40	41
RLY 5	23 134A	B48	42
RLY 6 & 7	40 98A	B49	43
RLY 6	76A 135	B53	46