

U. S. DEPARTMENT OF THE INTERIOR

GEOLOGICAL SURVEY

RDARL4, a FORTRAN interface for transferring
chemical analytical data from an Applied Research Laboratories
electron microprobe to a PDP-11 computing system

by

J. Stephen Huebner

Open-File Report 83-713

This report is preliminary and has not been reviewed for conformity with U.S. Geological Survey editorial standards. Any use of trade names is for descriptive purposes only and does not imply endorsement by the USGS.

CONTENTS

Page

Introduction -----	1
Hardware considerations -----	1
Description of program -----	3
Operator interaction during execution -----	7
Compiling the program -----	8
References cited -----	9
Appendix I. Description of auxillary mineral programs -----	10
Appendix II. Listing of RDARL4 -----	13
Appendix III. Examples of output:	
a. Hardcopy from line printer -----	29
b. Format of disk file input for FORTRAN program MINCLC -----	30
c. Format of disc file input for FORTRAN program PXVOID -----	30
d. Format of disk file input for FORTRAN program RECAMP -----	31
e. Format of disk file input for BASIC program OXIDES -----	31
f. Description of direct access file (virtual array) -----	32
g. Format of BASIC sequential file for statistics -----	32
h. Format of disk file input for BASIC mineral formula programs -----	33
i. Format of disk file for analyses of "known-unknowns" -----	34
j. Format of BASIC sequential file for archiving data -----	35
Appendix IV. On-line instructions for using RDARL4 -----	36

TABLE

Table 1. Structure of file of analyses created by microprobe -----	38
--	----

INTRODUCTION

Microprobe analysis at the U.S. Geological Survey does not end when a complete chemical analysis is obtained by the operator. Rather, an increasing number of USGS operators want to reconstitute their chemical analyses as atomic mineral formulas, norms of endmember mineral compositions, or compositional projections. In the future, operators will want to select analyses to construct data files for plotting and for statistical analyses. The traditional method of manually transcribing analyses from microprobe printout to another computing system is laborious and prone to transcription errors. RDARL4 is a machine interface between data files created by the microprobe and the BASIC and FORTRAN programs available in Reston for recalculating and plotting microprobe data.

The electron microprobe used in the USGS's Reston Microprobe Facility is a model SEMQ (Scanning Electron Microprobe Quantometer) manufactured by Applied Research Corporation (ARL) of Sunland, California. This instrument collects chemical data for up to 9 elements simultaneously or up to 14 elements asynchronously. Chemical analyses are reduced on line and may be stored on disk if the operator so chooses. Operators may obtain 100 to 170 analyses per 4 hour operating shift. Despite the high rate of productivity, demand for access to the microprobe is great. The efficient use of both the microprobe instrumentation and the large volumes of data obtained has been a major concern of the Reston Microprobe Facility staff. Use of the program RDARL4 permits operators to transfer large volumes of data to other computers quickly, accurately, and without competing for valuable instrument time.

HARDWARE AND SOFTWARE CONSIDERATIONS

The ARL-SEMQ microprobe automated control system is built around the PDP-11 series of computers and peripherals manufactured by the Digital Equipment Corporation (DEC), specifically the LSI-11/23 processor and RX-02 floppy disk controller and drive. The PDP-11 series of equipment is commonly used for control of major laboratory instrumentation and is available "off-line" elsewhere in the microprobe facility. DEC makes available software based upon its RT-11 operating system. This software permits data to be stored in words that are 1, 2, 4, or 8 bytes in length. The microprobe operating programs are written in Applied Research Laboratories Extended Basic (ARLEB) language and use assembly language subroutines and functions that can be called from ARLEB. Although the ARL system utilizes a conventional DEC processor and diskettes, ARL's software for disk storage uses a different disk directory and writes only 6-byte data words (called "entries"), a word length that is not supported by DEC. Stored on a RX-02 disk, an ARLEB data record is 42 6-byte entries long. Each entry contains one numeric value (or up to six ASCII characters).

A compact ARLEB disk storage routine was written in ARLEB language by James J. McGee of the USGS and inserted into the microprobe control program. This routine stores all the analytical information that most USGS operators want to retain. The first record of the data storage file contains instrumental parameters that are common to all microprobe analyses that follow in the file. Each subsequent record contains information for one analysis: title, stage position coordinates, and chemical analysis in terms of weight percent and cation formula unit (Table 1). Disk storage files are created during a microprobe operating session; after the operating shift the operator removes the data disk and either stores it for later use or goes to an off-line computer to retrieve the data.

Off line computing equipment utilizes the DEC LSI-11/23 processor, RX-02 disk controller and disk drives, a video terminal, and a line printer. The off-line equipment is completely compatible with the computing equipment that controls the microprobe. This interchangeability of hardware means that the only interfacing problems involve using the DEC RT-11 family of software to access a disk written under ARLEB. Fortunately, ARL provides "RT-LIB", a library of compiled MACRO (DEC assembly language) program modules that can be called from a FORTRAN program running under the RT-11 operating system whenever it is necessary to access the ARLEB disk. Once microprobe data has been read into an RT-11 based system, it is possible to create new disk files of any of a variety of file types: sequential (ASCII character) with either FORTRAN or BASIC carriage control and direct access (binary). This enables a USGS scientist to use any of a number of different programs (Appendix I) to recalculate the analytical data.

Program RDARL4 will read a data file that contains 100 or fewer analyses of 1 to 14 elements each. It will run under the RT-11 operating system in a PDP-11 minicomputer that has 64 kbytes of lower program memory and 18.2 kbytes of virtual memory for numeric and character arrays. The version described here is a general "brute-force" program designed to run quickly and to offer many options which can be executed without rereading the ARLEB data disc. The program could be easily modified to run in a smaller memory by decreasing the number of elements allowed per analysis, limiting the number of analyses that can be processed at one time, decreasing the number of output options, and/or by making uniform the data input requirements of the programs listed in Appendix I. Although the FORTRAN portion of the code could be easily modified to run under another operating system, numerous calls are made to modules of the "RT-LIB" library which, because it is written in PDP-11 assembly language, is not readily transportable to systems that do not use the PDP-11 family of processors. Conversion to a non-DEC system would be difficult and timeconsuming.

DESCRIPTION OF PROGRAM RDARL4

The program (Appendix II) has four major parts: definition and initialization of arrays and variables; reading of the ARLEB data disk; creation of output files; and subroutines to report disk errors and to handle commonly used interactions with the operator.

Program RDARL4 uses buffer arrays (DATA, IBUF) to hold information from a single record of the ARLEB data disc file from the time the record is read until its contents can be transferred to larger arrays that store all the titles, weight percent values, cation numbers, and stage coordinates, etc., that appear in a disk file (up to 101 records). These storage arrays have mnemonic names (for instance, ITITLE, WTPCT, PEROXY, CATION) and will be discussed in conjunction with program execution, below. Also important are the byte arrays (such as BPXVD and BOXIDE) that are used to store element names in the sequence that they will be written out for a given output option. These byte arrays are filled with the necessary ASCII characters using DATA statements.

Initial program execution involves keyboard input of the name of the ARLEB data file. The program then calls FIND, an RT-LIB module that accesses the ARLEB disk directory and returns the position, length (ISIZE), and ARLEB file-type (ITYPE, BFKIND), and an error code (IERR). The next call to RT-LIB, OPEN, opens the designated file for reading and returns only an error code. Each record is then read twice because there is no way to read a partial record, yet each record contains both ASCII and numeric information which must be handled differently. First, FDGET is used to retrieve ASCII information by returning the contents of the record, byte-by-byte, in an integer buffer array IBUF. Next, FFETCH reads the ARLEB record as a block of 42 6-byte numeric entries, converts each entry to a 4-byte floating point format, and returns the converted contents in the floating point buffer array DATA.

The first record of the ARLEB data file contains operating parameters that are common to all subsequent analyses, including the analyst's name, date of analysis, names of elements analyzed (in the sequence in which concentration data are to be stored), a 4-character abbreviation for the standard used for each element, and the X-ray line used to analyze each element. The complete list of contents of this first record is given in Table 1. FORTRAN code to select ASCII and numeric information for the first record from the buffer arrays IBUF and DATA and to place that information in more permanent variables and storage arrays in memory is used only once during each program execution.

Each subsequent record contains information that pertains to a single analysis. The program contains a loop to read and store these records, one-by-one. Exit from the loop is determined either by reaching the last block of the ARLEB file (ISIZE) or by encountering an analysis with a weight percent sum equal to zero. (The test for this condition is included because an ARLEB disk file is not necessarily completely filled with analyses.)

The information retrieved and stored is also listed in Table 1 and for each analysis includes a title, weight percent values, cation and anion numbers, and a flag (IAVLFG) which indicates analysis type and number of points averaged.

Execution of RDARL4 next gives the operator the choice of 11 output formats that may be used in almost any combination. The only restrictions are that reviewing the analyses on the video terminal, if that option is chosen, must be first and that sending the analyses to a lineprinter must precede creation of output files on disk. Simple reordering of blocks of code (and some renumbering of statements) would remove this constraint on operations.

The video terminal option is useful for operators who want to review the analyses in a data file, perhaps because they do not remember exactly the contents of a particular disk file. The format of the display is compact and includes all analyses in sequence. The lines displayed exceed 80 characters in length; operators who want to use this option may first want to set the video terminal for a 130 character line length. The only information displayed from the first record is the element names. The remaining information stored on the first record can be examined with the lineprinter (see below).

The line printer option gives a hard copy of the contents of the first record (operating parameters common to all analyses), then lets the operator specify the particular analyses to be printed (subroutine SELECT). Analyses need not be specified if only a record of the microprobe operating parameters is desired. An example of the printout is included in Appendix IIIa.

The following 9 options create disk files that can be used as input to supporting programs for mineral formula calculations, plotting, and other uses. (A tenth option, EXIT, terminates program execution.) Existing programs for which 5 of the 9 options were specifically prepared are briefly summarized in Appendix I. Data input to most of these programs is in the form of sequential files so that they may be viewed (on a video terminal or printer) and edited. The choice of sequential files also permits an operator to create an input file with an editor, without using RDARL4. To facilitate editing operations, these sequential files include the symbols of elements analyzed, whether or not an auxiliary program actually uses the element symbols (as in Program MINCLC) or merely reads and ignores them (as in the Programs which use disk files created by Options 2, 3, 4, and 7).

All auxiliary mineral formula programs except MINCLC require concentration data that is predetermined with respect to the nature of required elements and the sequence in which they appear. The microprobe operating system has no constraints on the number, nature, or sequence of elements analyzed and stored in an ARLEB data file. The reading and storage of data by RDARL4, as previously described, preserves the original number, nature, and sequence of elements in the byte array BELNAM. The new

sequence in which element data (one datum is two ASCII character element symbols followed by a floating point numeric concentration value) are to be put out is determined by comparing BELNAM, byte-by-byte, with an array such as BOXIDE, BPCVD, BCAMP, or BBASIC that contains the desired element symbols in the correct sequence:

```

DO 965, I=1,13
IFLAG=0
DO 965, J=1,NELMTS
IF(BELNAM(1,J).NE.BOXIDE(1,I)) GOTO 963
IF(BELNAM(2,J).NE.BOXIDE(2,I)) GOTO 963
KCOUNT(I)=J
IFLAG=1
963 IF(IFLAG.EQ.1) GOTO 965
KCOUNT(I)=16
965 CONTINUE

```

Whenever two pairs of bytes match (the element symbols are the same), the position of the desired element is stored in array KCOUNT. If no match is found (microprobe analysis did not include an element expected by the auxillary program), the array element is set to 16. Thus KCOUNT contains the sequence in which element concentrations are to be written out from WTPCT and CATIONS. A value of 16 will cause a WTPCT value of 0.00 and a CATION concentration of 0.000 to be written to the disk file. Array KCOUNT is used in implied DO LOOPS to control the actual sequence in which data are written from BELNAM and WTPCT or CATIONS to disk:

```

968 WRITE(2,971) N,IAVFLG(N),(ITITLE(L,N),L=1,24)
DO 970, I=1,13
WRITE(2,972)(BOXIDE(J,I),J=1,2),WTPCT(KCOUNT(I),N)
970 CONTINUE

```

The routine for constructing an input file for the general mineral formula program MINCLC is different. MINCLC can be configured to operate upon different sets of elements; at present the set consists of elements with atomic numbers 1 to 30 (hydrogen through zinc), although it could be changed or enlarged. During execution, MINCLC identifies the elements to which each concentration datum pertains, regardless of the number of the elements or their sequence in the input file. Program RDARL4 does not, therefore, contain code to write out elements in a particular sequence. Rather, the code is designed to write out only a number of data fields corresponding to the number of elements that were actually analyzed; records that contain fewer than 8 elements contain "\$\$" to indicate the last element-concentration pair of an analysis. RDARL4 also writes an initial record that tells MINCLC what to do ("****RDOXIONS" plus the number of anions and cations in the formula unit) and a title (from ITITLE) that precedes each analysis. Here another advantage of creating an output file that can be edited becomes apparent: if the operator includes more than one kind of mineral in a file created for MINCLC, he may want to change the formula basis by inserting additional control records (indicated by the initial "****") in that file.

Four disk file options were not designed to create an input file for a mineral recalculation program. Option 5 (BASIC virtual array) creates a direct access file for programs that have not yet been written. Both weight percent values and cation numbers are put out to the disk file. Option 6 creates a sequential file that can be used as input to the MINC BASIC statistical programs provided by DEC. Either weight percent values or cation numbers are put out. Option 8 (filing of "known-unknown" analyses) was written so that all analyses of well-characterized material (used to gage the quality of microprobe operation) can be filed for later use. Two uses are anticipated: detect and examine changes in microprobe performance with time; and evaluate the suitability of particular standards for particular elements and mineral groups. Option 8 stores operating parameters and weight percent values, but not cation numbers. Option 9 creates a sequential file that can be edited, searched, and concatenated using a keypad editor or BASIC language routines written by Elaine McGee of the USGS. In this manner, an operator can sort analyses into groups for plotting or archiving.

Four segments of code have been placed in subroutines because they are either unlikely or very likely to be called into use. Errors in reading the ARLEB disk or writing the DEC disk are rare, but if an error is encountered, the main program calls subroutine ERRORS, prints a brief description of the particular disk error encountered, and either halts execution (ARLEB disk error) or returns to the main program (DEC disk error). The other subroutines are used frequently.

Subroutine SETUP reminds the operator of the name of the ARLEB file by displaying it on the videoscreen, then asks the operator to provide a name for an output file that is to be written on disk. Because analyses from a microprobe operating shift can be used many different ways, operators may want to retain the ARLEB file name and append a file type that indicates the intended use of the output file. Thus, the ARLEB input filename "D1027" could be used in output files for several mineral recalculation programs, for example, D1027.MNC, D1027.PXV, and D1027.FEL.

Subroutine SELECT asks the operator to select the analyses that will be sent to the line printer or written to a disk file. The operator has several options: averaged analyses of unknowns; all the known-unknown analyses; or analyses, individually or in groups, designated by the sequence number (NANAL) in which they appeared in the ARLEB file and are, during execution, stored in arrays in memory. Each analysis is associated with a variable (IAVFLG) that indicates how many points were averaged to produce the analysis. Furthermore, if that analysis is of a "known-unknown", that is, an analysis of a well characterized material such as a microprobe standard that is used to gage overall system performance, the value 80 has been added to the number of points averaged to indicate that the analysis is of a "known-unknown". Thus averages of unknowns are those analyses with IAVFLG greater than 1 but less than 80 and "known-unknowns" are analyses with IAVFLG greater than or equal to 80.

Subroutine CHOICE gives the operator the possibility of adding analyses to the current output file; that file is not closed until the operator responds "N" to the question posed by CHOICE.

OPERATOR INTERACTION DURING EXECUTION

Operators should have some understanding of the flow or sequence of routines before attempting to run RDARL4. In Reston, the disk which contains the program contains a file of instructions (Appendix IV) that are printed at the video-terminal whenever the computer is started. The actual program begins by instructing the operator to place an ARLEB disk in drive #1, the right-hand drive on an RX-02 disk unit. After reading the file, the operator is first asked if the microprobe analyses are to be viewed on the video terminal; next, whether or not a printout is desired; and finally, whether or not a disk output file is to be written. If a disk file is to be created, the operator is given the choice of 9 file options or an exit from the program. Each of the disk options ends with an opportunity to select another option so that more than one kind of disk file can be produced without rereading the ARLEB disk. If another option is not selected, the operator has an opportunity to loop back to the beginning of the program to start again with a new ARLEB data file.

During execution, RDARL4 provides the operator with numerous prompts for information. These prompts are designed to be self explanatory, particularly if the operator has run the program before. However, operators will want to prepare in advance to respond to questions that involve more than a Y (yes), N (no), or C (continue) response. These questions concern:

Name of ARLEB data file (5 characters maximum)

Use to be made of the disk output file, if any (option number)

Output filename in format DYn:xxxxxx.YYY
where n=0 or 1 (left or right disc drive, respectively),
and YYY represents a file type (default DAT). Subroutine SETUP.

Cation to anion ratio (Option 1 - MINCLC output file) in the floating point format (F12.7, F12.7). For example, "4.0, 6.0" is an acceptable response

Weight % or cations (Option 6)

Number of analyses to be put out (Option 6)

Choice of individual analyses, averages only, or both
(Option 8 - Extraction of "known-unknown" analyses)

Choice of selecting averaged unknown analyses, all "known-unknown"
analyses, or designated analyses for most options. (Subroutine Select)

If the operator intends to use RDARL4 to create an output file on a disk in drive #1, that disk must have already been formatted when it is exchanged for the ARLEB disk during execution.

COMPILING THE PROGRAM

RDARL4 was compiled without overlays under DEC's RT-11 single job monitor that had been patched to load at address 28000 (to enable support for virtual memory arrays). The program is long; as listed in Appendix II the "no line-numbers" option of the compiler must be used. A command file for compiling and linking is given below. RTLIB and RX02 are contained in an ARL software kit, V.1.2-99.

```
DATE
TIME
FORTRAN/EXTEND/LIST:TT:/SHOW:0/NOLINENUMBERS/OBJECT:DY1: RDARL4
RLINK
RDARL4=DY1:RDARL4,DY0:RTLIB,RX02
C
TIME
```

References

- Smyth, Joseph R. (1980) Cation vacancies and the crystal chemistry of breakdown reactions in kimberlitic omphacites. *American Mineralogist*, vol. 65, p. 1185-1191.
- Stormer, J.C., Jr. (1983) The effects of recalculation estimates of temperature and oxygen fugacity from analyses of multicomponent oxides. *Amer. Mineral.*, vol 68, 586-594.

Appendix I: Brief descriptions of auxilliary mineral formula programs

In this appendix, references to variable and array names will be capitalized as they are in Program RDARL4, Appendix II. Examples of printable input files are given in Appendix III; the input files are the disk files generated by RDARL4.

Option 1. MINCLC is a FORTRAN program for converting a weight percent analysis to an oxygen-based formula unit. The program will attempt to achieve a specified cation-to-anion ratio by adjusting the proportions of the various valence states of the multivalent elements Ti, Cr, Mn, and Fe, if present. If the specified cation:anion ratio is 4:6, the program will calculate an idealized pyroxene site occupancy and pyroxene formula norm. It will also accept anions which substitute for oxygen. MINCLC was written in 1977 for punched card input by Phelps Freeborn, then of the USGS. When it was converted to run on an RT-11 based system by Huebner, January 25, 1982, the options were simplified and data input was made to be from a disk file, but the data format was left unchanged:

```
****RDOXIONS      OXYgens      CATions
NANALysis, TITLE: ITITLE
up to eight pairs of element symbols (BELNAM) and WTPCT
ditto; input terminated by $$
NANAL (etc.)
```

The first byte in the input file is a non-printing FORTRAN carriage control character, ASCII line feed.

Option 2. PXVOID is a FORTRAN program that recalculates pyroxene analyses to a pyroxene formula unit based upon six oxygens. Unlike MINCLC, PXVOID considers the possibility that deviations from the ideal 4:6 cation:oxygen ratio of pyroxenes are real and caused by nonstoichiometry. As originally written by Smyth (1980), data was input from the keyboard during execution. When converted to run under RT-11 by Huebner, March 7, 1982, data input was changed to be from disk file. The file format is similar to MINCLC with the following exceptions. The control record beginning with **** is replaced by an integer, 1 to indicate oxide weight percent input or 0 to indicate element weight percent (the corresponding RDARL4 variable is NTYPE). Element symbols are written in the required sequence from array BPXVD. The first byte in the input file is the non-printing FORTRAN carriage control character, ASCII line feed.

Appendix I (continued)

Option 3. RECAMP is a FORTRAN program that recalculates an amphibole analysis to several different formula units, using different assumptions. This program was originally coded by F. Spear and K. Kimball of the Massachusetts Institute of Technology, November 10, 1981, and requested that data be input from the keyboard. For USGS use, the program was modified by J.S. Huebner, March 3, 1982, to accept data from a disk file. The disk input file is particularly simple:

```
NANALysis   ITITLE
  9 pairs of element symbol (BCAMP) and oxide WTPCT.
```

Because the order of the elements is fixed, it is not necessary to include element symbols in the list. They were retained to facilitate examination and editing of the input file. The first byte of the input file is the non-printing FORTRAN carriage control character, ASCII line feed.

Option 4. OXCALC (following Stormer, 1983) is a BASIC program that converts microprobe chemical analyses of spinel and ilmenite to formula units and mole fractions of ulvospinel and magnetite in preparation for estimates of the temperature and pressure of last equilibration. OXCALC was modified by J.C. Stormer, February 18, 1983, for input of data from disk. Program RDARL4 contains a format statement that suppresses the initial non-printing FORTRAN carriage control character so that the file can be read under the BASIC operating system:

```
OPEN(UNIT=2,NAME=FILSPC,CARRIAGECONTROL='FORTRAN')
WRITE(2,966) NTYPE
966  FORMAT('+',I1)
```

The analyses which follow consist of 14 records each:

```
NANALysis, ITITLE
  13 pairs of element symbol (BOXIDE) and WTPCT
```

Option 7. AMPHI, MICA, CORD, FELDS, CHLOR, and GARN, are BASIC language programs written by Marta Kempa Flohr of the USGS to calculate mineral formulas, including formal site occupancies, for amphiboles, micas, cordierites, feldspars, chlorites, or garnets, respectively. These programs require oxide or element weight percent concentration data for 18 elements in the order in which the element symbols appear in byte array BBASIC. Program RDARL4 writes a sequential file without the initial byte of FORTRAN carriage control. The contents of the file are as follows:

```
NTYPE      (0 = element, 1 = oxide weight percent)
BBASIC (18 records, each containing an element symbol)
NANAL, ITITLE for first analysis in file
```

IQFILE (name of standard file), PEROXY (number of oxygens in cation's formula unit)

WTPCT value for silicon (first element), CATION value for silicon

. ,
. ,
. ,
. ,
. ,

WTPCT value for chlorine (last element), CATION value for chlorine
NANAL, ITITLE for second analysis in file

Option 9. This option creates sequential files on disk that can be concatenated and edited. Analytical data is read out in the order in which it was stored on the ARLEB disk. The initial byte of FORTRAN carriage control is suppressed so that the file created with this option can be searched and concatenated with a BASIC language program or an editor.

NELMTS, ***, PEROXY

BELNAM elements analyzed, listed one element-per-record

NANAL, ITITLE, PEROXY, IAVFLG

WTPCT, CATION first element is silicon

. ,
. ,
. ,
. ,
. ,

WTPCT, CATION last element is chlorine

WTPCT (15), CATION (15) weight percent and cation sums

NANAL, ITITLE

The file does not contain an initial FORTRAN carriage control character.

Appendix II. Listing of Program RDARL4

C RDARL4.FOR PERMITS AN RT-11 BASED SYSTEM TO READ AN ARL DATA DISK.
 C RDARL.OBJ LINKS TO ARL'S RTLIB, A COLLECTION OF COMPILED MACRO UTILITIES
 C THAT PERFORM VARIOUS DISK OPERATIONS AND CONVERT NUMERIC DATA FROM THE ARL
 C 3-WORD TO THE DEC 2-WORD FP FORMAT WHILE PRESERVING ASCII INFORMATION.
 C
 C WRITTEN BY J S HUEBNER AS RDARL1.FOR MARCH, 1982 (NUMERIC ONLY) AND
 C AS RDARL2.FOR ON JULY 13, 1982 (NUMERIC & ASCII)
 C OUTPUT OPTIONS BY HUEBNER, WRITTEN AS NEEDED. LAST EDIT 06/15/83
 C

C ARLEB USES A BLOCK OF 42 6-BYTE ENTRIES. FORMAT OF DATA FILED BY THE ARLEB
 C STORAGE ROUTINE WRITTEN BY MCGEE FOR THE USGS RESTON PROBE FACILITY:
 C

C FIRST RECORD CONTAINS OPERATING PARAMETERS:

ENTRY NUMBER	CONTENTS	RDARL4 VARIABLE NAME
1 - 4	OPERATORS' NAME	IOPNAM
5 - 7	ARL SYSTEM DATE	IARDT
8	NAME OF *QFILE	IQFILE
9	NAME OF DATA OUTPUT FILE	IARLOT
10	DATA REDUCTION FLAG	NREDCT
11	NUMBER OF ELEMENTS ANALYZED	NELMTS
12	OXIDE OR ELEMENT WT % FLAG	NTYPE
13	TEMPORARY COUNTER USED BY ARL	NCOUNT
14	FUTURE USE	NOUSED
15 - 28	ELEMENT NAME,X-RAY LINE	BELNAM,BXRLIN
29 - 42	STANDARD NAME	BSTAND

C SUBSEQUENT RECORDS CONTAIN ONE ANALYSIS EACH

1 - 8	TITLE	ITITLE
9	STAGE POSITION X	X
10	STAGE POSITION Y	Y
11	AVERAGING & KNOWN-UNKNOWN FLAG	IAVFLG
12 - 25	WEIGHT PERCENT	WTPCT
26	TOTAL WEIGHT PERCENT	WTPCT
27 - 40	CATION NUMBERS PER FORMULA UNIT	CATION
41	SUM OF CATIONS	CATION
42	NUMBER ANIONS PER FORMULA UNIT	PEROXY

```

DIMENSION DATA(42)          !buffers numeric records from FETCH
DIMENSION IBUF(128)         !buffers ASCII records from FDGET
INTEGER IOPNAM(12),IQFILE(3),IARLOT(3),KCOUNT(14),UPPER,NAME(3)
INTEGER*4 ITIME
BYTE BFILNM(15),BFKIND(18),BSTAND(4,14),BBUF(256)
BYTE BELNAM(2,14),BXRLIN(2,14),FILSPC(15),BFLAG1,BFLAG2,BFLAG3
BYTE BRPLY0,BPXVD(2,13),BCAMP(2,9),BBASIC(2,18),BOXIDE(2,13)
VIRTUAL ITITLE(24,100),WTPCT(16,100),CATION(16,100)
VIRTUAL NANAL(100),IAVFLG(100),PEROXY(100),X(100),Y(100)
EQUIVALENCE (NAME(1),BFILNM(5)),(BBUF(1),IBUF(1))
COMMON BFLAG1,BFLAG2,BFLAG3
DATA BFILNM/'D','Y','1',' ',' ','D','E','F','A','L','T',6*0/
DATA BFKIND/'B','A','S','I','C',' ',' ','S','Y','S','T','E','M',
&'D','A','T','A',' ',' '/
DATA BPXVD/'S','I',' ',' ',' ','A','L','T','I','C','R','F','C',
&'F','E','M','G','M','N','C','A','L','I','N','A','K',' '/
DATA BCAMP/'N','A','M','G','A','L','S','I','K',' ','C','A',
&'T','I','M','N','F','E'/
DATA BOXIDE/'S','I','T','I','A','L','V',' ',' ','C','R','F','C',
&'F','E','M','N','M','G','C','A','Z','N','N','I','N','B'/

```

Appendix II. Listing of Program RDARL4. (continued)

```

DATA BBASIC/'S','I','T','I','A','L','F','C','C','R','F','E',
&'Z','N','M','N','N','I','M','G','C','A','S','R','B','A',
&'N','A','K',' ','P',' ','F',' ','C','L'/
C
C INITIALIZE VARIABLES AND GET SYSTEM DATE AND TIME
C
TYPE*, ' '
TYPE*, 'PROGRAM RDARL4 PERMITS RT-11 TO READ ARLEB DATA FILES.'
TYPE*, 'IT WAS WRITTEN BY J S HUEBNER FOR THE RESTON PROBE FACILITY.'
TYPE*, 'THIS VERSION WAS COMPILED 06/15/83.'
CALL IDATE(IDATE1, IDATE2, IDATE3)
IF(IDATE2, EQ, 0) STOP 'SET SYSTEM DATE AND TIME, THEN BEGIN AGAIN'
CALL GTIM(ITIME)
CALL CVTTIM(ITIME, IHOURS, IMIN, ISEC, ITICK)
199 BFILNM(15)=0
IDRIVE=1
200 IERR=0
IREC=0
TYPE*, ' '
TYPE*, 'PLACE DISK WITH ARLEB DATA FILE IN DRIVE DY1:'
TYPE*, ' '
TYPE*, 'ENTER THE ARLEB FILENAME'
TYPE*, '(ARLEB accepts no more than 5 characters for the filename.)'
ACCEPT 210, NAME
210 FORMAT(3A2)
C
C FIND AND OPEN THE ARLEB OUTPUT DATA FILE
C
TYPE*, ' '
TYPE*, 'The disk access and data transfer will begin.'
ICHAN=IGETC()
TYPE*, 'Channel', ICHAN, ' allocated by the computer.'
CALL FIND (NAME, IDRIVE, IADD, ISIZE, ITYPE, IERR)
IF(IERR, NE, -1) CALL ERRORS(IERR, IREC)
IF(ITYPE-2) 310, 320, 330
310 N=1
GOTO 340
320 N=7
GOTO 340
330 N=13
340 TYPE 350, NAME, (BFKIND(I), I=N, N+5), IADD, ISIZE
350 FORMAT(1H0, 3A2, ' is a ', 6A1, ' file whose first record is ', /, 1H ,
&'block ', I4, ' and whose length is ', I3, ' blocks of 128 words.')
CALL OPEN (NAME, IDRIVE, IERR)
IF(IERR, NE, -1) CALL ERRORS(IERR, IREC)
C
C READ THE ARL OUTPUT FILE, BLOCK BY BLOCK, REFORMAT EACH RECORD (BLOCK)
C FOR STORAGE IN ARRAYS, AND DISPLAY HEADER INFORMATION ON VIDEO SCREEN
C
IREC=1 !read $QSET and $PSET information
CALL FDGET(IBUF, IADD, IDRIVE, IERR) !first without format conversion
IF(IERR, NE, -1) CALL ERRORS(IERR, IREC)
DO 125 I=1, 256
IF(BBUF(I), EQ, 0) BBUF(I)='040
125 CONTINUE
CALL FFETCH (IREC, IDRIVE, DATA, IERR) !with 3 to 2 word format convers.

```



```

IF(IERR.NE.-1) CALL ERRORS(IERR,IREC) !assign it to variables, arrays:
DO 130 I=1,12
IOPNAM(I)=IBUF(I)
130 CONTINUE
IARDT2=DATA(5)
IARDT1=DATA(6)
IARDT3=DATA(7)
DO 134 I=1,3
IQFILE(I)=IBUF(I+21)
IARLOT(I)=IBUF(I+24)
134 CONTINUE
NREDCT=DATA(10)
NELMTS=DATA(11)
NTYPE=DATA(12)
NCOUNT=DATA(13)
NOUSED=DATA(14)
DO 150 J=1,14
DO 145 I=1,2
BELNAM(I,J)=BBUF(6*J+78+I)
BXRLIN(I,J)=BBUF(6*J+81+I)
145 CONTINUE
DO 150 I=1,4
BSTAND(I,J)=BBUF(6*J+163+I)
150 CONTINUE
TYPE 152,(IOPNAM(M),M=1,12)
152 FORMAT(1H0,'THE MICROPROBE OPERATOR WAS ',12A2)
TYPE 153,IARDT2,IARDT1,IARDT3
153 FORMAT(1H,'DATE ANALYSES MADE: ',I2,'-',I2,'-19',I2)
TYPE*, ' '
155 FORMAT(9H ELEMENT ,I3,' is ',2A1,1X,2A1,' with standard ',4A1)
DO 160 K=1,NELMTS
TYPE 155,K,(BELNAM(I,K),I=1,2),(BXRLIN(I,K),I=1,2),(BSTAND(I,K),I=1,4)
160 CONTINUE
TYPE*, ' '
TYPE*, 'I AM STILL READING THE ARLEB DATA DISC - PLEASE BE PATIENT!'
TYPE*, ' '
DO 175 IREC=2,ISIZE !read $ANLZ output data
KOUNTR=IREC-1 !count # of analysis
CALL FDGET(IBUF,IADD+IREC-1,IDRIVE,IERR) !no format conversion
IF(IERR.NE.-1) CALL ERRORS(IERR,IREC)
DO 165 N=1,256
IF(BBUF(N).EQ.0) BBUF(N)='040
165 CONTINUE
CALL FFETCH (IREC,IDRIVE,DATA,IERR) !3 to 2 fmt conversion
IF(IERR.NE.-1) CALL ERRORS(IERR,IREC)
IF(DATA(26).EQ.0) KOUNTR=KOUNTR-1 !last+1 has wt%=0.00
IF(DATA(26).EQ.0.) TYPE*, 'LAST ANALYSIS IS #',KOUNTR
IF(DATA(26).EQ.0.) GOTO 176
NANAL(KOUNTR)=KOUNTR
DO 172 L=1,24
ITITLE(L,KOUNTR)=IBUF(L)
172 CONTINUE
X(KOUNTR)=DATA(9)
Y(KOUNTR)=DATA(10)
IAVFLG(KOUNTR)=DATA(11)
DO 174 M=1,15

```

```

      WTPCT(M,KOUNTR)=DATA(M+11)
      CATION(M,KOUNTR)=DATA(M+26)
174  CONTINUE
      PEROXY(KOUNTR)=DATA(42)
175  CONTINUE
176  CONTINUE
C
C REPORT END OF READ AND STORE-IN-ARRAY OPERATION
C
      LBLOCK=IREC+IADD-1
      TYPE 380
380  FORMAT(1H0,'          I JUST FINISHED READING THE ARLEB DISC.')
```

WRITE(7,381) LBLOCK

```

381  FORMAT(1H0,'          The last block read has number ',I4)
      TYPE*, ' '
382  TYPE*, '***** REMOVE ARLEB DISC FROM DRIVE DY1: *****'
      TYPE*, ' '
      TYPE*, '          THEN TYPE (C) TO CONTINUE'
      ACCEPT 701, IFLAG
      IF(IFLAG.NE.'C') GOTO 382
C
C VIEW DATA ON CRT
C
401  TYPE*, 'EXAMINE ANALYSES WITH VIDEO TERMINAL (Y or N)?'
      ACCEPT 701, IRPLY0
      IF(IRPLY0.NE.'Y'.AND.IRPLY0.NE.'N') GOTO 401
      IF(IRPLY0.NE.'Y') GOTO 703
      DO 410, I=1,KOUNTR
      NUMAVG=IAVFLG(I)
      IF(IAVFLG(I).GE.80) NUMAVG=NUMAVG-80
      TYPE 732, NANAL(I), (ITITLE(J,I), J=1,24)
      IF((IAVFLG(I).NE.0).AND.(IAVFLG(I).NE.80)) TYPE 445, NUMAVG
      TYPE 733, ((BELNAM(J,K), J=1,2), K=1,14)
      TYPE 734, (WTPCT(J,I), J=1,15)
      TYPE 735, (CATION(J,I), J=1,15), PEROXY(I)
445  FORMAT(1H , 'AVERAGE OF ', I2)
410  CONTINUE
C
C PRINT OUT DATA STORED IN ARRAYS
C
701  FORMAT(A1)
702  FORMAT(I3, I3)
      TYPE*, ' '
703  TYPE*, 'PRINT OUT HARD COPY OF DATA? (respond Y or N)'
      ACCEPT 701, IRPLY1
      IF(IRPLY1.NE.'Y'.AND.IRPLY1.NE.'N') GOTO 703
      IF(IRPLY1.NE.'Y') GOTO 800
      PRINT 710
710  FORMAT(1H , 'RDARL4 ENABLES RT-11 TO ACCESS CHEMICAL ANALYSES
& STORED IN ARLEB FORMAT BY THE RESTON ARL/SEM0 MICROPROBE')
      PRINT 711
711  FORMAT(1H0, 'IT WAS WRITTEN BY J S HUEBNER FOR THE RESTON
& MICROPROBE FACILITY AND COMPLETED JUNE 15, 1983')
      PRINT 712, IDATE1, IDATE2, IDATE3, IHOURS, IMIN
712  FORMAT(30H RT-11 RUN DATE AND TIME WERE , I2, '-', I2, '-19', I2,
& ' AT ', I2, ':', I2)
```

Appendix II. Listing of Program RDARL4. (continued)

```

PRINT 721,(IOPNAM(I),I=1,12),IARDT2,IARDT1,IARDT3
721  FORMAT(22HOMICROPROBE OPERATOR: ,12A2,25X,'MICROPROBE ANALYSIS
& DATE ',I2,'-',I2,'-19',I2)
PRINT 722,(IQFILE(I),I=1,3),(IARLOT(I),I=1,3) !intesar arrays
722  FORMAT(60H0THE TITLES OF THE *QSET INPUT AND $ANALZ OUTPUT
&FILES WERE ,3A2,' and ',3A2)
PRINT 723
723  FORMAT(1H0,' ELEMENT X-LINE STANDARD')
DO 724 J=1,NELMTS
PRINT 725,(BELNAM(I,J),I=1,2),(BXRLIN(I,J),I=1,2),(BSTAND(I,J),I=1,4)
724  CONTINUE
725  FORMAT(1H ,5X,2A1,10X,2A1,10X,4A1)
PRINT*, ' '
IF(NTYPE.EQ.0) PRINT*, 'ANALYSES STORED AS ELEMENT WEIGHT PERCENT'
IF(NTYPE.EQ.1) PRINT*, 'ANALYSES STORED AS OXIDE WEIGHT PERCENT'
PRINT*, ' '
IF(NREDCT.EQ.0) PRINT*, 'REDUCTION METHOD NOT SPECIFIED'
IF(NREDCT.EQ.0) GOTO 750
IF(NREDCT.GT.5) PRINT*, 'REDUCTION CODE =',NREDCT
IF(NREDCT.GT.5) GOTO 750
GOTO (726,727,728,729,730) NREDCT
726  PRINT*, 'DATA REDUCED WITH MAGIC-IV ROUTINE OF J W COLBY'
GOTO 750
727  PRINT*, 'DATA REDUCTION BY METHOD OF BENCE & ALBEE'
GOTO 750
728  PRINT*, 'DATA REDUCTION BY ARLs G-FACTOR METHOD'
GOTO 750
729  PRINT*, 'DATA NOT REDUCED BEYOND K-RATIO CALCULATION'
GOTO 750
730  PRINT*, 'DATA REDUCTION BY METHOD OF ZIEBOLD & OGILVIE'
PRINT*, ' '
REWIND 6
750  PRINT*, ' '
CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
IF(INDFLG.EQ.1) PRINT 755
IF(JAVFLG.EQ.1) PRINT 756
IF(KNWFLG.EQ.1) PRINT 757
IF(INDFLG.NE.1) LOWER=1
IF(INDFLG.NE.1) UPPER=KOUNTR
IF(LOWER.EQ.0.AND.UPPER.EQ.0) GOTO 731
DO 731 I=LOWER,UPPER
IF(INDFLG.EQ.1) GOTO 751
IF(JAVFLG.EQ.1.AND.IAVFLG(I).GT.1.AND.IAVFLG(I).LT.80) GOTO 751
IF(KNWFLG.EQ.1.AND.IAVFLG(I).GE.80) GOTO 751
GOTO 731
751  PRINT 732,NANAL(I),(ITITLE(J,I),J=1,24)
IF(X(I).NE.0) PRINT 740,X(I),Y(I)
PRINT 733,((BELNAM(J,K),J=1,2),K=1,14)
PRINT 734,(WTPCT(J,I),J=1,15)
NUMAVG=IAVFLG(I)
IF(IAVFLG(I).GE.80) NUMAVG=NUMAVG-80
IF((IAVFLG(I).NE.0).AND.(IAVFLG(I).NE.80)) PRINT 745,NUMAVG
PRINT 735,(CATION(J,I),J=1,15),PEROXY(I)
731  CONTINUE
REWIND 6
732  FORMAT(1H0,I3,' TITLE: ',24A2)

```

Appendix II. Listing of Program RDARL4. (continued)

```

740  FORMAT(1H+,60X,'X = ',F7.1,5X,'Y = ',F7.1)
733  FORMAT(5H      ,2X,14(2A1,5X),'TOTAL')
734  FORMAT(5H WT Z,14(F6.3,1X),F7.3)
745  FORMAT(1H+,124X,'AV OF ',I2)
735  FORMAT(5H CAT ,14(F6.3,1X),1X,F6.3,1X,'PER ',F7.3,' ANIONS')
755  FORMAT(1H0,'INDIVIDUALLY SELECTED ANALYSES OF ALL TYPES:')
756  FORMAT(1H0,'AVERAGES OF UNKNOWN ANALYSES ONLY:')
757  FORMAT(1H0,'ALL KNOWN-UNKNOWN ANALYSES:')
      CALL CHOICE(ANSWER)
      IF(ANSWER.EQ.'Y') GOTO 750

C
C  SELECT FORMAT FOR STORAGE OF DATA IN DISC FILE
C
800  TYPE*, 'DISC OUTPUT? (reply Y or N)'
      ACCEPT 701,IRPLY4
      IF(IRPLY4.EQ.'N') GOTO 1999
      TYPE*, ' '
      TYPE*, 'PLACE OUTPUT DISC in DRIVE 1'
801  TYPE*, ' '
      BFLAG1='040
      BFLAG2='040
      BFLAG3='040
      TYPE*, 'SELECT OUTPUT FILE OPTION FOR ONE OF THE FOLLOWING:'
      TYPE*, ' 1 MINCLC          5 BASIC VIRTUAL ARRAY'
      TYPE*, ' 2 PXVOID          6 BASIC SEQUENTIAL FILE FOR REGRESSIONS'
      TYPE*, ' 3 RECAP          7 BASIC MINERAL FORMULA RECALCULATIONS'
      TYPE*, ' 4 OXIDES          8 EXTRACT ANALYSES OF KNOWN-UNKNOWN'
      TYPE*, '                    9 BASIC SEQUENTIAL FILE FOR DATA STORAGE'
      TYPE*, '                   10 EXIT FROM PROGRAM'
      TYPE*, ' '
      ACCEPT 802,IRPLY5
802  FORMAT(I3)
      IF(IRPLY5.GT.10) TYPE*, 'REPLY OUT OF BOUNDS. TRY AGAIN.'
      IF(IRPLY5.GT.10) GOTO 801
      GOTO (805,840,860,960,880,900,920,940,1100,1999),IRPLY5

C
C  OPTION 1, MINCLC, COMPLETED JULY 14, 1982
C
805  TYPE*, 'YOU ELECTED A MINCLC FILE. NOW YOU MUST ENTER THE IDEAL'
      TYPE*, 'NUMBER OF CATIONS PER NUMBER OF OXYGENS (f12.7,f12.7):'
      READ(5,806,ERR=803) CAT,OXY
      IF(CAT.LE.0.01) GOTO 803
      IF(OXY.LE.0.01) GOTO 803
      GOTO 804
803  IERR=12
      CALL ERRORS(IERR,IREC)
      GOTO 805
806  FORMAT(2(F12.7))
804  CALL SETUP(IARLOT,FILSPC)
      OPEN(UNIT=2,NAME=FILSPC,ERR=807,CARRIAGECONTROL='LIST')
      GOTO 809
807  IERR=11
      CALL ERRORS(IERR,IREC)
      GOTO 804
809  WRITE (2,808,ERR=807) OXY,CAT
808  FORMAT(12H***RDOXIONS,38X,F10.0,F10.0)

```

Appendix II. Listing of Program RDARL4. (continued)

```

810  CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
      IF(INDFLG.NE.1) LOWER=1
      IF(INDFLG.NE.1) UPPER=KOUNTR
      DO 825 K=LOWER,UPPER
      IF(INDFLG.EQ.1) GOTO 811
      IF(JAVFLG.EQ.1.AND.IAVFLG(K).GT.1.AND.IAVFLG(K).LT.80) GOTO 811
      IF(KNWFLG.EQ.1.AND.IAVFLG(K).GE.80) GOTO 811
      GOTO 825
811  WRITE(2,826) K,(ITITLE(J,K),J=1,24)
      IF(NELMTS.GT.8) GOTO 818
      GOTO(813,814,815,816,817,818,819) NELMTS-2
813  WRITE(2,833)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=1,3)
      GOTO 825
814  WRITE(2,834)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=1,4)
      GOTO 825
815  WRITE(2,835)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=1,5)
      GOTO 825
816  WRITE(2,836)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=1,6)
      GOTO 825
817  WRITE(2,837)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=1,7)
      GOTO 825
818  WRITE(2,838)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=1,8)
      GOTO (825,819,820,821,822,823,824) NELMTS-7
819  WRITE(2,831)(BELNAM(I,9),I=1,2),(WTPCT(9,K))
      GOTO 825
820  WRITE(2,832)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=9,10)
      GOTO 825
821  WRITE(2,833)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=9,11)
      GOTO 825
822  WRITE(2,834)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=9,12)
      GOTO 825
823  WRITE(2,835)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=9,13)
      GOTO 825
824  WRITE(2,836)((BELNAM(I,J),I=1,2),WTPCT(J,K)),J=9,14)
825  CONTINUE
      CALL CHOICE(ANSWER)
      IF(ANSWER.EQ.'Y') GOTO 810
826  FORMAT(I3,' TITLE: '24A2)
831  FORMAT(2A1,F6.2,'$$')
832  FORMAT(2(2A1,F6.2),'$$')
833  FORMAT(3(2A1,F6.2),'$$')
834  FORMAT(4(2A1,F6.2),'$$')
835  FORMAT(5(2A1,F6.2),'$$')
836  FORMAT(6(2A1,F6.2),'$$')
837  FORMAT(7(2A1,F6.2),'$$')
838  FORMAT(8(2A1,F6.2))
      CLOSE(UNIT=2)
      GOTO 990
C
C  OPTION 2, PXVOID, COMPLETED AUGUST 18, 1982
C
840  DO 841, N=1,KOUNTR
      WTPCT(16,N)=0.0 !for elems not analyzed
841  CONTINUE
      TYPE*, 'You elected a Fortran sequential file for PXVOID'
      TYPE*, 'Your data must be in OXIDE WEIGHT PERCENT'

```

Appendix II. Listing of Program RDARL4. (continued)

```

TYPE*, ' '
DO 846, I=1,12                !sort element names
IFLAG=0
DO 846, J=1,NELMTS
IF(BELNAM(1,J).NE.BPXVD(1,I)) GOTO 845
IF(BELNAM(2,J).NE.BPXVD(2,I)) GOTO 845
KCOUNT(I)=J                !stores sequence of WTPCT
IFLAG=1
845 IF(IFLAG.EQ.1) GOTO 846
KCOUNT(I)=16
846 CONTINUE
TYPE 847,(KCOUNT(I),I=1,13)
847 FORMAT(27H SEQUENCE ELEMENTS PUT OUT:,13(2X,I2))
TYPE*, 'A VALUE OF 16 INDICATES ELEMENT NOT ANALYZED - 0.00%'
842 CALL SETUP(IARLOT,FILSPC)
OPEN(UNIT=2,NAME=FILSPC,ERR=843,CARRIAGECONTROL='LIST')
GOTO 844
843 IERR=11
CALL ERRORS(IERR,IREC)
GOTO 842
844 WRITE(2,855,ERR=843) NTYPE
848 CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
IF(INDFLG.NE.1) LOWER=1
IF(INDFLG.NE.1) UPPER=KOUNTR
DO 850 N=LOWER,UPPER
IF(INDFLG.EQ.1) GOTO 849
IF(JAVFLG.EQ.1.AND.IAVFLG(N).GT.1.AND.IAVFLG(N).LT.80) GOTO 849
IF(KNWFLG.EQ.1.AND.IAVFLG(N).GE.80) GOTO 849
GOTO 850
849 WRITE(2,856) N,(ITITLE(L,N),L=1,24)
WRITE(2,858) (BPXVD(J,1),J=1,2),WTPCT(KCOUNT(1),N),
&((BPXVD(J,I),J=1,2),WTPCT(KCOUNT(I),N),I=3,7)
WRITE(2,858) ((BPXVD(J,I),J=1,2),WTPCT(KCOUNT(I),N),I=8,13)
850 CONTINUE
855 FORMAT(I1)
856 FORMAT(I3,1X,24A2,20X)
858 FORMAT(6(2A1,F6.2))
CALL CHOICE(ANSWER)
IF(ANSWER.EQ.'Y') GOTO 848
CLOSE(UNIT=2)
GOTO 990

C
C OPTION 3, RECAMP, COMPLETED AUGUST 19, 1982
C
860 DO 861, N=1,KOUNTR
WTPCT(16,N)=0.0
861 CONTINUE
TYPE*, 'You elected a Fortran sequential file for RECAMP'
TYPE*, ' '
DO 866, I=1,9
IFLAG=0
DO 866, J=1,NELMTS
IF(BELNAM(1,J).NE.BCAMP(1,I)) GOTO 865
IF(BELNAM(2,J).NE.BCAMP(2,I)) GOTO 865
KCOUNT(I)=J
IFLAG=1

```

Appendix II. Listing of Program RDARL4. (continued)

```

865  IF(IFLAG,EQ.1) GOTO 866
      KCOUNT(I)=16
866  CONTINUE
862  CALL SETUP(IARLOT,FILSPC)
      OPEN(UNIT=2,NAME=FILSPC,ERR=863,CARRIAGE CONTROL='LIST')
      GOTO 867
863  IERR=11
      CALL ERRORS(IERR,IREC)
      BFLAG1='N'
      GOTO 862
867  CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
      IF(INDFLG.NE.1) LOWER=1
      IF(INDFLG.NE.1) UPPER=KOUNTR
      DO 870, N=LOWER,UPPER
      IF(INDFLG.EQ.1) GOTO 868
      IF(JAVFLG.EQ.1.AND.IAVFLG(N).GT.1.AND.IAVFLG(N).LT.80) GOTO 868
      IF(KNWFLG.EQ.1.AND.IAVFLG(N).GE.80) GOTO 868
      GOTO 870
868  WRITE(2,871,ERR=863) N,(ITITLE(L,N),L=1,18)
      WRITE(2,872) ((BCAMP(J,I),J=1,2),WTPCT(KCOUNT(I),N),I=1,9)
870  CONTINUE
871  FORMAT(I3,X,18A2)
872  FORMAT(9(2A1,F6.2))
      CALL CHOICE(ANSWER)
      IF(ANSWER.EQ.'Y') GOTO 867
      CLOSE(UNIT=2)
      GOTO 990

C
C  OPTION 5, DIRECT ACCESS FILE, COMPLETED JULY 6, 1982
C
880  TYPE*, 'YOU CHOSE A VIRTUAL (DIRECT ACCESS) DISC ARRAY FILE. ALL DATA'
      TYPE*, 'ARE PUT OUT; SELECTION OF DATA IS NOT YET POSSIBLE (8/31/82).'
      TYPE 881,(IARLOT(I),I=1,3)
881  FORMAT(1H0,'THE ARLEB OUTPUT FILENAME WAS ',3A2)
      TYPE*, 'ENTER BASIC VIRTUAL ARRAY FILESPEC AS DY???????????'
      TYPE*, ' '
      LENGTH=4*(NELMTS+1)+3 !length of each record in 2-byte words
      CALL ASSIGN(2,' ',-1)
      DEFINE FILE 2 (KOUNTR,LENGTH,U,NREC)
      NREC=1
      DO 890 I=1,KOUNTR
      WRITE(2'NREC) KOUNTR,NANAL(I),NELMTS,(WTPCT(J,I),J=1,NELMTS),
&WTPCT(15,I),(CATION(J,I),J=1,NELMTS),CATION(15,I)
890  CONTINUE
      CLOSE(UNIT=2)
      TYPE*,NREC,' RECORDS WERE PUT OUT TO DISC.'
      TYPE*, 'EACH RECORD CONSISTS OF',LENGTH,' 2-BYTE WORDS.'
      GOTO 990

C
C  OPTION 6, BASIC REGRESSION FILE, COMPLETED AUGUST 3, 1982
C
900  TYPE*, 'YOU SELECTED A MINC BASIC SEQUENTIAL FILE FORMATTED FOR THE'
      TYPE*, 'MINC CORRELATION MATRIX AND REGRESSIONS. THESE PROGRAMS ARE'
      TYPE*, 'LIMITED TO 75 ANALYSES. YOU WILL NEED TO HAVE SELECTED YOUR'
      TYPE*, 'ANALYSES IN ADVANCE AND KNOW THE TOTAL NUMBER TO BE WRITTEN'
      TYPE*, ' '

```

Appendix II. Listing of Program RDARL4. (continued)

```

TYPE*, 'IF YOU ARE NOT PREPARED TO PROCEED, TYPE ^C (control C)'
TYPE*, ' '
TYPE*, 'DO YOU WANT TO STORE WT % (1) OR CATIONS (2)?'
901  FORMAT(I3)
ACCEPT 901,IRPLY7
912  CALL SETUP(IARLOT,FILSPC)
OPEN(UNIT=2,NAME=FILSPC,ERR=913,CARRIAGECONTROL='FORTRAN')
GOTO 914
913  IERR=11
905  CALL ERRORS(IERR,IERR)
GOTO 912
914  NUMBER=0
TYPE*, 'HOW MANY ANALYSES SELECTED?'
ACCEPT 901,IRPLY8
WRITE(2,902,ERR=913) IRPLY8,NELMTS
902  FORMAT('+',I3,',',I2)
904  TYPE*, 'NOW ENTER THE FIRST AND LAST ANALYSIS NUMBERS IN A GROUP'
TYPE*, 'OF CONSECUTIVE ANALYSES (I2,I2)'
906  FORMAT(I3,I3)
ACCEPT 906,K,L
IF(IRPLY7.EQ.2) GOTO 910      !branch to write out cations
DO 907 I=K,L
903  FORMAT(1X,I3)
WRITE(2,903,ERR=919) I
NUMBER=NUMBER+1
DO 907 J=1,NELMTS
WRITE(2,909,ERR=919) WTPCT(J,I)
907  CONTINUE
IF(NUMBER.GE.IRPLY8) GOTO 916
GOTO 904
909  FORMAT(1X,F6.2)
GOTO 916
910  CONTINUE
911  FORMAT(1X,F6.3)
DO 915 I=K,L
WRITE(2,903,ERR=919) I
NUMBER=NUMBER+1
DO 915 J=1,NELMTS
WRITE(2,911,ERR=919) CATION(J,I)
915  CONTINUE
IF(NUMBER.GE.IRPLY8) GOTO 916
GOTO 904
916  CLOSE(UNIT=2)
TYPE*, ' '
TYPE*, 'USING THIS OPTION,'
TYPE*, 'YOU HAVE CREATED A FILE OF ',NUMBER,' ANALYSES.'
TYPE*, 'DO YOU WANT ANOTHER DISC FILE USING THIS OPTION (Y or N)?'
ACCEPT 701,IRPLY9
IF(IRPLY9.EQ.'Y') GOTO 900
GOTO 990
919  IERR=13
GOTO 905
C
C   OPTION 7, BASIC FORMULA RECALCULATIONS, COMPLETED JUNE 08, 1983
C
920  TYPE*, 'YOU CHOSE TO CREATE A BASIC SEQUENTIAL FILE FOR MINERAL

```


Appendix II. Listing of Program RDARL4. (continued)

```

& FORMULA CALCULTIONS.'
DO 921,N=1,KOUNTR
WTPCT(16,N)=0.0
CATION(16,N)=0.0
921 CONTINUE
DO 925,I=1,18
IFLAG=0
DO 925,J=1,NELMTS
IF(BELNAM(1,J).NE.BBASIC(1,I)) GOTO 923
IF(BELNAM(2,J).NE.BBASIC(2,I)) GOTO 923
KCOUNT(I)=J
IFLAG=1
923 IF(IFLAG.EQ.1) GOTO 925
KCOUNT(I)=16
925 CONTINUE
935 CALL SETUP(IARLOT,FILSPC)
OPEN(UNIT=2,NAME=FILSPC,ERR=936,CARRIAGECONTROL='FORTRAN')
GOTO 937
936 IERR=11
CALL ERRORS(IERR,IERR)
GOTO 935
937 WRITE(2,926,ERR=936) NTYPE
DO 922,L=1,18
WRITE(2,924) (BBASIC(K,L),K=1,2)
922 CONTINUE
924 FORMAT(1X,2A1)
926 FORMAT(1H+,I1)
927 CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
IF(INDFLG.NE.1) LOWER=1
IF(INDFLG.NE.1) UPPER=KOUNTR
DO 929, N=LOWER,UPPER
IF(INDFLG.EQ.1) GOTO 928
IF(JAVFLG.EQ.1.AND.IAVFLG(N).GT.1.AND.IAVFLG(N).LT.80) GOTO 928
IF(KNWFLG.EQ.1.AND.IAVFLG(N).GE.80) GOTO 928
GOTO 929
928 WRITE(2,930) N,(ITITLE(L,N),L=1,24)
WRITE(2,931) (IQFILE(I),I=1,3),PEROXY(N)
DO 929, I=1,18
WRITE(2,932) WTPCT(KCOUNT(I),N),CATION(KCOUNT(I),N)
929 CONTINUE
930 FORMAT(1X,I3,',',',',24A2)
931 FORMAT(1X,3A2,',',',',F6.3)
932 FORMAT(1X,F6.2,',',',',F6.3)
CALL CHOICE(ANSWER)
IF(ANSWER.EQ.'Y') GOTO 927
CLOSE(UNIT=2)
GOTO 990
C
C OPTION 8, EXTRACT AND FILE KNOWN-UNKNOWN ANALYSES, COMPLETED 10/7/82
C
940 TYPE*, ' '
TYPE*, 'MINC WILL FILE ON DISC ONLY THE ANALYSES OF KNOWN-UNKNOWN.'
L=0
M=0
NUMBER=0
DO 941, I=1,KOUNTR

```

Appendix II. Listing of Program RDARL4. (continued)

```

      IF(IAVFLG(I).EQ.80) L=L+1
      IF(IAVFLG(I).GT.80) M=M+1
941  CONTINUE
      IF(M+L.EQ.0) GOTO 959
      TYPE 942,L,M
942  FORMAT(1H0,'MINC FOUND ',I3,' INDIVIDUAL AND ',I3,' AVERAGE
& ANALYSES OF KNOWN-UNKNOWN')
      TYPE*, ' '
      TYPE*, 'DO YOU WANT TO FILE (1)INDIVIDUAL, (2)AVERAGE, (3)BOTH,'
      TYPE*, 'OR (4)NEITHER KIND OF ANALYSES ? (respond 1, 2, 3, or 4)'
      ACCEPT 958, IRPLY0
958  FORMAT(I1)
      IF(IRPLY0.EQ.4) GOTO 990
943  CALL SETUP(IARLOT,FILSPC)
      OPEN(UNIT=2,NAME=FILSPC,TYPE='NEW',ERR=948,CARRIAGECONTROL='FORTRAN')
      GOTO 949
948  IERR=11
      CALL ERRORS(IERR,IRES)
      GOTO 943
949  DO 950,K=1,KOUNTR
      IF(IRPLY0.EQ.1.AND.IAVFLG(K).EQ.80) GOTO 945
      IF(IRPLY0.EQ.1.AND.IAVFLG(K).EQ.81) GOTO 945
      IF(IRPLY0.EQ.2.AND.IAVFLG(K).GT.81) GOTO 945
      IF(IRPLY0.EQ.3.AND.IAVFLG(K).GE.80) GOTO 945
      GOTO 950
945  NUMBER=NUMBER+1
      IF(NUMBER.GE.1) GOTO 946
      WRITE(2,944,ERR=948)K,(IOPNAM(I),I=1,12),IARDT2,IARDT1,IARDT3
      GOTO 947
946  WRITE(2,953)K,(IOPNAM(I),I=1,12),IARDT2,IARDT1,IARDT3
947  WRITE(2,954)(IQFILE(I),I=1,3),(IARLOT(I),I=1,3),NELMTS,NREDCT
      NUMAVG=IAVFLG(K)-80
      WRITE(2,955)(ITITLE(I,K),I=1,24),NUMAVG
      WRITE(2,956)((BELNAM(I,J),I=1,2),(WTPCT(J,K)),J=1,NELMTS)
      WRITE(2,952)WTPCT(15,K)
      WRITE(2,957)((BSTAND(I,J),I=1,4),J=1,NELMTS)
950  CONTINUE
      CLOSE(UNIT=2)
      TYPE*, ' '
      TYPE*,NUMBER,' ANALYSES WERE JUST FILED'
944  FORMAT(2H+#,I4,' OPERATOR: ',12A2,2X,'DATE: ',I2,'-',I2,'-19',I2)
953  FORMAT(2H +#,I4,' OPERATOR: ',12A2,2X,'DATE: ',I2,'-',I2,'-19',I2)
954  FORMAT(8H *QSET=,3A2,' DATA STORED IN ',3A2,2X,'NELMTS=',I2,3X,
&'NREDCT=',I1)
955  FORMAT(2H ,24A2,2X,'AVERAGE OF ',I2)
956  FORMAT(2H ,14(2A1,F6.2))
952  FORMAT(2H , 'SUM=',F6.2)
957  FORMAT(2H ,4A1,13(1X,4A1))
      GOTO 990
959  TYPE*, 'NONE FOUND'
      GOTO 990
C
C  OPTION 4, OXIDES, COMPLETED AUGUST 31, 1982
C
960  TYPE*, 'YOU CHOSE THE FORMAT FOR OXIDE FORMULAS AND GEOTHERMOMERTY'
      DO 961, N=1,KOUNTR

```

Appendix II. Listing of Program RDARL4. (continued)

```

WTPCT(16,N)=0.0
961 CONTINUE
DO 965, I=1,13
IFLAG=0
DO 965, J=1,NELMTS
IF(BELNAM(1,J).NE.BOXIDE(1,I)) GOTO 963
IF(BELNAM(2,J).NE.BOXIDE(2,I)) GOTO 963
KCOUNT(I)=J !sequence elements put out
IFLAG=1
963 IF(IFLAG.EQ.1) GOTO 965
KCOUNT(I)=16 !if elem not present, wtX=0.0
965 CONTINUE
975 CALL SETUP(IARLOT,FILSPC)
OPEN(UNIT=2,NAME=FILSPC,ERR=976,CARRIAGECONTROL='FORTRAN')
GOTO 977
976 IERR=11
CALL ERRORS(IERR,IERR)
GOTO 975
977 WRITE(2,966,ERR=976) NTYPE
966 FORMAT('+',I1)
967 CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
IF(INDFLG.NE.1) LOWER=1
IF(INDFLG.NE.1) UPPER=KOUNTR
DO 970, N=LOWER,UPPER
IF(INDFLG.EQ.1) GOTO 968
IF(JAVFLG.EQ.1.AND.IAVFLG(N).GT.1.AND.IAVFLG(N).LT.80) GOTO 968
IF(KNWFLG.EQ.1.AND.IAVFLG(N).GE.80) GOTO 968
GOTO 970
968 WRITE(2,971) N,IAVFLG(N),(ITITLE(L,N),L=1,24)
DO 970, I=1,13
WRITE(2,972)(BOXIDE(J,I),J=1,2),WTPCT(KCOUNT(I),N)
970 CONTINUE
971 FORMAT(I3,',',I2,',',24A2)
972 FORMAT(1X,2A1,',',F6.2)
CALL CHOICE(ANSWER)
IF(ANSWER.EQ.'Y') GOTO 967
CLOSE(UNIT=2)
GOTO 990

C
C OPTION 9. BASIC SEQUENTIAL FILE THAT CAN BE CONCATENATED,
C SEARCHED, OR EDITED. USE FOR ARCHIVING GROUPS OF ANALYSES.
C
1100 TYPE*, 'YOU ELECTED A DISK FILE OF THE ELEMENTS ANALYZED AND IN THE'
TYPE*, ' SAME ORDER IN WHICH THEY APPEARED IN THE ARL PRINTOUT. THIS'
TYPE*, ' FORMAT CAN BE SEARCHED, CONCATENATED, AND EDITED.'
1101 CALL SETUP(IARLOT,FILSPC)
OPEN(UNIT=2,NAME=FILSPC,ERR=1102,CARRIAGECONTROL='FORTRAN')
GOTO 1103
1102 IERR=11
CALL ERRORS(IERR,IERR)
GOTO 1101
1103 WRITE(2,1105,ERR=1102) NELMTS,PEROXY(1)
1105 FORMAT(1H+,I3,',',***,',',F6.3)
DO 1110, J=1,NELMTS
WRITE(2,1115) (BELNAM(I,J),I=1,2)
1110 CONTINUE

```

Appendix II. Listing of Program RDARL4. (continued)

```

1115  FORMAT(1X,2A1)
      WRITE(2,1116)
1116  FORMAT(6H TOTAL)
1120  CALL SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
      IF(INDFLG.NE.1) LOWER=1
      IF(INDFLG.NE.1) UPPER=KOUNTR
      DO 1130, N=LOWER,UPPER
      IF(INDFLG.EQ.1) GOTO 1128
      IF(JAVFLG.EQ.1.AND.IAVFLG(N).GT.1.AND.IAVFLG(N).LT.80) GOTO 1128
      IF(KNWFLG.EQ.1.AND.IAVFLG(N).GT.80) GOTO 1128
      GOTO 1130
1128  WRITE(2,1131) NANAL(N),(ITITLE(L,N),L=1,24),PEROXY(N),IAVFLG(N)
      DO 1129,I=1,NELMTS
      WRITE(2,1132) WTPCT(I,N),CATION(I,N)
1129  CONTINUE
      WRITE(2,1132) WTPCT(15,N),CATION(15,N)
1130  CONTINUE
1131  FORMAT(1X,I3,' ',',',24A2,' ',',',F6.3,' ',',',I3)
1132  FORMAT(1X,F6.2,' ',',',F6.3)
      CALL CHOICE(ANSWER)
      IF(ANSWER.EQ.'Y') GOTO 1120
      CLOSE(UNIT=2)
      GOTO 990

C
C  BRANCH TO ANOTHER OPTION OR EXIT
C
990  TYPE*, ' '
      TYPE*, 'ANOTHER OUTPUT FILE OPTION? (Y or N)'
      ACCEPT 701, IRPLY6
      IF(IRPLY6.EQ.'Y') GOTO 801
1999  TYPE*, 'READ A NEW ARLEB DATA FILE? (Y or N)'
      ACCEPT 701, IRPLY6
      IF(IRPLY6.EQ.'Y') GOTO 199
      STOP 'GOODBYE. HAVE A NICE DAY.'
      END

C
C--SUBROUTINE ERRORS IDENTIFIES THE DISK ERROR--
C
      SUBROUTINE ERRORS(IERR,IREC)
      IF(IERR.GE.11) GOTO 5          !ARLEB errors <11; RT-11 errors >10
      TYPE*, 'ERROR AT RECORD ',IREC,' ANALYSIS ',IREC-1          !ARLEB only
5     GOTO (10,20,30,40,50,60,70,80,90,100,110,120,130),IERR+1
10    STOP 'ILLEGAL DISC ADDRESS. BEGIN AGAIN.'
20    STOP 'BAD READ/WRITE. DEFECTIVE DISK? BEGIN AGAIN.'
30    STOP 'PARITY ERROR; DEFECTIVE DATA TRANSMISSION. BEGIN AGAIN.'
40    STOP 'LIBRARY FULL (NO ROOM FOR FILE HEADER).'
50    STOP 'DISC FULL (NOT ENOUGH FREE SPACE).'
60    STOP 'NON-EXISTANT DATA RECORD REQUESTED. CHOOSE NEW NAME, BEGIN AGAIN.'
70    STOP 'DRIVE NOT READY. CHECK AND CORRECT, THEN BEGIN AGAIN.'
80    STOP 'NON-EXISTANT DRIVE REQUESTED.'
90    TYPE*, 'FILE NOT FOUND or WRONG FILE TYPE or INCORRECT DISK.'
      STOP 'CHANGE TO ARLEB DISK or CHOOSE NEW NAME. PLEASE BEGIN AGAIN.'
100   STOP 'BAD FETCH OR FILE. DENSITY INCORRECT? PLEASE BEGIN AGAIN.'
110   TYPE*, 'DISC IS NOT A VALID ARL BASIC SYSTEM (BOOTABLE) DISC.'
      STOP 'PLACE ANOTHER DISK IN THE RIGHT-HAND DRIVE AND BEGIN AGAIN.'
120   TYPE*, 'ERROR WHILE TRYING TO OPEN AN OUTPUT FILE. MAKE SURE THAT A'

```

Appendix II. Listing of Program RDARL4. (continued)

```

TYPE*, 'PROPERLY FORMATTED DISK WITH FREE SPACE IS PRESENT AND THAT THE'
TYPE*, 'DISK DRIVE IS READY. CAREFULLY REENTER THE NAME OF THE OUTPUT'
TYPE*, 'FILE, FOLLOWING THE SUGGESTED SYNTAX. IF YOU DO NOT SUCCEED,'
TYPE*, 'TERMINATE THE PROGRAM WITH A CONTROL C.'
TYPE*, ' '
130 TYPE*, 'UNACCEPTABLE VALUES. CHECK SYNTAX, PLEASE TRY AGAIN.'
TYPE*, ' '
140 TYPE*, 'ERROR WHILE WRITING. REOPEN THE OUTPUT FILE.'
TYPE*, ' '
RETURN
END

C
C--SUBROUTINE SETUP ACCEPTS FILENAME
C
SUBROUTINE SETUP(IARLOT,FILSPC)
INTEGER IARLOT(3)
BYTE FILSPC(15)
TYPE 200,(IARLOT(I),I=1,3)
200 FORMAT(29H THE ARL OUTPUT FILENAME WAS ,3A2,/)
210 TYPE 220
220 FORMAT(52H ENTER NAME OF FILE TO BE CREATED AS DY?:?????.???)
ACCEPT 230,FILSPC
230 FORMAT(15A1)
IF(FILSPC(1).EQ.' ') GOTO 210
FILSPC(15)=0
TYPE*, ' '
RETURN
END

C
C--SUBROUTINE SELECT LETS YOU SELECT ANALYSES TO BE PUT OUT
C
SUBROUTINE SELECT(LOWER,UPPER,INDFLG,JAVFLG,KNWFLG)
INTEGER UPPER
COMMON BFLAG1,BFLAG2,BFLAG3
INDFLG=0
JAVFLG=0
KNWFLG=0
LOWER=0
UPPER=0
IF(BFLAG1.EQ.'Y') GOTO 323
TYPE 300
300 FORMAT(1H0,'CHOOSE ANALYSES TO BE PUT OUT.')
TYPE 301
301 FORMAT(1H0,'YOU WILL BE ASKED TO CHOSE BETWEEN AVERAGED UNKNOWN,')
TYPE 302
302 FORMAT(1H , 'ALL KNOWN-UNKNOWN, AND THEN INDIVIDUAL ANALYSES.')
305 TYPE 310
310 FORMAT(1H0,'ALL AVERAGES OF UNKNOWN ANALYSES (Y or N)?')
ACCEPT 320, BFLAG1
320 FORMAT(A1)
IF(BFLAG1.NE.'Y'.AND.BFLAG1.NE.'N') GOTO 305
IF(BFLAG1.EQ.'Y') JAVFLG=1
IF(BFLAG1.EQ.'Y') GOTO 360
323 IF(BFLAG2.EQ.'Y') GOTO 330
324 TYPE 325
325 FORMAT(1H0,'ALL THE KNOWN-UNKNOWN (Y or N)?')

```

Appendix II. Listing of Program RDARL4. (continued)

```

ACCEPT 320, BFLAG2
IF(BFLAG2.NE.'Y'.AND.BFLAG2.NE.'N') GOTO 324
IF(BFLAG2.EQ.'Y') KNWFLG=1
IF(BFLAG2.EQ.'Y') GOTO 360
330 TYPE 340
340 FORMAT(1H0,'Enter FIRST and LAST analysis numbers in a group ',/,
&1H ', 'of CONSECUTIVE analyses (I3,I3). Just type <RET> to exit.')
```

```

ACCEPT 350,LOWER,UPPER
350 FORMAT(I4,I3)
IF(LOWER.EQ.0.OR.UPPER.EQ.0) GOTO 355
IF(LOWER.GT.100.OR.LOWER.LT.0) TYPE*, 'First number too large or
& small. Try again.'
IF(UPPER.GT.100.OR.UPPER.LT.0) TYPE*, 'Second number too large or
& small. Try again.'
IF((UPPER-LOWER).LT.0) TYPE*, 'Second number too small. Try again.'
IF((UPPER-LOWER).LT.0) GOTO 330
IF(LOWER.GT.100.OR.LOWER.LT.0) GOTO 330
IF(UPPER.GT.100.OR.UPPER.LT.0) GOTO 330
INDFLG=1
GOTO 360
355 TYPE 370
370 FORMAT(1H0,'YOU SELECTED NOTHING OR TYPED A ZERO. TRY AGAIN (Y or N)?')
ACCEPT 320, BFLAG3
IF(BFLAG3.EQ.'Y') GOTO 330
360 RETURN
END

C
C--SUBROUTINE CHOICE ENABLES A SECOND SELECTION
C
SUBROUTINE CHOICE(ANSWER)
409 TYPE*, 'SELECT MORE ANALYSES FOR THIS OUTPUT FILE (Y or N)?'
ACCEPT 410,ANSWER
410 FORMAT(A1)
IF(ANSWER.NE.'Y'.AND.ANSWER.NE.'N') GOTO 409
RETURN
END
```

Appendix III. Examples of Output created by Program RDARL4.

a. Hardcopy printout.

RDARL4 ENABLES RT-11 TO ACCESS CHEMICAL ANALYSES STORED IN ASCII FORMAT BY THE RESTON ARI/SEMIO MICROPROBE
 IT WAS WRITTEN BY J S HUERNER FOR THE RESTON MICROPROBE FACILITY AND COMPLETED JUNE 15, 1983
 RT-11 RUN DATE AND TIME WERE 6-28-1983 AT 15:28
 MICROPROBE OPERATOR: HUERNER - HD & GNT MICROPROBE ANALYSIS DATE 10-29-1982
 THE TITLES OF THE #QSET INPUT AND #ANALZ OUTPUT FILES WERE #GAPY and D1027

ELEMENT	X-LINE	STANDARD
FE	KA	FXHH
SI	KA	FXAD
HG	KA	FXAD
CA	KA	FXAD
AL	KA	FSTA
NA	KA	FSTA
TI	KA	OXIL
CR	KA	OXTR
MN	KA	DLST

ANALYSES STORED AS OXIDE WEIGHT PERCENT
 DATA REDUCTION BY METHOD OF HENCE & ALBEE

AVERAGES OF UNKNOWN ANALYSES ONLY:

19 TITLE: HD BAILLAL LAKE USNM R3005															TOTAL	AV OF
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN							
WT %	3.405	54.037	16.697	25.389	0.864	0.261	0.005	0.018	0.334	0.674	0.000	0.000	0.000	0.000	101.009	4
CAT	0.103	1.964	0.904	0.989	0.037	0.018	0.000	0.001	0.010	0.059	0.000	0.000	0.000	0.000	4.027	PER 6.000 ANIONS

16 TITLE: HD DUCKTOWN USNM 119197															TOTAL	AV OF
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN							
WT %	5.693	53.549	15.044	24.184	0.399	0.495	0.028	0.014	1.176	0.674	0.000	0.000	0.000	0.000	100.582	5
CAT	0.176	1.978	0.828	0.957	0.017	0.035	0.001	0.000	0.037	0.037	0.000	0.000	0.000	0.000	4.030	PER 5.000 ANIONS

24 TITLE: H UVAROVITE USNM 139715 ZONER															TOTAL	AV OF
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN							
WT %	2.189	38.217	0.637	33.043	7.937	0.143	0.043	18.580	0.836	0.674	0.000	0.000	0.000	0.000	101.625	9
CAT	0.072	1.513	0.038	1.401	0.370	0.011	0.001	0.581	0.028	0.021	0.000	0.000	0.000	0.000	4.016	PER 6.000 ANIONS

29 TITLE: KAK GNT															TOTAL	AV OF
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN							
WT %	210.617	42.296	19.185	5.111	23.571	0.159	0.409	0.081	0.487	0.674	0.000	0.000	0.000	0.000	101.915	4
CAT	0.312	1.489	1.006	0.193	0.978	0.011	0.011	0.002	0.015	0.079	0.000	0.000	0.000	0.000	4.016	PER 6.000 ANIONS

INDIVIDUALLY SELECTED ANALYSES OF ALL TYPES:

1 TITLE: HD TEST															X = 8893.0	Y = 37153.5	TOTAL
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN								
WT %	223.134	48.486	4.257	20.762	0.083	0.474	0.027	0.002	2.542	0.674	0.000	0.000	0.000	0.000	99.767		
CAT	0.786	1.969	0.258	0.904	0.004	0.037	0.001	0.000	0.087	0.353	0.000	0.000	0.000	0.000	4.046	PER 6.000 ANIONS	

2 TITLE: HD TEST															X = 8893.0	Y = 37153.5	TOTAL
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN								
WT %	222.820	49.072	4.206	20.753	0.499	0.469	0.005	0.005	2.530	0.674	0.000	0.000	0.000	0.000	100.359		
CAT	0.767	1.973	0.252	0.894	0.024	0.037	0.000	0.000	0.086	0.224	0.000	0.000	0.000	0.000	4.033	PER 6.000 ANIONS	

3 TITLE: HD TEST															X = 8893.0	Y = 37153.5	TOTAL
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN								
WT %	222.979	49.986	4.111	20.852	0.112	0.454	0.000	0.033	2.627	0.674	0.000	0.000	0.000	0.000	101.154		
CAT	0.766	1.992	0.244	0.890	0.005	0.035	0.000	0.001	0.089	0.164	0.000	0.000	0.000	0.000	4.022	PER 6.000 ANIONS	

4 TITLE: HD TEST															X = 8893.0	Y = 37153.5	TOTAL
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN								
WT %	223.137	49.733	3.881	20.754	0.141	0.452	0.008	0.000	2.622	0.674	0.000	0.000	0.000	0.000	100.727		
CAT	0.775	1.993	0.232	0.891	0.007	0.035	0.000	0.000	0.089	0.128	0.000	0.000	0.000	0.000	4.021	PER 6.000 ANIONS	

5 TITLE: HD M12330 AV OF 4															TOTAL	AV OF
WT %	FE	SI	HG	CA	AL	NA	TI	CR	MN							
WT %	223.017	49.320	4.114	20.780	0.209	0.462	0.009	0.008	2.580	0.674	0.000	0.000	0.000	0.000	100.500	4
CAT	0.774	1.982	0.246	0.895	0.010	0.036	0.000	0.000	0.088	0.105	0.000	0.000	0.000	0.000	4.031	PER 6.000 ANIONS

Appendix III. Examples of Output created by Program RDARL4.

b. Printout of sequential disk file for auxillary program MINCLC

```

****RDXIONS
10 TITLE: HD BAILAL LAKE USNM R3005
FE 3.41SI 54.04MG 16.70CA 25.39AL 0.86NA 0.26TI 0.00CR 0.02
MN 0.33$$
16 TITLE: HD DUCKTOWN USNM 119197
FE 5.69SI 53.55MG 15.04CA 24.18AL 0.40NA 0.50TI 0.03CR 0.01
MN 1.18$$
24 TITLE: H UVAROVITE USNM 139715 ZONED
FE 2.19SI 38.22MG 0.64CA 33.04AL 7.94NA 0.14TI 0.04CR 18.58
MN 0.84$$
29 TITLE: KAK GNT
FE 10.62SI 42.30MG 19.18CA 5.11AL 23.57NA 0.16TI 0.41CR 0.08
MN 0.49$$
1 TITLE: HD TEST
FE 23.13SI 48.49MG 4.26CA 20.76AL 0.08NA 0.47TI 0.03CR 0.00
MN 2.54$$

```

6.

4.

c. Printout of sequential disk file for auxillary Program PXVOID

```

1
10 HD BAILAL LAKE USNM R3005
SI 54.04AL 0.86TI 0.00CR 0.02FC 0.00FE 3.41
MG 16.70MN 0.33CA 25.39LI 0.00NA 0.26K 0.00
16 HD DUCKTOWN USNM 119197
SI 53.55AL 0.40TI 0.03CR 0.01FC 0.00FE 5.69
MG 15.04MN 1.18CA 24.18LI 0.00NA 0.50K 0.00
24 H UVAROVITE USNM 139715 ZONED
SI 38.22AL 7.94TI 0.04CR 18.58FC 0.00FE 2.19
MG 0.64MN 0.84CA 33.04LI 0.00NA 0.14K 0.00
29 KAK GNT
SI 42.30AL 23.57TI 0.41CR 0.08FC 0.00FE 10.62
MG 19.18MN 0.49CA 5.11LI 0.00NA 0.16K 0.00
5 HD M12330 AV OF 4
SI 49.32AL 0.21TI 0.01CR 0.01FC 0.00FE 23.02
MG 4.11MN 2.58CA 20.78LI 0.00NA 0.46K 0.00
29 KAK GNT
SI 42.30AL 23.57TI 0.41CR 0.08FC 0.00FE 10.62
MG 19.18MN 0.49CA 5.11LI 0.00NA 0.16K 0.00

```


Appendix III. Examples of Output created by Program RDARL4.

d. Printout of sequential disk file for auxillary program RECAMP

```

1 HD TEST
NA 0.47MG 4.26AL 0.08SI 48.49K 0.00CA 20.76TI 0.03MN 2.54FE 23.13
2 HD TEST
NA 0.47MG 4.21AL 0.50SI 49.07K 0.00CA 20.75TI 0.01MN 2.53FE 22.82
3 HD TEST
NA 0.45MG 4.11AL 0.11SI 49.99K 0.00CA 20.85TI 0.00MN 2.63FE 22.98
4 HD TEST
NA 0.45MG 3.88AL 0.14SI 49.73K 0.00CA 20.75TI 0.01MN 2.62FE 23.14
5 HD M12330 AV OF 4
NA 0.46MG 4.11AL 0.21SI 49.32K 0.00CA 20.78TI 0.01MN 2.58FE 23.02
25 KAK GNT
NA 0.15MG 19.06AL 23.61SI 42.38K 0.00CA 4.99TI 0.48MN 0.50FE 10.63
1 HD TEST
NA 0.47MG 4.26AL 0.08SI 48.49K 0.00CA 20.76TI 0.03MN 2.54FE 23.13

```

e. Printout of sequential disk file for auxillary Program OXCALC

```

1
10, 4,HD BAILAL LAKE USNM R3005
SI, 54.04
TI, 0.00
AL, 0.86
V , 0.00
CR, 0.02
FC, 0.00
FE, 3.41
MN, 0.33
MG, 16.70
CA, 25.39
ZN, 0.00
NI, 0.00
NB, 0.00
16, 5,HD DUCKTOWN USNM 119197
SI, 53.55
TI, 0.03
AL, 0.40
V , 0.00
CR, 0.01
FC, 0.00
FE, 5.69
MN, 1.18
MG, 15.04
CA, 24.18
ZN, 0.00
NI, 0.00
NB, 0.00
24, 9,H UVAROVITE USNM 139715 ZONED

```

Appendix IIIf. Description of BASIC virtual array file.

The direct access file or BASIC virtual array created by Option 5 is written in internal (binary) computer format; its contents cannot readily be examined or edited. The record length in 2-byte words is

$$\text{LENGTH} = 4 * (\text{NELMTS} + 1) + 3$$

Each record includes the number of analyses in the file (KOUNTR) and the number of elements analyzed (NELMTS). Thus, the length of the file can be computed once the first record is read using a default record size of 63, the maximum size that can be written. Data for each record is contained in the following sequence:

KOUNTR
NANAL (I)
NELMTS
WTPCT up to 14 value depending upon NELMTS
WTPCT total of weight present values
CATION up to 14 values
CATION sum of cation numbers

g. BASIC Sequential File for statistics

5, 9
1
23.13
48.49
4.26
20.76
0.08
0.47
0.03
0.00
2.54
2
22.82
49.07
4.21
20.75
0.50
0.47
0.01
0.00
2.53
3
22.98
49.99
4.11
20.85
0.11
0.45
0.00
0.03
2.63
4
23.14

Appendix IIIh. BASIC Sequential File or auxillary programs used to calculate mineral formula of amphiboles, micas, cordierite, feldspars, chlorite, and garnets.

```

1
SI
TI
AL
FC
CR
FE
ZN
MN
NI
MG
CA
SR
BA
NA
K
P
F
CL
10,HD BAILAL LAKE USNM R3005
*GAPY , 6.000
54.04, 1.964
0.00, 0.000
0.86, 0.037
0.00, 0.000
0.02, 0.001
3.41, 0.103
0.00, 0.000
0.33, 0.010
0.00, 0.000
16.70, 0.904
25.39, 0.989
0.00, 0.000
0.00, 0.000
0.26, 0.018
0.00, 0.000
0.00, 0.000
0.00, 0.000
0.00, 0.000
16,HD DUCKTOWN USNM 119197
*GAPY , 6.000
53.55, 1.978
0.03, 0.001
0.40, 0.017
0.00, 0.000
0.01, 0.000
5.69, 0.176
0.00, 0.000
1.18, 0.037
0.00, 0.000
15.04, 0.828
24.18, 0.957

```

Appendix IIIi. Printout of disk file for storing "known-unknown"
microprobe analyses with Option 8 of Program RDARL4.

```

# 1 OPERATOR: HUEBNER -= HD & GNT          DATE: 10-29-1982
*QSET=*GAFY DATA STORED IN D1027  NELMTS= 9  NREDCT=2
HD TEST                                AVERAGE OF 0
FE 23.13SI 48.49MG  4.26CA 20.76AL  0.08NA  0.47TI  0.03CR  0.00MN  2.54
SUM= 99.77
PXHD PXAD PXAD PXAD FSTA FSTA OXIL OXTB OLST
# 2 OPERATOR: HUEBNER -= HD & GNT          DATE: 10-29-1982
*QSET=*GAFY DATA STORED IN D1027  NELMTS= 9  NREDCT=2
HD TEST                                AVERAGE OF 0
FE 22.82SI 49.07MG  4.21CA 20.75AL  0.50NA  0.47TI  0.01CR  0.00MN  2.53
SUM=100.36
PXHD PXAD PXAD PXAD FSTA FSTA OXIL OXTB OLST
# 3 OPERATOR: HUEBNER -= HD & GNT          DATE: 10-29-1982
*QSET=*GAFY DATA STORED IN D1027  NELMTS= 9  NREDCT=2
HD TEST                                AVERAGE OF 0
FE 22.98SI 49.99MG  4.11CA 20.85AL  0.11NA  0.45TI  0.00CR  0.03MN  2.63
SUM=101.15
PXHD PXAD PXAD PXAD FSTA FSTA OXIL OXTB OLST
# 4 OPERATOR: HUEBNER -= HD & GNT          DATE: 10-29-1982
*QSET=*GAFY DATA STORED IN D1027  NELMTS= 9  NREDCT=2
HD TEST                                AVERAGE OF 0
FE 23.14SI 49.73MG  3.88CA 20.75AL  0.14NA  0.45TI  0.01CR  0.00MN  2.62
SUM=100.73
PXHD PXAD PXAD PXAD FSTA FSTA OXIL OXTB OLST
# 5 OPERATOR: HUEBNER -= HD & GNT          DATE: 10-29-1982
*QSET=*GAFY DATA STORED IN D1027  NELMTS= 9  NREDCT=2
HD M12330 AV OF 4                      AVERAGE OF 4
FE 23.02SI 49.32MG  4.11CA 20.78AL  0.21NA  0.46TI  0.01CR  0.01MN  2.58
SUM=100.50
PXHD PXAD PXAD PXAD FSTA FSTA OXIL OXTB OLST
# 25 OPERATOR: HUEBNER -= HD & GNT        DATE: 10-29-1982
*QSET=*GAFY DATA STORED IN D1027  NELMTS= 9  NREDCT=2
KAK GNT                                AVERAGE OF 0
FE 10.63SI 42.38MG 19.06CA  4.99AL 23.61NA  0.15TI  0.48CR  0.10MN  0.50
SUM=101.89
PXHD PXAD PXAD PXAD FSTA FSTA OXIL OXTB OLST

```

Appendix IIIj: Printout of sequential file for archiving, concatenating,
and editing analyses

```

    9,***, 6.000
FE
SI
MG
CA
AL
NA
TI
CR
MN
TOTAL
 10,HD BAILAL LAKE USNM R3005           , 6.000,  4
   3.41, 0.103
  54.04, 1.964
   16.70, 0.904
   25.39, 0.989
    0.86, 0.037
    0.26, 0.018
    0.00, 0.000
    0.02, 0.001
    0.33, 0.010
 101.01, 4.027
  16,HD DUCKTOWN USNM 119197           , 6.000,  5
   5.69, 0.176
  53.55, 1.978
   15.04, 0.828
   24.18, 0.957
    0.40, 0.017
    0.50, 0.035
    0.03, 0.001
    0.01, 0.000
    1.18, 0.037
 100.58, 4.030
  24,H UVAROVITE USNM 139715 ZONED     , 6.000,  9
   2.19, 0.072
  38.22, 1.513
   0.64, 0.038
  33.04, 1.401
   7.94, 0.370
   0.14, 0.011
   0.04, 0.001
  18.58, 0.581
   0.84, 0.028
 101.63, 4.016
  29,KAK GNT                             , 6.000,  4
   10.62, 0.312
   42.30, 1.489
   19.18, 1.006
    5.11, 0.193
   23.57, 0.978

```

Appendix IV. Printout of instructions contained in Reston Disk File

(Press the NO SCROLL key at the lower left of the keyboard to halt the motion of the display; press the key again to advance the display.)

WELCOME TO RDARL4!

RDARL4 is a Fortran program that enables the MINC computer to read the data files written by the ARLEB operating system in the ARL-SEMQ probe. For input you must use a disk that contains data written in the abbreviated storage format coded by Jim McGee that appears in an SEMQ operating program such as \$ANBA. (The lengthy storage format provided by ARL is not compatible with RDARL4!) RDARL4 will coach you by asking a number of questions; your responses will be alphabetic characters (Y, N, or C), integer numbers (e.g., 12 calls for a one- or two-digit number without decimal point), rarely a floating point number (6.0 or 4.0 are acceptable for F12.7 format), and file names. For disk output, you will need a formatted MINC disk (see the INITIALIZE command, MINC Programming reference, volume 3). You will also need to have made some decisions in advance. Specifically, if you elect to create an output file on disk, you will be asked what use will be made of the disk output (RDARL4 takes care of preparing the correct format for that use). All options will ask you to select the analyses that you wish included in the disk output file. In addition, the option that writes a disk file for MINC regressions needs to know how many analyses you will select. Your current choices for output are:

Option#

- 1 HARD-COPY PRINTOUT
- 1 MINCLC FORTRAN INPUT DISK FILE
- 2 PXVOID " " "
- 3 RECOMP " " "
- 5 MINC-BASIC VIRTUAL ARRAY ON DISK
- MINC BASIC SEQUENTIAL FILES ON DISK for:
- 6 REGRESSION PROGRAMS on the MINC APPLICATIONS PACKAGE disk)
- 7 SILICATE FORMULA RECALCULATION
- 4 OXIDE FORMULA RECALCULATION
- 8 FORTRAN SEQUENTIAL DISK FILE OF KNOWN-UNKNOWN ANALYSES
- 9 BASIC SEQUENTIAL DISK FILE FOR EDITING, CONCATENATING, ARCHIVING

Appendix IV. Printout of instructions contained in Reston Disk File (continued)

If you plan to select more than one kind of output file, you will need to specify different output file names. One way to do this is to use a file name that embeds a mnemonic abbreviation (MIN, PXV, REC, OXD, etc.) for the kind of file. For example, DY1:DA709.MIN would be a file created for the mineral formula program MINCLC.

Now set the computer system's current date and time by using the following format after the monitor's prompt (.):

.DATE 19-AUG-82<return>

.TIME 13:49<return>

If you intend to display the data on this video screen, you can prevent "wraparound" by pressing the keys SETUP, 9, SETUP to change the number of characters per line from 78 to 130. However, in so doing you will erase the current screen display, so read the next three lines first!

To run the program, type RUN RDARL4<return>.

To interrupt the program, type ^C twice. (^C means "control C")

Have Fun!

Steve Huebner

June 10, 1983

Table 1. Structure of a file of chemical analyses created by microprobe

First Record Contains Parameters Common to All Analyses

<u>Entry Number</u>	<u>Contents</u>	<u>RDARL4 Variable or Array Name</u>
1 - 4	Operator's Name	IOPNAM
5 - 7	ARL System Date	IARDT
8	Name of *QFILE	IQFILE
9	Name of Data Output File	IARLOT
10	Data Reduction Flag	NREDCT
11	Number of Elements Analyzed	NELMTS
12	Oxide or Element Wt. % Flag	NTYPE
13	Temporary Counter Used by ARL	NCOUNT
14	Future Use	NOUSED
15 - 28	Element Names, X-ray Lines (up to 14)	BELNAM,BXRLIN
29 - 42	Standard Names (up to 14)	BSTAND

Subsequent Records Contain One Analysis Each

<u>Entry Number</u>	<u>Contents</u>	<u>RDARL4 Array Name</u>
1 - 8	Title	ITITLE
9	Stage Position X	X
10	Stage Position Y	Y
11	Averaging & Known-Unknown Flag	IAVFLG
12 - 25	Weight Percent Values (14)	WTPCT
26	Total Weight Percent	WTPCT
27 - 40	Cation Numbers (14)	CATION
41	Sum of Cations	CATION
42	Number Anions Per Formula Unit	PEROXY