

U. S. DEPARTMENT OF INTERIOR
GEOLOGICAL SURVEY

Coal Petrographic Laboratory
Procedures and Safety Manual

by

T. A. Moore and R. W. Stanton

Open-file Report 85-20

This report is preliminary and has not been reviewed for conformity with USGS editorial standards and stratigraphic nomenclature.

PREFACE

This manual has a dual purpose: 1) to provide to new employees, interested colleagues, the industry, and the public a guide to the methods and procedures used in sample preparation and analysis, within a coal petrographic laboratory, and 2) to act as a reference for safety and consistency in sample preparation in the laboratory. Precise methodology ensures reliable quality control and provides maximum safety potential. In writing this handbook we have tried to set procedures that will allow employees to work under a minimum of supervision.

Many of the sample preparation procedures used in this manual follow the American Society for Testing Materials' (ASTM) recommendations. The appropriate ASTM standard number is noted at the beginning of these procedures. References are cited for sample preparation methods taken from the literature. After each procedure, necessary safety precautions are listed; these **must** be followed each time a particular technique is performed.

Although this manual is as up-to-date as possible, new techniques for sample preparation and analysis are constantly being improved and better materials for use in preparations are constantly being developed. Therefore, updating of this manual will be periodically necessary.

TABLE OF CONTENTS

	Page
PREFACE	i
LIST OF FIGURES	iii
INTRODUCTION	1
Purpose of Laboratory	1
Types of Samples	1
General Rules	2
LABORATORY ORGANIZATION	3
Sample Flow	3
Record Keeping	3
SAMPLE HANDLING	7
Sample Splitting and Contamination	7
Importance of Clear Labeling	7
SAMPLE PREPARATION, ANALYSIS, AND SAFETY PROCEDURES	8
Sample Crushing and Grinding	8
Jaw Crusher	8
Large High-Speed Rotary Mill Grinder	10
Small High-Speed Rotary Mill Grinder	12
Casting of Pellet and Block Samples	15
Pellet Making — Anthracite and Bituminous samples.	15
Pellet Making — Lignite and Subbituminous samples.	17
Microblock Casting	19
Microsample Casting	21
Hand Sample Casting	23
Vacuum Impregnating	25
Polishing and Etching of Pellets, Microblocks, Microsamples and Hand Samples	27
Pellet, Microblock and Microsample Polishing	27
Hand Sample Polishing	30
Etching of Polished Pellets and Other Mounts	32
Concentration and Preparation of Dispersed Organic Matter	34
Petrographic and Other Physical Analyses	40
Density Determination	40
Low Temperature Ashing (LTA)	43
High Temperature Ashing (HTA)	46
Volatile Matter Determination	48
X-Ray Radiography of Drill Core Samples	50
Vitrinite Reflectance	53
Maceral Analysis	60
EMERGENCY PROCEDURES AND PRECAUTIONS	66
Eye Wash Station	66
Fire Extinguisher	66
First Aid	66
Acid Spills	66
ACKNOWLEDGMENTS	66
CITED AND SELECTED REFERENCES	67

LIST OF FIGURES

	Page
Figure 1: Flow diagram of sample processing	4
Figure 2: Flow diagram of record keeping scheme	5
Figure 3: Sample status sheet	6
Figure 4: Diagram of microscope set up	54
Figure 5: Calibration matrix	56
Figure 6: Classification sheet for Anthracite and Bituminous Coals	62
Figure 7: Classification Sheet for Lignite and Subbituminous Coals	64

1. INTRODUCTION

1.1 Purpose of Laboratory

The purpose of this laboratory is to prepare coal and associated rock samples for physical and optical analyses. Physical analyses performed in our laboratory consist of chemistry determinations, low temperature ashing (LTA), high temperature ashing (HTA), volatile matter determinations, and X-ray radiography. Optical analyses consist of vitrinite reflectance and spectral fluorescence (measures of coal diagenesis) measurements and maceral (organic) analysis. Subsplittings of samples are also sent to the Branch of Analytical Laboratories (USGS) for chemical and mineralogical analyses.

From the results of these analyses, coal samples can be characterized by their composition, rank, and washability potential. Coal bed variability can also be determined when a suite of samples is analyzed from one coal bed. The detailed study of coal beds is important for both utilization and geologic studies. In addition to being of value in mine planning, coal research gives clues to the geologic processes that formed the coal bed.

1.2 Types of Samples

Coal is the major rock type processed in the laboratory. Occasionally, other rocks that are associated with coal beds are prepared (i.e. claystone, shale, limestone, sandstone). Peat samples are also periodically processed and analyzed.

When coal is received as a bulk sample in the laboratory, it is usually in one of two forms: (1) core sample or (2) coal from a channel sample. A core sample arrives as an intact, oriented unit that readily allows subsampling and description. A channel sample, usually collected from a fresh mine face, arrives as broken and unoriented pieces in a bag. Channel samples collected from outcrops can be weathered, greatly limiting the types of analyses that can be performed reliably.

Although the samples may differ in mode of collection, sample preparations for each type are similar.

1.3 General Safety Rules

- 1. Safety glasses are mandatory for any work performed in the laboratory.**
- 2. Contact lenses are not permitted when using potting resins, acids, or organic solvents.**
- 3. Any visitors and minors may be present only with the permission of the lab supervisor, and must wear safety glasses at all times.**
- 4. The last employee to leave for the day should lock the doors, check status of ovens, furnaces, and all other electrical devices.**
- 5. All reagents should be clearly labeled and stored in designated places.**
- 6. All heavy items and chemicals should be stored below counter level.**
- 7. Fume masks must be worn when pouring resins.**
- 8. Rubber or disposable plastic gloves must be worn when using acetone, methanol, acids, toluene, resins, or colloidal silica polishing compound.**
- 9. No open-toed shoes may be worn in the laboratory or grinding room.**
- 10. The following must be used under the fume hood:**
 - a. Acetone (used for cleaning molds)**
 - b. Methanol (used for cleaning molds)**
 - c. Toluene (used for adding to epoxy and removing lap cloths)**
 - d. Epoxy-Resins (used in preparing pellets)**
 - e. Shellac (used for spraying core)**
 - f. Polyester Resins (used for potting block and core samples)**
 - g. Acids (used for special preparations)**

N.B.: Acids must not be used or stored with any organic compounds (items a-f).

2 LABORATORY ORGANIZATION

2.1 Sample Flow

Our laboratory follows ASTM standard D2013 (ASTM, 1983a) for the correct procedures of crushing, grinding, sieving, and splitting coal samples. The steps in processing of an idealized sample are illustrated in Figure 1. When a bulk sample is received in the laboratory, it is first crushed to minus 8 mesh (<2.4 mm), then split into three subsamples; one-third is sent to the Branch of Analytical Laboratories (USGS) for mineralogical and chemical analyses; one-third is retained as a laboratory split that is processed further by us; and the remaining third is kept as a storage split. When insufficient sample (less than 100 g) is available for three splits, the storage split is omitted.

Using a large high-speed rotary mill grinder with a 2 mm sieve screen, the laboratory split is crushed to less than 20 mesh (850 um). Three subsplits obtained by using a riffle sampler (opening > 2.5 mm), are as follows: (1) a split for pellet making, (2) a split for density determination, and (3) a split for storage. Again, if insufficient sample is available for all three splits, the storage split is not taken.

After a density determination has been performed on the density subsample, it is ground to less than 60 mesh (250 um) using a small high-speed rotary mill grinder. Two subsplits are made from the ground sample: one for a high temperature ash (HTA) determination and the other for a low temperature ash (LTA) determination. Because little sample is needed for the HTA, only a very small subsplit is taken (about 2). After the LTA has been performed, the remaining unashed sample is retained as a storage split. The ashed remains are held for X-ray mineralogy.

To summarize, a sample progresses through three basic preparation steps: (1) preliminary crushing and splitting at 8 mesh, (2) crushing to 20 mesh for the petrographic splits, and (3) 60 mesh fine-grinding of the density subsplit for ash determinations. This preparation procedure gives representative subsamples with the least amount of contamination.

2.2 Record Keeping

Laboratory organization not only entails the physical handling the samples, but also the maintenance of sample status records. It is just as important to know a sample's stage of preparation as it is to know where it is physically located. Because our laboratory, at any particular time, may have 200 to 300 samples at nine or ten different stages of processing, written ledgers are used to note sample progression.

In our laboratory, eight notebooks are kept (figure 2). Six of these notebooks are used primarily for on-the-spot data entry for specific procedures. These are the Density, HTA, LTA, and Volatile Matter data entry books and the Pellet and Microblock data books. A sample's progress in these and all other procedures, is noted in a seventh notebook, the Sample Status Book. The sample status sheets (figure 3) were designed to logically follow the flow of sample preparation. The eighth notebook, the Values Book, compiles the final data from the LTA, HTA, Density, and Volatile Matter determinations, thus allowing easy access and comparison of data for any sample.




FIGURE 1: Flow diagram of idealized sample preparation




FIGURE 2: Flow diagram of record keeping

[illegible]

6

3 SAMPLE HANDLING

3.1 Sample Splitting and Contamination

All final interpretations of data depend on the quality of sample preparation in the laboratory. Ideally, any subsplit that is analyzed will be representative of the original bulk sample. To assure that the sample is representative, mechanical riffles should be used. Particle size must also be considered when subsplits are small. Improper splitting gives unrepresentative splits and leads to erroneous results.

Sample contamination is another common problem within the laboratory. Grinding machine, splitting riffles, mortar and pestles, sieves, and other equipment that come into contact with the coal or rock samples must be thoroughly cleaned between samples (this means using acetone and disposable wipes). In addition, laboratory counter tops, fume hoods, and the general working area should be kept as clean as possible to prevent dust contamination of samples.

3.2 Importance of Clear Labeling

Personal notation does not work within the laboratory. Each sample bag and container must be neatly and clearly marked as to its contents so others may understand it. A mislabeled sample can have disastrous effects on interpretations. Any sample that is spilled, dropped, or lost should be immediately reported to a supervisor (i.e. **NO** sample is better than a mismarked or unrepresentative sample). If any question arises about a sample number, it should be noted on the bag or container and a supervisor should be informed.

4 SAMPLE PREPARATION, ANALYSIS, AND SAFETY PROCEDURES

4.1 Sample crushing and grinding (ASTM 2013, 1983a)

4.1A Jaw Crusher (8 mesh grinding)

Equipment	jaw crusher large riffler (3/8"-1/2" chutes) acetone disposable wipes plastic bags permanent marker pen Sample Status Book
Safety Equipment	dust respirator safety glasses/goggles ear mufflers rubber gloves closed-toe shoes laboratory coat and/or apron

4.1A.1 Preliminary

- a. Thoroughly clean crushing plates of jaw crusher with acetone and disposable wipes.
- b. Check minimum gap between plates to insure the spacing is 1/4".
- c. Place cardboard shelf under crushing plates.
- d. Clean and place plastic sample catcher bin under plates and on top of cardboard shelf.
- e. Turn on hood vacuum.

4.1A.2 Crushing

- a. Turn on jaw crusher.
- b. Place bulk sample between crushing plates. Crush until all sample has passed or the bin becomes full.
- c. Turn off machine and open crushing plates.
- d. Sweep all residue sample off the plates and lower platform into bin.
- d. Thoroughly clean plates and platform with acetone and high pressure air.

4.1A.3 Splitting

- a. Split sample into thirds using large riffler. Mark the bags "Analytical labs", "Petrography", and "Storage". If there is insufficient sample for the 3 splits omit the storage split.
- b. Clean machine with acetone.

4.1A.4 Record Keeping

Record all grinding and splitting in the Sample Status Book.

4.1A.5 Safety Procedures

a. Using the Jaw Crusher

- dust respirator, ear protection, and safety glasses/goggles must be worn at all times.
- laboratory coat and/or apron must be worn
- no open-toed shoes or sandals may be worn
- keep fingers and clothing away from crushing plates

b. Acetone

- avoid breathing fumes
- wear protective gloves
- avoid skin and eye contact

4.1B Large High-Speed Rotary Mill Grinder (for 20 mesh grinding)

Equipment large high-speed rotary mill grinder
disposable wipes
mortar and pestle
permanent marker pen
riffler (medium, 1/4" chutes)
20 mesh sieve
2 oz. jar(s)
4 oz. jar(s)
razor knife
Sample Status Book

Safety Equipment dust respirator
safety glasses/goggles
ear mufflers
protective gloves
closed-toe shoes
laboratory coat and/or apron

4.1B.1 Preliminary

- a. Check that the proper sizing screen is in the machine (2 mm screen for 20 mesh).
- b. **Unplug machine** and clean machine thoroughly using high pressure air, acetone, and disposable wipes.
- c. Clean the metal sample-catching cup.

4.1B.2 Grinding

- a. Close mill door and plug in machine.
- b. Close slide between overhead sample holder and grinding chamber.
- c. Place sample into overhead sample holder (do not fill over 3/4 full), proceed through step 6, then close slider and pour in the remaining sample.
- d. Put cap on sample holder.
- e. Turn on machine.
- f. Open slider.
- g. Let sample grind completely.
- h. Turn off machine and listen for "click" that indicates the blades have stopped.
9. Open mill door.
- h. Check to see that all of the sample has passed through the screen (if not, close door and continue to grind until all of the sample passes through the screen).

- i. Shake spring-loaded bottom of mill to assure that all of the sample has fallen into metal cup.
- j. Remove metal cup from under the mill and place sample into a 20 mesh sieve.
- k. Clean machine and metal cup thoroughly with high pressure air and acetone.

4.1B.3 Sieving

- a. Place crushed sample into 20 mesh sieve (step 12 under 4.1b.2 above).
- b. Gently shake sample through screen.
- c. Place remaining sample which does not pass through the screen into a mortar and pestle.
- d. Grind by hand until all of the sample passes through the 20 mesh screen.
- e. Clean sieving equipment and mortar and pestle with high pressure air and acetone.

4.1B.4 Splitting

- a. Use medium-size riffler (>1/4" riffle chutes) and obtain successive subsplits:
 - Split for Density (use 4 oz. jar and make a subsplit of about 20 g).
 - Split for Pellet (use 2 oz. jar, for a 60-g subsplit).
- b. Bag and heatseal remainder and mark as Storage split -20 mesh.
- c. Clean riffler equipment using high pressure air.

4.1B.5 Record Keeping

Record all grinding and splits in Sample Status Book.

4.1B.6 Safety Procedures

- a. **Using the large High-Speed Rotary Knife Grinder**
 - dust respirator, ear protection, and safety glasses/goggles must be worn at all times
 - laboratory coat and/or apron should be worn
 - no open-toed shoes or sandals may be worn
 - do not ever open grinder door while blades are still moving
 - always unplug machine before cleaning inside of the chamber

2. Acetone

- avoid breathing fumes
- avoid skin and eye contact
- wear protective gloves

4.1C Small High-Speed Rotary Mill Grinder (for 60 mesh grinding)

Equipment	small high-speed rotary mill grinder belt set for high speed grinding (pulley ratio 2x motor) 20 and 40 mesh delivery tubes collection jar with special lid oven, capable of maintaining 105°C acetone disposable wipes 60 mesh sieve with bottom pan plastic bags and heatsealer vials permanent marker pens microsplitter (1/8" chutes) brush Sample Status Book
Safety Equipment	dust respirator safety glasses/goggles closed-toe shoes rubber gloves laboratory coat and/or apron

4.1C.1 Preliminary

- a. Preheat oven to 60°C for subbituminous coals or 105°C for higher rank coals.
- b. Place opened sample jars (density jars) in oven for one hour.
- c. Clean high speed rotary mill and other equipment thoroughly with high pressure air and a brush.
- d. Thoroughly clean 40 mesh delivery tube and collection jar.

4.1C.2 Grinding and Sieving

- a. Turn on grinder.
- b. Place dried sample in hopper letting the inverted sample jar act as a cover.
- c. Let machine grind for about 45 seconds.
- d. Turn off machine and remove the collection jar.
- e. Place the sample from the collection jar into a 60 mesh sieve.

- f. Place collection jar back into grinder and turn on grinder again.
- g. Sieve sample.
- h. Place the portion that does not pass sieve back into grinder (do this while grinder is still on).
- i. Continue to process until all of the sample passes the 60 mesh sieve.
- j. Turn off machine and with the sieve placed below the discharge tube remove and clean glass plate and discharge tube.
- k. Clean machine thoroughly with acetone and air.

4.1C.3 Splitting

- a. Use small-sized riffler (microsplitter, 1/8" riffle chutes).
- b. obtain a successive split for high temperature ash (HTA) (use small vials to split off about 2-3 grams).
- c. Mark vial HTA split -60 mesh.
- d. Place remainder of sample in plastic bag.
- e. Heatseal and mark bag LTA split -60 mesh.
- f. Clean riffle equipment using high pressure air.

4.1C.4 Record Keeping

Enter all grinding and splitting information in Sample Status Book.

4.1C.5 Suggestions

- a. It is important to thoroughly clean the mesh on the delivery tube. If there are problems cleaning the mesh, place tube into an ultrasonic cleaner. Dry tube by placing in oven.
- b. If the mesh of the delivery tubes consistently clogs, try heating the tube before grinding.
- c. Grinder and tube should be cleaned with acetone often when grinding shaley and dirty coals.
- e. Rotary blade gap should be 0.015 inch.

4.1C.6 Safety Procedures

a. Using the small High-Speed Rotary Mill Grinder

- use in vented hood
- dust respirator and safety glasses must be worn at all times
- laboratory coat and/or apron should be worn
- no open-toed shoes or sandals may be worn
- do not remove glass plate while blades are moving
- do not insert fingers into chamber for any reason unless machine is unplugged
- after adjusting blades always check to see if rotor will turn 360° before switching machine on.

b. Acetone

- avoid breathing fumes
- avoid skin and eye contact
- protective gloves must be worn

4.2 Casting of pellet and block samples

4.2A Pellet Making

Equipment stainless steel pellet molds (I.D.=1")
with 1/4" plugs and 1" plugs
epoxy-resin
wooden stirring sticks
paper cups
plastic mixing rod
marking pen
heating oven
masking tape
release agent (fluorocarbon)
hydraulic press
analytical balance (triple beam)
Pellet Book
Sample Status Book

Safety Equipment rubber gloves
rubber apron or laboratory coat
organic vapor respirator
safety glasses/goggles
fume hood

4.2a.1 Procedure for Anthracite and Bituminous Coal (ASTM D2797, 1983b)

a. Preliminary

1. Label 3 molds with sample number.
2. Label 2 paper cups with sample number.
3. Spray release agent (fluorocarbon) on molds of the 2 cups.
4. Pour equal (approximately 30 g) amounts of coal into each of the 2 cups.

b. Mixing and Pouring of Epoxy (performed in fume hood)

1. Mix epoxy-resin in correct proportions (e.g. 4 parts resin (20 g) to one part activator (5 g) (by weight)).
2. Add 2-3 ml of toluene to epoxy (this will decrease the viscosity of the epoxy).
3. Mix epoxy with coal in cup until coal is "just wet."
4. Place plug into the bottom of the mold.
5. Place equal amounts of coal/epoxy mixture into each of the 4 molds.

c. Pressing and Compacting the Mixture

1. Place pellet mold under hydraulic press and pump to 7000 p.s.i. and release; continue to pump and release until dial stabilizes (approximately 8-10 times).
2. Pump to 7000 p.s.i. again and let set for at least 3 min. under pressure.
3. Remove, repeat for each mold.
4. Put aside to cure overnight.

d. Popping and Labeling of Pellets

1. Extract pellet from mold by using the side plunger on the hydraulic press.
2. Persist in removal of all plugs.
3. Replace pellet into the same mold, remembering to space the bottom 1/2" (i.e. use two thin plugs for proper spacing).
4. Place appropriate typed color-coded label on top surface of pellet.
5. Pour fresh, clear, epoxy-resin to the rim of the mold.
6. Let cure overnight (for fast curing, place in oven for one hour 60°C).
7. Pop pellet out using hydraulic press (to avoid damaging the pellet, place 1/4" plug on coal surface and press out on plug).

e. Cleaning of Molds (performed in fume hood)

1. Place empty, used molds and plugs in an acetone bath overnight.
2. Scrub molds with steel wool in fume hood
3. Rinse with acetone.
4. Let dry.
5. Store the clean, dried molds and plugs away from dirt and dust.

f. Record Keeping

1. Record in Pellet Book the date that pellet was made and labeled.
2. Record in Sample Status Book that pellet was made.

4.2A.2 Procedure for Lignite and Subbituminous Coals

a. Preliminary

1. Dry sample in oven in a petri dish for at least 2 hours at 60°C.
2. Label 3 molds with sample number and place in oven with sample for 2 or more hours.
3. Label 2 paper cups with the sample number.

b. Mixing and Pouring of Epoxy

1. Mix epoxy-resin in correct proportions (e.g. 4 parts resin to one part activator (by weight) (about 20 g to 5 g).
2. Add 2-3 ml of toluene to epoxy of dry, (to reduce viscosity).
3. Pour approximately equal amounts of dry, heated coal into each of the 2 cups (total 60 g).
4. Mix epoxy into coal until the mixture appears "gummy" or "oozy".
5. Remove molds and plugs from oven.
6. Spray release agent (fluorocarbon) on molds and plugs.
7. Place equal amount of "gummy" coal mixture into each of the 4 molds.

c. Pressing and Compacting the Mixture

1. Place pallet mold under hydraulic press and pump to 5000 p.s.i. and release; pump again until dial stabilizes (approximately 8-10 times).
2. Pump again to 5000 p.s.i. and leave for 2-3 min. under pressure.
3. Remove, repeat for each mold.
4. Put aside to cure overnight.

d. Popping and Labeling of Pellets

Same as for anthracite and bituminous coals

e. Cleaning of Molds

Same as for anthracite and bituminous coals

f. Record Keeping

1. Record date pellet was made and labeled in Pellet Book.
2. Record in Sample Status Book that pellet was made.

4.2a.3 Safety Procedures

1. Handling Epoxy

- organic vapor respirator, safety glasses/goggles, and rubber gloves must be worn
- do not wear contact lenses
- avoid breathing fumes, use under fume hood
- laboratory coat and/or apron should be worn
- avoid skin and eye contact

2. Hydraulic press

- safety glasses/goggles should be worn
- avoid overloading press
- properly secure mold under press
- make sure press features are tightened

C. Cleaning Molds – use of acetone

- rubber gloves must be worn
- safety glasses/goggles should be worn
- use acetone under fume hood, avoid breathing fumes
- use hand lotion or cream on hands to prevent excessive drying

4.2B Microblock Casting

Equipment stainless steel pellet molds (I.D.=1")
with 1/4" plugs
rock saw
epoxy-resin
wooden stirring sticks
paper cups
marker pen
masking tape
release agent (fluorocarbon)
analytical scale (triple beam)
Sample Status Book
Microblock Book

Safety Equipment rubber gloves
laboratory apron or coat
organic vapor respirator
safety glasses/goggles
fume hood

4.2B.1 Cutting

- a. If sample is oriented, clearly mark orientation.
- b. Using a water lubricated rock saw (or band saw) cut a cube the correct size to fit into a stainless steel pellet mold (I.D. = 1").
- c. Let sample dry at room temperature for 8 hours or put into oven at 60° for one hour.

4.2B.2 Casting

- a. Label mold with sample number.
- b. Place two 1/4" plugs in bottom of mold and tape into place.
- c. Spray release agent (fluorocarbon) into mold.
- d. Mix a small amount of epoxy-resin in the correct proportion (4 parts Armstrong C-4 to one part activator "D" (by weight)).
- e. Place sample into mold and pour a few drops of epoxy into the bottom.
- f. Let set up (about 2 hours--this will prevent "floating" of the sample).
- g. Place label on sample in mold.
- h. Mix and pour fresh epoxy-resin to the rim of the mold.
- i. Let cure overnight or in an oven at 60°C for one hour.

4.2B.3 Record Keeping

- a. Record in Microblock Book date sample was prepared.
- b. Record in Sample Status Book that a microblock was made.

4.2B.4 Safety Procedures

a. Rock Saw or Band Saw

- safety glasses/goggles and ear protection must be worn
- laboratory coat and/or apron should be worn
- make sure all saw features are tightened
- keep fingers and clothing away from blade and moving parts
- never force blade to cut too fast
- be sure the proper blade for the job is used

b. Handling Epoxy

- organic vapor respirator, safety glasses/goggles, and rubber gloves must be worn
- avoid breathing fumes, use under fume hood
- avoid skin and eye contact
- laboratory coat and/or apron should be worn
- do not wear contact lenses

4.2C Microsample Casting

Equipment stainless steel mold (I.D. = 1") with
plugs
epoxy-resin
wooden stirring sticks
paper cups
release agent (fluorocarbon)
hand drill or drill press
Sample Status Book
Pellet Book

Safety Equipment organic vapor respirator
rubber gloves
laboratory apron or coat
safety glasses/goggles
fume hood

4.2C.1 Preliminary – Making Epoxy "Blank"

- a. Make epoxy-resin blank by pouring epoxy mix into a stainless steel mold spaced at the bottom with two 1/4" plugs.
- b. Let cure for one day or one hour in an oven at 60°C.
- c. Pop out epoxy blank.
- d. Polish one side of epoxy blank to at least 5 micrometer polish.
- e. Drill a small, shallow hole of appropriate size (depending on the amount and the dimensions of the sample grain(s),

4.2C.2 Mixing and Casting

- a. Mix crushed sample with epoxy either in a paper cup or in the prepared, drilled hole (depending on amount of sample).
- b. Make sure to fill hole with mixture level with the surface of the pellet and tap pellet mold on counter surface gently to bring air bubbles to surface.
- c. Wipe away any excess epoxy from the surface.
- d. Let cure one day or 1–2 hours in an oven at 60°C.

4.2C.3 Record Keeping

- a. Record in Pellet Book that sample was prepared.
- b. Record in Sample Status Book that a microsample was made.

4.2C.4 Safety Procedures

a. Handling Epoxy

- organic vapor respirator, safety glasses/goggles, and rubber gloves must be worn
- avoid breathing fumes, use under fume hood
- avoid skin and eye contact
- laboratory coat and/or apron should be worn
- do not wear contact lenses

b. Using Drill

- safety glasses/goggles must be worn
- keep fingers and clothing away from drill bit and moving parts

4.2D Hand Sample Casting

Equipment polyester or epoxy resin and hardener
cake/pie pans(s)
release agent (silicone)
paper cups
dye (optional)

Safety Equipment organic vapor respirator
rubber gloves
safety glasses/goggles
laboratory apron or coat
fume hood

4.2D.1 Cutting

- a. Coal block or core piece should be taped with masking tape to prevent breakage.
- b. Sample should be cut to produce desired sample size.
- c. Cut sample should be dried in oven (60°C) for 2 hours

4.2D.2 Casting

- a. Spray tray with silicone release agent.
- b. Mix polyester (or epoxy resin) in appropriate proportions.
- c. Add dye, if desired.
- d. Place sample, cut-face down, in pan and pour plastic. (Note: depending on the rank of the sample, it may be best to pour plastic over all or only part of the block.
- e. Let sample cure for a day or more at room temperature inside a fume hood.
- f. After removing sample from pan, excess sample and plastic can be trimmed off edges and surface (opposite cut face). A rim can be made by using masking tape around edges and a label and another layer of clear resin can be applied to the first cut face side.

4.2D.3 Suggestions

- A. If the polyester or epoxy-resin does not completely setup after a day, place in an oven for an hour or two at 60°C. Remove and let cool; plastic should have hardened.
- B. Addition of some dyes may retard hardening process; to compensate, add a higher proportion of hardener.

4.2D.4 Safety Procedures

1. Rock Saw or Band Saw

- safety glasses/goggles must be worn
- laboratory coat and/or apron should be worn
- make sure all saw features are tightened
- keep fingers and clothing away from blade and moving parts
- never force blade to cut too fast
- make sure proper blade is used

2. Handling of Plastic-Resin

- organic vapor respirator, safety glasses/goggles and rubber gloves must be worn
- avoid breathing fumes, use under fume hood
- avoid skin and eye contact
- laboratory coat and/or apron should be worn
- do not wear contact lenses

4.2E Vacuum Impregnating

Equipment vacuum bell jar
vacuum pump
stainless steel mold (I.D. = 1") and
1/4" plugs
cake/pie pan(s)
polyester or epoxy-resin

Safety Equipment organic vapor respirator
safety glasses/goggles
rubber gloves
laboratory apron or coat
fume hood

[It is often desirable to vacuum-impregnate coals that are low-rank, "dirty", or are cores that will be slabbed (cut into 1/4" thickness) for X-ray radiography. Vacuum impregnating causes pores in the coal to fill with resin; the more resin that penetrates the coal, the less likely that a sample will crack or flake. Vacuum impregnating may be used with core samples, hand samples, microblocks, and crushed coal.]

4.2E.1 Procedure for Hand Sample, Core Sample and Microblock Impregnation

- a. Place cut and prepared sample into tray (for hand samples and core sample) or stainless steel mold (for microblock).
- b. Mix correct amount of resin (either polyester or epoxy). Pour mixture over sample but do not fill to capacity.
- c. Put sample into vacuum chamber (Note: If sample tray covers evacuating/venting hole, raise it up on blocks or place on a platform).
- d. Place vacuum bell jar over sample.
- e. Turn on vacuum pump. Check seal around vacuum bell jar. If there is no seal, check that the bell jar is seated securely on the rubber mat.
- f. Watch sample and vent air into chamber if "boiling" over of resin occurs.
- g. Evacuate for 5-10 minutes.
- h. Turn off vacuum pump and vent air into chamber slowly.
- i. Remove sample to allow to let it cure overnight under a fume hood.

4.2E.2 Procedure for Crushed Samples

- a. Mix -20 mesh, dried, crushed sample with epoxy-resin in a paper cup (see section 4.2a). Use copious amounts of epoxy.

- b. Place sample into a stainless steel mold (with one 1/4" bottom plug) and pour a small amount of epoxy on top.
- c. Put sample and open-topped mold into vacuum chamber.
- d. Follow steps 5 through 8 of preceding procedure.
- e. Remove sample and place top plug into mold.
- f. Using the hydraulic press, place 3000 p.s.i. on sample. Release immediately and repeat 3-4 times.
- g. Let cure overnight. If the sample still flakes or cracks, try using more epoxy and longer evacuating time and less pressure with fewer releases).

4.2E.3 Safety Procedures

a. Handling of Polyester or Epoxy-Resin

- organic vapor respirator, safety glasses/goggles, and rubber gloves must be worn
- avoid breathing fumes, use under fume hood
- avoid skin and eye contact
- laboratory coat and/or apron should be worn

b. Vacuum Bell Jar and Pump

- check oil level in pump
- remain in room while evacuating
- vent air into bell jar slowly

4.3 Polishing of pellets, microblocks, microsamples, and hand samples

4.3A Pellet, Microblock, and Microsample Polishing

(Modified from Cole and Berry, 1965; ASTM D2797, 1983b)

Equipment	automatic polishing machines short-napped polishing cloth 600 & 400 grit paper, 8" diam. 15 micrometers diamond wheel 1.0 micrometers aluminum oxide polishing compound 0.06 micrometer colloidal silica polishing suspension sample holder sample leveler ultrasonic cleaner Sample Status Book Pellet or Microblock Book
Safety Equipment	safety glasses/goggles laboratory apron or coat

4.3A.1 Polishing Pellet Tops

- a. Mount 6 pellets into pellet holder with the tops down (Note: If polishing less than 6 pellets, insert epoxy 'blanks' into remaining holes).
- b. Level and tighten.
- c. Place on 15-micrometer diamond wheel or 400 grit paper for 1-2 min. or until surface is completely flat. Use a steady stream of water. For proper sequence of turning on and off of the lapidary wheel see suggestion B.
- d. Rinse in tap water.
- e. Relevel pellets in holder.
- e. Rough polish for 2 min. on short-napped cloth with 1.0 micrometer aluminum oxide (50:1 mixture of distilled water and polishing compound).
- f. Rinse with tap water.

4.3A.2 Polishing of Sample Bottoms

- a. After polishing the tops of 6 pellets, relevel with the bottoms down.
- b. Grind for 2 min. on a 15-micrometer diamond wheel. Start with an initial pressure of 30 pounds and increase by 5 pound increments until a brownish slurry occurs. Use a slow stream of tap water while the wheel is turning.

- c. Relevel.
- d. Rinse with tap water. Clean ultrasonically for 2 min. then rinse.
- e. Final grind on 600 grit paper for 2 min. Start with an initial pressure of 30 pounds and increase by 5 pound increments until brownish slurry occurs. Use a slow stream of tap water.
- f. Rinse with tap water. Clean ultrasonically for 2 min. then rinse.
- g. Rough polish for 2-3 min. on short-napped cloth with 1.0 micrometer aluminum oxide (1:1 mixture of polishing compound and distilled water). Start with an initial pressure of 30 pounds and increase by increments of 5 pounds until a brownish slurry appears.
- h. Without changing the pressure or stopping the machine, flush with distilled water for 1 1/2 min. to remove particles and excess polishing compound.
- i. Without releasing pressure, turn off motor.
- j. Rinse with distilled water. Clean ultrasonically for 2 min. then rinse.
- k. Final polish for 2 min. on short-napped cloth with 0.06 micrometer colloidal silica. Start with an initial pressure of 30 pounds and increase by increments of 5 pound until slurry occurs.
- l. Without changing the pressure or stopping the machine, flush with distilled water for 1 1/2 min.
- m. Turn off motor without releasing pressure.
- n. Rinse in distilled water. Clean ultrasonically (2 min). Rinse again.
- o. Blow dry and place pellets in desiccator.

4.3A.3 Record Keeping

- a. Record in Pellet Book or Microblock Book that the sample was polished.
- b. Note in Sample Status Book that sample was polished.

4.3A.4 Suggestions

- A. Grinding times may vary for polishing of microsamples depending on thickness of sample.
- B. Sequence for turning on and off polishing laps are as follows:
 - a. turn on water with nothing on the wheel
 - b. turn on motor and run hand over the lap to remove any dirt or excess grit
 - c. turn off motor
 - d. place pellet holder under automatic polishing unit and secure

- e. turn on water
- f. place a little pressure on pellet holder then turn on motor
- g. place appropriate pressure on pellet holder
- h. polish for the recommended time
- i. turn off motor
- j. release pressure
- k. turn off water.

For hard coals, increased polishing times may be needed.

If microprobe or SEM analyses are to be conducted on the Sample dionized water must be substituted for tap water in all steps and new cloths must be put on laps.

4.3A.5 Safety Procedures

A. Automatic Polishing Machine

- safety glasses/goggles should be worn
- keep fingers and clothing away from moving parts, do not insert fingers under machine when in motion
- make sure automatic polishing machine is securely seated and the shaft is locked into position before turning on wheel.

B. Ultrasonic Cleaner

- keep fingers out of cleaner while it is on

C. Colloidal Silica Polishing Suspension

- rubber gloves and safety glasses/goggles must be worn
- avoid skin and eye contact (basic solution)
- after use, lap must be kept moist to prevent drying of polishing suspension and creating hazardous dust and seal in plastic bag in drawer.

4.3B Hand Sample Polishing

Equipment 12" diameter lapidary wheel 60, 120,
240, 320, 600 silicon carbide grit
paper and/or powders
5.0, 3.0, 0.05 micrometers aluminum
oxide polishing compounds
band saw or rock saw

Safety Equipment safety glasses/goggles
laboratory apron or coat

4.3B.1 Polishing cut surface

- a. Remove excess plastic from side of sample to be polished using a band saw or rock saw (be sure to make a straight, even cut).
- b. Use 60 grit silicon carbide powder or grit paper to remove any saw marks or plastic from face of sample. Grind until flat, uniform surface exists (use a slurry of water of grit).
- c. Clean ultrasonically for 2 min. (place sample face down but raised up off the bottom).
- d. Use 120 grit paper or powder and grind surface until all scratches from the previous step have been removed and a uniform surface exists.
- e. Clean ultrasonically for 2 min.
- f. Continue grinding through 240, 320, and 600 grit, removing all scratches from the previous step and clean ultrasonically after each step.
- g. Rough polish on short napped cloth with 5.0 micrometer aluminum oxide.
- h. Rinse. Clean ultrasonically (2 min.). Rinse.
- i. Polish on short napped cloth with 0.3 micrometer aluminum oxide.
- j. Rinse. Clean ultrasonically (2 min.). Rinse.
- k. Final polish on short napped cloth with 0.05 micrometer aluminum oxide or 0.06 colloidal silica polishing suspension.
- l. Rinse. Clean ultrasonically (2 min.). Rinse and blow dry.

4.3B.2 Suggestions

- a. Grinding time at each step will vary considerably between samples; degree of hardness and size of sample controls the amount of time needed.
- b. When polishing, a clockwise rotation (opposite that of the turning lapidary wheel) around the wheel will give the best polishing results.

4.3B.3 Safety Procedures

a. Band Saw or Rock Saw

- safety glasses/goggles and ear protection must be worn
- laboratory coat and/or apron should be worn
- make sure all saw features are secure and tight
- keep fingers and clothing away from blade and moving parts
- never force blade to cut too fast
- make sure proper blade is used

b. Lapidary Wheels

- never place fingers under wheel while moving

c. Ultrasonic Cleaner

- keep fingers out of cleaner when on

d. Colloidal Silica Polishing Suspension

- rubber gloves and safety glasses/goggles must be worn
- avoid skin and eye contact
- after use, lap must be thoroughly rinsed of all polishing suspension

4.3C Etching of Polished Pellets and Other Mounts

(Modified from Stach, p. 338, 1982)

Equipment 50, 100, 500 ml beakers
glass stirring rods
bunsen burner or hot plate
watch glass(es) and/or petri dish(es)
10ml graduated cylinder
sulfuric acid (H_2SO_4 concentrated)
potassium permanganate (KMnO_4)
sodium sulfite (Na_2SO_3)

4.3C.1 Preliminary

- a. Mix etching solution of 100 ml water + 25g KMnO_4
+5ml H_2SO_4 (concentrated).

NOTE: always add acid to water

- b. Mix rinsing solution of 100ml water + 25g Na_2SO_3
+ 5ml H_2SO_4

- c. Stir rinsing solution until all Na_2SO_3
has dissolved.

- d. Heat the etching solution in a water bath until all the KMnO_4 has
dissolved.

4.3C.2 Etching

- a. If a dilution of the etching solution is necessary: (i.e. for low
rank coals)
-- dilute to the desired concentration
-- reheat to boiling for one min. then continue with the etching
procedure.
- b. Immediately pour part of the heated etching solution into a watch
glass.
- c. Submerge face of polished pellet (or other mount) into etching
solution (see table below for recommended etching times).
- d. Remove mount and rinse immediately with water to remove excess
etching solution (1-2 sec.).
- e. Submerge pellet (or other mount) into rinsing solution for one
min. or until purplish stain has been completely removed.
- f. Clean ultrasonically with water for one min.
- g. Blow dry with air immediately.

4.3C.3 Suggestions

A. Recommended etching times for different ranks of coals

Rank	Rel. % Sol.	Approx. Time (sec.)	Example of Coal bed
Low Vol. Bit.	100	55-60	Pocahontas
Med. Vol. Bit.	100	45-55	Upper Freeport
Hi Vol. Bit.	100	20-30	Waynesburg
Subbituminous	75	4-5	Powder River coals

Note: These are only recommended times; different coals beds may require variations in time and solution concentrations.

4.3C.4 Safety Procedures

e. Handling of chemicals

- all chemicals must be used under a fume hood
- heavy rubber gloves, laboratory coat and rubber apron must be worn when handling chemicals and when etching is performed
- safety glasses/goggles and organic vapor respirator must be worn
- avoid skin and eye contact with all chemicals H_2SO_4 will cause severe burns)
- when mixing solutions always add acid to water

b. Use of Bunsen Burner or Hot Plate

- use under fume hood
- keep all chemicals and flammables away from heat
- place burner or hot plate on a level surface

c. Ultrasonic Cleaner

- when the ultrasonic cleaner is working keep fingers out

4.3D Concentration and Preparation of Dispersed Organic Matter (Modified from Barrows and others, 1979)

Equipment vacuum oven
mortar and pestle
acid resistant polyethylene or teflon
containers
hot plate or bunsen burner
500 ml beaker(s)
Hydrochloric Acid, 10% (HCl)
Sodium Carbonate (Na_2CO_3)
Zinc Bromide (ZnBr_2)
Hydrofluoric Acid, 48% (HF)
centrifuge tubes
centrifuge
glass bottles

Safety Equipment organic vapor respirator
safety glasses/goggles
heavy rubber gloves
laboratory coat
rubber apron
fume hood

4.3D.1 Preliminary

- a. Crush sample through jaw crusher (-8 mesh) [see section 4.1 for proper procedure].
- b. Place sample in polyethylene or teflon container.

4.3D.2 Chemical Disaggregation

BEFORE USING, READ FOLLOWING MATERIAL SAFETY DATA SHEETS ON PAGES 36 – 39...

- a. Treat sample with HCl for 1–3 days to dissolve carbonate minerals.
- b. Place sample in acid resistant polyethylene or teflon container resting in a water bath at 50°C under a fume hood.
- c. Decant and rinse with water (3 rinses).
- d. Treat sample with 48% HF for 1–3 days to dissolve silicate minerals.
- e. Decant and add Na_2CO_3 solution to neutralize unspent acid. Rinse with water (7 rinses).
- f. Dissolve HF reaction products with 10% HCl (one day).

4.3D.3 Mechanical Separation

- a. Transfer sample to centrifuge tubes.
- b. Centrifuge and decant water.
- c. Float kerogen fraction in ZnBr_2 solution (s.g.=2.0), centrifuging for a short time.

- d. Decant ZnBr_2 solution and kerogen.
- e. Discard heavy fraction (unreacted shale, pyrite, and other acid insoluble minerals).
- f. Dilute ZnBr_2 with water.
- g. Centrifuge to sink kerogen and decant.
- h. Rinse kerogen with water.
- 9. Centrifuge and decant water (3 rinses).

4.3D.4 Pellet Preparation

- a. Transfer sample to glass bottles and dry in a vacuum oven at 50°C (one day).
- b. Scrape kerogen from bottles and carefully grind in mortar and pestle.
- c. If sample size is large, follow procedure 4.2a for making pellets and then polish according to procedure 4.3a.
- d. If sample size is small follow procedure 4.2c for microsample mounting then polish according to procedure 4.3a.

4.3D.5 Safety Procedures

- a. **Using Hydrofluoric Acid (HF)**
BEFORE USING, READ FOLLOWING MATERIAL SAFETY DATA SHEETS ON PAGES 36 – 39...
 - **WARNING:** Hydrofluoric acid is an extremely hazardous chemical at any concentration or temperature. Physical contact will cause severe burns and inhalation of the vapor will irritate the respiratory system. Use all necessary measures to prevent physical contact with the acid. Consider any accident with hydrofluoric acid a serious accident. **LOCATION OF CALCIUM GLUTONATE GEL SHOULD BE KNOWN IN ADVANCE OF ANY AND ALL USE OF HF ACID.**
 - Heavy rubber gloves, safety glasses/goggles, laboratory coat and rubber apron must be worn at all times.
 - Organic vapor respirator must be worn
 - HF must be kept and used under a fume hood. Use hood's glass door in halfway closed position .
 - do not allow HF to stand in glass containers—it is very corrosive of glass.

Should physical contact with HF result, affected area must be flushed with water immediately for 15 minutes and help called for. Calcium glutonate gel should be applied and massaged into affected area. Victim should then be taken to nearest hospital Emergency Room for treatment.

- c. **Hydrochloric Acid and Zinc Bromide**
 - heavy rubber gloves, safety glasses/goggles, laboratory coat and rubber apron must be worn
 - organic vapor respirator must be worn
 - avoid skin and eye contact
 - must be used under a fume hood



National Safety Council

RESEARCH AND DEVELOPMENT SECTION

SAFETY AND HEALTH INFORMATION SHEET

Title: HYDROFLUORIC ACID EXPOSURE

Hydrofluoric acid (hydrogen fluoride, HF) and substances that can form HF as a result of hydrolysis by moisture are toxic on inhalation and severely corrosive to the skin, eyes and mucous membranes in both the vapor and the liquid states. Concentrated aqueous solutions may cause immediate pain upon contact with the skin. Contact with dilute solutions will result in delayed pain up to 24 hours later. HF penetrates the skin barrier and continues to cause damage until proper treatment is completed. HF forming substances include:

1. Hydrogen Fluoride (HF)
Boron Trifluoride (BF₃)
Phosphorous Pentafluoride (PF₅)

Form copious HF fumes in moist air; visible; extremely irritating; corrosive; early warning of presence.

2. Bromine Pentafluoride (BrF₅)
Bromine Trifluoride (BrF₃)
Iodine Pentafluoride (IF₅)

Form colored vapors from halogen but not HF type fumes; visible; irritating; corrosive; warning of presence.

3. Carbonyl Fluoride (COF₂)
Phosphorous Trifluoride (PF₃)
Fluorine (F₂)
Chlorine Trifluoride (ClF₃)
Sulfur Tetrafluoride (SF₄)
Silicone Tetrafluoride (SiF₄)

Little if any color; generally bad odors; irritating; corrosive; give little warning; can oxidize and hydrolyze in the body.

Precautionary Measures

Because of the extremely hazardous nature of HF the following measures are recommended to protect employees from exposure:

1. Training in the proper use of HF or HF-producing materials should be given before any exposure. The use of HF or HF-producing materials should be limited to laboratories with total exhaust ventilation systems. All work should be accomplished within a fume hood or otherwise be protected by

The information and recommendations contained in this publication have been compiled from sources believed to be reliable and to represent the best current opinion on the subject. No warranty, guarantee, or representation is made by the National Safety Council as to the absolute correctness or sufficiency of any representation contained in this and other publications, and the National Safety Council assumes no responsibility in connection therewith; nor can it be assumed that all acceptable safety measures are contained in this (and other publications), or that other or additional measures may not be required under particular or exceptional conditions or circumstances.

an exhausted enclosure unless prohibited by size of apparatus. In the latter case, local exhaust pick-ups should be located in areas where accidental release of materials is most likely to occur. Shields or see-through plastic wraps are also recommended for likely splash points (valves, fittings, etc.).

2. Quantities of HF or HF-producing materials on hand in the work area should be kept to an absolute minimum consistent with project needs.
3. Protective clothing/devices should be used. Chemical splash goggles are mandatory for any use. Face shields and impervious clothing (gloves, aprons, coats, etc.) should also be used depending on the degree of risk to exposure. Respiratory equipment should be considered for escape/rescue purposes.
4. A means for continuous flushing of the eyes and body with water must be provided. A hand-held aerated emergency spray nozzle that is connected to a continuous supply of water by means of a short (5 ft) (1.52m) length of hose is recommended within the lab. A safety shower should be accessible, preferably within 25 feet (7.6m) of the work area.
5. Persons working with HF or HF-producing materials should not work alone, i.e., when not within audio or visual range of another individual for more than a few minutes at a time.
6. Persons working with or supervising work that involves HF or HF-producing materials should request a review of the project by safety/medical staff before work begins. All project participants should be fully familiar with the hazards of HF; review of the manufacturer's Material Safety Data Sheet is recommended.

Emergency Procedures

Whenever anyone has been seriously exposed to HF or HF-producing materials, emergency action must be initiated immediately. Any person performing first aid must avoid contaminating themselves when assisting an exposed individual. First aid treatment must be administered by the injured coworkers at once as follows.

Liquid or Vapor Contact With Skin

Individuals who have had contact with HF or HF-producing materials must be showered immediately at the work site using available emergency shower facilities. Contaminated clothing should be removed as rapidly as possible. It is essential that all exposed areas of the body be washed with copious amounts of water for at least 15 minutes.



National Safety Council

RESEARCH AND DEVELOPMENT SECTION

SAFETY AND HEALTH INFORMATION SHEET

Title: HYDROFLUORIC ACID EXPOSURE

Immediately begin soaking the area in an iced solution of 0.13 percent benzalkonium chloride ("Zepharin Chloride"). If the burned area cannot be immersed in the solution, use saturated compresses changed every two minutes. Call a physician and continue soaking until further instruction. If no physician is available, continue soaking for two hours and then apply a generous quantity of calcium gluconate gel.

Contact With Eyes

If liquid HF has entered the eyes or if the eyes have been exposed to high concentrations of vapors, they should be flushed with large quantities of clean water for 15 minutes or more using emergency facilities at the accident site. The eyelids should be held apart during irrigation to ensure contact of water with all tissues on the surface of the eye and lids. A physician, preferably an eye specialist, should be called in attendance at the first possible moment. DO NOT USE ZEPHARIN CHLORIDE OR CALCIUM GLUCONATE GEL IN THE EYES.

Inhalation

Any person who has inhaled HF or HF-producing materials must be moved immediately to an area with clean atmosphere. He should remain quiet, preferably lying down, and be kept warm and comfortable. If breathing has stopped, begin artificial respiration (mouth-to-mouth resuscitation). Oxygen rich air (approximately 40 percent oxygen) should be started as soon as possible. Prompt professional attention is required. Under no circumstances should a patient be permitted to return to work or to go home until examined and discharged by a physician who is aware of the nature of the exposure.

Ingestion

Encourage the exposed person to drink a large quantity of water without delay. Then, administer milk or two ounces of Milk of Magnesia.

Storage

Only quantities of HF or HF-producing materials which are needed for a one day use, should be in the work area. Additional supplies should be kept in a chemical storage room.

Spill, Leak, or Disposal

Notify safety personnel, provide adequate ventilation, and remove ignition sources since hydrogen may be generated by reactions with metals. Use protective clothing and equipment. HF vapor should be passed through a packed tower scrubber. Spills should be covered with lime to form a slurry. Do not flush to sewers or waterways.

The information and recommendations contained in this publication have been compiled from sources believed to be reliable and to represent the best current opinion on the subject. No warranty, guarantee or representation is made by the National Safety Council as to the absolute correctness or sufficiency of any representation contained in this and other publications, and the National Safety Council assumes no responsibility in connection therewith; nor can it be assumed that all acceptable safety measures are contained in this (and other publications), or that other or additional measures may not be required under particular or exceptional conditions or circumstances.

Serial no. RSD-83-4

Date: September, 1983

The neutralized slurry can be scraped up for disposal in an approved hazardous waste landfill. Liquid wastes may be neutralized in an impoundment or open tank with lime in a remote location away from buildings and people. Follow Federal, State, and local regulations.

References

"Occupational Health Guidelines for Hydrogen Fluoride," September, 1978.

NIOSH/OSHA Occupational Health Guidelines for Chemical Hazards, DHHS (NIOSH) Publication No. 81-123.

U.S. Governmental Printing Office, Washington, DC 20402.

4.4 Petrographic and other physical analyses

4.4A Density Determination (Stanton, 1982)

Equipment gas comparison pycnometer
analytical balance (sensitivity of
0.01g, with a 200 g capacity)
desiccator and drying agent
oven, capable of obtaining and holding
105°C
Helium
gas drying tube, drying agent
vacuum pump
disposable wipes
acetone
Density Book
Values Book
Sample Status Book

4.4A.1 Preliminary

- a. Obtain prepared density split (20 – 30 grams, minimum 5 grams at 20 mesh).
- b. Place sample into petri dish and label.
- c. Dry in oven for at least 4 hours at 105 – 110°C.
- d. After drying, place in desiccator and permit to cool to room temperature before making volume determination.
- e. Turn on balance and allow it to warm up for about 20 min.
- f. Helium gas source for pycnometer should be set for 3 p.s.i. Gas should be run through lines and drying cylinder before use each day to purge system
- g. A zero measurement with 4 separate determinations must be made at the beginning of each day. The average zero correction that results will be checked by one determination after each sample. The procedure for a zero measurement check is identical to that for a sample except the sample cup is empty. Also, after the initial zero check, a volume determination should be made on one of the two known volume ball standards.

4.4A.2 Density Determination

- a. Clean sample cup with acetone.
- b. Tare cup on balance.
- c. Place all of the prepared density split in cup.
- d. Weigh cup and sample and record.

- e. Make sure reference wheel on pycnometer is at clockwise stop.
- f. Rotate measuring wheel to estimated volume of sample.
- g. Firmly lock sample cup into unit.
- h. Open purge valve and coupling valve.
- i. Turn on vacuum pump.
- j. Open vacuum valve for 60 sec. then close valve and turn off pump, in that order (Note: with large and very powdery samples (shales) it may be necessary to have vacuum introduced gradually to avoid material getting sucked up into machine).
- k. Open gas valve for 10 sec. then close.
- l. Open vent valve for 10 sec. then close.
- m. Rotate reference wheel counterclockwise to stop.
- n. Rotate measuring wheel counterclockwise past the starting number by 5 cubic centimeters.
- o. Rotate measuring wheel clockwise to starting number.
- p. Wait 10 sec.
- q. Close coupling valve.
- r. Turn both wheels (reference and measuring) clockwise simultaneously so that pointer is kept just left of zero on scale. Continue until reference wheel stop is reached.
- s. Rotate measuring wheel counterclockwise to bring pointer to zero.
- t. Open coupling valve. If no significant movement of needle occurs then record value to nearest 0.01 cm³.
- u. Open purge valve and remove sample cup.
- v. Repeat until 3 values are obtained from one sample.

4.4A.3 Record Keeping

- a. Record humidity, temperature, zero correction, and date for each sample in the Density Book.
- b. Record all volume measurements in the Density Book.
- c. Record final average corrected density measurement in the Values Book.
- d. Record that a density determination was made in the Sample Status Book.

4.4A.4 Suggestions

- a. The zero correction factor is determined by finding the number that must be added to the zero check to equal zero. This number is then added to the average volume of sample to equal corrected volume.
- b. For the daily zero measurement check, for values should be obtained consecutively that are within 0.03 cm³. For zero measurement checks, after each sample, the value should be with 0-.03 cm³ of the average daily zero check.
- c. This machine requires patience, consistency and practice by the operator to give repeatable results. For this reason, best results are obtained from one operator doing large groups of samples.
- d. Some possible sources of error are:
 - leakage around cup lip
 - sample absorption emittance of gases (incomplete drying of sample)
 - above average room humidity and temperature
 - incomplete evacuation
 - interruption or slowing of procedures after coupling valve is closed
 - sticky zero pointer

4.4A.5 Safety Procedures

- a. High Pressure Gas Cylinder
 - never bump or knock valve stem
 - always turn off gas at valve stem when finished
 - when transporting cylinders the safety cap must be over the stem
 - cylinder has to be securely anchored to wall with chains
- b. Acetone
 - avoid breathing fumes
 - always wear rubber gloves
 - avoid skin and eye contact

4.5B Low Temperature Ashing (LTA)

Equipment plasma asher
petri dished
analytical balance accurate to at
least 0.0001g
mortar and pestle
100 mesh sieve screen
oven, capable of heating to and
maintaining 105^oC
desiccator with drying agent
oxygen tank
1 N ammonium acetate (pH7)
LTA Book
Values Book
Sample Status Book

Procedure of Anthracite and Bituminous Coals

4.5B.1 Preliminary

- a. Obtain the -60 mesh split marked or LTA.
- b. Pulverize sample with mortar and pestle to -100 mesh.
- c. Place in bottle and label.
- d. Zero check balance.
- e. Weigh out 3g of sample in a clean, preweighed petri dish (dishes are cleaned, oven dried, and desiccated at least 24 hours prior to use). Record weight in LTA Book.
- f. Dry sample in oven at 105^oC for one hour.
- g. Zero check balance.
- h. Weigh sample dish and record.
- i. Calculate % moisture and record in LTA Book.

4.5B.2 Ashing Procedure

- a. Place sample in asher chamber.
- b. Close chamber door.
- c. Close main door.
- d. To turn on LTA – push "on" button.
- e. Then turn on pump (switch box is on pump).
- f. Wait about 5–10 min. for pump down to evacuate to a 1mm vacuum.
- g. Push "RF" button.

- h. One 1mm is reached:
 - turn up RF power to 125 watts (or about 30 watts per chamber)
 - turn up oxygen flow to 1.5 cc

9 .Chambers should exhibit pale blue plasma.

4.5B.3 Stirring and Weighing Samples

- a. Turn "RF" power knob counterclockwise.
- b. Push "RF" button off.
- c. Open main door.
- d. When vacuum has been purged, open chamber doors.
- e. Remove and stir sample by gently shaking sample dish to expose the unoxidized grains (stirring should be done 3 times daily).
- f. Weighting of sample is done one a day.
- g. To weigh, remove and place glass over dish..
- h. Take into balance room and wait 5 min.
- i. Weigh, after first zero checking balance.
- j. Sample is done when the weight changes by less than 0.0001 gram or begins to gain weight (i.e. by absorbing moisture).
- k. Either remove completed sample or note sample weight in the LTA Book.

4.5B.4 Record Keeping

- a. Record % moisture in LTA Book.
- b. Record final ash % in LTA Book and in the Values Book.
- c. Note in Sample Status Book that LTA has been done.

Procedure for Lignite and Subbituminous Coals

(Modified from Given and Yarzab, 1978)

4.5B.1. Preliminary

- a. Obtain the -60 mesh split marked or LTA.
- b. Further pulverize sample to -100 mesh with mortar and pestle.
- c. Place in bottle and label.
- d. Zero check balance.
- e. Weigh out 3g of sample in a clean, preweighed petri dish (dishes are cleaned, oven dried, and desiccated at least 24 hours prior to use).

- f. Dry sample in oven for one hour at 105°C.
- g. Zero check balance.
- h. Weigh sample and dish and calculate moisture.
- i. Mix and stir sample for 1-24 hours (depending on coal) with 1 N ammonium acetate (pH7).
- j. Reweigh to determine weight loss.

4.5B.2 Ashing Procedure

Same as for Anthracite and Bituminous coals

4.5B.3 Stirring and Weighing Sample

Same as for Anthracite and Bituminous coals

4.5B.4 Record Keeping

- a. Record % moisture in LTA Book.
- b. Record final ash % in LTA Book and the Values Book.
- c. Note in Sample Status Book that LTA has been done.

4.5B.5 Safety Procedures

a. Oxygen Tank

- secure oxygen tank with a chain and tank holder to a wall
- never increase the flow of oxygen over 2.0 cc

b. 1 N Ammonium Acetate

- avoid skin and eye contact
- rubber gloves must be used
- organic vapor respirator must be used

4.5C High Temperature Ashing (HTA) (ASTM D3174, 1983c)

Equipment electric muffle furnace
17 ml porcelain crucibles
analytical balance with sensitivity of
0.0001g
desiccator with drying agent
oven capable of heating to and
maintaining 105°C
spatula
tongs
HTA Book
Values Book
Sample Status Book

Safety Equipment safety glasses/goggles
heat resistant gloves
laboratory coat

4.5C.1 Preliminary

- a. Obtain the -60 mesh split marked for HTA.
- b. Place samples into petri dishes and then into an oven for one hour at 105°C.
- c. Place clean crucibles (17ml) into oven with samples.
- d. After an hour remove sample and crucibles; place into desiccator and let cool to room temperature.
- d. Zero check balance.
- e. Weight crucible and record in HTA Data Book as dish weight (D).
- f. Place about a gram of sample into crucible and record as dish and sample weight (DS).

4.5C.2 Ashing

- a. Place sample into cold furnace (i.e. room temperature, about 25°C).
- b. Turn on fan for furnace.
- c. Heat furnace to 500°C within first hour (8°C per min.)
- d. heat furnace to 750°C in the second hour (4°C per min).
- e. Maintain temperature at 750°C for 60 min).
- f. Turn off furnace.
- g. Let furnace cool to 400°C.

- h. Remove sample and place in desiccator (Note: make sure nothing in desiccator is flammable).
- i. Let cool to room temperature (about 8–12 min.).
- j. Zero check balance.
- k. Weigh crucible and ash, then record in ledger as dish and ash weight (DA).
- l. Calculate ash: $\{DA-D/DS-D\} \times 100 = \% \text{ Ash.}$

4.5C.3 Record Keeping

- a. Record in HTA Book % Ash.
- b. Record Final % in Values Book under HTA %.
- c. Note in Sample Status Book that HTA was performed.

4.5C.4 Suggestions

a. Input for Programmable Furnace

- Initial rate (1st hour) $8^{\circ}\text{C}/\text{min.}$
- Transition temperature 500°C
- Hold time 0 min.
- Final rate (2nd hour) $4^{\circ}\text{C}/\text{min.}$
- Final temperature 750°C
- Final hold time 60 min.

b. Repeatability (%)

- In Lab:
 - No carbonates 0.2
 - Carbonates present 0.3
 - Coals with more than 12% ash
 - with carbonate and pyrite 0.5
- Between Labs:
 - No carbonates 0.3
 - Carbonates present 0.5
 - Coals with more than 12% ash
 - with carbonate and pyrite 1.0

4.5C.5 Safety Procedures

Operating Furnace

- safety glasses/goggles must be worn
- laboratory coat should be worn
- when removing sample, heat resistant gloves must be worn

4.5D Volatile Matter Determination (ASTM D3173, 1983d; D3175, 1983e)

Equipment platinum/fused quartz crucibles
analytical balance, sensitivity at
least 0.0001g
electric muffle furnace
tongs
oven capable of reaching and
maintaining 105°C
hydrofluoric acid, (HF), 10%
Volatile Matter Book

Safety Equipment heat resistant gloves
laboratory apron and coat
safety glasses/goggles
fume hood
heavy rubber gloves
organic vapor respirator

4.5D.1 Preliminary – Moisture Determination (ASTM D3173, 1983d)

- a. Zero check balance.
- b. Weigh out approximately one gram of -60 mesh coal.
- c. Place sample in a preweighed, dried, and cleaned platinum or fused quartz crucible [Note: If using the fused quartz crucible NEVER touch with bare hands, ALWAYS use tongs or lens paper.
- d. Place in heating oven for one hour at 105°C.
- e. Reweigh and determine % moisture using the following equation:
$$\% \text{ Moisture} = \{(A-B)/A\} \times 100$$

where: A = grams of sample used, and
B = grams of sample after heating

4.5D.2 Volatile Matter Determination (ASTM D3175, 1983e)

- a. Heat muffle furnace to 950°C + 20°C [Note: It is best to heat furnace to near 970°C to compensate for rapid heat loss when opening furnace door].
- b. Place crucible with top in place into heated furnace.
- c. After 2-3 min. inspect crucible to verify that the lid is still properly seated (if necessary, reseal the lid).
- d. After heating for exactly 7 min., remove crucible from furnace and without disturbing lid, allow to cool in desiccator.
- e. Zero check analytical balance.
- f. Weigh as soon as cold

- g. Determine % volatile matter as follows:
 $\% \text{ volatile matter} = C - D$
 where: C = % weight loss and
 D = % moisture

4.5.3 Cleaning of Fused Quartz Crucibles

BEFORE USING HF... READ PAGES 36 - 39.

- a. Prepare a solution of 10 % hydrofluoric acid. (Note: Hydrofluoric acid is Extremely Dangerous, use every precaution to assure safety as outlined on page 29. If you have never used this acid before, seek help initially).
2. Rinse crucible in hood in dilute HF to remove carbon deposits.
3. Continue rinsing for 3-5 min. (or until clean) followed by a thorough rinsing with deionized water and complete drying.

4.5.4 Record Keeping

Record volatile matter determination in Volatile Matter Book and in Values Book.

4.5D.5 Suggestions

- a. Washed and rinsed crucibles must be dried as quickly as possible. Water droplets left on the crucible's surface will accumulate atmospheric contaminants that are likely to cause localized devitrification.
- b. Temperature regulation during the procedure is critical; maintain it as exactly as possible.
- c. In-lab repeatability should be as follows:

Anthracite	0.3%
Semianthracite and Bituminous	D.5%
Subbituminous	0.7%
Lignite and Peat	1.0%

4.5D.6 Safety Procedures

- a. Handling Hydrofluoric Acid
 - **WARNING:** Hydrofluoric acid is an extremely hazardous chemical at any concentration or temperature. Physical contact will cause serious burns and inhalation of the vapors will irritate the respiratory system. Use all necessary measures to prevent physical contact with the acid. Consider any accident with hydrofluoric acid a serious accident. -READ PAGES 36 - 39.
 - heavy rubber gloves, safety glasses/goggles, laboratory coat, and rubber apron must be worn at all times
 - organic vapor respirator must be worn
 - must be used under a fume hood
- B. Operating Furnace
 - safety glasses/goggles must be worn
 - laboratory coat and/or apron should be worn
 - when removing samples, asbestos gloves must be worn

4.4E X-Ray Radiography

Equipment	X-ray radiography unit band saw or rock saw masking tape [wide] marker pen X-ray film (e.g. instant positive transparent radiographic 8 x 10 film) instant developing unit 13" x 11" plastic core tube holder 13" PVC (3" diam.) or acrylic (2 1/4" diam.) core tubes
Safety Equipment	safety glasses/goggles ear protection laboratory coat and/or apron radiation badge dust respirator

4.4E.1 Preliminary

- a. If the core is not received in plastic tubes, carefully place it into appropriate sized tubes.
- b. For 2 1/4" core, place into split acrylic tubes. Seal ends with plastic caps and mark orientation on tube (i.e. the "up" arrow).
- c. For 3" core place into a half tube of PVC and tape core into place; note orientation on tape and core.

4.4E.2 X-Ray Radiography

- a. Place closed 2 1/4" core tubes or taped 3" half tubes into plastic core tube holder.
- b. Tape two pieces of instant X-radiograph film together so they measure 15 1/4" x 11 1/2".
- c. Place film in chamber at the appropriate level.
- d. Position plastic tube holder with tubes over X-radiograph film.
- e. Secure all portals to chamber area with lead partition.
- f. Close main door to chamber area and place lead partition around unit.
- g. Set unit's controls (see table below).
- h. Step back from X-ray machine until exposure has been made, remove tubes.
- i. When exposure has been made, remove tubes.
- j. Process film through instant processor according to instructions.

4.4E.3 Core Cutting

- a. For 2 1/4" cores, remove end caps and half of the split tube; tape exposed core and note orientation on the taped portion.
- b. Cut the 2 1/4" taped cores using band saw.
- c. Cut the 3" cores (already taped) using band saw.
- d. Annotate radiograph using observations that can be seen both on the cut surface and the radiograph.
- e. If samples are taken, note sample intervals on radiographs.

4.4E.4 Suggestions

- a. Sequence for turning on X-ray radiograph unit

- 1) turn on plug switch near socket
- 2) turn key to "on" position
- 3) select MA setting
- 4) slowly turn KV knob to desired setting
- 5) to turn off unit, first, slowly turn KV knob to zero
- 6) turn key to "off" position
- 7) turn off plug switch.

- b. X-Ray Exposure Table:

coal ¹ contrast typed	core diameter	KV	MA	exposure time (min.)	source to sample distance (in.)
Med High	3"	60	3	24	40
Med High	2 1/4"	60	3	12	40
Med High	slabbed (1/4")	60	3	2	16
Low	slabbed (1/4")	40	3	12-15	16

¹Refers to amount of contrast between the "lighter" (less dense) and "darker" (more dense) components in the coal.

4.4E.5 Safety Procedures

a. X-Ray Unit

- X-ray badge and ring must be worn at all times
- all portals must be sealed off before turning on X-rays
- lead shield must be in front of door before turning on X-rays
- never open door while X-rays are on
- never stand close to machine while X-rays are on
- have machine periodically checked for radiation leaks

b. Band Saw or Rock Saw

- dust respirator and safety glasses and ear protection must be worn
- laboratory coat and/or apron should be worn
- keep fingers and clothing away from blade and other moving parts
- never force the sample to feed too fast
- ear protection should be worn with the band saw
- make sure the proper blade is being used

4.4F Vitrinite Reflectance (ASTM D2798, 1983f)

Equipment petrographic microscope with reflected light capabilities, capable of 500X magnification, photomultiplier with peak picker, and rotating circular stage
immersion oil, nondrying, refractive index 1.515 – 1.519
sample-leveling press
modeling clay
glass reflectance calibration standards
glass slide (2"X1")

4.4F.1 Setting-up the Microscope

[Numbers refer to figure 4]

- a. Turn on:
 - large digital voltmeter (1)
 - small digital voltmeter (2)
 - lamp power supply (3)
 - high voltage power supply (4)
 - stable power supply (5)
- b. Increase Current to:
 - lamp power supply (to highest setting) (6)
 - stable power supply (to 8.25 amps) (7)
- c. Switch on peak picker (8)
- s. Check object lens type (50X.85p) (9).
- f. Push out half-stop (10).
- g. Push in slide and check that number "8" is displayed (11).
- h. Pull out MPV 2 (12).
- i. Switch condenser out (line up knob with reference line) (13).
- j. Check to see that spectral wedge is out (rotate dial to red dot) (14).
- k. Place switch towards illustrated hand (15).
- l. Place switch down for through-the-eye-piece-viewing during metering (16).
- m. Push red lighted button to engage high voltage power supply (17).
- n. Place glass reflectance calibration standard under scope and pick a reference glass using the coordinates for calibration matrix.




FIGURE 4: Diagram of microscope set-up (refer to text for explanation)





FIGURE 4: (cont'd)

GLASS *	1	2	3	4	5	6
1	<u>0.300</u>	0.508	0.911	1.024	1.351	1.649
2	0.306	<u>0.496</u>	0.906	1.022	1.360	1.654
3	0.304	0.504	<u>0.921</u>	1.007	1.363	1.650
4	0.307	0.508	0.924	<u>1.039</u>	1.386	1.670
5	0.324	0.554	0.928	1.049	<u>1.413</u>	1.703
6	0.307	0.509	0.914	1.028	1.361	<u>1.661</u>
Calibration Values	0.306	0.509	0.923	1.020	1.372	1.656
Stage Coordinates	4x116	4x109	4x100	13x116	13x109	13x100
<div> <div>DATE _____</div> <div>BY _____</div> </div>						

FIGURE 5: Example of a calibration matrix

- o. Engage measuring square (18).
- p. while watching the objective and calibration standard, focus down until the two become very close.
- q. While looking through the eye piece, focus up until measuring square is sharp and defined.
- r. Check volt meters to see that both read zero and adjust dark current if necessary (19).
- s. Engage filter S542-21 (20).
- t. Engage photometer by pushing foot switch--red light should turn on (21).
- u. Let machine warm-up for 20-30 min.

4.4F.2 Centering and Calibrating Instrument

- a. Turn off overhead lights.
- b. After machine has warmed-up, level a sample and focus as before.
- c. To assure sharp focus of measuring square, focus on sample and engage focusing square. Both sample and measuring square should be sharply in focus. If not adjust the focusing square (22).
- d. To center objective,
 - place cross-hair over a small object (e.g. pyrite, micrinite)
 - rotate stage and note path of object
 - center accordingly.
- e. To calibrate photometer,
 - place glass calibration standard under scope and focus as before
 - move to a reference glass using coordinates from matrix
 - clear peak picker by pressing in track button (light will blink) then release
 - rotate stage
 - if voltage does not agree with matrix-calibration value, adjust voltage (23)
 - clear peak piker and rotate stage again; continue this until photometer is calibrated to reference glass
 - move to another reference glass and calibrate as before
 - return to first reference glass and check reading; if into adjustment is necessary, instrument is calibrated; if adjustment is needed, then proceed as before until both glasses are calibrated.

4.4F.2 Measurement of Vitrinite Reflectance

- a. Place leveled sample under microscope and focus as before.
- b. Systematically transects the pellet measuring only the structureless vitrinites (telocollinite) that are scratch-free and pyrite-free.
- c. Whenever a particle of telocollinite is found, center on a scratchless area.
- d. Focus sharply and engage measuring square (mirror housing lever).
- e. Clear peak picker.
- f. Rotate stage (by pressing foot pedal that engages motor).
- g. Record maximum reflectance.
- h. Continue measurements until 25 counts are recorded.
- i. Check calibration.
- j. Measure another 25 particles on this mount; then check calibration.
- k. Measure 50 counts on next pellet (a total of 100 counts for one sample).

4.4F.3 Preparation of Calibration Matrix

- a. Set up microscope according to steps 1 to 20 in A and 1 to 4 of B.
- b. Place glass calibration standard on stage.
- c. Center objective onto the first glass.
- d. Record position coordinates on stage.
- e. Clear peak picker.
- f. Rotate stage.
- g. Recore value—adjust voltage until calculated value of glass No. 1 is read as peak.
- h. Read and record both position and peak of all other glasses without any adjustment of voltage.
- i. Enter each value into row 1 of matrix (figure 5).
- j. Return to coordinates of glass No. 2 and adjust voltage to read calculated value of glass No. 2.
- k. Move stage to coordinates of each glass and measure the maximum peak on each value on next row of matrix (figure 5).

- l. Repeat steps 10 and 11 until matrix is completed.
- m. Throw out values associated with highly inconsistent readings (greater than 0.2% difference) (such as glass 5 in figure 5).
- n. Average all other readings except the adjusted reading (underlined in fig. 5).
- o. The average will be the calibration value to be used at the specified coordinate of each glass.

4.4G Maceral Analysis (ASTM D2799, 1983g; ICCP, 1963)

Equipment petrographic microscope with reflected
light capabilities, 500X with eye
piece with cross hair and mechanical
stage point counter
immersion oil, nondrying, refractive
index 1.515 to 1.519
sample leveling press
modeling clay
2" X 1" glass slides

4.4G.1 Setting-up the Microscope

(Numbers refer to figure 4)

- a. Turn on:
 lamp power supply (3)
 mercury vapor lamp (24).
 (for blue-light [fluorescence] analysis only)
- b. Increase current to lamp power supply (to highest setting) (6).
- c. Ignite mercury vapor lamp (25).
- d. Check object lens type (50X.85p) (9).
- e. Engage half stop (10).
- f. Switch condenser out (line up knob with reference line) (13).
- g. Push on MPV 2 (12).
- h. Level sample using modeling clay and sample leveling press.
- i. Secure mechanical stage on platform.
- j. Place immersion oil on pellet face and place under objective on the mechanical stage.
- k. While watching the objective and pellet, focus down until the two become quite close.
- l. While looking through the eye piece, focus up until the image is sharp and clear.
- m. For blue-light [fluorescence] analysis engage:
 yellow filtration (26)
 blue filtration (27)
 fluorescing light source (28)

4.4G.2 Anthracite and Bituminous Maceral Analyses

1. For anthracite and bituminous coals use the classification sheet displayed in figure 6.
2. Perform 500 point counts using white-light and classification in section A of figure 6.
3. If fluorescing macerals are present, perform 500 counts on the sample pellet using the blue-light classification (section B of figure 6).
4. Finally, perform 500 counts on the etched pellet ("c" pellet) using white-light and the classification for etched macerals (section C of figure 6).
5. Recalculate and combine the analyses and record in section D of figure 7.

4.4G.3 Lignite and Subbituminous Maceral Analyses

- a. For lignite and subbituminous coals use the classification sheet displayed in figure 7.
- b. Perform 500 point counts using white-light and the classification in section A. of figure 7.
- c. Perform 500 counts on the same pellet using the fluorescing light source and the classification in section B of figure 7.
- d. finally, perform 500 counts on the etched pellet ("c" pellet) using white light and classification for etched macerals (section C. of figure 7).
- e. Recalculate and combine the analyses and record in section D of figure 7.

4.4G.4 Safety aspects

Never, ever, look directly at the blue light illumination without the barrier filter in place.

U.S.GEOLOGICAL SURVEY COAL LAB
 *****BITUMINOUS COAL DATA SHEET*****
 A. WHITE LIGHT ANALYSIS

MACERAL	(A1)	(A2)	(A3) Totals (A1+A2)	
Vitrinite----				
Vitrodet.-----				(A4)_____
Sporinite----				(tot. Vitrinite)
Cutinite-----				
Resinite-----				(A5)_____
"Resodet."-----				(tot. Liptinite)
Alginite-----				
Liptodet.-----				
Fusinite -----				
Semifusinite----				
Micrinite-----				(A6)-----
Macrinite-----				(tot. Inert.)
Inertodet.-----				
Sclerotinite----				

D. SUMMARY OF DATA

0.001(B6 - A6) = D1
 0.001(____ - ____ = ____

WHITE LIGHT TOTALS	MACERAL VARIETY	CRYPTO C5 X D1 =	REVISED %	REPT %	
	Telinite	X			
	Corpo in telo	X			Recal.
	Desmocollinite	X			- From
	Detrocollinite	X			Analysis
	Corpo in desmo	X			C
	Pcoll in desmo	X			
	Sporinite	/			
	Cutinite	/			
	Resinite	/			
	"Resodet."	/			From
	Exsudatinite	/			- Analysis
	Fluorinite	/			B
	Bituminite	/			
	Alginite	/			
	Liptodetrinite	/			
	Fusinite	/			
	Semifusinite	/			
	Micrinite	/			From
	Macrinite	/			- Analysis
	Inertodet.	/			A
	Sclerotinite	/			

PELLET ANALYSTS Lta _____ SET _____
 Hta _____
 A. by _____ Density _____ SAMPLE NO. _____
 Pyr S _____
 B. by _____ Tot S _____ DATE _____

FIGURE 6: Maceral classification sheet for high-rank coals

C. USGS DATA SHEET FOR ETCHED PELLETS

MACERAL (crypto)	(C1) by	(C2) by	(C3)Total [1+2]	(C5)=[C3-C4]*100 tot ea/tot vit.
Telinite				
Corpo in telinite				
Desmocollinite				
Detrocollinite				
Corpo in desmo				
Poricoll.in desmo				
all others				XXXXXXXXXXXXX

(C4)Total Vitrinite Counts_____

B. BLUE LIGHT ANALYSIS				
MACERAL	(B1)	(B2)	B3=B1+B2	comments
Sporinite				
Cutinite				
Resinite				
"Resodetrinite"				
Exsudatinite				
Fluorinite				
Bituminite				
Alginite				
Liptodetrinite				
Vit + INERT.			(B6)	

(B4)TOTAL LIPTINITE CTS_____

SAMPLE NO. _____

FIGURE 6: (cont'd)

U.S.GEOLOGICAL SURVEY COAL LAB
 ***** LOW RANK COAL DATA SHEET *****
 A. WHITE LIGHT ANALYSIS

MACERAL	(A1)	(A2)	(A3) Totals (A1+A2)	
Huminite-----				(A4)-
Sporinite-----				(tot. huminite)
Cutinite-----				
Resinite-----				
"Resodet."-----				
Suberinite-----				(A5)-
Alginite-----				(tot Liptinite)
Liptodet.-----				
Fusinite-----				
Semifusinite-----				
Micrinite-----				(A6)-
Macrinite-----				(tot. Inert.)
Inertodet.-----				
Sclerotinite-----				

D. SUMMARY OF DATA

0.001(B6 - A6) = D1
 0.001(____-____) = ____

WHITE LIGHT TOTALS	MACERAL VARIETY	CRYPTO C5 X D1 =	REVISED %	REPT %	
	Humotelinite	X			
	Eucorpo/telo	X			
	Poricorpo/telo	X			Recal.
	Eugelinite	X			From
	Porigelinite	X			- Analysis
	Detrogelinite	X			C
	Eucorpo/gelo	X			
	Poricorpo/gelo	X			
	Sporinite	////			
	Cutinite	////			
	Resinite	////			
	"Resodet."	////			
	Suberinite	////			From
	Exsudatinite	////			- Analysis
	Fluorinite	////			B
	Bituminite	////			
	Alginite	////			
	Liptodetrinite	////			
	Fusinite	////			
	Semifusinite	////			
	Micrinite	////			From
	Macrinite	////			- Analysis
	Inertodet.	////			A
	Sclerotinite	////			

PELLET ANALYSTS Lta _____ SET _____
 Hta _____
 A. by _____ Density _____ SAMPLE NO. _____
 Pyr S _____
 B. by _____ Tot S _____ DATE _____

FIGURE 7: Maceral classification sheet for low-rank coals

C. USGS DATA SHEET FOR ETCHED PELLETS

MACERAL (crypto)	(C1) by	(C2) by	(C3)Total [1+2]	(C5) = [C3/C4]*100 tot ea/tot hum.
Humotelinite				
Eucorpo in telo				
Poricorpo in telo				
Eugelinite				
Porigelinite				
Detrogelinite				
Eucorpo/gelo				
Poricorpo/gelo				
all others				XXXXXXXXXXXX

(C4)Total Huminite Counts_____

B. BLUE LIGHT ANALYSIS

MACERAL	(B1)	(B2)	B3=B1+B2	comments
Sporinite				
Cutinite				
Resinite				
"Resodetrinite"				
Suberinite				
Exsudatinite				
Fluorinite				
Bituminite				
Alginite				
Liptodetrinite				
HUM + INERT.			(B6)	

(B4)TOTAL LIPTINITE CTS_____

SAMPLE NUMBER _____

FIGURE 7: (cont'd)

5 EMERGENCY PROCEDURES AND PRECAUTIONS

5.1 Eye Wash Station

The eye wash station should be positioned in the laboratory free of obstructions. The water should be changed at least once every 6–12 months. Be acquainted with the correct operation of the eye wash station before an emergency.

5.2 Fire Extinguisher

The fire extinguisher should be located free of obstructions within the laboratory. The correct type of extinguisher (CO₂, water air pressure, dry chemical, etc.) should be used depending on type of hazard that exist in the laboratory. Before using dial X 7222, if no answer, 9-911 to report fire then extinguish.

5.3 First Aid

A first aid box consisting of band-aids, antiseptic, bandages, etc., should be placed in an easily accessible place in the laboratory. Any accident should always be reported to a supervisor regardless of how serious.

5.4 Acid Spills

In case of acid spills the emergency number is X 7222 and should be called IMMEDIATELY—DO NOT ATTEMPT TO CLEAN UP YOURSELF.

6 ACKNOWLEDGMENTS

Many people have helped to develop and improve the procedures in this manual. We wish to thank particularly Derek Widmayer, Leslie Ruppert, Peter Swallow, and Philippa Benson who contributed to the refinement of these sample preparation procedures. Also, our thanks to Brenda Pierce who helped in the drafting and reviewing of this manual. Through the years other industry and university personnel have through conversation and their gracious hospitality influenced many of the details in our procedures. These individuals include: Dr. J. J. Renton, Mr. Ralph Gray, Dr. Alan Davis, Dr. William Spackman, Dr. Jack Crelling and Mr. William Grady.

- Allen, T., 1974, Particle size measurement, Chapman and Hall, London, England, 455p.
- American Society for Testing and Materials, 1983a, Preparing coal samples for analysis, ASTM designation D2013: 1983 Annual Book of ASTM Standards, section 5, vol. 05.05, Philadelphia, p. 295-309.
- , 1983b, Standard method of preparing coal and coke, from coal, ASTM designation D3174: 1983 Annual Book of ASTM Standards, section 5, vol. 05.05, Philadelphia, p. 375-379.
- , 1983c, Ash in the analysis sample of coal and coke from coal, ASTM designation D3174: 1983 Annual Book of ASTM Standards, section 5, vol. 05.05, Philadelphia, p. 398-401.
- , 1983d, Moisture in the analysis sample of coal and coke, ASTM designation D3173: 1983 Annual Book of ASTM Standards, section 5, vol. 05.05, Philadelphia, p. 395-397.
- , 1983e, Volatile matter in the analysis sample of coal and coke, ASTM designation D3175: 1983 Annual Book of ASTM Standards, section 5, vol. 05.05. Philadelphia, p. 402-407
- , 1983f, Microscopical determination of the reflectance of the organic components in a polished specimen of coal, ASTM designation D2798: 1983 Annual book of ASTM Standards, section 5, vol. 05.05, Philadelphia, p. 380-383.
- , 1983g, Microscopical determination of volume percent of physical components of coal, ASTM designation D2799: 1983 Annual Book of ASTM Standards, section 5, vol. 05.05, Philadelphia, p. 384-385.
- Barrows. M. H., Cluff, R. M., and Harvey, R. D., 1979, Petrology and maturation of dispersed organic matter in the New Albany Shale Group, in Proceedings 3rd Eastern Gas Shales Symposium, Hilma Barlow, ed., Morgantown Energy Technology Center and Eastern Section AAPG, META/SP-79-6.
- Cole, D. L., and Berry, W. F., 1965, Preparation and polishing of coal and coke for petrographic analysis: Bituminous Coal Research, Inc., Monroeville, Pa., 38p.
- ICCP (International Committee of Coal Petrology), 1963, International handbook of coal petrology: Centre National de al Recherche Scientifique, Paris, France.
- Given, P. H., and Yarab, R. F., 1978, Analysis of the organic substance of coal: problems posed by the presence of mineral matter, in: C. Karr, ed., Analytical Methods of Coal and Coal Products, Academic Press, New York, v. 2, p. 3-41.
- Hawley, G. G., 1977, The condensed chemical dictionary, 9th edit: Van Nostrand Reinhold Company, New York, N.Y., 957p.

Jones, R. H., 1970, A bore core radiography unit, Journal of Physics E., Science Instr., v. 3, p. 363-366.

Laboratory Safety Seminar, 1982, National Hazards Control Institute, Alpha, N.J., 140 p.

Laboratory Safety, 1979, J.T. Baker Chemical Company, Phillipsburg, N.J.

Laboratory Waste Disposal Manual, 1969, Manufacturing Chemists Association, Washington, D.C., 170 p.

National Safety Council, 1983, Hydrofluoric acid exposure: Safety and Health Information Sheet, RSD-83-4, 4p.

Renton, J.J. 1977, Some practical considerations in X-radiography, West Virginia Geol. and Econ. Survey, MERC/SP/77-7, 27p.

Safe Handling of Compressed Gases in Laboratory and Plant, Matheson Gas Products, 13 p.

Stach, E., Machowsky, M-Th., Teichmuller, M., Taylor, G.H., Chandra, D., and Teichmuller, R., 1982, Stach's Textbook of coal petrology: Berlin, Gebruder Borntraeger, 428 p.

Stanton, R.W., 1982, Determination of true density of pulverized coal samples, U.S. Geological Survey Open-File Report 82-835.

Stanton, R.W., Moore, T.A., Ruppert, L.F., and Cecil, C.B., 1983, X-ray radiographic description of coal core as an aid for estimating coal quality 1983 GAS Book of Abstracts with Programs, p. 695.

Stockham, J.D., and Fochtman, E.G., 1977, Particle size analysis: Ann Arbor Science Publishers, Inc., Ann Arbor, Mich., 140 p.

Survey Safety Handbook, 1979, U.S. Geological Survey, Reston, Va.