

U.S. DEPARTMENT OF THE INTERIOR

U.S. GEOLOGICAL SURVEY

MASSPEC: A PC program to control and to process  
data from an automated mass-spectrometer

by

Robert C. Bigelow, Robert B. Vaughn, and Stanley E. Church<sup>1</sup>

Open-File Report 99-161

This report is preliminary and has not been reviewed for conformity with U.S. Geological Survey editorial standards. Any use of trade, product, or firm names is for descriptive purposes only and does not imply endorsement by the U.S. Government.

<sup>1</sup>Denver, Colorado

# **MASSPEC: A PC program to control and to process data from an automated mass-spectrometer**

*by Robert C. Bigelow, Robert B. Vaughn, and Stanley E. Church*

## **I. INTRODUCTION:**

The Geological Survey has used a single-filament, solid-source Mass Spectrometer (MS) for research on isotopic mass ratios in earth materials. The National Bureau of Standards (NBS) designed and built the spectrometer as a single-collector instrument that emphasizes measurement precision for scientific research rather than high-speed data acquisition. A custom-programmed personal computer controls the magnetic field of the spectrometer and acquires relative isotopic abundance data from the MS. A program, called MASSPEC, automates the collection of isotopic data and computes isotopic ratios subject to operator-specified statistical parameters.

The spectrometer came equipped with a Hewlett-Packard HP85 computer/controller that served adequately for many years, but is now obsolete. NBS supplied an HP85 program in a dialect of HP BASIC that one of us, (SC), modified to improve its ability to filter obviously erroneous data, but it was still necessary to import the data into a spreadsheet and reduce it manually to remove the effects of noise and sample instability. Since an MS-DOS Personal Computer (PC) was already employed to run the spreadsheet, we decided to replace the HP85 with a PC that was programmed to run the MS. The objectives of the MASSPEC program were to replace the functionality of the HP85 controller and to provide a higher degree of data reduction and more flexible control over measurement parameters than was previously possible. Some HP-85 BASIC program concepts were retained but most of the program could not be translated for the PC since it used HP-85 concepts and system features that are difficult to duplicate on a PC. This report describes the MASSPEC program.


## **II. MASS SPECTROMETER SYSTEM OVERVIEW:**

A small amount of sample material is placed on a filament which is then mounted in the source of the MS. Heating the filament drives off sample ions which are accelerated through a fixed voltage and magnetically deflected to generate an ion current in a Faraday cup<sup>1</sup>. Adjusting the magnetic field of the MS tunes the machine to measure ions of a given atomic mass and charge, usually a singly ionized isotope. A programmable<sup>2</sup> Varian FR-41 Gaussmeter measures and controls the magnetic field. The ion collector current develops a measurable voltage across a known resistance ( $10^9$ ,  $10^{10}$ , or  $10^{11}$  ohms) at the input of an NBS-designed electrometer. An HP 3456A Digital Voltmeter (DVM) integrates the output voltage from the electrometer and converts it to digital form with 6 ½ digit resolution.

The entire system can be controlled and run manually; however, four devices, the Gaussmeter, the DVM, the electrometer<sup>3</sup> and a printer also can be controlled through GPIB (IEEE-488) interfaces. The PC communicates with those instruments via a National Instruments GPIB-PCIIA board, an IEEE-488 I/O interface. The block diagram (Fig 1) of the MS control system diagrams the relationships among these automated parts of the system. The computer sends GPIB commands to the Gaussmeter to change the magnetic field and then waits a predetermined amount of time to allow the magnetic field to stabilize before taking readings from the DVM. The electrometer converts the ion signal from the MS to a voltage in the range 0 to 10 Volts. The DVM determines the voltage from the electrometer by averaging over a period of either 10 or 100 power line cycles (60 Hz. in N. America, 50 Hz. in Europe) and broadcasts its reading over the GPIB interface. The advantage of this technique is that the timing accuracy for the readings is very precise; However, the computer must re-synchronize

with the DVM each time it accepts readings from it. After a full suite of readings is recorded by the PC, the results are stored on disk and transmitted to a printer.

### III. MEASUREMENT DESIGN:



**Figure1: Mass Spectrometer Layout**


The objective is to determine the ratios of isotopes in the samples with as high a degree of precision as can be achieved. Large numbers of readings are taken and then strenuously filtered to remove questionable readings. The philosophy behind the MASSPEC program is that it is better to drop a good (if doubtful) reading than to risk including bad ones in the computations. However, all readings are preserved so that the investigator can override the program's judgement if he so chooses.

In an ideal MS, all isotopes would be measured simultaneously; however, that is not possible with our present instrument. Instead, the magnetic field is tuned to measure each isotope in turn, and that sequence of measurements is repeated for a predetermined number of times. The data must be interpolated to infer simultaneous data values. Systematic noise, emission irregularities and dropping signal levels are all limiting factors of accuracy in these measurements. The heated filament drives off ionized sample material and the rate of ion production is roughly proportional to the remaining amount of sample. Hence, the signal strength (ions per second) usually decreases gradually over time as the

sample depletes. Sometimes abrupt shifts in signal levels are observed that are attributed to sputtering, changes of sample distribution on the filament and other irregularities. Much of the MASSPEC program is specifically designed to reject spurious or questionable data caused by those irregularities.

#### IV. MEASUREMENT DETAILS:

As illustrated in Figure 2, peak signals from each isotopic species appear at well defined




**Figure 2:** Ion current as a function of magnetic field. Points "r", "s" and "s'" are reference and two sample isotopic peaks, respectively. Background is measured near points "b" on at least one side of each isotopic peak.

magnetic field settings and do not overlap. Both the present and previous programs determine background by averaging measurements at a fixed offset (points b, Figure 2) to each side of the isotope peaks. This technique works very well when isotope peaks are cleanly separated with no overlap or interference because the background changes very slowly with magnetic field. MASSPEC also permits background points to be set individually in Gauss units near each isotopic peak

and can determine background from a single background point per isotope if desired.


A typical determination of isotopic ratios consists of a repetitive series of measurements on each isotope. The fundamental building block of these measurements is the **set**, a fixed number



**Figure 3:** A set is a series of measurements of one isotope, a sequence is one set of each isotope measured serially, beginning with the reference isotope. Only the mean and standard deviation of each set is recorded, subject to filtering


(usually 10) of measurements over a period of several seconds at a single magnetic field setting (figure 3). A **good set** is one that has at least a predefined minimum number of valid data points as determined by filtering, (described below). Conversely, a **bad set**, is one that has too few valid data points after filtration. After a set is filtered, the mean and standard deviation of the set are stored, as are the number of valid data points remaining in the set. Individual data measurements comprising a set

are not retained. Measurements proceed by tuning the mass spectrometer to one isotope after another and making a set of measurements at each tuning. One isotopic species is designated as a reference (r) so the other isotopes (s, s' ...) may be ratioed against it. A **sequence** begins with a reference isotope set, followed by one set of each of the other isotopes (Figure 3). This procedure is repeated until a predetermined number of good sequences (ie containing no bad sets) have been measured. An additional constraint is that the good sequences must all be within "groups". As shown in Figure 4, a **group** consists of an uninterrupted number of good sequences plus a following reference set. The additional reference set is required to allow interpolation when computing ratios. A bad set (e.g.: set  $s_7$  in Figure 4) forces the program to start a new group beginning with a new reference set. Since the number of groups is a measure of the number of bad sets, the program aborts if the group number gets too large.



**Figure 4:** Plot illustrating the relationship between sequences and groups. Each bar represents one set. The weight,  $W$  is a fraction of the period,  $t$ , between the beginnings of successive sequences and is used to compute interpolations. The weight shown is for the first sample ion; the second (and succeeding ) ions each have their own weights. Note that the bad set,  $s_7$ , does not belong to a sequence or group

Both MASSPEC and the HP85 program filter out defective measurements by comparing each measurement to the mean, and rejecting those that differ by more than a certain multiple of the Standard Deviation (SD). This is referred to as Rule 1 in Figure 5. Rule 1 is sufficient if most of the measurements are well behaved with only a few outliers (bad data points) and the signals are well above the noise. However, if the data are quite noisy, the standard deviation itself may be too large to be much help in determining what is or is not acceptable data. To eliminate the effects of very noisy data, especially at small signal levels, MASSPEC also compares the measured Relative Standard Deviation (RSD) to a predefined limit. If the RSD is larger than the limiting value, the measured SD is replaced for filtering purposes by the limiting RSD multiplied by the measured mean value and the data comparison is made as before. The advantage to this scheme, called Rule 2, is that sets with very noisy or abruptly changing measurements are detected, filtered, and rejected if necessary. Figure 5 diagrams the effects of both data filtering rules in terms of the SD and RSD. After filtration, MASSPEC



**Figure 5:** The effect of two rules limiting the allowable values of the variance of the means. Rule 1 limits the variance at large values of means to fixed multiples of SD. At small mean values, Rule 2 limits the allowable variance to no more than a fixed RSD

re-computes the mean average of the remaining data and reports the recalculated mean, the number of valid data points, and the standard deviation. The same filtering system is used at several levels in MASSPEC. It is applied to background sets, isotopic sets, and to ratios of sets. The desired output from MASSPEC is the computed ratio of each sample (ie: non-reference) isotope to the reference isotope along with the relevant statistics. Along with the rest of the data, MASSPEC records the time each set is measured. This makes it possible, when forming ratios, to correct for the fact that sample and reference isotopes are measured at slightly different times. We assume linear interpolation is sufficient to provide the correction. The usual procedure is to interpolate between reference isotope sets with time weighting ( $W$ , in figure 4) to derive a "reference" at the time the sample isotope set was measured. Optionally, time weighted interpolation can be applied to sample isotope sets and ratios formed with an intervening reference set; this is referred to in the program listing as an "inverse sense" ratio. For each sample isotope to reference isotope pair, the mean of all ratios is formed and each ratio is compared to it. Rules 1 and 2 are applied and outlier ratios are thrown out. The mean ratio is re-computed from the remaining ratios along with its standard deviation and reported.

It is possible to do many different experiments using MASSPEC without modifying the program itself. Experiments are defined for MASSPEC by entries in an auxiliary file called ISOPARAM.REC which contains the parameters for each experiment. The parameters include the gauss settings for each isotope peak; their associated background settings; and the filter parameters for each stage of measurement (See Appendix A). ISOPARAM.REC is an ASCII file; any parameter may be changed and new experiments may be added (up to a limit of 20) by editing the file with a text editor.

## V. EXAMPLE:

As an example of how the parameters affect the program, consider the excerpt from a typical ISOPARAM.REC file listed in Appendix A. For convenience, Table 1 reproduces part of ISOPARAM.REC. Suppose we want to measure  $^{206}\text{Pb}$ ,  $^{207}\text{Pb}$  and  $^{208}\text{Pb}$  using  $^{206}\text{Pb}$  as the reference.

The initial entries for those three lead isotopes define the gauss settings for each and where the background measurements are made. MASSPEC takes "NumReads" (10) readings for each set and computes the mean (M) and standard deviation (SD) for each isotope. The program invokes Rule 2 by comparing the SD to a "limiting" SD formed by multiplying IsoMaxDev by M and uses the lesser of those quantities when applying Rule 1. Any reading whose absolute difference from the mean exceeds IsoSigmaTol (1.5) times the SD violates Rule 1 and is rejected (again, see Figure 5).

206	' Iso(1).Mass	INTEGER	Mass Units Iso(1)=reference isotope
6703.0	' Iso(1).Gauss	SINGLE	Gaussmeter setting for isotope peak
1.0	' Iso(1).Range	SINGLE	Voltage range for Chart Recorder
... (similar entries for isotopes 2 and 3)...			
10	' NumReads	INTEGER	Readings per set for each isotope
7	' IsoGoodReads	INTEGER	Min number of readings in a Good Isotope Set
.003	' IsoMaxDev	SINGLE	Maximum allowed Rel. Std. Deviation for an Isotope set (Rule 2)
1.5	' IsoSigmaTol	SINGLE	Tolerance: Multiplies Std. Deviation to form limit for outliers (Rule 1)

**Table 1:** A partial, somewhat reformatted list from Appendix A ISOPARAM.REC of parameters for one isotope.

Once outlier data points have been eliminated, MASSPEC checks to be sure that IsoGoodReads (7) good readings remain to form a valid measurement and if so recalculates the mean, M, with those data points. While it would be possible to continue to re-compute SDS, throw out additional points and compute a new means until they converged, our experience with this program and its predecessor indicates that one cycle of corrections is adequate to yield good statistics when the data are good. Most of the remaining parameters apply to Baselines and Ratios and are applied in the same way. as the Isotope parameters. The Delay parameter specifies the time for the magnetic field to settle whenever the MS is re-tuned; usually before beginning each set of measurements.

## VI. PROGRAM OVERVIEW:

The MASSPEC program is written in Microsoft Quick BASIC version 4.5. Initially, we had hoped to translate the NBS HP85 BASIC program as a basis for MASSPEC, but the required changes were too great for that to be a productive approach. Nevertheless, many of the concepts of the NBS program were retained and incorporated in MASSPEC. The MASSPEC Program listing is in Appendix B (q.v.). The program is available on diskette.

The main part of the program is divided into 10 sections with 12 subprograms and functions.

- A. The Declaration section has a \$INCLUDE 'QBDECL.BAS' statement that invokes the library for National Instruments IEEE488 PC card. This statement is required so that the program can invoke functions supplied by National Instruments. In the setup phase, the program initializes variables and types. In the 'Open Files' section, the program enters a large loop at the 'SetupFiles:' label that terminates in a 'SELECT CASE' statement near the end of the main program. The program asks for a sample name which becomes the name for 3 files; the backup to the report file <Samplename>.Bnn, a log file

- <Samplename>.LOG and a <Samplename>.Rnn file that is intended to be used with a spreadsheet, although this has not been tested. The LOG file is intended to record everything that happens in case the program aborts or for some other reason is unable to finish. The .Rnn file is written to a floppy disk so that it will be an automatic backup. If the same sample is re-measured, the logfile is appended and the other two files are created with incrementing extensions (Rnn, Bnn) so that data is not overwritten. This creates a lot of housekeeping but it means that, in general, data is not lost. After several optional questions about measurement parameters and documentation, the operator is given a chance to revise his entries and the program enters the "GET Parameters" section.
- B. The ISOPARAMS.REC file is scanned twice. The first time, every experiment definition header is read in and served up in a menu. Once the operator chooses the desired definition the program re-reads ISOPARAMS.REC and imports the chosen parameters from it.
  - C. In the Change Parameters section, two menu options allow some of the parameters to be overridden. The RANGE question scales the output to the chart recorder but does not affect the data seen by the computer. The NUMBER OF CYCLES refers to the integration time of the DVM. 100 is the usual period and should give the highest accuracy but the measurements complete much faster at 10 cycle integration times. Sometimes the sample on the filament depletes quickly or the data are extremely regular and it is necessary or desirable to integrate for only 10 line cycles. The program now goes fully automatic and does not require further operator intervention until the end.
  - D. The Baseline is measured at one point or preferably two bracketing points near the isotopic peak. This determines systematic errors, particularly the positive offset (usually about 500 mV) that is deliberately introduced to avoid non-linearities in the MS system electronics. The program will terminate if there are more than BaseMaxBadSets, because that would indicate an extraordinarily unstable system. Good baseline information is summarized to the screen.
  - E. As the program moves into the measurement phase several variables are set up that keep their meanings throughout the program. The basic measurement is a set of readings that are averaged and checked for bad points as explained previously. The variable I indicates the isotope. The first isotope (I=1) is always the reference isotope. The variable J counts how many usable sets of measurement have been collected irrespective of which isotope. The variable N counts how many complete cycles through sequence of I values have taken place. N is the sequence number which indicates how much data will be available to compute ratios later. Sequences always start with a reference set. The group number, G, groups valid sequences together and indicates the number of restarts because of defective data as previously explained. If there are rules violations, the group number increments and if it becomes too large the program prints an error message and quits. Some debug statements, mostly PRINT USING, have been left in the program. They can be invoked by uncommenting them.
  - F. Once sufficient good data (as defined by N) has been collected the program enters the



computation stage, Analyze Ratios. First, reference sets are interpolated then ratios are formed between the sample and reference data. In the normal ratio procedure, bracketing reference data are interpolated to derive a value that represents the reference isotope signal at the time a sample isotope is measured. In the optional "reverse sense" ratioing method, bracketing sample data are interpolated to infer values at the time a reference signal is measured. Once all the ratios are formed, they are averaged and checked against their standard deviation in the same way that sets were. Ratios have their own limiting parameters in ISOPARAMS.REC which often are tighter than the parameters for the individual isotopes.

- G. Once the computations are finished, the data are written out to the ASCII delimited <Samplename>.Rnn file for possible use in a spreadsheet, are printed out, and the <Samplename>.Bnn files are written out. Both files are closed. The Data are printed to the screen and to the printer in the Report section, using routines from the Display GoSubs section that follows the Loop section
- H. In the Loop Logic section of the main program, the operator is given choices as to whether to repeat the measurement exactly, change sample IDs only and keep the parameters, change the parameters for the same sample, change everything or quit. Based on his decision, the program branches as indicated in the CASE statement. If the measurement is partly or fully repeated, some parameters are retained to simplify the operators task.
- I. The Display Gosubs and Error Messages sections follow the main body of code. The Display GOSUBS generate formatted output to the screen, the printer, and the Log file.

### III. SUBPROGRAMS AND FUNCTIONS

There are both logical and hardware-specific subprograms for the voltmeter and for the gaussmeter. This arrangement of subprograms has two major advantages. First, it insures that a minimum of reprogramming will need to be done should a different IEEE PC interface board, Digital VoltMeter (DVM) or Gaussmeter be substituted for the present units in the future. Second, it allowed us to substitute dummy routines for the IEEE-specific routines and test the program on a separate PC, which avoided tying up the MS system. The Measure Subprogram performs statistical calculations on sets of measurements whereas the GetDVM subroutine has the IEEE-specific calls for reading the DVM. Likewise the Move Subprogram has the logic and the time delays for the Gaussmeter but the OutputGauss Subprogram has the IEEE specific calls. The Move Subprogram produces a rotating or whirling line display that is intended to reassure the operator that the system really is functioning while it moves to a new Gauss setting; however, it does not work well with some PCs. OutputGauss has an algorithm for converting 5 digit Gaussmeter settings to the binary coded decimal form required by the Gaussmeter along with the range data coded in for the electrometer. The forms of GetDVM and OutputGauss were based in large part on examples in National Instruments documentation. The conversion of data to Binary Coded Decimal plus Range control follows NBS documentation supplied with the electrometer and the NBS-modified Varian Gaussmeter. The Range data is sent to the Gaussmeter which then passes it to the electrometer.

The Menu subprogram provides a moderately flexible means of presenting up to 20 multiple choice menu items in a vertical format. Keystrokes on the grey vertical arrow keys (or on their numeric

keypad equivalents with the "Num Lock" key off) highlight the selections; hitting "Enter" selects the choice. The menu choices and their consequences are in arrays which are usually at the beginning of the program. The Interpol8 Subprogram defines the interpolation routines. The weighting function is DelT (for Delta-T), functionally the same as W, above. Ratio Statistics Subprogram actually computes the ratio statistics. It is possible to have near cancellation in the SigSq# variable and it can even compute to a negative. If it is less than  $10^{-15}$ , it taken as zero. There are several very small Functions or Subprograms that are mostly self-explanatory. The Getxxx functions get different variable types from ISOPARAM.REC

## VIII. ENDNOTES AND REFERENCES

1. The theory and operation of mass spectrometers is covered in, for example, Hamilton, E.I. and Aherns, L.H., 1965, *Applied Geochemistry*: London and New York, Academic Press, Solid sources are discussed on pp 17-23. Lead isotopes measurements are discussed in Chapter 9, pp 167-203
2. The Varian Gaussmeter and the electrometer were modified to be controlled over the GPIB bus at by Schiedler at NBS (unpublished report). The NBS report refers to the electrometer as a "parametric amplifier"( PA). Our present instrument has newer electronics but otherwise is similar to one described in Shields, W.R., ed., 1966, NBS Technical Note 277; *Analytical Mass Spectrometry Section: Instrumentation and Procedures for Isotopic analysis*, 25 July 1966: [Washington, D. C.], National Bureau of Standards.
3. The electrometer provides a fixed range output to drive the DVM and a range-switched output to drive a chart recorder. The range settings are derived from GPIB signals sent to the Gaussmeter. Since the output to the DVM is fixed, range settings should not affect the data seen by the PC. The NBS program indicates such an effect may exist; if so, it has been too small for us to observe and we have ignored it.

## APPENDIX A: ISOPARAM.REC

Each ISOPARAM.REC file must start with a file date header as shown. MASSPEC checks the revision information in ISOPARAM.REC to insure that it has the correct format and sequences of parameters. Whenever the form of ISOPARAM.REC changes, a new header must be defined and code in MASSPEC changed. This insures that incompatible ISOPARAM.REC and MASSPEC versions cannot be used together. The parameters are listed in groups that define a measurement sequence in MASSPEC. Each parameter appears on its own line and is contained in the first 12 characters of the line; the remainder of the line is comment. All entries in a group of parameters must be defined. Each experiment definition starts with a header, such as "PB206-208", that identifies it as a menu choice in the program. Definitions must be separated from each other by at least one blank line.

The header is followed by the parameters for each isotope. For each isotope, the mass in AMU is followed by the Gauss setting for the isotopic peak and the Gauss settings at which the Baseline(s) near it will be measured. A negative BaseGauss2 reading indicates it is not taken. After each isotope is defined, definitions similar to Table 1 appear for the Isotopes, the Baselines, and the Ratios. The last two entries are the maximum allowable number of bad baseline sets and the time delay between sets to let the electrometer (PA) settle the magnetic field settle after it is retuned. If more than 3 (typically) baseline readings in a row are noisy or bad, the measurement should be aborted. The optimum time delay is determined by trial and error as a compromise between affecting the first readings in a set and slowing the measurements so much that drift becomes unacceptable

### Sample ISOPARAM.REC file

```
' Last revised 11 Mar 93 for Rev 6.7
' The first non-blank line must be the revision line (11Mar93)
' add BaseGauss1 & 2; remove GaussOffset from records (11Mar93)
' This file holds parameters for the Masspec program
' up to 20 different parameter sets can be defined
' Sets must be separated by at least one line that is blank
' in the first 12 characters (these lines qualify).
' Pb 206, 207, 208 default sequence

PB 206-208
3
206
6703.0
1.0
6711.0
6695.0
207
6719.3
1.0
6711.3
6727.3
208
6736.0
3.0
6728.0
6744.0

' IsoName  STRING * 10      Name of Sequence of Isotopes
' NumIsos  INTEGER          No. of Isotopes
' Iso(1).Mass  INTEGER      Mass Units  Iso(1)=reference isotope
' Iso(1).Gauss  SINGLE      Gaussmeter setting for isotope peak
' Iso(1).Range  SINGLE      Voltage range for Chart Recorder
' BaseGauss1  SINGLE      Baseline Gauss; 1 required
' BaseGauss2  SINGLE      2nd optional; ignore if 0 or negative.
' Iso(2).Mass  INTEGER      Mass Units
' Iso(2).Gauss  SINGLE      Gauss
' Iso(2).Range  SINGLE      Volts
' BaseGauss1  SINGLE      Baseline Gauss reading, one required
' BaseGauss2  SINGLE      2nd optional; ignore if 0 or negative.
' Iso(3).Mass  INTEGER      Mass Unit
' Iso(3).Gauss  SINGLE      Gauss
' Iso(3).Range  SINGLE      Volts
' BaseGauss1  SINGLE      Baseline Gauss reading, one required
' BaseGauss2  SINGLE      2nd optional; ignore if 0 or negative.
```

10	' NumReads INTEGER	Isotope; No. of readings per set
7	' IsoGoodReads INTEGER	Min No of reads in a Good Isotope Set
10	' BaseNumReads INTEGER	No. of sets read for each isotope
7	' BaseGoodReads INTEGER	Min no.reads in Good Baseline set
10	' NumRatios INTEGER	Minimum number of valid ratios
.003	' RatioMaxDev SINGLE	Max.Rel.Std.Dev for Ratios
1.5	' RatioSigmaTol SINGLE	Tolerance * Ratio Std.Dev=Reject
.003	' IsoMaxDev SINGLE	Max.Rel.Std.Dev for isotopes
1.5	' IsoSigmaTol SINGLE	Tolerance * Std.Dev.= Reject level
.003	' BaseMaxDev SINGLE	Max Rel.Std.Dev. for BaseLines
1.5	' BaseSigmaTol SINGLE	Tolerance * Base.Std.Dev.=Rej.Level
3	' BaseMaxBadSets INTEGER	Maximum allowed bad baseline sets
2	' Delay INTEGER	Min. time to let PA settle

Mandatory: at least one blank line between definitions

	' Pb 204, 206 default sequence	
PB 204-206	' IsoName STRING * 10	Name of Sequence of Isotopes
2	' NumIsos INTEGER	No. of Isotopes
206	' Iso(1).Mass INTEGER	MU Iso(1) is reference isotope
6703.0	' Iso(1).Gauss SINGLE	Gaussmeter reading for isotope
1	' Iso(1).Range SINGLE	Volts
6711.0	' BaseGauss1 SINGLE	Baseline Gauss reading, one required
6695.0	' BaseGauss2 SINGLE	2nd optional; ignore if 0 or negative.
204	' Iso(2).Mass INTEGER	Mass Units
6670.0	' Iso(2).Gauss SINGLE	Gauss
0.1	' Iso(2).Range SINGLE	Volts
6662.0	' BaseGauss1 SINGLE	Baseline Gauss reading, one required
6678.0	' BaseGauss2 SINGLE	2nd optional; ignore if 0 or negative.
10	' NumReads INTEGER	readings per set per isotope
7	' IsoGoodReads INTEGER	Min reads in a Good Isotope Set
10	' BaseNumReads INTEGER	No. of sets read for each isotope
7	' BaseGoodReads INTEGER	Min no.reads in Good Baseline set
10	' NumRatios INTEGER	Minimum number of valid ratios
.003	' RatioMaxDev SINGLE	Max.Rel.Std.Dev for Ratios
1.5	' RatioSigmaTol SINGLE	Tolerance * Ratio Std.Dev=Reject
.003	' IsoMaxDev SINGLE	Max.Rel.Std.Dev for isotopes
1.5	' IsoSigmaTol SINGLE	Tolerance * Std.Dev.= Reject level
.003	' BaseMaxDev SINGLE	Max Rel.Std.Dev. for BaseLines
1.5	' BaseSigmaTol SINGLE	Tolerance * Base.Std.Dev.=Rej.Level
3	' BaseMaxBadSets INTEGER	Maximum allowed bad baseline sets
2	' Delay INTEGER	Min. time to let PA settle

up to 18 more definitions may be included; formats must conform to foregoing examples

## APPENDIX B: PROGRAM 'MASSPEC' LISTING

Note: a '\' at the end of a line indicates continuation onto the next line for listing purposes. Continuation lines are not supported by Quick Basic 4.5. To compile, use IEEEv9x5.BAS or edit this program to remove continuation lines and this heading. Required files and programs for compilation are: Microsoft Quick Basic 4.5, QBDECL.BAS from National Instruments Inc., IEEEv9x5, and ISOPARAM.REC. The PC must have a National Instruments IEEE interface installed that is supported by QBDECL.BAS. Possibly useful Debug statements have been left in place but are commented out by default.

```
' Program MASSPEC    R.C.Bigelow 23Mar99 Additional Comments
' Last major revision is IEEEv9x5.BAS 8Mar96
' Make Exe from the compiler menu will not recognize the need for V or X compiler
' options so compile will fail. Run compiler from command line with options:
'   C:\QB45> BC IEEEv9x5.BAS /E/O/V/C:512/X
DEFINT A-Z
' $INCLUDE: 'QBDECL.BAS'           'National Instruments IEEE488 software
' $DYNAMIC
DECLARE FUNCTION GetAns$ (Prompt$, Default$)
DECLARE FUNCTION GetStr12$ ()
DECLARE FUNCTION GetSingle! ()
DECLARE FUNCTION GetInt% ()
DECLARE FUNCTION Interpol8# (M() AS ANY, SampIndex%, RefIndex%, MeasIndex%)
DECLARE FUNCTION MaxOf# (A#, B#)
DECLARE FUNCTION YNans% (message$)
DECLARE SUB GetInfo (OpName$, FilAmps$, FilVolts$, SubmitName$, Comment$)
DECLARE SUB GetSampName (SampName$)
DECLARE SUB MoveTo (Gauss!, Range%, Delay)
DECLARE SUB Measure (SigmaTol!, MaxStDev!, Reads, Mean#, Sigma#, Time!, GoodReads)
DECLARE SUB RatioStats (Ratio() AS DOUBLE, Results() AS ANY, I%, J%, T%)
DECLARE SUB ReadDVM (VmRead#)
DECLARE SUB Statistics (VALUE AS DOUBLE, Sum AS DOUBLE, SqSum AS DOUBLE, Number%)
DECLARE SUB OutputGauss (Gauss!, Range%)
DECLARE SUB VertMenu (Header$, Menu$(), Size, RowNum, Choice)
'-----
'Comments
' Program IEEE.BAS gathers data from an NBS designed Mass spectrometer.
' Program controls a Varian Gaussmeter (Gauss), an NBS designed
' Parametric Amplifier (PA) attached to the Gaussmeter, and
' a Hewlett Packard High Model 3456A Digital Voltmeter (DVM) via a
' National Instruments IEEE-488 controller in an IBM-type 386 PC.
' DEFINITIONS
' Sequence: an ordered list of isotopes to be measured, beginning with the
'           Reference Isotope followed by one or more Sample Isotopes.
' Set:      a series of (usually 10) consecutive measurements of a single
'           isotope. Mean and standard deviation (Sigma) are reported;
'           individual measurements are not.
' Tolerance Measurements are rejected if their difference from the mean is
'           greater than the Tolerance (usually 1.5 to 3) times Sigma
' MaxDev    Maximum value a base or isotope Sigma can have. Helps detect
'           bad sets( see below)
' For Isotope Measurements:
' Good Set  A Set with no less than 'IsoGoodReads' measurements (usually 7)
'           and a Sigma less than < 'IsoMaxDev'. The Mean and Sigma for all
'           Good Sets are reported and used to compute Ratios if possible.
' Bad Set:  Violates Good Set rules and causes a break in measurements.
'           Increments the 'Group' number for subsequent measurements and
'           forces immediate re-measure of Ref. Isotope Set. Bad Set
'           data is reported but is not used to compute ratios.
' Group:    A collection of good data sets, unbroken by a bad data set.
```

```

'           Begins with a Ref. Isotope set. Within a group, only sets between
'           the initial and final Ref. Isotope set are used to construct
'           Ratios. All sets in a Group have the same 'Group'number.
' Ratios:    For each Sample isotope, an inferred Reference Isotope value is
'           determined by interpolation between bracketing Ref. Isotope
'           measurements. The interpolation is time-weighted to yield
'           an apparent value for the Ref. Isotope at the time the Sample
'           Isotope was measured. The ratio is the mean value of the Sample
'           Isotope set divided by the apparent Reference Isotope value.
'           Constructed for each sequence within a group for which there is
'           a Good Set of initial and final reference isotope measurements.
' For Baseline measurements:
' Good Set  A Set with no less than 'BaseGoodReads' measurements (usually 7)
'           and a Rel. Std. Deviation < 'BaseMaxDev'. The mean and standard
'           deviation for all Good Sets are reported and the average of
'           'BaseNumReads' (usually 10) mean values becomes the Iso.Baseline
'           parameter for each isotope
' Bad Set:  Violates Good Set rules and increments 'BadReadFlag'; which, if
'           it exceeds BaseMaxBadSets (usually 3), aborts the measurement.
' End Comments
' -----
' Setup

CONST FALSE = 0, TRUE = NOT FALSE      'TRUE = -1 in Quick BASIC 4.0 & up
TYPE IsoParams
    Mass      AS INTEGER                ' Isotopic Mass, from Parameters file
    Gauss     AS SINGLE                 ' Gaussmeter reading for above
    Range     AS SINGLE                 ' PA range from Params. File; modifiable
    BaseGauss1 AS SINGLE                ' Baseline Gaussmeter reading, one required
    BaseGauss2 AS SINGLE                ' 2nd optional; ignore if 0 or negative.
END TYPE

TYPE MeasRecord      ' Measured data
    Mass      AS INTEGER                ' As in IsoParams; links the two arrays
    Avrg      AS DOUBLE                 ' Av. Volt. from 1 set of (usually) 10 samples
    StDv      AS DOUBLE                 ' Standard Deviation for Avrg
    Tick      AS SINGLE                 ' Time (Sec.) set was measured
    Valid     AS INTEGER                ' Number of valid readings
END TYPE

' Subscripts are Sample #, Isotope #, and ratio type
REDIM Ratio(30, 3, 2) AS DOUBLE

TYPE ResultRecord
    SumRatio   AS DOUBLE                ' Sum of ratios
    SqSumRatio AS DOUBLE                ' Squares of sums of ratios
    NumRatios  AS INTEGER                ' Counts all ratios of a given type
    RatioAv    AS DOUBLE                ' Average of Ratios
    RatioDev   AS DOUBLE                ' Std.Dev. of ratios
    GoodRatios AS INTEGER                ' Number of Ratios within Tolerance
END TYPE

' The Results Subscripts are Isotope # and Ratio Type where Ratio Type is:
' 0 to indicate combined results for both kind of ratios
' 1 for "normal" Ratio of Sample vs time-weighted, interpolated Reference
' 2 for "reverse sense" Ratio of weighted, interpolated Sample vs. Reference
REDIM Results(3, 2) AS ResultRecord

```

```

'Entries in the IsoParams.Rec (Params) file are as follows:      ---
DIM IsoName          AS STRING * 10 ' Name of this Sequence of Isotopes
DIM NumIsos          AS INTEGER    ' No. of Isotopes,
REDIM Iso(3)         AS IsoParams  ' Holds params for each isotope; REDIMed
                                   ' later to actual size. Iso(1) is usual-
                                   ' ly the reference isotope

DIM NumReads         AS INTEGER    ' Number of readings per set per isotope
DIM IsoGoodReads     AS INTEGER    ' Min no of reads in a Good Isotope Set
DIM BaseNumReads     AS INTEGER    ' Number of Base sets to read per isotope
DIM BaseGoodReads    AS INTEGER    ' Min no of reads in a Good Baseline set
DIM NumRatios        AS INTEGER    ' Minimum number of valid ratios
DIM RatioMaxDev      AS SINGLE     ' Max. Relative Std.Dev for Ratios
DIM RatioSigmaTol    AS SINGLE     ' Tolerance for Ratio Std.Dev
DIM IsoMaxDev        AS SINGLE     ' Max. Permissible Std.Dev
DIM IsoSigmaTol      AS SINGLE     ' Tolerance for Isotope Std. Dev.
DIM BaseMaxDev       AS SINGLE     ' Max Rel. Stand. Dev. in BaseLine set
DIM BaseSigmaTol     AS SINGLE     ' Tolerance for Base Std. Dev.
DIM BaseMaxBadSets   AS INTEGER    ' Maximum allowed bad baseline sets
DIM Delay            AS INTEGER    ' Min. time to let PA settle

SampName$ = "": FilAmps$ = "": FilVolts$ = "
OpName$ = "": SubmitName$ = "

' Scratch variables for measurement calls and computing ratios
DIM Mean AS DOUBLE, Sigma AS DOUBLE, Sum AS DOUBLE, SqSum AS DOUBLE
DIM Time AS SINGLE

'RefIso is the index of the Reference Isotope (conventionally 1)
RefIso = 1

DIM RangeMenu(7) AS STRING
RangeMenu(0) = ".003 Volts": RangeMenu(1) = "0.01 Volts"
RangeMenu(2) = ".03 Volts": RangeMenu(3) = "0.1 Volts"
RangeMenu(4) = "0.3 Volts": RangeMenu(5) = "1.0 Volts"
RangeMenu(6) = "3.0 Volts": RangeMenu(7) = "10 Volts"

DIM RevSenseMenu(0 TO 1) AS STRING
RevSenseMenu(0) = "Omit Reverse Sense Ratios"
RevSenseMenu(1) = "Include Reverse Sense Ratios"

DIM IntegrateMenu(0 TO 1) AS STRING
IntegrateMenu(0) = "Integrate for 10 Line Cycles"
IntegrateMenu(1) = "Integrate for 100 Line Cycles"

RepeatHeader$ = "Choose one of the following:"
DIM RepeatMenu(4) AS STRING
RepeatMenu(0) = "New Sample, New parameters"
RepeatMenu(1) = "New Sample, Same parameters"
RepeatMenu(2) = "Same Sample, New Paramaters"
RepeatMenu(3) = "Redo this measurement (no changes)"
RepeatMenu(4) = "Quit"

REDIM RatioType$(3)
RatioType$(0) = "Combined Data"
RatioType$(1) = "Samp/Interp-Ref"
RatioType$(2) = "Interp-Samp/Ref"
RatioType$(3) = "All types"
COLOR 7, 1
'End Setup
'-----

```



```

'-----
'Open Files
GetInitialData:
CLS
CALL GetSampName(SampName$)
CALL GetInfo(OpName$, FilAmps$, FilVolts$, SubmitName$, Comment$)
SameParams = FALSE
SetupFiles:
COLOR 7, 1: CLS : TickMark$ = TIME$
CLOSE
OPEN "Cons:" FOR OUTPUT AS #1
OPEN "LPT1:" FOR OUTPUT AS #2 'for parallel port printer
' OPEN "LPT2" FOR OUTPUT AS #2 'for HP-IB printer; Note no colon

' Trick to make Rext$ display two digits (eg: .R00). start at 100 then use
' RIGHT$(STR$(..),2)) to lop off the 1; always leaves 2 digits in the string.
FOR Rnn = 100 TO 199
    CLOSE #5
    Rext$ = ".R" + RIGHT$(STR$(Rnn), 2) 'part of trick.
    RecFile$ = "C:" + SampName$ + Rext$ 'File #5 ASCII, Comma delimit file
    OPEN RecFile$ FOR APPEND AS #5
    IF LOF(5) = 0 THEN EXIT FOR ' preserves Rnn
NEXT Rnn
IF Rnn > 199 THEN 'This means something is probably very wrong!
    PRINT "There are now 99 "; SampName$; ".Rnn files! This program"
    PRINT "cannot make any more Record files with this Sample ID."
    INPUT "(hit 'ENTER'; Program will stop)", X$
    CLOSE : STOP
END IF

OpenDataFile:
PRINT " The Record file for this sample is: "; RecFile$: PRINT
PRINT " Please be sure a formatted diskette with room for data is"
INPUT " in Drive B: then hit <Enter>.", X$: PRINT : PRINT
CLOSE #4
Bnn = Rnn
Bext$ = ".B" + RIGHT$(STR$(Bnn), 2)
Datafile$ = "B:\\" + SampName$ + Bext$ 'File #4 Backup for Record file
ON ERROR GOTO BadOpenDataFile
OPEN Datafile$ FOR APPEND AS #4
ON ERROR GOTO 0
IF LOF(4) > 0 THEN 'This means something is probably very wrong!
    CLOSE : CLS : KILL RecFile$: COLOR 14, 4
    PRINT "There is already a file named "; Datafile$; " but it does"
    PRINT "not correspond to the new Record file, "; RecFile$
    PRINT "Please rename "; Datafile$; ", move it to other storage "
    INPUT "and/or replace the Diskette.(hit 'ENTER'; Program will stop)", X$
    COLOR 7, 1
    STOP
END IF

PRINT " The Data file for this sample is: "; Datafile$: PRINT
CLOSE #3
Logfile$ = "C:" + SampName$ + ".LOG" 'File #3 records everything
OPEN Logfile$ FOR APPEND AS #3
PRINT #3, "-----Begin Log For "; SampName$; " -----"
PRINT " The Log file for this sample is: "; Logfile$
IF LOF(3) > 0 THEN
    PRINT " That Logfile already exists. Your current session will be appended to it"
END IF: PRINT

```

```

PRINT
PRINT "      Sample Name: "; SampName$; TAB(36);
PRINT " Time: "; TickMark$; TAB(52);
PRINT "Date: "; DATE$
PRINT "  Filament Volts: "; FilVolts$; TAB(40);
PRINT "  Filament Amps: "; FilAmps$
PRINT " Operators Name: "; OpName$; TAB(40);
PRINT "Submitter's Name: "; SubmitName$
PRINT "      Comments: "; Comment$

LOCATE 24, 5
IF NOT YNans("      Is this all correct?") THEN GOTO GetInitialData

'End Open files
'-----
'GET Parameters

IF SameParams = TRUE THEN GOTO MeasureBase      'Keep same parameters

ON ERROR GOTO ParamFileErr
OPEN "IsoParam.Rec" FOR INPUT AS #6
ON ERROR GOTO 0      'turn off error routine
REDIM NameList(20) AS STRING
'look for string
IsoParamRev$ = "Last revised 11 Mar 93 for Rev 6.7"

DO 'Find first non-blank line
    LINE INPUT #6, Line1$
LOOP UNTIL LEN(Line1$) > 0
IF INSTR(Line1$, IsoParamRev$) = 0 THEN
    CLOSE : CLS : COLOR 14, 4
    PRINT "Wrong revision of ISOPARAMS.REC. The revision is:"
    PRINT MID$(Line1$, 14)
    PRINT "it should be:"
    PRINT IsoParamRev$
    INPUT "(hit 'ENTER'; Program will stop)", X$
    COLOR 7, 1: STOP
END IF

I = 0: Flag = FALSE
DO
    NameList(I) = GetStr12$
    IF (NameList(I) = "") THEN
        Flag = TRUE      ' Blanks separate Param records
    ELSEIF (Flag = TRUE) THEN
        I = I + 1: Flag = FALSE
    ELSE
        Flag = FALSE
    END IF
LOOP UNTIL EOF(6)
NameList(I) = ""      ' Last one is invalid
MaxI = I - 1      ' Subscript for last Valid NameList entry
CLS
Choice = 0
Header$ = "Select the desired Isotope Range"
CALL VertMenu(Header$, NameList(), MaxI, 2, Choice)
CLOSE #6

NewParam:
OPEN "IsoParam.Rec" FOR INPUT AS #6      ' Again
DO      ' skip unwanted entries

```

```

LOOP UNTIL GetStr12$ = NameList(Choice)
IsoName = NameList(Choice)
NumIsos = GetInt
REDIM Iso(NumIsos) AS IsoParams
FOR I = 1 TO NumIsos
Iso(I).Mass = GetInt: Iso(I).Gauss = GetSingle: Iso(I).Range = GetSingle
Iso(I).BaseGauss1 = GetSingle: Iso(I).BaseGauss2 = GetSingle:
NEXT I
NumReads = GetInt
IsoGoodReads = GetInt
BaseNumReads = GetInt
BaseGoodReads = GetInt
NumRatios = GetInt
RatioMaxDev = GetSingle
RatioSigmaTol = GetSingle
IsoMaxDev = GetSingle
IsoSigmaTol = GetSingle
BaseMaxDev = GetSingle
BaseSigmaTol = GetSingle
BaseMaxBadSets = GetInt
Delay = GetInt
'End Get Parameters
'-----
'Change Parameters
ChangeParams:
REDIM RangeNo(NumIsos) AS INTEGER
FOR I = 1 TO NumIsos
RangeNo(I) = 0
ON ERROR GOTO fini
DO UNTIL (ABS(VAL(RangeMenu(RangeNo(I))) - Iso(I).Range) < .00001)
RangeNo(I) = RangeNo(I) + 1
IF RangeNo(I) > 7 THEN ERROR 248
LOOP
ON ERROR GOTO 0
RowNum = 2
COLOR 7, 1: CLS
Header$ = "Select Parametric Amplifier Range for Mass "
Header$ = Header$ + STR$(Iso(I).Mass) + ":"
CALL VertMenu(Header$, RangeMenu(), 7, RowNum, RangeNo(I))
Iso(I).Range = VAL(RangeMenu(RangeNo(I)))
NEXT I
CALL VertMenu("", RevSenseMenu(), 1, 12, RevSense)
Choice = 0
CALL VertMenu("", IntegrateMenu(), 1, 16, Choice)
TMax = 1 + RevSense ' =1 for normal ratios, =2 to include RevSense
IntegTime = 10 ^ (Choice + 1) 'returns 10 for 10 Hz and 100 for 100 Hz
'End Change Parameters
'-----
'Measure Baseline
MeasureBase:
ErrFlag = FALSE 'Set for too many bad readings (BadReadFlag or G too large)
P = 3
GOSUB DisplayHeader

COLOR 7, 1: CLS
CALL ReadDVM(-100#)
' Debug --- CALL ReadDVM(-10#)
Time0! = TIMER

REDIM Baseline(NumIsos, 2) AS MeasRecord
REDIM NumBases(NumIsos) AS INTEGER

```

```

PRINT
PRINT "      Base Sigma Tol: "; BaseSigmaTol; "      Base Max Dev: "; BaseMaxDev
PRINT "Base Gauss  Mean      Sigma      MaxRejLev.  RejLevel  Good Readings"
FOR I = 1 TO NumIsos      'Isotope Index as usual
' (Iso(I).BaseGauss2 > 0) Evaluates as TRUE (-1) or FALSE (0)
  NumBases(I) = 1 - (Iso(I).BaseGauss2 > 0)
  SqSum# = 0
  FOR L = 1 TO NumBases(I)
    SELECT CASE L
      CASE 1
        BaseGauss! = Iso(I).BaseGauss1      ' required baseline reading
      CASE 2
        BaseGauss! = Iso(I).BaseGauss2      ' optional baseline reading
    END SELECT

    CALL MoveTo(BaseGauss!, RangeNo(I), Delay)
    BadReadFlag = 0      ' keeps track of bad readings
    PRINT #3, "Base Gauss: "; BaseGauss!; "      Base Sigma Tol: "; BaseSigmaTol;
    PRINT #3, "      Base Max Dev: "; BaseMaxDev
    PRINT #3, "Mean      Sigma      MaxRejLev.  RejLevel  Good Readings"
    DO      ' get at least GoodReads good baseline readings without
      ' more than BaseMaxBadSets in a row bad readings
      IF BadReadFlag >= BaseMaxBadSets THEN
        ErrFlag = TRUE
        PRINT BadReadFlag; " bad Baseline readings in a Row. Bailing out..."
        GOTO WriteRecordFiles
      END IF

      CALL Measure(BaseSigmaTol, BaseMaxDev, BaseNumReads, Mean#, Sigma#, T!, GoodReads)
      RejectLevel# = BaseSigmaTol! * Sigma#
      MaxRejLevel# = BaseMaxDev! * Mean#
      PRINT USING "#####.##      "; BaseGauss!;
      PRINT USING "#.#####      "; Mean#; Sigma#; MaxRejLevel#; RejectLevel#;
      PRINT USING "#####      "; GoodReads
      PRINT #3, USING "#.#####      "; Mean#; Sigma#; MaxRejLevel#; RejectLevel#;
      PRINT #3, USING "#####      "; GoodReads
      BadReadFlag = BadReadFlag + 1
    LOOP UNTIL GoodReads >= BaseGoodReads
    PRINT
    PRINT #3,
    Baseline(I, L).Avrg = Mean#
    Baseline(I, L).StDv = Sigma#
    Baseline(I, L).Valid = GoodReads
    Baseline(I, 0).Avrg = Baseline(I, 0).Avrg + Mean# * GoodReads
    Baseline(I, 0).Valid = Baseline(I, 0).Valid + GoodReads
    SqSum# = SqSum# + Mean# ^ 2
  NEXT L
  Baseline(I, 0).Avrg = Baseline(I, 0).Avrg / Baseline(I, 0).Valid
NEXT I

FOR P = 1 TO 3 STEP 2
  GOSUB DisplayBaseData
NEXT P
' Debug -- CLOSE : INPUT "halt", X$
'End Measure Baseline
'-----
'Measure Samples
' Meas holds all data from one full set of ratio measurements
' change @ ver 7x5: now terminate after a full set; previously, only a
' single reference was done in the last set. Change accommodates
' reverse sense ratios
REDIM Meas(2 * NumRatios, NumIsos) AS MeasRecord

```

```

REDIM Group(2 * NumRatios + 1) 'Allow plenty of room
G = 1      'Counts Groups of good data sets; bad data set forces new Group
N = 0      'Number of ratios completed so far
J = 0      'index to measurement array

CALL ReadDVM(-IntegTime * 1#)

'Debug--
' LPRINT : LPRINT "SetNum      Group      FullSets  ";
' FOR I = 1 TO NumIsos
'   LPRINT USING "Mass ###      "; Iso(I).Mass;
' NEXT I: LPRINT
CLS
MeasFmt1$ = "###  ###      ###  #.#####  #.#####  ##"
MeasFmt2$ = "          ###  #.#####  #.#####  ##"
PRINT "Measured Values:"
PRINT : PRINT "Set  Group  Mass  Mean      Sigma      GoodReads"
PRINT #3, "Measured Values:"
PRINT #3, : PRINT #3, "Set  Group  Mass  Mean      Sigma      GoodReads"

DO ' Loop on J and N values
'
' Debug-- LPRINT USING "###      "; J; G; N;

FOR I = 1 TO NumIsos      ' Denotes isotope in sequence
CALL MoveTo(Iso(I).Gauss, RangeNo(I), Delay)
CALL Measure(IsoSigmaTol, IsoMaxDev, NumReads, Mean, Sigma, Time!, GoodReads)
Mean = Mean - Baseline(I, 0).Avrg
' Keep all data in LOG file, even bad data! Note no CR (so error messages
' Can be printed on the same line )
IF I = 1 THEN
PRINT USING MeasFmt1$; J; G; Iso(I).Mass; Mean; Sigma; GoodReads;
PRINT #3, USING MeasFmt1$; J; G; Iso(I).Mass; Mean; Sigma; GoodReads;
ELSE
PRINT USING MeasFmt2$; Iso(I).Mass; Mean; Sigma; GoodReads;
PRINT #3, USING MeasFmt2$; Iso(I).Mass; Mean; Sigma; GoodReads;
END IF

' Debug --
' LPRINT USING "### "; I; J;
' LPRINT USING "#.##### "; Mean; Sigma; GoodReads

IF (ABS(Mean) < IsoSigmaTol * Sigma) THEN
Mean = 0!      'forces bad read
PRINT " ERR: |Mean| < "; IsoSigmaTol; "* Sigma"
PRINT #3, " ERR: |Mean| < "; IsoSigmaTol; "* Sigma"
ELSEIF (GoodReads < IsoGoodReads) THEN 'Not enough good reads this time
Mean = 0!      'forces bad read
PRINT " ERR: Too Few GoodReads"
PRINT #3, " ERR: Too Few GoodReads"
ELSEIF (IsoSigmaTol * Sigma > IsoMaxDev * Mean) THEN 'Sigma too big
SigDivMean$ = STR$(CSNG(Sigma / Mean))
Mean = 0!      'forces bad read
PRINT " ERR: RSD"; SigDivMean$; " too big"
PRINT #3, " ERR: RSD"; SigDivMean$; " too big"
ELSE ' need CR at end of printed line for ok data
PRINT
PRINT #3,
END IF

```

```

IF (Mean = 0!) THEN
    G = G + 1                                ' start a new Group
    IF G > NumReads THEN
        PRINT " Too many bad data sets.  Bailing out"
        ErrFlag = TRUE
        GOTO WriteRecordFiles
    END IF
    N = N - 1                                ' don't count partial sets of data
    EXIT FOR                                  ' back to reference isotope
END IF

Meas(J, I).Mass = Iso(I).Mass
Meas(J, I).Avrg = Mean
Meas(J, I).StDv = Sigma
Meas(J, I).Tick = Time! - Time0!
Meas(J, I).Valid = GoodReads
Group(J) = G
IF (J = 2 * NumRatios) THEN GOTO ExceedsArray

NEXT I:
    PRINT : PRINT #3,
' Debug -- LPRINT

' If the Ref. Isotope is zero, there are no valid readings for the current
' value of J, and we do not store a row of nulls
    IF Meas(J, 1).Avrg > 0! THEN J = J + 1

    LastSample = J                          ' net total of J + 1 rows of information with
                                            ' NumRatios good sets so far (keep updating)

    N = N + 1                                ' counts ratios
    IF (N = NumRatios) THEN EXIT DO 'Normal Exit
LOOP    ' terminates on previous EXIT DO

    CALL MoveTo(Iso(1).Gauss, RangeNo(1), Delay)    'Move back to Ref

P = 3
GOSUB DisplayIsoMeas

'End Measure Samples
'-----
'Analyze Ratios
' Subscripts in Ratio are sample#, Isotope #, and ratio type
' Results record includes three kinds of ratios: "normal" and "reverse
' sense" interpolated ratios (Subscripts 1 and 2) and the combined results
' for both kinds of ratios (Subscript 0)
' RefIso is the index of the Reference Isotope (conventionally 1)

REDIM Ratio(LastSample + 1, NumIsos, 2) AS DOUBLE
REDIM Results(NumIsos, 2) AS ResultRecord
FOR I = 1 TO NumIsos
    IF I <> RefIso THEN ' Form Ratios for Sample Isotopes only
        FOR J = 0 TO LastSample - 1 ' Index to Meas() Array
' "Normal" ratio; interpolate between bracketing Refs to time of Sample
            Test! = Meas(J, I).Avrg * Meas(J, RefIso).Avrg * Meas(J + 1, RefIso).Avrg
            IF (Test!) > 0 THEN ' no zero elements in Test!
                InterpRef# = Interpol8#(Meas(), I, RefIso, J)
                Ratio(J, I, 1) = Meas(J, I).Avrg / InterpRef#
                CALL RatioStats(Ratio(), Results(), I, J, 1)
            END IF 'Test!
        
```

```

' "Reverse Sense" ratio; interpolate between bracketing Sampls. to time of Ref
  IF RevSense THEN
    Test! = Meas(J, I).Avrg * Meas(J + 1, I).Avrg * Meas(J + 1, RefIso).Avrg
    IF (Test!) > 0 THEN      ' no zero elements in Test!
      InterpAv# = Interpol8#(Meas(), RefIso, I, J + 1)
      Ratio(J + 1, I, 2) = InterpAv# / Meas(J + 1, RefIso).Avrg
      CALL RatioStats(Ratio(), Results(), I, J + 1, 2)
    END IF 'Test!
  END IF 'RevSense
NEXT J
' Now filter out of spec ratios
  FOR T = 1 TO TMax      'TMax = 2 => RevSense is TRUE. TMax = 1 otherwise
  WHILE Results(I, T).NumRatios <> Results(I, T).GoodRatios
    Results(I, T).NumRatios = Results(I, T).GoodRatios 'save GoodRatios
  'Activate this code to compare kinds of ratios separately
  '   Mean = Results(I, T).RatioAv
  '   Tol! = RatioSigmaTol * Results(I, T).RatioDev
  ' Activate this code to use the combined ratios for comparison Also
  ' generates Normal Ratio data if Reverse Sense isn't used.
    Mean = Results(I, 0).RatioAv
    Tol! = RatioSigmaTol * Results(I, 0).RatioDev

    IF Tol! > RatioMaxDev * Mean THEN
      Tol! = RatioMaxDev * Mean
      BEEP: PRINT "ERR: Ratio Deviation clamped to RatioMaxDev * Mean"
      PRINT #3, "Ratio Tol! ="; Tol!; " Results("; I; ", 0).RatioDev = ";
      PRINT #3, Results(I, 0).RatioDev
    END IF

    FOR J = 0 TO LastSample
      IF Ratio(J, I, T) > 0 AND (ABS(Ratio(J, I, T) - Mean) > Tol!) THEN
        Ratio(J, I, T) = -Ratio(J, I, T)
        CALL RatioStats(Ratio(), Results(), I, J, T)
      END IF
    NEXT J
  WEND
NEXT T
END IF 'I<>RefIso
NEXT I

P = 3
  GOSUB DisplayRatios
  GOSUB DisplayAnalysis

'End Analyze Ratios
'-----
'Begin Write Recordfiles
' Write files as comma delimited ASCII. These files are formatted and
' ordered for use with a spreadsheet. The information is the order: header
' information, followed by results, calculated ratios, Isotope measurements
' and baseline measurements. Nine or more significant figures are kept for
' double precision variables.

WriteRecordFiles:
  ON ERROR GOTO DiskWriteError
  FOR P = 4 TO 5
    WRITE #P, SampName$, DATE$, TickMark$, FilVolts$, FilAmps$, IsoName
    WRITE #P, STR$(IntegTime) + " Hz", SubmitName$, OpName$
    WRITE #P, Comment$
  
```

```

FOR I = 1 TO NumIsos
    PRINT #P, USING " ###_, "; Iso(I).Mass;
NEXT I: PRINT #P, USING " ###"; Iso(RefIso).Mass
PRINT #P,
Fs$ = CHR$(47)

IF ErrFlag THEN GOTO RawData' write computations only for complete runs
FOR I = 1 TO NumIsos
    IF I <> RefIso THEN
        WRITE #P, Iso(I).Mass, Fs$, Iso(RefIso).Mass
        FOR T = 0 TO 2
            PRINT #P, USING "\          \"; RatioType$(T);
            PRINT #P, USING ", ###.#####"; Results(I, T).RatioAv;
            PRINT #P, USING ", #.#####"; Results(I, T).RatioDev;
            PRINT #P, USING ", ###"; Results(I, T).GoodRatios
        NEXT T: WRITE #P,
    END IF
NEXT I

FOR I = 1 TO NumIsos
    IF I <> RefIso THEN
        WRITE #P, Iso(I).Mass, Fs$, Iso(RefIso).Mass
        FOR J = 0 TO LastSample
            IF Ratio(J, I, 1) <> 0 OR Ratio(J, I, 2) <> 0 THEN
                PRINT #P, USING " ###_, ##"; J; Group(J);
                FOR T = 1 TO 2
                    IF Ratio(J, I, T) <> 0 THEN
                        PRINT #P, USING ", ###.#####"; Ratio(J, I, T);
                    ELSE
                        PRINT #P, ",    0.0          ";
                    END IF
                NEXT T: PRINT #P,
            END IF 'Ratios <> 0
        NEXT J: PRINT #P,
    END IF
NEXT I
PRINT #P,
RawData:
FOR I = 1 TO NumIsos
    PRINT #P, USING "###_, #####.#"; Iso(I).Mass; Iso(I).Gauss
    FOR J = 0 TO LastSample
        IF Meas(J, I).Avrg > 0 THEN
            PRINT #P, USING " ###_, ##"; J; Group(J);
            PRINT #P, USING ", ###.#####"; Meas(J, I).Avrg;
            PRINT #P, USING ", #.#####"; Meas(J, I).StDv;
            PRINT #P, USING ", ###"; Meas(J, I).Valid;
            DelTime! = MaxOf(Meas(J, I).Tick - Meas(J + (J > 0), I).Tick, 0)
            IF (Meas(J + (J > 0), I).Tick = 0) THEN DelTime! = 0
            PRINT #P, USING ", ###.#####"; DelTime!;
            PRINT #P,
        END IF
    NEXT J: PRINT #P,
NEXT I

FOR I = 1 TO NumIsos
    PRINT #P, USING " ###, "; Iso(I).Mass
    FOR L = 1 TO NumBases(I)
        IF L = 1 THEN PRINT #P, USING " #####.#, "; Iso(I).BaseGauss1;
        IF L = 2 THEN PRINT #P, USING " #####.#, "; Iso(I).BaseGauss2;
        PRINT #P, USING "#.#####", "; Baseline(I, L).Avrg;
        PRINT #P, USING "#.#####", "; Baseline(I, L).StDv;
    
```



```

        PRINT #P, USING "###"; Baseline(I, L).Valid;
        PRINT #P,
    NEXT L
    PRINT #P, USING "          , #.#####,"; Baseline(I, 0).Avrg;
    PRINT #P, USING "          , ###"; Baseline(I, 0).Valid
    PRINT #P,
NEXT I: PRINT #P, CHR$(12)
NEXT P
ON ERROR GOTO 0
'End Write Recordfiles
'-----
'Report
CLS :

FOR P = 1 TO 2

    GOSUB DisplayHeader

    GOSUB DisplayBaseData

    GOSUB DisplayIsoMeas

    IF ErrFlag = FALSE THEN 'Only if the measurement completed properly
        GOSUB DisplayRatios
        GOSUB DisplayAnalysis
    END IF
NEXT P

PRINT #2, CHR$(12); : PRINT #3, CHR$(12);
PRINT #1, CHR$(13):

LOCATE 25, 2: COLOR 15, 1
INPUT ; "Hit 'ENTER' to continue ", X$: COLOR 7, 1

'End Report
'-----
'Begin LoopLogic
CLS :
' size = 4 and RowNum = 2 (rows from top)
RepeatFlag = 0 'Default choice to be highlighted
CALL VertMenu(RepeatHeader$, RepeatMenu(), 4, 2, RepeatFlag)
COLOR 7, 1: CLS
SameParams = FALSE
SELECT CASE RepeatFlag
    CASE 0
        ' New Sample, Change parameters
        SampName$ = ""
        ' Activate the following line to pre-clear all information
        '   FilAmps$ = "": FilVolts$ = "": OpName$ = "": SubmitName$ = ""
        CALL GetSampName(SampName$)
        CALL GetInfo(OpName$, FilAmps$, FilVolts$, SubmitName$, Comment$)
    CASE 1
        ' New Sample, Same parameters
        SampName$ = ""
        CALL GetSampName(SampName$)
        SameParams = TRUE
    CASE 2
        ' Same Sample, Change Paramaters

        ' Activate the following line to pre-clear all information
        '   FilAmps$ = "": FilVolts$ = "": OpName$ = "": SubmitName$ = ""

        CALL GetInfo(OpName$, FilAmps$, FilVolts$, SubmitName$, Comment$)

```

```

CASE 3          ' Redo this measurement (no changes)
    SameParams = TRUE
CASE ELSE      ' Quit
    STOP
END SELECT

GOTO SetupFiles

PRINT " At 'Begin Error messages'. Should not ever get here!": GOTO fini

'End LoopLogic
'-----
'Begin Display GoSubs
DisplayHeader:
    PRINT #P,
    PRINT #P, "      Sample Name: "; SampName$; TAB(36);
    PRINT #P, " Time: "; TickMark$; TAB(52);
    PRINT #P, "Date: "; DATE$
    PRINT #P, " Filament Volts: "; FilVolts$; TAB(40);
    PRINT #P, " Filament Amps: "; FilAmps$
    PRINT #P, "Integration Time: "; STR$(IntegTime); " Line Cycles"
    PRINT #P, " Operators Name: "; OpName$; TAB(40);
    PRINT #P, "Submitter's Name: "; SubmitName$
    PRINT #P, "      Comments: "; Comment$
    PRINT #P,
RETURN

DisplayBaseData:
    PRINT #P,
    PRINT #P, "Baseline Measurements:"
    FOR I = 1 TO NumIsos
        PRINT #P, USING " Mass: ### "; Iso(I).Mass
        PRINT #P, "Base Gauss ";
        PRINT #P, "Baseline      Sigma      Reads ";
        PRINT #P,
        FOR L = 1 TO NumBases(I)
            IF L = 1 THEN PRINT #P, USING " #####.## "; Iso(I).BaseGauss1;
            IF L = 2 THEN PRINT #P, USING " #####.## "; Iso(I).BaseGauss2;
            PRINT #P, USING "##### "; Baseline(I, L).Avrg;
            PRINT #P, USING "##### "; Baseline(I, L).StdV;
            PRINT #P, USING "###"; Baseline(I, L).Valid;
            PRINT #P,
        NEXT L
        PRINT #P, USING "      Average: ##### "; Baseline(I, 0).Avrg;
        PRINT #P, USING "      #####"; Baseline(I, 0).Valid;
        PRINT #P, : PRINT #P,
    NEXT I
RETURN

DisplayIsoMeas:
    PRINT #P, "Isotope Measurements:"
    FOR I = 1 TO NumIsos
        PRINT #P, USING " Mass: ### Gauss: #####.##"; Iso(I).Mass; Iso(I).Gauss
        PRINT #P, " No. Group ";
        PRINT #P, "##### ##### ##### #####";
        PRINT #P, " Mean      Sigma      Rel.Std.Dev      Reads      Time Int.";
        PRINT #P,
    
```

```

FOR J = 0 TO LastSample
  IF Meas(J, I).Avrg > 0 THEN
    PRINT #P, USING "### ###      "; J; Group(J);
    PRINT #P, USING "###.#####      "; Meas(J, I).Avrg;
    PRINT #P, USING "###.#####      "; Meas(J, I).StDv;
    PRINT #P, USING "###.#####      "; Meas(J, I).StDv / Meas(J, I).Avrg;
    PRINT #P, USING "###      "; Meas(J, I).Valid;
    DelTime! = MaxOf(Meas(J, I).Tick - Meas(J + (J > 0), I).Tick, 0)
    IF (Meas(J + (J > 0), I).Tick = 0) THEN DelTime! = 0
    PRINT #P, USING "###.###"; DelTime!;
    PRINT #P,
  END IF
NEXT J: PRINT #P,
NEXT I
RETURN

DisplayRatios:
FOR I = 1 TO NumIsos
  IF I <> RefIso THEN
    PRINT #P, USING " Ratios: ###/###"; Iso(I).Mass; Iso(1).Mass
    '
    "### ###      ###.#####      ###.#####      "
    PRINT #P, " Set Group   Samp/Interp-Ref      ";
    IF RevSense THEN PRINT #P, "Interp-Samp/Ref ";
    PRINT #P,
    FOR J = 0 TO LastSample
      '
      print ratios only if one or more of them is <> 0
      IF Ratio(J, I, 1) <> 0 OR Ratio(J, I, 2) <> 0 THEN
        PRINT #P, USING "### ###      "; J; Group(J);
        FOR T = 1 TO TMax
          IF Ratio(J, I, T) <> 0 THEN
            PRINT #P, USING "#####.#####      "; Ratio(J, I, T);
          ELSE
            PRINT #P, "      -----      ";
          END IF
        NEXT T: PRINT #P,
      END IF 'Ratios <> 0
    NEXT J: PRINT #P,
  END IF
NEXT I
PRINT #P,
RETURN

DisplayAnalysis:
PRINT #P, "Analysis of Ratios  "
FOR I = 2 TO NumIsos
  PRINT #P, USING "Mass Ratio: ###/###"; Iso(I).Mass; Iso(1).Mass;
  PRINT #P, "      Mean      (Inverse)      ";
  PRINT #P, " Sigma      Rel.Std Dev      Valid"
  ' (N < 3) evaluates to -1 if it is TRUE and to 0 if it is FALSE
  FOR TT = 0 TO 2 * (-(RevSense <> 0)) ' Either Samp/Interp-Ref only or
    ' all three
    T = (TT + 1) * -(TT < 2) ' T sequence 1, 2, 0 or just 1
    PRINT #P, RatioType$(T);
    PRINT #P, USING "#####.#####      "; Results(I, T).RatioAv;
    IF Results(I, T).RatioAv > 0 THEN
      PRINT #P, USING "(#####.#####) "; 1 / (Results(I, T).RatioAv);
      PRINT #P, USING "#####.#####      "; Results(I, T).RatioDev;
      PRINT #P, USING "###.#####      "; Results(I, T).RatioDev / Results(I, T).RatioAv;
    END IF
  NEXT TT
NEXT I

```

```

        ELSE
            PRINT #P, "    ----    ";
            PRINT #P, "    ----    ";
            PRINT #P, "    ----    ";
        END IF
        PRINT #P, USING "###"; Results(I, T).GoodRatios
    NEXT TT: PRINT #P,
NEXT I
RETURN

'End Display GoSubs
'-----
'Begin ERROR messages
ExceedsArray:
    BEEP
    PRINT " Can't store more than "; 2 * NumRatios; " measurements.  Bailing out..."
    RESUME fini

ParamFileErr:
    BEEP: PRINT "Can't find disk file ISOPARAM.REC.  Bailing out..."
    RESUME fini

DiskWriteError:
    PRINT "Error on Writing to disk: "; ERR
    IF P = 4 THEN PRINT "Hard Disk Not responding or full"
    IF P = 5 THEN PRINT "Diskette Not ready or full"
    RESUME fini

fini:
' Debug--
COLOR 14, 4: INPUT ; "    Hit Return to Stop or ^C to Break", X$: COLOR 7, 1
CLOSE : END

BadOpenDataFile:
PRINT "Open error = "; ERR
SELECT CASE ERR
    CASE 71
        PRINT "Diskette drive B: is not ready; either there is no disk ";
        PRINT "or it isn't properly inserted."
    CASE ELSE
        PRINT "There is some kind of error with Diskette Drive B:"
END SELECT
PRINT "Hit any key to continue": X$ = INPUT$(1)
RESUME OpenDataFile

```

## SUBROUTINES AND FUNCTIONS:

```

REM $STATIC
FUNCTION GetAns$ (Prompt$, Default$)
' If there is no default, an answer is not required.  Prompt$ and Default$
' added in Version 6x4
'RCB 10Mar93
COLOR 7, 1
PRINT Prompt$;
IF LEN(Default$) > 0 THEN 'answer is required or you get the default
    PRINT "(Default: "; Default$; ") ";
    INPUT ; "", Entry$: PRINT
    IF LEN(Entry$) = 0 THEN
        Entry$ = Default$
    END IF
ELSE 'Entry optional
    INPUT ; "", Entry$: PRINT
END IF
GetAns$ = Entry$
END FUNCTION

SUB GetInfo (OpName$, FilAmps$, FilVolts$, SubmitName$, Comment$)
    PRINT " Please Enter:"
    OpName$ = GetAns$(" Your initials or last name: ", OpName$)
    FilAmps$ = GetAns$(" The Filament Amperage: ", FilAmps$)
    FilVolts$ = GetAns$(" The Filament Voltage: ", FilVolts$)
    SubmitName$ = GetAns$(" The Submitters Name: ", SubmitName$)
    Comment$ = GetAns$("Comment: ", Comment$)
    PRINT :
END SUB

FUNCTION GetInt
GetInt = CINT(VAL(GetStr12$))
END FUNCTION

SUB GetSampName (SampName$)

GetSampName:
    SampMsg$ = " Please enter a unique Sample ID ( up to 8 Char): "
    DO
        SampName$ = GetAns$(SampMsg$, SampName$)
    LOOP UNTIL LEN(SampName$) > 0
    SampName$ = LEFT$(SampName$, 8)
    PRINT
END SUB

DEFSNG A-Z
FUNCTION GetSingle
GetSingle = VAL(GetStr12$)
END FUNCTION

DEFINT A-Z
FUNCTION GetStr12$
LINE INPUT #6, Lin$
X$ = LTRIM$(LEFT$(Lin$, 12))
GetStr12$ = X$
END FUNCTION

```

```

FUNCTION Interpol8# (M() AS MeasRecord, SampIndex%, RefIndex%, MeasIndex%)
' M will be a 2 dimensional array (typically Meas(J,I) ). J is the index
' to the Meas array. SampIndex and RefIndex are the sample and reference
' isotope indices (typically RefIndex = 1 and SampIndex = 2 or 3).
' The measurements are taken in the sequence Ref.. Samp.. Ref.. Samp....
' Delt is a weighting function based on the times of measurement for a
' Sample Isotope and its bracketing Reference Isotopes. It determines
' Interpol8, the linearly interpolated value for the Ref. Isotope at the
' time the Sample Isotope was measured. Note that calling this function with
' Sample and Reference roles reversed finds an interpolated value of the
' sample at the time the reference was measured
,
DIM Delt AS SINGLE
SI = SampIndex%: R = RefIndex%:
JS = MeasIndex% ' Set index for the sample
JR = MeasIndex% ' Set index for the reference
' The "JR" subscript for the reference must be adjusted to
' bracket the sample in the time sequence. Time(JR) < Time(JS) < Time(JR+1)
DO UNTIL M(JR, R).Tick < M(JS, SI).Tick
    JR = JR + 1
LOOP
Delt = (M(JS, SI).Tick - M(JR, R).Tick) / (M(JR + 1, R).Tick - M(JR, R).Tick)
Interpol8# = (M(JR + 1, R).Avrg * Delt + M(JR, R).Avrg * (1 - Delt))
' LPRINT " Samp   Ref   JSamp  JRef   Times/Averages @ R Interp and R+1"
' LPRINT USING " ##      "; SI; R; JS; JR;
' LPRINT USING " #.#####      "; M(JR, R).Tick; M(JS, SI).Tick; M(JR + 1, R).Tick
' LPRINT USING " ##      "; SI; R; JS; JR;
' LPRINT USING " #.#####      "; M(JR, R).Avrg; (M(JR + 1, R).Avrg * Delt + M(JR, R).Avrg
* (1 - Delt)); M(JR + 1, R).Avrg
' LPRINT
END FUNCTION

FUNCTION MaxOf# (A#, B#)
IF A# >= B# THEN
    MaxOf = A#
ELSE
    MaxOf = B#
END IF
END FUNCTION

SUB Measure (SigmaTol!, MaxStdDev!, Reads, Mean#, Sigma#, Time!, GoodReads)
' Calling routine supplies SigmaTol! MaxStdDev!, Reads
' Measure returns mean#, Sigma#, Time! and Good Readings
' RejLevel is a multiple of Sigma# but limited to a maximum fraction of Mean#
DIM Sum AS DOUBLE, SumOfSqs AS DOUBLE
REDIM Sample(Reads) AS DOUBLE

CALL ReadDVM(Mean#) 'Dummy read to synchronize with DVM
FOR J = 1 TO Reads
    CALL ReadDVM(Sample(J))
    Sum = Sum + Sample(J)
    SumOfSqs = SumOfSqs + (Sample(J)) ^ 2
NEXT J
Time! = TIMER
Mean# = Sum / Reads
Sigma# = SQR(SumOfSqs / Reads - (Mean# ^ 2))
RejectLevel# = SigmaTol! * Sigma# ' a multiple of Sigma
MaxRejLevel# = MaxStdDev! * Mean# ' a fraction of Mean

```

```

IF RejectLevel# > MaxRejLevel# THEN
'   PRINT "MEASURE: WARNING! Sigma exceeds ";SigmaTol!; * Max. Reject Level"
'   PRINT #3, "MEASURE: WARNING! Sigma exceeds ";SigmaTol!; * Max. Reject Level"
    RejectLevel# = MaxRejLevel#
END IF
GoodReads = Reads
FOR J = 1 TO Reads
    IF (ABS(Sample(J) - Mean#) > RejectLevel#) THEN
        Sum = Sum - Sample(J)
        GoodReads = GoodReads - 1
    END IF
NEXT J
IF GoodReads > 0 THEN
    Mean# = Sum / GoodReads
ELSE
    Mean# = 0:
END IF
EXIT SUB
END SUB

```

```

SUB MoveTo (Gauss!, Range, Delay)
'Rev to return cursor to its col,row on input. all printing is on row 25
'no CRs allowed so screen will not roll up unless input row is 25, then
'it rolls up
Finish! = TIMER + Delay
CALL OutputGauss(Gauss!, 7)          ' Range 7 to avoid saturating the PA
C = CSRLIN
R = POS(0)
DO WHILE R > 24
    PRINT : R = POS(0)
LOOP
LOCATE 25, 1: COLOR 7, 1: PRINT SPACE$(79);
LOCATE 25, 1: PRINT "      Changing the Magnet      ";
P = POS(0) - 1: COLOR 1, 7
DO UNTIL (Finish! < TIMER)
    LOCATE 25, P: PRINT "|"; : LOCATE 25, P: PRINT "|";
    LOCATE 25, P: PRINT "/"; : LOCATE 25, P: PRINT "/";
    LOCATE 25, P: PRINT "-"; : LOCATE 25, P: PRINT "-";
    LOCATE 25, P: PRINT "\"; : LOCATE 25, P: PRINT "\";
LOOP
COLOR 7, 1
CALL OutputGauss(Gauss!, Range)      ' Changes range only
LOCATE 25, 1: COLOR 7, 1: PRINT SPACE$(79);
LOCATE C, R
END SUB

```

```

SUB OutputGauss (Gauss!, Range)
' Contains IEEE-488 interface specific calls
' Transmits 4 bytes to GaussDevice; computation of array G() follows
' NBS Program. Output is in BCD nibbles except for first byte in sequence:
' Braces indicate G() array values as 16 bit integers
' [0 Range], [1st BCD 2nd BCD], [3rd BCD 4th BCD], [5th BCD LOAD]
' Range is 1 to 8 (not 0 to 7 as in the rest of the program)
' LOAD is the load control byte, 14; 10 Kilogauss overrange is not used
' To use it increase LOAD to 15.

```

```

STATIC InitFlag AS INTEGER, OldGauss AS SINGLE, OldRange AS INTEGER
STATIC GaussDevice AS INTEGER
IF Gauss! = OldGauss AND OldRange = Range THEN ' we are already there
    EXIT SUB
END IF

```

```

OldGauss = Gauss!                                ' otherwise save the new ones
OldRange = Range
REDIM G(4) AS INTEGER

InitIt:
  IF InitFlag = 0 THEN
    GaussName$ = "GAUSSM"
    ' GaussDevice = 8 is the usual IEEE-488 address of the Gaussmeter
    CALL IBFIND(GaussName$, GaussDevice)
    IF GaussDevice < 0 THEN
      PRINT "Can't Find Gaussmeter on IEEE-488 Buss. Bailing out..."
      ON ERROR GOTO fini
      ERROR 254
    END IF
    MASK = &H8000                                ' error condition
    CALL IBTRAP(MASK, 3)                          ' sets up applications monitor for mode 3
                                                    ' error trapping
    CALL IBCLR(GaussDevice)                       ' clear Gaussmeter
    InitFlag = 1                                  ' should not need to do this again
  END IF                                          ' now we are ready

IF Gauss! < .1 OR Gauss! >= 9999.95 OR Range < 1 OR Range > 7 THEN
  ON ERROR GOTO fini
  PRINT "Error in OutputGauss: Gauss = "; Gauss!; " Range = "; Range
  ERROR 253
END IF

' Form a string Orggggg0 where r is the output range and n is a digit from
' the gauss readings, without decimal point. Convert string to 8 BCD digits
' (4 bytes). Then add hex 'E', Load Control trigger for the Gauss controller

Gauss$ = STR$(FIX(Range * 100000! + 10 * Gauss!)) + "0"

FOR I = 1 TO 7 STEP 2
  G(I \ 2) = &H10 * VAL(MID$(Gauss$, I, 1)) + VAL(MID$(Gauss$, I + 1, 1))
NEXT I
G(3) = G(3) + &HE 'adds trigger character
Gauss$ = ""
FOR I = 0 TO 3
  Gauss$ = Gauss$ + CHR$(G(I))
NEXT I
CALL IBWRT(GaussDevice, Gauss$)                  ' Send the bytes to the Gaussmeter
END SUB

SUB RatioStats (Ratio() AS DOUBLE, Res() AS ResultRecord, I, J, T)
'
' All definitions are as in main program. Logic is:
'   Sum = Sum + Value
'   SqSum = SqSum + Value ^ 2
'   Number = Number + 1
'   Average = Sum/Number
'   StdDev = SquareRoot(SqSum/Number - Average^2)
' Note that the combined stats (X=0) and the individual stats (X = T)
' are updated during each call
IF Ratio(J, I, T) = 0 THEN
  LPRINT " DEBUG RatioStats: was called with a zero value for Ratio"
  LPRINT " Set Index "; J; " Isotope Index "; I; " and Ratio type "; T
  STOP ' this shoudn't happen
END IF

```



```

FOR X = 0 TO T STEP T
  Res(I, X).SumRatio = Res(I, X).SumRatio + Ratio(J, I, T)
  IF Ratio(J, I, T) > 0 THEN
    Res(I, X).SqSumRatio = Res(I, X).SqSumRatio + Ratio(J, I, T) ^ 2
    Res(I, X).GoodRatios = Res(I, X).GoodRatios + 1
  ELSE
    Res(I, X).SqSumRatio = Res(I, X).SqSumRatio - Ratio(J, I, T) ^ 2
    Res(I, X).GoodRatios = Res(I, X).GoodRatios - 1
  END IF
  IF Res(I, X).GoodRatios = 0 THEN
    LPRINT "DEBUG RatioStats: Number of Ratios is zero for:"
    LPRINT " Set Index "; J; " Isotope Index "; I; " and Ratio type "; X
    STOP ' this shoudn't happen
  END IF
  Res(I, X).RatioAv = Res(I, X).SumRatio / Res(I, X).GoodRatios
  SigSq# = Res(I, X).SqSumRatio / Res(I, X).GoodRatios
  SigSq# = SigSq# - Res(I, X).RatioAv ^ 2
  IF SigSq# < .0000000000000001# THEN SigSq# = 0
  Res(I, X).RatioDev = SQR(SigSq#)
NEXT X
,
END SUB

```

```

SUB ReadDVM (VmRead#)
' Last change 3May93 RCB: use VmRead# as a flag on initialization so 10 or
' 100 Hz can be chosen on the fly. (10 is default). If the routine is
' called without previous initialization, it defaults
  STATIC InitFlag AS INTEGER
  STATIC DvmDevice AS INTEGER
  DIM DvmInput AS STRING * 14
  IF VmRead# < -1# THEN InitFlag = 0 ' Force reinitialization
InitDVM:
  IF InitFlag = 0 THEN
    DVMname$ = "HP-DVM"
'   DVMdevice = 22 '22 is the usual IEEE-488 address of the DVM
    CALL IBFIND(DVMname$, DvmDevice)
    IF DvmDevice < 0 THEN
      PRINT "Can't Find VoltMeter on IEEE-488 Buss. Bailing out..."
      ON ERROR GOTO fini
      ERROR 254
    END IF
    IF CINT(ABS(VmRead#)) = 100 THEN
      DVMcommand$ = "H6STG100STIF1R4Z0" 'Setup command from prev. program
    ELSE
      DVMcommand$ = "H6STGW10STIF1R4Z0" 'Setup command for testing
    END IF
    Stat = ILWRT(DvmDevice, DVMcommand$, 17) ' Send command to DVM
'   CONST EEND=    &H2000 ' EOI or EOS detected
'   CONST CMPL=    &H100  ' I/O completed
    IF (Stat <> CMPL) THEN
      PRINT "Bad Status ("; HEX$(Stat); ") on HP DVM. Bailing out....."
      ON ERROR GOTO fini
      ERROR 251
    END IF
    InitFlag = 1 ' should not need to do this again
  END IF ' now we are ready
  Stat = ILRD(DvmDevice, DvmInput, 14)

```

```

IF (Stat <> EEND + CMPL) THEN
    PRINT "DVM Read error(s); Status is: "; HEX$(Stat); " Err code: "; IBERR;
    PRINT "  Bailing out....."
    ON ERROR GOTO fini
    ERROR 250
ELSE
    VmRead# = VAL(DvmInput)
END IF
END SUB

SUB VertMenu (Header$, Menu$(), Size, RowNum, Choice)
' VGA Colors
' Required Inputs: Header$, Menu$(), Size = index of largest valid entry
' in Menu$() and RowNum (Rows from the top of screen) Choice is default
' value to be highlighted on input.
' Menu subscripts start at 0. Screen is not cleared.
' Menu choices must fit one to a line; up to 20 choices allowed (full screen)
' Returns to caller: Choice (Menu subscript from 0 to Size)
' If key is an up arrow or a down arrow Choice changes and Key$ is reassigned
' The Key$ variable is 2 bytes for arrow keys. The value in Key$ is
' actually a null plus either "H" or "P". Capslock is ignored.
' Numlock and Shift affect the numeric keypad as usual.

Upkey$ = CHR$(0) + "H": Downkey$ = CHR$(0) + "P" ' Actual values
up$ = CHR$(24): dn$ = CHR$(25) ' these print as up and down arrows
CR$ = CHR$(13)
COLOR 15, 2
LOCATE 25, 5: COLOR 15, 6
PRINT "Select your entry with the "; up$; dn$; " keys, then hit 'Enter'";
COLOR 15, 2
LOCATE RowNum, 5
PRINT Header$: RowNum = RowNum + 2
LOCATE RowNum, 5
FOR I = 0 TO Size
    LOCATE , 5: PRINT Menu$(I)
NEXT I
LOCATE RowNum + Choice, 5
' Choice = 0
GOSUB Hilight

WHILE (1)
    DO: Key$ = INKEY$: LOOP WHILE Key$ = ""
    IF (Key$ = Upkey$) THEN
        GOSUB Normal
        Choice = Choice - 1
        IF Choice = -1 THEN Choice = Size
        GOSUB Hilight
    ELSEIF (Key$ = Downkey$) THEN
        GOSUB Normal
        Choice = (Choice + 1) MOD (Size + 1)
        GOSUB Hilight
    ELSEIF Key$ = CR$ THEN
        COLOR 7, 1
        GOTO ArrowFini
    ELSE
        LOCATE 24, 5: COLOR 14, 4
        BEEP:
    PRINT "Sorry, I don't understand that! Hit 'Q' to quit; any other key to continue";
    X$ = INPUT$(1)

```

```

        IF UCASE$(X$) = "Q" THEN
            ON ERROR GOTO fini
            ERROR 245
        END IF
        LOCATE 24, 5: COLOR 15, 1: PRINT SPACE$(75); : LOCATE 24, 5: COLOR 15, 6
        PRINT "Select your entry with the "; up$; dn$; " keys, then hit 'Enter'";
        GOSUB Hilight
    END IF
WEND
Hilight:
    COLOR 15, 6: LOCATE RowNum + Choice, 5
    PRINT Menu$(Choice);
RETURN
Normal:
    COLOR 15, 2: LOCATE RowNum + Choice, 5
    PRINT Menu$(Choice);
RETURN
ArrowFini:
RowNum = RowNum + Size + 2
'COLOR 7, 1: CLS
END SUB

FUNCTION YNans (message$)
' this routine returns -1 (true) for yes and 0 (false) for no. The
' <cr> option returns -1 RCB 28Feb92
PRINT
DO
.....    PRINT message$; " ( y/n or <cr> ): ";
    COLOR 1, 7: PRINT "Y"; : COLOR 7, 1: LOCATE , POS(0) - 1
    A$ = INPUT$(1): PRINT A$;
    IF A$ = CHR$(13) THEN A$ = "Y"
    IF UCASE$(A$) = "Y" THEN
        Ans = -1
    ELSEIF UCASE$(A$) = "N" THEN
        Ans = 0
    ELSEIF UCASE$(A$) = "Q" THEN
        ON ERROR GOTO fini
        PRINT : ERROR 245
    ELSE
        Lp% = CSRLIN: GOSUB WipeIt: LOCATE 23, 5
        TRef! = TIMER + 4: COLOR 14, 4
        PRINT "Sorry, I don't understand that! Please answer "; : LOCATE 24, 5
        PRINT "'y' or <cr> for 'yes'; 'n' for 'no' or 'q' to quit";
        DO: LOOP WHILE TIMER < TRef!
        GOSUB WipeIt: LOCATE Lp%, 1
        Ans = 1
    END IF
LOOP UNTIL (Ans < 1)
Lp% = CSRLIN
YNans = Ans
LOCATE 24, 5: PRINT SPACE$(50);
LOCATE Lp% + 1, 1
EXIT FUNCTION

WipeIt:
    COLOR 7, 1: FOR W = 23 TO 24
        LOCATE W, 1: PRINT SPACE$(80); : NEXT W
RETURN

END FUNCTION

```