

U.S. Department of the Interior
U.S. Geological Survey

**Methods of Analysis by the U.S. Geological Survey
Organic Geochemistry Research Group—
Determination of Selected Herbicides and Their
Degradation Products in Water Using Solid-Phase
Extraction and Gas Chromatography/Mass
Spectrometry**

By J.L. KISH, E.M. THURMAN, E.A. SCRIBNER, and L.R. ZIMMERMAN

Open-File Report 00-385

Lawrence, Kansas
2000

U.S. Department of the Interior

Bruce Babbitt, Secretary

U.S. Geological Survey

Charles G. Groat, Director

The use of firm, trade, and brand names in this report is for identification purposes only and does not constitute endorsement by the U.S. Geological Survey.

For additional information write to:

District Chief
U.S. Geological Survey
4821 Quail Crest Place
Lawrence, KS 66049-3839

Copies of this report can be purchased from:

U.S. Geological Survey
Information Services
Building 810
Box 25286, Federal Center
Denver, CO 80225-0286

CONTENTS

Abstract	1
Introduction	1
Determination of Herbicides and Degradation Products in Water Using Solid-Phase Extraction and Gas Chromatography/Mass Spectrometry	2
Method of Analysis	2
Scope and Application	2
Summary of Method	2
Interferences	3
Apparatus and Instrumentation	3
Reagents and Consumable Materials	4
Sampling Methods, Sample-Collection Equipment, and Cleaning Procedures	4
Standards	4
Gas Chromatograph/Mass Spectrometry Performance	5
Evaluation of Gas Chromatograph Performance	5
Evaluation of Mass Spectrometer Performance	5
Calibration	6
Procedure	6
Calculation of Results	8
Qualitative Identification	8
Quantitation	9
Reporting of Results	9
Method Performance	9
Conclusions	10
References Cited	13

FIGURES

1. Chromatogram showing total ions of eight herbicides and five degradation products in 1.0-nanogram-per-microliter standard solution	5
2. Graphs showing examples of high, intermediate, and low calibration curves for dimethenamid	8

TABLES

1. Herbicides and degradation products suitable for determination using method described, with registry numbers and molecular weights	3
2. Retention times and relative retention times of selected herbicides, degradation products, surrogate compounds, and internal standard analyzed using method described	7
3. Mean recoveries and relative standard deviations of selected herbicides and degradation products in reagent-water samples using gas chromatography/mass spectrometry	11
4. Mean recoveries and relative standard deviations of selected herbicides and degradation products in surface-water samples using gas chromatography/mass spectrometry	11
5. Mean recoveries and relative standard deviations of selected herbicides and degradation products in ground-water samples using gas chromatography/mass spectrometry	12
6. Mean observed concentrations, standard deviations, and method detection limits of selected herbicides and degradation products in reagent-water samples spiked at 0.05 microgram per liter	12
7. Absolute recoveries for selected herbicides and degradation products spiked at concentrations of 1.0 microgram per liter	13

CONVERSION FACTORS, MISCELLANEOUS ABBREVIATIONS, AND ABBREVIATED WATER-QUALITY UNITS

Conversion Factors

Multiply	By	To obtain
liter (L)	33.82	ounce
gram (g)	0.002205	pound
kilopascal	0.1450	pound per square inch
meter (m)	3.281	foot

Temperature can be converted to degrees Celsius ($^{\circ}\text{C}$) or degrees Fahrenheit ($^{\circ}\text{F}$) by the equations:

$$^{\circ}\text{C} = 5/9 (^{\circ}\text{F} - 32)$$

$$^{\circ}\text{F} = 9/5 (^{\circ}\text{C}) + 32.$$

Miscellaneous Abbreviations

cubic centimeter (cm^3)
inside diameter (i.d.)
mass to charge (m/z)
micrometer (μm)
milliabsorbance units (mAU)
milligram (mg)
millimeter (mm)
millimole (mM)
milliseconds (ms)
minute (min)
nanogram (ng)
volt (V)

Abbreviated Water-Quality Units

liter per minute (L/min)
microgram per liter ($\mu\text{g/L}$)
microliter (μL)
milligram per milliliter (mg/mL)
milliliter (mL)
milliliter per minute (mL/min)
nanogram per microliter ($\text{ng}/\mu\text{L}$)

Methods of Analysis by the U.S. Geological Survey Organic Geochemistry Research Group—Determination of Selected Herbicides and Their Degradation Products in Water Using Solid-Phase Extraction and Gas Chromatography/Mass Spectrometry

By J.L. Kish¹, E.M. Thurman¹, E.A. Scribner¹, and L.R. Zimmerman²

Abstract

A method for the extraction and analysis of eight herbicides and five degradation products using solid-phase extraction from natural water samples followed by gas chromatography/mass spectrometry is presented in this report. This method was developed for dimethenamid; flufenacet; fluometuron and its degradation products, demethylfluometuron (DMFM), 3-(trifluoromethyl)phenylurea (TFMPU), 3-(trifluoromethyl)aniline (TFMA); molinate; norflurazon and its degradation product, demethylnorflurazon; pendamethalin; the degradation product of prometryn, deisopropylprometryn; propanil; and trifluralin. The eight herbicides are used primarily in the southern United States where cotton, rice, and soybeans are produced. The exceptions are dimethenamid and flufenacet, which are used on corn in the Midwest.

Water samples received by the U.S. Geological Survey's Organic Geochemistry Research Group in Lawrence, Kansas, are filtered to remove suspended particulate matter and then passed through disposable solid-phase extraction columns containing octadecyl-bonded porous

silica (C-18) to extract the compounds. The herbicides and their degradation products are removed from the column by ethyl acetate elution. The eluate is evaporated under nitrogen, and components then are separated, identified, and quantified by injecting an aliquot of the concentrated extract into a high-resolution, fused-silica capillary column of a gas chromatograph/mass spectrometer under selected-ion mode.

Method detection limits ranged from 0.02 to 0.05 $\mu\text{g/L}$ for all compounds with the exception of TFMPU, which has a method detection limit of 0.32 $\mu\text{g/L}$. The mean absolute recovery is 107 percent. This method for the determination of herbicides and their degradation products is valuable for acquiring information about water quality and compound fate and transport in water.

INTRODUCTION

This report describes a method that uses solid-phase extraction (SPE) followed by gas chromatography/mass spectrometry (GC/MS) for the analysis of six herbicides and five degradation products, which are used primarily in the southern United States to enhance cotton, rice, and soybean production, and for two herbicides used in corn-growing areas of the Midwest. This method was developed by the U.S. Geological Survey (USGS) Organic Geochemistry Research Group in Lawrence, Kansas (Thurman and others,

¹U.S. Geological Survey, Lawrence, Kansas.

²University of Kansas Center for Research, Inc., Lawrence, Kansas.

1992; Meyer and others, 1993; Thurman and Mills, 1998; Zimmerman and Thurman, 1999; Thurman and others, 2000).

Reconnaissance studies in the Midwest have shown widespread detection of herbicides. Approximately three-fourths of all pre-emergent herbicides in the United States are applied to row crops in a 10-State area of the midwestern United States where herbicides frequently are detected in surface water (Thurman and others, 1991; Gianessi and Puffer, 1995). Because many herbicides and their metabolites are water soluble, they may leach into ground water (Hallberg, 1989; Thurman and others, 1991; Kolpin and others, 1993) as well as transported in surface runoff (Wauchope, 1978; Leonard, 1988).

Equally important to water quality is the application of herbicides to cotton and rice in the southern United States. Cotton and rice receive three to five times more herbicides per acre than do corn or soybeans. Cotton-growing areas of the United States extend from the East Coast (The Carolinas) to the Mississippi River Delta, the Texas High Plains, and the arid deserts of the Southwest (Arizona and California). These areas of the country have different climate, precipitation, and soil types, which result in different weed and insect pressures, as well as different runoff potentials; therefore, leaching patterns often are different. Because of these considerations, the types and amounts of herbicides applied may vary considerably throughout cotton-growing areas (Coupe and others, 1998; Thurman and others, 2000).

The analytical method described in this report was developed by the USGS to determine concentrations of the following herbicides and their degradation products: dimethenamid; flufenacet; fluometuron and its degradation products, demethylfluometuron (DMFM), 3-(trifluoromethyl)phenylurea (TFMPU), 3-(trifluoromethyl)aniline (TFMA); molinate; norflurazon and its degradation product, demethylnorflurazon; pendamethalin; the degradation product of prometryn, deisopropylprometryn; propanil; and trifluralin. The GC/MS method of analysis described in this report has been assigned the method code "O-2132-99." This unique code represents the automated method of analysis for organic compounds as it is described in this report and can be used to identify the method.

This report provides a detailed description of the method, including the apparatus, reagents, instrument calibration, and the SPE procedure required for sample analysis. Method detection limits, mean extraction

recoveries, and relative standard deviations for the GC/MS methods also are presented.

DETERMINATION OF HERBICIDES AND DEGRADATION PRODUCTS IN WATER USING SOLID-PHASE EXTRACTION AND GAS CHROMATOGRAPHY/MASS SPECTROMETRY

Method of Analysis (O-2132-99)

Scope and Application

The method described in this report is suitable for the determination of low concentrations (in micrograms per liter) of selected cotton, rice, soybean, and corn herbicides and their degradation products in natural water samples. Registry numbers and molecular weights are shown in table 1 for each herbicide and degradation product. This method is applicable to herbicides and their degradation products that are (1) efficiently partitioned from the water phase onto an octadecyl (C-18) silica phase that is chemically bonded to a solid silica matrix and (2) sufficiently volatile and thermally stable for gas chromatography. Suspended particulate matter is removed from the samples by filtration, so this method is suitable only for dissolved-phase herbicides and their degradation products.

Herbicides were selected for analysis because of their extensive use in the United States and their importance to studies being conducted by the USGS. The calibration range for the method is equivalent to concentrations from 0.05 to 5.0 $\mu\text{g/L}$ without dilution.

Summary of Method

Natural water samples are filtered at the collection site using glass-fiber filters with a 0.7- μm nominal pore diameter to remove suspended particulate matter. In the laboratory, filtered water samples are passed through a preconditioned C-18 column. The adsorbed compounds are removed from the C-18 with ethyl acetate. The eluate is evaporated further under nitrogen. The sample components are separated, identified, and quantified by injecting an aliquot of the concentrated extract into a high-resolution, fused-silica capillary column of a GC/MS system under selected-ion mode

Table 1. Herbicides and degradation products suitable for determination using method described, with registry numbers and molecular weights

[CAS, Chemical Abstract Service; DP, degradation product; AMID, amide; PU, phenylurea; TC, thiocarbamate; PDZ, pyridazinone; DNA, dinitroaniline; --, not available]

Herbicide or degradation product	Class	CAS registry number	Molecular weight (atomic mass units)
Deisopropylprometryn	DP	--	199
Demethylfluometuron (DMFM)	DP	--	218
Demethylnorflurazon	DP	--	289
Dimethenamid	AMID	87674-68-8	275
Flufenacet	AMID	142459-58-3	363
Fluometuron	PU	2164-17-2	232
Molinate	TC	2212-67-1	187
Norflurazon	PDZ	27314-13-2	303
Pendimethalin	DNA	40487-42-1	281
Propanil	AMID	709-98-8	218
3-(trifluoromethyl)aniline (TFMA)	DP	--	161
3-(trifluoromethyl)phenylurea (TFMPU)	DP	--	204
Trifluralin	DNA	1582-09-08	335

(SIM). Compounds eluting from the GC column are identified by comparing their measured ions and retention times to reference ions and retention times obtained by the measurement of control standards under the same conditions used for the water samples. The concentration of each identified compound is measured by relating the MS response of the quantitation ion produced by that compound to the MS response of the quantitation ion produced by a compound that is used as an internal standard. Surrogate compounds, whose concentrations are known in every sample, are measured with the same calibration procedure.

Interferences

Organic compounds having identical mass characteristic ions and GC retention times to those of the herbicides and their degradation products of interest may interfere.

Apparatus and Instrumentation

• *Analytical balances*—Capable of accurately weighing 0.0100 g ± 0.0001 g.

- *Autopipettes*—10-, 100-, and 200- μ L, variable-volume autopipettes with disposable tips (Rainin, Woburn, MA, or equivalent).
- *Tekmar six-position AutoTrace*—Automated SPE workstation (Tekmar-Dohrmann, Cincinnati, OH).
- *Extraction software*—Tekmar AutoTrace Extraction software, version 1.33 (Tekmar-Dohrmann, Cincinnati, OH).
- *Automated solvent evaporator*—The heat-bath temperature needs to be maintained at 45 °C and the nitrogen gas pressure at 103 kilopascals (Turbovap LV or equivalent, Zymark Inc., Hopkinton, MA).
- *Fused-silica capillary column*—Cross-linked methyl siloxane capillary column (12 m x 0.2 mm i.d., 0.33- μ m film thickness) (HP Ultra 1 or equivalent, Hewlett Packard, Wilmington, DE).
- *GC/MS benchtop system*—Hewlett Packard (Wilmington, DE) model 5890 series II Plus or equivalent GC with autoinjector connected to a Hewlett Packard model 5970 or equivalent MS detector.
- *GC conditions*—Oven, 60 °C (hold 1 min), then ramp to 200 °C at 6 °C/min, then 30 °C/min to 250 °C, and hold for 4 min;

injection port, 210 °C; carrier gas, helium; injection volume, 2 µL, splitless injection.

- MS conditions—Multiplier, 400 over autotune; detector, 280 °C; dwell time, 25 ms; mass ions monitored are listed in table 2 (see section on “Calibration”).
- *Moisture sieve and oxygen scrubber for carrier gas.*
- *Data system*—Computer and printer compatible with the GC/MS system used.
- *Software*—HP DOS ChemStation Software, 1030A version C (Hewlett Packard, Wilmington, DE) is used to acquire and store data and for peak integration.

Reagents and Consumable Materials

- *Sample bottles*—Baked 125-mL amber glass bottles (Boston round) with Teflon-lined lids.
- *Sample filters*—A 0.70-µm glass-fiber filter (Gilsen, Middleton, WI).
- *Reagent water*—Generated by purification through activated charcoal filtration and deionization with a high-purity, mixed-bed resin, followed by another activated charcoal filtration, and finally distillation in an autostill (Wheaton or equivalent, Millville, NJ).
- *Analytical standards*—Standards of the triazine and chloroacetanilide compounds.
- *SPE columns*—C-18 SEP-PAK Plus, containing 360 mg of 40-µm C-18 bonded-silica packing (Waters, Milford, MA).
- *Disposable snap-cap finish centrifuge tubes*—10 mL (Kimble or equivalent, Vineland, NJ).
- *Solvents*—Ethyl acetate [American Chemical Society (ACS) grade] and methanol [high-performance liquid chromatography grade (HPLC)].
- *Gas for evaporation*—Nitrogen, ultrapure grade.
- *Pasteur pipettes*—(Kimble or equivalent, Vineland, NJ).
- *0.1-mL autosampler vials*—Plastic vials with glass cone inserts (Wheaton, Millville, NJ).
- *GC carrier gas*—Helium, ultrapure grade.

Sampling Methods, Sample-Collection Equipment, and Cleaning Procedures

Following USGS protocol, surface-water samples are collected with a depth-integrating technique at three or more locations across each stream (Ward and Harr, 1990). The water samples from each site are composited in a single glass container or Teflon bottle.

Samples are withdrawn from the compositing container and filtered through a 0.70-µm glass-fiber filter using a peristaltic pump. Filters are leached with about 200 mL of sample prior to filtration of sample. The filtered material for analysis is collected in baked 125-mL amber glass bottles with Teflon-lined lids. Samples are chilled immediately and shipped to the laboratory within 3 days of collection. At the laboratory, samples are logged in, assigned identification numbers, and refrigerated at 4 °C until extracted and analyzed.

Standards

- *Stock standard solutions*—Obtain herbicides, degradation products, internal standard, and surrogate standards as pure materials from commercial vendors or chemical manufacturers. If pure materials are obtained, prepare solutions of 1 mg/mL by accurately weighing, to the nearest 0.1 mg, 50 mg of the pure material in a 50-mL volumetric flask and dilute with methanol. Transfer the stock solutions to clean 2-mL vials and store in a freezer. These solutions are stable for about 24 months.
- *Primary fortification standard*—A solution containing all the compounds for analysis at a known concentration. This standard is used in the preparation of control samples. Prepare a 1.23-ng/µL concentration, primary fortification standard solution by combining appropriate volumes of the individual stock solutions in a 1-L volumetric flask. Use adjustable autopipette to dispense an appropriate volume and dilute with methanol. Store at less than 4 °C. This solution is stable for about 24 months.
- *Internal standard solution*—Herbicides are ratioed against an internal standard to determine concentration. Prepare a solution of phenanthrene-*d*₁₀ in ethyl acetate at a concentration of 0.2 ng/µL. The internal standard may be purchased as a 100-µg/mL solution in methylene chloride. Dilute 800 µL in 4 L of ethyl acetate. Store at less than 4 °C. This solution is stable for about 12 months.
- *Surrogate solution*—A secondary compound used to determine the final concentration. The surrogate has the same function as the internal standard except it is spiked into the water. Prepare solutions of atrazine-*d*₅ and terbuthylazine in methanol at concentrations of 1.23 ng/µL. If

pure materials are obtained, prepare solutions of 1 mg/mL by accurately weighing, to the nearest 0.001 mg, 50 mg of the pure material in a 50-mL volumetric flask and dilute with methanol. Transfer the stock solutions to clean vials and store in a freezer. This solution is stable for about 24 months.

- **Calibration solution**—Prepare a calibration solution using the primary fortification standard in water that contains all target herbicides and metabolites. This solution is spiked at concentrations from 0.05 to 5.0 µg/L (0.05, 0.1, 0.2, 0.5, 1.0, 2.0, and 5.0 µg/L) and the surrogate standards, atrazine-d₅ and terbuthylazine, at a constant concentration of 1 µg/L.

Gas Chromatography/Mass Spectrometry Performance

Evaluation of Gas Chromatograph Performance

Gas chromatograph performance is evaluated by peak shape, internal standard response, and by comparison of response factors relative to response

factors obtained using a new capillary column and freshly prepared calibration solutions (a standard curve). An example of the separation and peak shape of cotton, rice, soybean, and corn herbicides and degradation products and internal standards is shown in a total ion chromatogram of a 1.0-µg/L standard solution (fig. 1).

If peak shape deteriorates or if response factors fail to meet the calibration criteria, the injection liner is changed, or maintenance on the capillary column is performed to bring the gas chromatograph into compliance. Part of the inlet end on the capillary column may be removed to restore performance. Specifically, poor peak shape and a loss of response for the herbicides and degradation products susceptible to loss on injection indicate a need for immediate action.

Evaluation of Mass Spectrometer Performance

Mass spectrometer performance is evaluated by assessing isotopic ratios, contamination, electron multiplier sensitivity, and abundance.




Figure 1. Chromatogram showing total ions of eight herbicides and five degradation products in 1.0-nanogram-per-microliter standard solution. Retention times shown above each peak correspond to compounds listed in table 2.

- Tune the mass spectrometer before each GC/MS sample set (approximately 43 injections or three extraction sample sets) using the procedure and software supplied by the manufacturer. Parameters in the tuning software are set to give ± 0.15 atomic mass unit resolution at masses 69, 219, and 502 in the spectrum of perfluorotributylamine (PFTBA). With the resolution of the 69 ion at 100-percent abundance, the mass 219 ion should be 35 ± 20 percent, and the mass 502 ion should be more than 3 percent relative abundance; however, their masses may vary depending on the mass spectrometer used. Check mass assignments to ensure accuracy to ± 0.15 atomic mass units and that mass peak widths measured at one-half the peak height range from about 0.50 to 0.60 atomic mass units.
- Also, during the tuning of the mass spectrometer, check the mass spectrometer for the presence of excessive water and air, which indicate leaks in the vacuum. If detected, locate and fix leaks.
- Initially adjust the electron multiplier of the mass spectrometer to ensure that the established reporting level for each selected compound can be achieved.
- Initial calibration data are entered into a computer spreadsheet (Microsoft Excel, Microsoft, Inc., Seattle, WA), and ratios are calculated for each quantitation ion relative to the internal standard (phenanthrene- d_{10}). Graphs are made from the GC/MS data by plotting the correlation curve with the phenanthrene- d_{10} ratios of a single ion on the x axis and the concentrations of the standards used on the y axis. Three graphs are made for each ion, one with concentrations ranging from 0 to 0.20 $\mu\text{g/L}$, another with concentrations ranging from 0.20 to 2.0 $\mu\text{g/L}$ (fig. 2), and the final curve ranging from 0.2 to 5.0 $\mu\text{g/L}$. This gives an intermediate and a low curve to keep the response linear. The final curve is a quadratic curve that is used to give high-end results. The low curve is plotted with one point at 0. The spreadsheet determines slopes, y intercepts, and correlation coefficient values (r^2) for the graphs.
- Initial calibration data acquired using a new capillary column and fresh calibration solutions are acceptable if the correlation coefficient (r^2) value for all curves is greater than or equal to 0.99 for all compounds.
- Subsequent daily response factors calculated for the majority of compounds need to agree within ± 20 percent of the mean response factor for the compounds analyzed. A response factor is equal to the area of the quantitation ion for the selected compound or surrogate divided by the area of the quantitation ion for the internal standard. Analyze at least two calibration standards with each sample set, one high calibration standard ranging from 0.5 to 2.0 $\mu\text{g/L}$ and one low standard ranging from 0.05 to 0.20 $\mu\text{g/L}$ to verify instrument response in each range.

Calibration

- Acquire initial calibration data by using a new capillary column and freshly prepared calibration solutions. Use these data in the subsequent evaluation of GC/MS performance.
- Acquire data for each calibration solution by injecting 2 μL of each solution into the GC/MS according to the conditions already described. Calculate the relative retention time (RRT_c) for each selected compound and the surrogate compounds in the calibration solution or in a sample as follows:

$$RRT_c = \frac{RT_c}{RT_i}, \quad (1)$$

where RT_c = uncorrected retention time of the quantitation ion of the selected compound or surrogate compound, and

RT_i = uncorrected retention time of the quantitation ion of the internal standard (phenanthrene- d_{10}).

See table 2 for an example of retention times and relative retention times.

Procedure

- *Sample preparation*—In the automation of sample extraction, the Autotrace workstation is used (Tekmar-Dohrman, Cincinnati, OH). Should an environmental sample contain less than 123 mL, distilled water is added to bring the volume to the required 123 mL. Any volume added is recorded. An extraction sample set consists of eight samples, one duplicate sample, two standard control samples (one high concentration and one low concentration), and one blank control sample. Each bottle is spiked with the surrogate standard solution, atrazine- d_5 and terbuthylazine, at a con-

Table 2. Retention times and relative retention times of selected herbicides, degradation products, surrogate compounds, and internal standard analyzed using method described

[min, minute; m/z, mass-to-charge ratio; --, not applicable]

Compound	Retention time (min)	Relative retention time (dimensionless)	Quantitation ion (m/z)	Confirmation ion 1 (m/z)	Confirmation ion 2 (m/z)
Herbicides or degradation products (in order of increasing retention time)					
3-(trifluoromethyl)aniline (TFMA)	4.73	0.25	161	142	114
Molinate	14.78	.77	126	55	187
Fluometuron	17.05	.89	72	232	--
Trifluralin	17.95	.93	264	306	--
3-(trifluoromethyl)phenylurea (TFMPU)	17.96	.93	161	204	142
Demethylfluometuron (DMFM)	18.40	.96	161	58	142
Deisopropylprometryn	20.26	1.05	184	199	157
Propanil	21.19	1.10	161	217	--
Dimethenamid	21.23	1.10	154	230	203
Flufenacet	23.43	1.22	151	211	123
Pendimethalin	24.20	1.26	252	281	162
Demethylnorflurazon	26.37	1.37	289	145	291
Norflurazon	26.60	1.38	303	145	102
Surrogate compounds					
Atrazine-d ₅	19.00	--	205	220	--
Terbutylazine	19.60	--	214	229	--
Internal standard					
Phenanthrene- d ₁₀	19.26	1.00	188	--	--

centration of 1.0 µg/L (100 µL of 1.23-ng/µL surrogate solution into 123 mL) with an autopipette.

- *Workstation preparation*—Before a sample set is extracted on the automated workstation, each port is flushed with 15 mL of methanol:water (4:1) and then again with distilled water. All SPE columns, test tubes, reagents, working solvents, surrogate spike, and samples then are loaded onto the instrument.
- *Conditioning the SPE columns*—The workstation conditions each SPE column by sequentially passing 1 mL methanol, 1 mL ethyl acetate, 1 mL methanol, and 3 mL distilled water through each column at a flow rate of 10 mL/min by positive pressure.
- *Loading the sample*—123 mL of sample is passed through the SPE column at a flow rate of 20 mL/min.
- *Eluting the SPE column*—Each SPE column is eluted with 4.0 mL ethyl acetate to remove the compounds at a flow rate of 2 mL/min.

- *Spiking of internal standard solution*—After all the samples in a set have been loaded and eluted, 500 µL of 0.2-ng/µL phenanthrene-d₁₀ solution is hand spiked into each eluate.
- *Separation of ethyl acetate and residual water*—Due to water in the SPE column that is eluted with the ethyl acetate, the ethyl acetate is transferred off of the aqueous phase into another tube. This is done manually using a pasteur pipette.
- *Evaporation*—The spiked eluate then is evaporated to approximately 75 µL under nitrogen in a water bath at 45 °C.
- *Transfer to vials*—Using a baked disposable Pasteur pipette, the eluate is withdrawn from the 10-mL glass centrifuge tube into a pipette, and transferred to an appropriately labeled GC autosampler vial containing a 0.1-mL insert for GC/MS analysis. The GC autosampler vial is capped and stored at less than 4 °C until analysis by GC/MS.




Figure 2. Examples of high, intermediate, and low calibration curves for dimethenamid. (A) low calibration curve for dimethenamid with a linear curve fit, (B) intermediate calibration curve for dimethenamid with a linear curve fit, and (C) high calibration curve for dimethenamid with a quadratic curve fit.

- *Sample analysis and data evaluation*—Ensure that GC/MS conditions for the analysis of the selected compounds in sample extracts are the same as those used in the analysis of the calibration solutions. Prior to the analysis of any sample extracts, ensure that the GC performance evaluation criteria have been met. Inject 2 μ L of the sample extract and acquire data using the GC/MS conditions described.

Calculation of Results

Qualitative Identification

- The expected retention time (*RT*) of the GC peak of the quantitation ion for the selected compound of interest needs to be within ± 6 seconds of the expected retention time that is based on the RRT_c obtained from the analysis of the internal stan-

dard. Calculate the expected retention time as follows:

$$RT = (RRT_c)(RT_i), \quad (2)$$

where RT = expected retention time of the selected compound or surrogate compound,

RRT_c = relative retention time of the selected compound or surrogate compound, and

RT_i = uncorrected retention time of the quantitation ion of the internal standard.

- Mass-spectral verification for each selected compound is done by comparing the relative integrated abundance values of the selected ions monitored with the relative integrated abundance values obtained from the standard control samples. The relative ratios of the ions need to be within ± 20 -percent of the relative ratios of those obtained in the absence of any obvious interferences. Slopes for compounds that interact with the GC inlet are modified to meet the ± 20 -percent criteria. As samples are analyzed, the heated inlet is coated with an involatile residue. Over time, this residue builds up and causes specific sorption of some analytes, which is an inhibitory factor. The charts monitor the surrogate-to-internal-standard ratio; if the ratio is not within 20 percent of its mean, then a new standard curve is analyzed. These modifications from 0 to 20 percent are referred to as the correction factor in equation 4.

Quantitation

- Calculate the volume of sample processed as follows:

$$DF = \left(\frac{123}{123 - V_{np}} \right) \left(\frac{123}{123 - V_a} \right), \quad (3)$$

where DF = dilution factor,

V_{np} = volume not pumped, in milliliters, and

V_a = volume added, in milliliters.

The dilution factor is incorporated into the calculation for determining final concentrations in samples.

- If a selected compound has passed the aforementioned qualitative identification criteria, calculate the concentration in the sample as follows:

$$C = \left[\left(\left(\frac{A_c}{A_i} \right) (m) + y \right) (DF) \right] x CF, \quad (4)$$

where C = concentration of the selected compound or surrogate compound in the sample, in micrograms per liter;

A_c = area of the quantitation ion for the selected compound or surrogate identified;

A_i = area of the quantitation ion for the internal standard;

m = slope of correlation curve between the selected compound and phenanthrene- d_{10} from the original calibration data;

y = y intercept of correlation curve between the selected compound and phenanthrene- d_{10} from the original calibration data;

DF = dilution factor as described in equation 3; and

CF = correction factor.

Reporting of Results

Concentrations of herbicides and degradation products are reported from 0.05 to 5.0 $\mu\text{g/L}$ without dilution. If the concentration is greater than 5.0 $\mu\text{g/L}$, the sample extract is diluted (volume increased to approximately 150 μL with eluting solvent) and re-analyzed. If the concentration is greater than 10 $\mu\text{g/L}$, the sample is re-extracted with a 1:10 dilution (sample:distilled water) and re-analyzed for those compounds with concentrations greater than 10 $\mu\text{g/L}$.

Method Performance

A reagent-water sample, a surface-water sample collected from Poison Creek in Valley County, Idaho, and a ground-water sample collected from a well in Valley County, Idaho, were used to test the performance of the GC/MS method. The surface- and ground-water samples were collected in 45-L carboys and were split into 123-mL samples. One set of seven samples was spiked with 0.2 $\mu\text{g/L}$ of each herbicide and degradation product of interest, and one set of samples was spiked with 1.0 $\mu\text{g/L}$ of each herbicide and degradation product of interest. In addition, unspiked samples of surface and ground water were extracted and analyzed to determine background

concentrations of the herbicides and degradation products. All subsamples were analyzed in one laboratory, the USGS Organic Geochemistry Research Laboratory in Lawrence, Kansas, using one GC/MS system. Each sample set was extracted and analyzed on different days intermittently between sample sets so that comparison of different matrices and concentrations included bias from day-to-day variation. Method recoveries from the analyses are listed in tables 3 through 5.

Mean recovery: Mean recovery in reagent-, surface-, and ground-water samples was determined by comparing the mean calculated concentration from seven replicate samples as shown in the "Quantitation" section to the spiked concentration (0.2 µg/L).

Corrections for background concentrations: Neither surface- nor ground-water samples required correction for background concentrations of compounds. All unfortified reagent-water samples also had no detections.

Method detection limits (MDL's): An MDL is defined as the minimum concentration of a compound that can be identified, measured, and reported with a 99-percent confidence that the compound concentration is greater than zero. MDL's were determined according to procedures outlined by the U.S. Environmental Protection Agency (1992). Seven replicate samples of reagent water fortified at 0.05 µg/L were analyzed to determine MDL's (table 6). Each sample was analyzed on different days during May and June 1998, so day-to-day variation was included.

The MDL was calculated using the following equation:

$$MDL = (S)(t_{(n-1, 1-\alpha=0.99)}), \quad (5)$$

where S = standard deviation of the replicate analyses, in micrograms per liter, at the spiked concentration;
 $t_{(n-1, 1-\alpha=0.99)}$ = Student's t value for the 99-percent confidence level with $n-1$ degrees of freedom (U.S. Environmental Protection Agency, 1992); and
 n = number of replicate analyses.

The estimated MDL's ranged from 0.02 to 0.32 µg/L (table 6). According to the U.S. Environmental Protection Agency (1992) procedure, the fortified concentrations should be no more than five times the estimated MDL. The fortified concentrations were within five times the MDL.

Standard deviation: Standard deviations for replicate concentrations in reagent-, surface-, and ground-water samples were determined for all compounds.

Absolute recovery: Absolute recovery of each compound was determined by comparing eight replicate samples processed using this procedure to solvent spiked with the compounds injected directly into the GC/MS at 1.0-µg/L concentrations. Compound quantitation-ion ratios to internal-standard target-ion ratios were compared. Absolute recoveries are listed in table 7. Absolute recovery is different from the mean recoveries listed in tables 4–6 in that mean recoveries are calculated from an initial calibration curve that is processed in the same manner as the samples, thus correcting for routine analyte losses. Absolute recoveries ranged from 74 to 123 percent, with a mean absolute recovery for all compounds of 107 percent (table 7).

CONCLUSIONS

This report presents a method for routine analysis of eight herbicides and five degradation products in natural water samples. From the data presented in this report, SPE and GC/MS combine to produce a sensitive and reliable method for the determination of low concentrations of the selected herbicides and degradation products in natural water samples.

Concentrations of herbicides and degradation products are reported from 0.05 to 5.0 µg/L without dilution. Method detection limits ranged from 0.02 to 0.05 µg/L for all compounds with the exception of TFMPU, which has a method detection limit of 0.32 µg/L. The mean absolute recovery for all compounds is 107 percent.

Table 3. Mean recoveries and relative standard deviations of selected herbicides and degradation products in reagent-water samples using gas chromatography/mass spectrometry

[µg/L, micrograms per liter; RSD, relative standard deviation]

Herbicide or degradation product	Reagent water			
	Seven replicate samples spiked at 0.2 µg/L		Seven replicate samples spiked at 1.0 µg/L	
	Mean recovery (percent)	RSD	Mean recovery (percent)	RSD
Deisopropylprometryn	120	0.12	103	0.10
Demethylfluometuron (DMFM)	111	.04	105	.15
Demethylnorflurazon	86	.04	79	.17
Dimethenamid	109	.06	98	.05
Flufenacet	108	.05	97	.12
Fluometuron	101	.02	92	.09
Molinate	101	.05	91	.13
Norflurazon	83	.05	78	.15
Pendimethalin	130	.08	98	.15
Propanil	111	.05	104	.13
3-(trifluoromethyl)aniline (TFMA)	90	.07	90	.17
3-(trifluoromethyl)phenylurea (TFMPU)	81	.11	75	.18
Trifluralin	113	.07	82	.11

Table 4. Mean recoveries and relative standard deviations of selected herbicides and degradation products in surface-water samples using gas chromatography/mass spectrometry

[µg/L, micrograms per liter; RSD, relative standard deviation]

Herbicide or degradation product	Surface water			
	Seven replicate samples spiked at 0.2 µg/L		Seven replicate samples spiked at 1.0 µg/L	
	Mean recovery (percent)	RSD	Mean recovery (percent)	RSD
Deisopropylprometryn	112	0.03	96	0.16
Demethylfluometuron (DMFM)	130	.06	125	.24
Demethylnorflurazon	95	.04	94	.18
Dimethenamid	109	.02	104	.08
Flufenacet	127	.03	112	.06
Fluometuron	109	.05	97	.14
Molinate	83	.02	87	.14
Norflurazon	90	.04	94	.19
Pendimethalin	90	.02	98	.08
Propanil	116	.04	111	.10
3-(trifluoromethyl)aniline (TFMA)	68	.02	80	.15
3-(trifluoromethyl)phenylurea (TFMPU)	95	.12	92	.37
Trifluralin	86	.01	78	.10

Table 5. Mean recoveries and relative standard deviations of selected herbicides and degradation products in ground-water samples using gas chromatography/mass spectrometry

[µg/L, micrograms per liter; RSD, relative standard deviation; --, not applicable]

Herbicide or degradation product	Ground water			
	Seven replicate samples spiked at 0.2 µg/L		Seven replicate samples spiked at 1.0 µg/L	
	Mean recovery (percent)	RSD	Mean recovery (percent)	RSD
Deisopropylprometryn	99	0.15	97	0.20
Demethylfluometuron (DMFM)	139	.15	147	.32
Demethylnorflurazon	63	.07	89	.19
Dimethenamid	88	.12	108	.12
Flufenacet	91	.11	112	.13
Fluometuron	100	.11	100	.17
Molinate	71	.09	88	.15
Norflurazon	63	.07	87	.18
Pendimethalin	90	.12	106	.18
Propanil	32	.07	123	.11
3-(trifluoromethyl)aniline (TFMA)	65	.05	79	.18
3-(trifluoromethyl)phenylurea (TFMPU)	--	--	114	.39
Trifluralin	90	.11	100	.19

Table 6. Mean observed concentrations, standard deviations, and method detection limits of selected herbicides and degradation products in reagent-water samples spiked at 0.05 microgram per liter

[µg/L, micrograms per liter; MDL, method detection limit]

Herbicide or degradation product	Mean observed concentration (µg/L)	Standard deviation (µg/L)	MDL (µg/L)
Deisopropylprometryn	0.05	0.01	0.02
Demethylfluometuron (DMFM)	.05	.02	.05
Demethylnorflurazon	.03	.01	.04
Dimethenamid	.04	.01	.03
Flufenacet	.05	.01	.04
Fluometuron	.04	.01	.03
Molinate	.03	.01	.03
Norflurazon	.03	.01	.03
Pendimethalin	.04	.01	.03
Propanil	.05	.01	.04
3-(trifluoromethyl)aniline (TFMA)	.05	.01	.03
3-(trifluoromethyl)phenylurea (TFMPU) ¹	.16	.11	.32
Trifluralin	.03	.01	.02

¹TFMPU was analyzed using a fortified concentration of 0.2 microgram per liter.

Table 7. Absolute recoveries for selected herbicides and degradation products spiked at concentrations of 1.0 microgram per liter

[%, percent]

Herbicide or degradation product	Absolute recovery (%)
Deisopropylprometryn	109
Demethylfluometuron (DMFM)	119
Demethylnorflurazon	127
Dimethenamid	107
Flufenacet	121
Fluometuron	116
Molinate	109
Norflurazon	123
Pendimethalin	78
Propanil	92
3-(trifluoromethyl)aniline (TFMA)	111
3-(trifluoromethyl)phenylurea (TFMPU)	102
Trifluralin	74
Mean absolute recovery	107

REFERENCES CITED

Coupe, R.H., Thurman, E.M., and Zimmerman, L.R., 1998, Relation of usage to the occurrence of cotton and rice herbicides in three streams of the Mississippi Delta: *Environmental Science & Technology*, v. 32, p. 3673–3680.

Gianessi, L.P., and Puffer, C.M., 1995, Herbicides use in the United States—national summary report, revised April 1991: Washington, D.C., Resources for the Future, 128 p.

Hallberg, G.R., 1989, Pesticide pollution in groundwater in the humid United States: *Agriculture, Ecosystems, and Environment*, v. 26, p. 299–367.

Kolpin, D.W., Burkart, M.R., and Thurman, E.M., 1993, Hydrogeologic, water-quality, and land-use data for the reconnaissance of herbicides and nitrate in near-surface aquifers of the midcontinental United States, 1991: U.S. Geological Survey Open-File Report 93–114, 61 p.

Leonard, R.A., 1988, Herbicides in surface waters, in Grover, R., ed., *Environmental chemistry of herbicides*, volume I: Boca Raton, Florida, CRC Press, p. 45–87.

Meyer, M.T., Mills, M.S., and Thurman, E.M., 1993, Automated solid-phase extraction of herbicides from water for gas chromatographic-mass spectrometric analysis: *Journal of Chromatography*, v. 629, p. 55–59.

Thurman, E.M., Bastian, K.C., and Mollhagan, Tony, 2000, Occurrence of cotton herbicides and insecticides in playa lakes of the High Plains of West Texas: *The Science of the Total Environment*, v. 248, p. 189–200.

Thurman, E.M., Goolsby, D.A., Meyer, M.T., and Kolpin, D.W., 1991, Herbicides in surface waters of the Midwestern United States—the effect of spring flush: *Environmental Science & Technology*, v. 25, p. 1794–1796.

Thurman, E.M., Goolsby, D.A., Meyer, M.T., Mills, M.S., Pomes, M.L., and Kolpin, D.W., 1992, A reconnaissance study of herbicides and their metabolites in surface water of the Midwestern United States using immunoassay and gas chromatography/mass spectrometry: *Environmental Science & Technology*, v. 26, p. 2440–2447.

Thurman, E.M., and Mills, M.S., 1998, *Solid-phase extraction*: New York, John Wiley & Sons, Inc., 344 p.

U.S. Environmental Protection Agency, 1992, Guidelines establishing test procedures for the analysis of pollutants (appendix B, part 136, Definition and procedures for the determination of the method detection limit): U.S. Code of Federal Regulations, Title 40, revised as of July 1, 1992, p. 565–567.

Ward, R.D., and Harr, C.A., 1990, Methods for collection and procession of surface-water and bed-material samples for physical and chemical analyses: U.S. Geological Survey Open-File Report 90–140, 71 p.

Wauchope, R.D., 1978, The pesticide content of surface water draining from agricultural fields—a review: *Journal of Environmental Quality*, v. 7, p. 459–472.

Zimmerman, L.R., and Thurman, E.M., 1999, Method of analysis by the U.S. Geological Survey Organic Geochemistry Research Group—determination of triazine and chloroacetanilide herbicides in water by solid-phase extraction and capillary-column gas chromatography/mass spectrometry with selected-ion monitoring: U.S. Geological Survey Open-File Report 98–634, 21 p.