



Chapter G

The determination of forty elements in geological and botanical samples by inductively coupled plasma-atomic emission spectrometry

By Paul H. Briggs

U.S. Geological Survey Open File Report 02-223-G

Analytical methods for chemical analysis of geologic and other materials, U.S. Geological Survey

Contents

Principle G-1

Interferences G-2

Scope G-1

Apparatus G-2

Reagents G-2

Safety precautions G-3

Procedure G-3

Digestion of samples G-3

ICP-AES analysis G-3

Calculation G-5

Assignment of uncertainty G-5

Bibliography G-18

Tables

1. Reporting limits for 40 elements by ICP-AES G-2
2. Operating conditions for determination of 40 elements by ICP-AES G-4
3. Analytical performance summary for 40 elements by ICP-AES G-5

The determination of forty elements in geological and botanical samples by inductively coupled plasma-atomic emission spectrometry

By Paul H. Briggs

Codes: E011

T01

Accepted: 10/06/93

Revised: 11/06/01

Principle

Forty major, minor, and trace elements are determined in geological and botanical materials by inductively coupled plasma-atomic emission spectrometry (ICP-AES). The sample is decomposed using a mixture of hydrochloric, nitric, perchloric, and hydrofluoric acids at low temperature (Crock and others, 1983). The digested sample is aspirated into the ICP-AES discharge where the elemental emission signals are measured simultaneously for the forty elements. Calibration is performed with digested rock reference materials and a series of multi-element solution standards (Lichte and others, 1987).

Interferences

ICP-AES interferences may result from spectral interferences, background shifts, and matrix effects (Thompson and Walsh, 1983). Inter-element correction factors and background corrections are applied using the proprietary data system software (Perkin-Elmer, 1999). Approximately 100 spectral interference corrections are performed on each sample. Further corrections are made when an element influences other elements beyond the "normal correction." It is common to not report an effected element due to the extraordinary interference of the affecting element. Matrix effects can generally be negated by proper matching of standard and sample.

Sample decomposition using this multi-acid digestion technique is suited to dissolve most rock types, soils, sediments and ashed botanical samples. As with any technique there are going to be exceptions. The method does not fully dissolve refractory or resistant minerals and some secondary minerals. Examples of incomplete digestion are as follows: Ba in barite, Cr in chromite, Ti in rutile, Sn in cassiterite, Al in corundum, and rare earth elements in a monazite. Samples that contain elements in concentrations beyond the linear working range where normally the element is a trace constituent, will have to be diluted (i.e., Mg in a dolomite, Pb in a galena, Zn in a sphalerite, Cu in a chalcopyrite). This dilution increases the lower reporting limits.

Scope

Analysis by ICP-AES for major, minor, and trace elements is useful for a variety of geochemical investigations. The lower and upper reporting limits used for this method are shown in table 1. Approximately 150 samples can be analyzed daily by the ICP-AES instrumentation.

Table 1.—Reporting limits for 40 elements by ICP-AES

<i>Element</i>	<i>Concentration range</i>		<i>Element</i>	<i>Concentration range</i>	
Aluminum, Al	0.005	50 %	Gallium, Ga	4	50,000 ppm
Calcium, Ca	0.005	50 %	Holmium, Ho	4	5,000 ppm
Iron, Fe	0.005	25 %	Lanthanum, La	2	50,000 ppm
Potassium, K	0.005	50 %	Lithium, Li	2	50,000 ppm
Magnesium, Mg	0.005	5 %	Manganese, Mn	4	50,000 ppm
Sodium, Na	0.005	50 %	Molybdenum, Mo	2	50,000 ppm
Phosphorous, P	0.005	50 %	Niobium, Nb	4	50,000 ppm
Titanium, Ti	0.005	25 %	Neodymium, Nd	4	50,000 ppm
Silver, Ag	2	10,000 ppm	Nickel, Ni	3	50,000 ppm
Arsenic, As	10	50,000 ppm	Lead, Pb	4	50,000 ppm
Gold, Au	8	50,000 ppm	Scandium, Sc	2	50,000 ppm
Barium, Ba	1	35,000 ppm	Tin, Sn	5	50,000 ppm
Beryllium, Be	1	5,000 ppm	Strontium, Sr	2	15,000 ppm
Bismuth, Bi	10	50,000 ppm	Tantalum, Ta	40	50,000 ppm
Cadmium, Cd	2	25,000 ppm	Thorium, Th	4	50,000 ppm
Cerium, Ce	4	50,000 ppm	Uranium, U	100	100,000 ppm
Cobalt, Co	2	25,000 ppm	Vanadium, V	2	30,000 ppm
Chromium, Cr	2	50,000 ppm	Yttrium, Y	2	25,000 ppm
Copper, Cu	2	15,000 ppm	Ytterbium, Yb	1	5,000 ppm
Europium, Eu	2	5,000 ppm	Zinc, Zn	2	15,000 ppm

Apparatus

- Perkin Elmer Optima 3000 simultaneous ICP-AES
- Hot plate with 50-position aluminum heating block
- 30-mL Teflon vessels with caps (Savillex)
- Acid dispensers (Labindustries)
- Repeating pipet (Eppendorf)
- Drying oven set at 95°C
- 13x100 mm disposable polypropylene test tubes with caps

Reagents

Hydrochloric acid, HCl reagent grade, 37 percent

Nitric acid, HNO₃ reagent grade, 70 percent

Hydrofluoric acid, HF reagent grade, 48 percent

Perchloric acid, HClO₄ reagent grade, 70 percent

Deionized water (DI)

One percent nitric acid solution: 10 mL 70 percent conc HNO₃ diluted in 1000 mL DI water

Aqua regia: three parts conc HCl and one part conc HNO₃; solution is not stable and must be prepared immediately before use

Lutetium internal standard (Lu): 500 µg Lu/mL, as Lu₂O₃ in 5 percent (v/v) HCl

Safety precautions

All laboratory personnel are required to wear safety glasses, rubber gloves, and lab coats when working in the laboratory. All sample digestions must be performed in a perchloric acid hood; the latter is washed down after each day's use. Refer to the *CHP* and *MSDS* for specific precautions, effects of overexposure, disposal, and first-aid treatment, for reagents used in the digestion procedure and operating the ICP-AES instrumentation. Calcium gluconate gel should be available in labs where HF is in use.

Procedure

Digestion of samples

1. Weigh 0.200 g sample into Teflon vessel.
2. Add 100 μ L Lu internal standard to each vessel with repeating pipet.
3. Rinse side walls of Teflon vessel with a minimum amount of DI water.
4. In the fume hood, slowly add 3 mL HCl and allow any reaction to subside.
5. Add 2 mL HNO₃, 1 mL HClO₄, and 2 mL HF. Place sample solution vessel on hot plate with aluminum heat block at a controlled temperature of 110°C in a perchloric acid fume hood.
6. Evaporate sample solution to complete dryness on hot plate (usually overnight).
7. Remove from hot plate, cool to touch and add 1 mL HClO₄ and 2 to 3 mL DI water.
8. Return to hot plate and evaporate to hard dryness. The temperature of the hot plate is increased to 160°C. This step usually takes a few hours.
9. Remove dried sample from hot plate and cool.
10. Add 1.0 mL *aqua regia* with repeating pipet and let react for 15 min.
11. Add 9.0 mL 1 percent HNO₃ and thread screw cap tightly on vessel. Place vessel in drying oven for 1 hour at a controlled temperature of 95°C.
12. Remove sample solution and cool. Transfer solution into labeled disposable polypropylene test tube and cap with test tube cap.
13. Analyze sample solution by ICP-AES.

ICP-AES analysis

The ICP-AES instrument is calibrated at the start of each day using a well characterized in-house geological reference material (Pikes Peak granite) and nine single and or multi-element solution standards. The major and trace elements are determined by comparing the element intensities obtained from the standards to those obtained from the samples. There are three method preparation blanks digested with each sample set. A blank subtraction is performed to negate the effect of the reagents. Table 2 shows instrumental operating conditions and element wavelengths for this method.

Table 2.—Operating conditions for determination of 40 elements by ICP-AES

RF power to the torch 1300W
 Plasma argon flow rate 15 L/min
 Auxiliary argon flow rate 0.5 L/min
 Nebulizer argon flow rate 0.75 L/min
 Sample pump rate 1 mL/min
 Observation height above
 load coil 15 mm
 Equilibration time 15s
 Auto integration time 0.5 to 20 s
 Nebulizer Modified Babington

<i>Element</i>	<i>Wavelength, nm</i>	<i>Element</i>	<i>Wavelength, nm</i>
Ag	328.0	Mg	279.0
Al	308.2	Mn	257.6
As	188.9	Mo	202.0
Au	242.7	Na	589.5
Ba	413.1	Nb	292.7
Be	313.0	Nd	406.1
Bi	223.0	Ni	231.6
Ca	430.2	P	213.6
Cd	226.5	Pb	220.3
Ce	413.7	Sc	424.6
Co	228.6	Sn	189.9
Cr	267.7	Sr	460.7
Cu	324.7	Ta	240.0
Eu	381.9	Th	350.9
Fe	273.9	Ti	223.0
Ga	294.3	U	385.9
Ho	345.6	V	292.4
K	766.4	Y	371.0
La	408.6	Yb	328.9
Li	670.7	Zn	213.8

Most elements in each sample data set are normalized (i.e. standardized) using well defined in-house reference materials (RM) that have undergone the sample digestion process. The normalized sample set is quality control checked by an independent, established RM that has undergone the sample digestion process also. The RM QC check is accepted if the recovery is within the upper and lower control limits of three times the standard deviation of the certified value, and the concentrations of the elements are >10 times the lower reporting limit. If the standardization is not accepted the sample set and in-house standards are redigested and reanalyzed.

Calculation

$$\text{Element concentration} = \text{IRU}/\text{IRS} \times \text{CONSTD} \times \text{WT SOLN}/\text{WT SAMPLE} + \text{IEC}$$

where

IRU = intensity of element/intensity of Lu
 IRS = intensity of calibration standard/intensity of Lu
 CONSTD = conc of calibration standard
 WT SOLN = weight of final solution
 WT SAMPLE = weight of sample
 IEC = interelement corrections

Assignment of Uncertainty

The analytical results for the selected reference materials, duplicate samples, and method blanks are summarized in table 3. Please note: Some *pv* data has been converted from the oxide using the conversion factors in appendix A, table A1. Proposed values taken from the following published reference materials compilations: GSD-3, GSD-5 AND AGV-1 (Govindaraju, 1994), SRM 2709, 2710 and 2711 (NIST Certificate of Analysis, 1992), TILL-4 (CCRMP, 1995), and SRM 1571 *pv* values converted to ash weight from ash% of 9.2 (NBS Certificate of Analysis, 1977) .

Table 3.—Analytical performance summary for 40 elements by ICP-AES See page ix of the introduction to this Methods Manual for an explanation of the abbreviations used in the analytical performance summary tables.

Reference	Description	n	mean	s	pv	% RSD	% R
Aluminum, Al (%)							
SRM 1571	orchard leaves	7	0.44	0.01	---	2	---
GSD-3	stream sediment	35	6.11	0.1	6.4	2	95
SRM 2710	soil	5	6.33	0.08	6.44	1	98
SRM 2711	soil	46	6.53	0.02	6.53	0.3	100
TILL-4	till	17	7.32	0.2	7.62	2	96
SRM 2709	soil	46	7.42	0.2	7.5	2	99
GSD-5	stream sediment	40	8.08	0.2	8.19	2	99
AGV-1	andesite	6	9.13	0.2	9.07	2	101
Arsenic, As (ppm)							
AGV-1	andesite	6	<10	---	0.88	---	---
GSD-3	stream sediment	35	17.6	3	17.6	19	100
SRM 2709	soil	46	17.7	4	17.7	25	100
GSD-5	stream sediment	40	76.4	5	75	7	102
SRM 1571	orchard leaves	7	104	7	109	7	95
SRM 2711	soil	46	106	5	105	5	101
TILL-4	till	17	112	11	111	10	101
SRM 2710	soil	5	634	27	626	4	101

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Barium, Ba (ppm)							
TILL-4	till	17	407	15	395	4	103
SRM 1571	orchard leaves	7	456	4	478	1	95
GSD-5	stream sediment	40	464	24	440	5	106
GSD-3	stream sediment	35	626	19	615	3	102
SRM 2710	soil	5	707	7	707	1	100
SRM 2711	soil	46	723	7	726	1	100
SRM 2709	soil	46	939	20	968	2	97
AGV-1	andesite	6	1210	35	1226	3	99
Beryllium, Be (ppm)							
GSD-3	stream sediment	35	1.2	0.07	1.5	6	81
AGV-1	andesite	6	2.0	0.08	2.1	4	95
SRM 1571	orchard leaves	7	<2	---	---	---	---
SRM 2711	soil	46	2.2	0.1	---	5	---
GSD-5	stream sediment	40	2.2	0.1	2.3	5	98
SRM 2710	soil	5	2.5	0.07	---	3	---
TILL-4	till	17	3.0	0.1	3.7	4	82
SRM 2709	soil	46	3.7	0.6	---	16	---
Bismuth, Bi (ppm)							
SRM 2709	soil	46	<10	---	---	---	---
SRM 2710	soil	5	<10	---	---	---	---
SRM 2711	soil	46	<10	---	---	---	---
TILL-4	till	17	<10	---	40	---	---
GSD-3	stream sediment	35	<10	---	0.79	---	---
AGV-1	andesite	6	<10	---	0.057	---	---
GSD-5	stream sediment	40	<10	---	2.4	---	---
SRM 1571	orchard leaves	7	<20	---	---	---	---
Cadmium, Cd (ppm)							
TILL-4	till	17	<2	---	---	---	---
SRM 2709	soil	46	<2	---	0.38	---	---
GSD-5	stream sediment	40	<2	---	0.82	---	---
AGV-1	andesite	6	<2	---	0.069	---	---
GSD-3	stream sediment	35	<2	---	0.1	---	---
SRM 1571	orchard leaves	7	<4	---	1.2	---	---
SRM 2710	soil	5	20.3	2	21.8	9	93
SRM 2711	soil	46	40.1	3	41.7	7	96

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Calcium, Ca (%)							
GSD-3	stream sediment	35	0.127	0.0050	0.18	4	71
TILL-4	till	17	0.880	0.02	0.89	2	99
SRM 2710	soil	5	1.27	0.006	1.25	0.5	102
SRM 2709	soil	46	1.93	0.05	1.89	3	102
SRM 2711	soil	46	2.87	0.05	2.88	2	100
AGV-1	andesite	6	3.54	0.04	3.53	1	100
GSD-5	stream sediment	40	3.86	0.08	3.79	2	102
SRM 1571	orchard leaves	7	22.5	0.6	22.72	3	99
Cerium, Ce (ppm)							
SRM 1571	orchard leaves	7	11.2	0.8	---	7	---
SRM 2709	soil	46	45.9	2	42	5	109
SRM 2710	soil	5	59.4	3	57	5	104
GSD-3	stream sediment	35	67.4	4	64	6	105
AGV-1	andesite	6	69.3	3	67	4	103
SRM 2711	soil	46	74.6	4	69	6	108
TILL-4	till	17	81.1	5	78	7	104
GSD-5	stream sediment	40	93.2	6	89	6	105
Chromium, Cr (ppm)							
AGV-1	andesite	6	3.3	1	10.1	44	33
SRM 1571	orchard leaves	7	29.7	3	28.3	11	105
TILL-4	till	17	40.9	3	53	8	77
SRM 2711	soil	46	44.3	2	47	5	94
SRM 2710	soil	5	55.7	32	39	58	143
GSD-5	stream sediment	40	71.6	3	70	4	102
GSD-3	stream sediment	35	91.9	4	87	5	106
SRM 2709	soil	46	126	5	130	4	97
Cobalt, Co (ppm)							
SRM 1571	orchard leaves	7	<2	---	2.17	---	---
SRM 2710	soil	5	2.2	5	10	213	22
TILL-4	till	17	8.2	2	8	29	102
SRM 2711	soil	46	9.5	2	10	19	95
SRM 2709	soil	46	14.0	2	13.4	14	104
GSD-3	stream sediment	35	15.7	3	11.7	19	134
AGV-1	andesite	6	17.1	3	15.3	19	112
GSD-5	stream sediment	40	20.2	2	18.9	11	107

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Copper Cu (ppm)							
SRM 2709	soil	46	32.7	1	34.6	3	94
AGV-1	andesite	6	57.5	1	60	2	96
SRM 2711	soil	46	115	2	114	2	101
GSD-5	stream sediment	40	141	4	137	3	103
SRM 1571	orchard leaves	7	145	9	130	6	112
GSD-3	stream sediment	35	178	6	177	3	100
TILL-4	till	17	246	10	237	4	104
SRM 2710	soil	5	2900	25	2950	1	98
Europium, Eu (ppm)							
SRM 2709	soil	46	<1	---	0.9	7	99
SRM 2711	soil	46	1.1	0.06	1.1	6	99
SRM 2710	soil	5	1.1	0.06	1	6	110
TILL-4	till	17	1.2	0.07	<1	6	---
GSD-3	stream sediment	35	1.2	0.08	1.3	7	89
GSD-5	stream sediment	40	1.2	0.08	1.4	7	89
AGV-1	andesite	6	1.7	0.1	1.64	6	107
SRM 1571	orchard leaves	7	<2	---	---	9	---
Galium, Ga (ppm)							
SRM 1571	orchard leaves	7	<8	---	0.87	---	---
GSD-3	stream sediment	35	14.8	1	15.9	8	93
SRM 2709	soil	46	15.2	1	14	8	108
SRM 2711	soil	46	15.4	1	15	6	103
TILL-4	till	17	16.8	1	---	6	---
GSD-5	stream sediment	40	20.5	2	20.3	8	101
AGV-1	andesite	6	20.6	0.7	20	3	103
SRM 2710	soil	5	38.3	1	34	4	113
Gold, Au (ppm)							
SRM 2709	soil	46	<8	---	---	---	---
SRM 2711	soil	46	<8	---	0.03	---	---
SRM 2710	soil	5	<8	---	0.6	---	---
GSD-3	stream sediment	35	<8	---	---	---	---
AGV-1	andesite	6	<8	---	---	---	---
TILL-4	till	17	<8	---	---	---	---
GSD-5	stream sediment	40	<8	---	---	---	---
SRM 1571	orchard leaves	7	<16	---	---	---	---

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Holmium, Ho (ppm)							
SRM 2709	soil	46	<4	---	0.54	---	---
SRM 2710	soil	5	<4	---	0.6	---	---
SRM 2711	soil	46	<4	---	---	---	---
TILL-4	till	17	<4	---	---	---	---
GSD-3	stream sediment	35	<4	---	0.9	---	---
AGV-1	andesite	6	<4	---	0.67	---	---
GSD-5	stream sediment	40	<4	---	1.1	---	---
SRM 1571	orchard leaves	7	<8	---	---	---	---
Iron, Fe (%)							
SRM 1571	orchard leaves	7	0.3	0.02	0.326	7	104
SRM 2711	soil	46	2.9	0.05	2.89	2	99
SRM 2710	soil	5	3.3	0.06	3.38	2	99
SRM 2709	soil	46	3.5	0.07	3.5	2	101
TILL-4	till	17	4.0	0.1	3.94	2	101
GSD-5	stream sediment	40	4.1	0.09	4.08	2	100
GSD-3	stream sediment	35	4.7	0.1	4.57	3	103
AGV-1	andesite	6	4.7	0.08	4.74	2	100
Lanthanum, La (ppm)							
SRM 1571	orchard leaves	7	14.4	0.9	---	6	---
SRM 2709	soil	46	24.5	1	23	5	106
SRM 2710	soil	5	32.1	2	34	5	94
GSD-3	stream sediment	35	37.9	3	39	7	97
SRM 2711	soil	46	40.6	3	40	6	101
AGV-1	andesite	6	42.3	2	38	4	111
TILL-4	till	17	43.6	3	41	8	106
GSD-5	stream sediment	40	44.0	3	46	7	96
Lead, Pb (ppm)							
SRM 2709	soil	46	16.9	2	18.9	12	89
AGV-1	andesite	6	36.1	1	36	4	100
GSD-3	stream sediment	35	40.5	3	40	7	101
TILL-4	till	17	48.1	3	50	7	96
GSD-5	stream sediment	40	112	7	112	6	100
SRM 1571	orchard leaves	7	516	56	489	11	105
SRM 2711	soil	46	1110	76	1162	7	95
SRM 2710	soil	5	5360	471	5532	9	97

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Lithium, Li (ppm)							
SRM 1571	orchard leaves	7	6.4	0.2	6.52	2	99
AGV-1	andesite	6	10.6	0.9	12	8	88
SRM 2711	soil	46	27.3	0.8	---	3	---
TILL-4	till	17	30.0	1	30	4	100
GSD-3	stream sediment	35	32.5	0.9	33	3	99
SRM 2710	soil	5	40.2	0.7	---	2	---
GSD-5	stream sediment	40	47.4	1	45	3	105
SRM 2709	soil	46	55.9	1	---	2	---
Magnesium, Mg (%)							
GSD-3	stream sediment	35	0.39	0.007	0.41	2	94
GSD-5	stream sediment	40	0.57	0.009	0.58	2	98
TILL-4	till	17	0.73	0.01	0.76	2	96
SRM 2710	soil	5	0.85	0.006	0.853	1	100
AGV-1	andesite	6	0.91	0.006	0.92	1	99
SRM 2711	soil	46	1.05	0.01	1.05	1	100
SRM 2709	soil	46	1.52	0.02	1.51	2	101
SRM 1571	orchard leaves	7	6.69	0.1	6.74	2	99
Manganese, Mn (ppm)							
GSD-3	stream sediment	35	412	10	387	2	106
TILL-4	till	17	480	7	465	1	103
SRM 2709	soil	46	548	11	538	2	102
SRM 2711	soil	46	638	9	638	1	100
AGV-1	andesite	6	734	12	774	2	95
SRM 1571	orchard leaves	7	977	19	989	2	99
GSD-5	stream sediment	40	1150	25	1162	2	99
SRM 2710	soil	5	9500	163	10100	2	94
Molybdenum, Mo (ppm)							
GSD-5	stream sediment	40	<2	---	1.2	---	---
SRM 2709	soil	46	<2	---	2	---	---
AGV-1	andesite	6	<2	---	2.7	---	---
SRM 2711	soil	46	<2	---	1.6	---	---
SRM 1571	orchard leaves	7	<4	---	3.26	---	---
TILL-4	till	17	14.4	0.8	16	5	90
SRM 2710	soil	5	16.8	2	19	12	88
GSD-3	stream sediment	35	91.5	2	92	2	100

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Neodymium, Nd (ppm)							
SRM 1571	orchard leaves	7	<8	---	---	---	---
SRM 2710	soil	5	18.3	1	---	5	---
SRM 2709	soil	46	19.3	1	---	5	---
GSD-3	stream sediment	35	26.6	1	30	5	89
AGV-1	andesite	6	31.0	1	33	5	94
SRM 2711	soil	46	31.6	2	31	5	102
TILL-4	till	17	32.3	2	30	7	108
GSD-5	stream sediment	40	33.1	2	35	6	95
Nickel, Ni (ppm)							
SRM 1571	orchard leaves	7	16.0	0.3	14.1	2	113
SRM 2710	soil	5	16.4	0.3	14.3	2	114
AGV-1	andesite	6	16.8	0.2	16	1	105
TILL-4	till	17	19.4	0.7	17	3	114
SRM 2711	soil	46	21.5	0.7	20.6	3	105
GSD-3	stream sediment	35	28.5	1	25.6	4	111
GSD-5	stream sediment	40	37.2	1	34	4	109
SRM 2709	soil	46	80.6	3	88	4	92
Niobium, Nb (ppm)							
SRM 2710	soil	5	<4	---	---	---	---
SRM 1571	orchard leaves	7	<8	---	---	---	---
GSD-3	stream sediment	35	16.1	6	16	39	101
SRM 2709	soil	46	17.0	6	---	38	---
AGV-1	andesite	6	20.5	3	15	15	137
GSD-5	stream sediment	40	20.8	7	19	36	110
SRM 2711	soil	46	22.6	7	---	33	---
TILL-4	till	17	23.4	11	15	46	156
Phosphorus, P (%)							
GSD-5	stream sediment	40	0.062	0.003	0.06	5	104
GSD-3	stream sediment	35	0.063	0.003	0.06	5	104
SRM 2709	soil	46	0.065	0.003	0.062	4	105
SRM 2711	soil	46	0.085	0.004	0.086	4	99
TILL-4	till	17	0.095	0.002	0.09	2	106
SRM 2710	soil	5	0.118	0.004	0.106	4	111
AGV-1	andesite	6	0.224	0.009	0.22	4	102
SRM 1571	orchard leaves	7	2.50	0.1	2.28	5	110

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Potassium, K (%)							
GSD-5	stream sediment	40	1.8	0.05	1.69	3	105
GSD-3	stream sediment	35	1.9	0.05	2.02	3	96
SRM 2709	soil	46	2.0	0.1	2.03	6	97
SRM 2710	soil	5	2.1	0.03	2.11	2	98
SRM 2711	soil	46	2.4	0.1	2.45	5	98
AGV-1	andesite	6	2.5	0.09	2.42	3	102
TILL-4	till	17	2.5	0.07	2.7	3	94
SRM 1571	orchard leaves	7	13.8	1	15.98	9	86
Scandium, Sc (ppm)							
SRM 1571	orchard leaves	7	<4	---	---	---	---
SRM 2710	soil	5	9.6	0.3	8.7	3	111
SRM 2711	soil	46	9.7	0.3	9	3	108
TILL-4	till	17	11.2	0.2	10	2	112
SRM 2709	soil	46	12.3	0.3	12	3	102
AGV-1	andesite	6	12.8	0.2	12.2	2	105
GSD-3	stream sediment	35	13.8	0.4	14.3	3	97
GSD-5	stream sediment	40	14.4	0.3	14.5	2	99
Silver, Ag (ppm)							
TILL-4	till	17	<2	---	---	---	---
SRM 2709	soil	46	<2	---	0.41	---	---
GSD-5	stream sediment	40	<2	---	0.36	---	---
AGV-1	andesite	6	<2	---	0.078	---	---
GSD-3	stream sediment	35	<2	---	0.59	---	---
SRM 1571	orchard leaves	7	<4	---	---	---	---
SRM 2711	soil	46	4.6	0.2	4.63	4	99
SRM 2710	soil	5	35.1	1	35.3	4	100
Sodium, Na (%)							
SRM 1571	orchard leaves	7	0.092	0.005	0.089	5	103
GSD-3	stream sediment	35	0.20	0.006	0.24	3	85
GSD-5	stream sediment	40	0.27	0.008	0.3	3	91
SRM 2710	soil	5	1.12	0.03	1.14	3	98
SRM 2709	soil	46	1.16	0.03	1.16	3	100
SRM 2711	soil	46	1.17	0.03	1.14	2	103
TILL-4	till	17	1.81	0.06	1.82	3	99
AGV-1	andesite	6	3.13	0.04	3.2	1	98

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Strontium, Sr (ppm)							
GSD-3	stream sediment	35	87.1	3	90	3	97
TILL-4	till	17	116	2	109	2	106
GSD-5	stream sediment	40	210	4	204	2	103
SRM 2709	soil	46	231	5	231	2	100
SRM 2711	soil	46	247	4	245.3	2	101
SRM 2710	soil	5	325	3	240	1	135
SRM 1571	orchard leaves	7	438	10	402	2	109
AGV-1	andesite	6	669	4	662	1	101
Tantalum, Ta (ppm)							
SRM 2709	soil	46	<40	---	---	---	---
SRM 2710	soil	5	<40	---	---	---	---
SRM 2711	soil	46	<40	---	---	---	---
TILL-4	till	17	<40	---	1.6	---	---
GSD-3	stream sediment	35	<40	---	1.08	---	---
AGV-1	andesite	6	<40	---	0.9	---	---
GSD-5	stream sediment	40	<40	---	1.4	---	---
SRM 1571	orchard leaves	7	<80	---	---	---	---
Thorium, Th (ppm)							
AGV-1	andesite	6	5.4	0.5	6.5	10	82
SRM 1571	orchard leaves	7	<8	---	0.7	---	---
GSD-3	stream sediment	35	9.0	0.7	9.2	8	98
SRM 2709	soil	46	11.5	0.9	10.8	8	107
SRM 2711	soil	46	14.2	1	13.6	8	104
SRM 2710	soil	5	14.4	2	13.5	11	107
GSD-5	stream sediment	40	15.7	1	15.2	7	104
TILL-4	till	17	17.4	2	17.4	11	100
Tin Sn (ppm)							
SRM 2711	soil	46	<5	---	---	---	---
SRM 2709	soil	46	<5	---	---	---	---
AGV-1	andesite	6	<5	---	4.2	---	---
GSD-3	stream sediment	35	<5	---	3.4	---	---
GSD-5	stream sediment	40	<5	---	4.6	---	---
SRM 2710	soil	5	6.8	2	---	26	---
TILL-4	till	17	9.8	3	---	27	---
SRM 1571	orchard leaves	7	<10	---	---	---	---

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Titanium, Ti (%)							
SRM 1571	orchard leaves	7	<.1	---	---	---	---
SRM 2710	soil	5	0.13	0.01	0.283	4	45
SRM 2711	soil	46	0.28	0.01	0.306	3	90
SRM 2709	soil	46	0.34	0.01	0.342	2	99
TILL-4	till	17	0.42	0.02	0.49	5	87
GSD-5	stream sediment	40	0.45	0.03	0.54	6	83
GSD-3	stream sediment	35	0.50	0.05	0.64	10	79
AGV-1	andesite	6	0.67	0.03	0.63	5	107
Uranium, U (ppm)							
SRM 2709	soil	46	<100	---	3	---	---
SRM 2710	soil	5	<100	---	25	---	---
SRM 2711	soil	46	<100	---	---	---	---
TILL-4	till	17	<100	---	5	---	---
GSD-3	stream sediment	35	<100	---	1.86	---	---
AGV-1	andesite	6	<100	---	1.92	---	---
GSD-5	stream sediment	40	<100	---	2.6	---	---
SRM 1571	orchard leaves	7	<200	---	0.32	---	---
Vanadium, V (ppm)							
SRM 1571	orchard leaves	7	5.6	0.2	---	3	---
TILL-4	till	17	69.1	1	67	1	103
SRM 2710	soil	5	73.8	1	76.6	2	96
SRM 2711	soil	46	81.2	1	81.6	2	100
GSD-5	stream sediment	40	105.6	2	109	2	97
SRM 2709	soil	46	112.8	1	112	1	101
GSD-3	stream sediment	35	119.3	2	120	2	99
AGV-1	andesite	6	122.1	0.8	121	1	101
Ytterbium, Yb (ppm)							
AGV-1	andesite	6	1.4	0.1	1.72	8	80
GSD-3	stream sediment	35	1.6	0.3	2.6	19	61
SRM 2709	soil	46	1.7	0.1	1.6	7	103
TILL-4	till	17	1.7	0.2	3.4	13	50
SRM 1571	orchard leaves	7	<2	---	---	---	---
GSD-5	stream sediment	40	2.1	0.2	2.9	10	72
SRM 2710	soil	5	2.2	0.08	1.3	4	172
SRM 2711	soil	46	2.7	0.2	2.7	6	100

**Table 3.—Analytical performance summary for 40 elements by ICP-AES---
Continued.**

Reference	Description	n	mean	s	pv	% RSD	% R
Yttrium, Y (ppm)							
SRM 1571	orchard leaves	7	5.3	0.5	---	9	---
GSD-3	stream sediment	35	16.1	2	22	12	73
SRM 2709	soil	46	16.4	0.7	18	4	91
TILL-4	till	17	17.9	0.8	33	4	54
AGV-1	andesite	6	19.3	0.8	20	4	97
GSD-5	stream sediment	40	20.4	1	26	6	79
SRM 2710	soil	5	20.9	0.7	23	3	91
SRM 2711	soil	46	26.4	1	25	5	106
Zinc, Zn (ppm)							
GSD-3	stream sediment	35	47.2	2	52	4	91
TILL-4	till	17	65.0	2	70	4	93
AGV-1	andesite	6	88.8	2	88	3	101
SRM 2709	soil	46	103	3	106	3	97
GSD-5	stream sediment	40	251	7	243	3	103
SRM 1571	orchard leaves	7	300	8	272	3	110
SRM 2711	soil	46	346	8	350.4	2	99
SRM 2710	soil	5	6540	405	6952	6	94

Table 3.—Continued--Duplicate samples results

Duplicate samples	k	n	mean	s	%RSD	min	max	<total	<pairs
Ag ppm	6	2	27.0	0.9	3	2	120	150	75
Al %	81	2	5.94	0.1	2	0.05	11	0	0
As ppm	29	2	127	4	3	10	2800	104	52
Au ppm	0	2	-	-	-	-	-	162	81
Ba ppm	81	2	625	14	2	24	1200	0	0
Be ppm	65	2	3.78	0.1	2	1	110	32	16
Bi ppm	0	2	-	-	-	-	-	162	81
Ca %	81	2	2.46	0.1	4	0.08	29	0	0
Cd ppm	26	2	28.7	1	4	2	190	110	55
Ce ppm	73	2	78.8	6	7	18	180	16	8
Co ppm	69	2	17.6	2	9	1	180	24	12
Cr ppm	80	2	181	7	4	1	8500	2	1
Cu ppm	81	2	70.5	4	6	4	1700	0	0
Eu ppm	2	2	2.10	0.02	0.9	2	2	158	79
Fe %	81	2	3.21	0.05	2	0.06	17	0	0
Ga ppm	70	2	16.0	0.5	3	5	26	22	11
Ho ppm	0	2	-	-	-	-	-	162	81
K %	81	2	3.45	0.2	5	0.1	25	0	0
La ppm	78	2	37.9	3	7	2	96	6	3
Li ppm	77	2	32.9	1	3	3	100	8	4
Mg %	81	2	1.24	0.02	2	0.09	5	0	0
Mn ppm	81	2	1069	51	5	49	7700	0	0
Mo ppm	25	2	25.6	2	6	2	210	112	56
Na %	79	2	1.01	0.03	3	0.02	3	4	2
Nb ppm	66	2	15.8	2	15	5	56	30	15
Nd ppm	74	2	34.0	2	7	8	83	14	7
Ni ppm	81	2	69.2	3	5	3	1700	0	0
P %	81	2	0.444	0.01	3	0.009	7	0	0
Pb ppm	73	2	497	38	8	6	9700	16	8
Sc ppm	69	2	11.5	0.2	2	2	25	24	12
Sn ppm	12	2	11.8	3	23	5	56	138	69
Sr ppm	81	2	217	3	2	26	910	0	0
Ta ppm	0	2	-	-	-	-	-	162	81
Th ppm	71	2	11.3	1	8	4	27	20	10
Ti %	74	2	0.304	0.01	4	0.05	0.69	14	7
U ppm	0	2	-	-	-	-	-	162	81
V ppm	79	2	74.5	1	1	3	160	4	2
Y ppm	78	2	18.1	0.6	3	2	95	6	3
Yb ppm	68	2	2.00	0.1	7	1	5	26	13
Zn ppm	81	2	583	24	4	8	6300	0	0

Table 3.–Continued--Method blank results 3s values are considered the lower limit of detection (LOD), and 5s values are considered the lower limits of determination (LLD)

Method blank	n	Mean	s	3s	5s
Ag ppm	30	-0.01	0.01	0.04	0.06
Al %	30	0.00006	0.00002	0.00005	0.0001
As ppm	30	-0.7	1	4	6
Au ppm	30	-0.1	0.06	0.2	0.3
Ba ppm	30	0.08	0.09	0.3	0.5
Be ppm	30	0.0005	0.001	0.002	0.003
Bi ppm	30	0.5	0.5	1	2
Ca %	30	-8E-06	0.0001	0.0003	0.0004
Cd ppm	30	0.005	0.02	0.05	0.08
Ce ppm	30	0.2	0.1	0.4	0.7
Co ppm	30	-0.2	0.1	0.3	0.5
Cr ppm	30	-0.1	0.2	0.6	0.9
Cu ppm	30	0.001	0.01	0.03	0.1
Eu ppm	30	0.01	0.003	0.01	0.01
Fe %	30	-0.00004	0.00001	0.00004	0.00006
Ga ppm	30	-0.06	0.08	0.2	0.4
Ho ppm	30	0.03	0.02	0.06	0.1
K %	30	0.00002	4E-05	0.0001	0.0002
La ppm	30	-0.02	0.03	0.08	0.1
Li ppm	30	0.0005	0.01	0.03	0.05
Mg %	30	-0.00002	1E-05	4E-05	7E-05
Mn ppm	30	0.0005	0.005	0.01	0.02
Mo ppm	30	-0.08	0.1	0.4	0.6
Na %	30	0.000002	6E-06	2E-05	3E-05
Nb ppm	30	-0.5	0.2	0.5	0.8
Nd ppm	30	0.05	0.07	0.2	0.4
Ni ppm	30	0.1	0.05	0.1	0.2
P %	30	0.00003	4E-05	0.0001	0.0002
Pb ppm	30	-0.3	0.3	1.0	2
Sc ppm	30	-0.005	0.004	0.01	0.02
Sn ppm	30	-0.2	0.5	1	2
Sr ppm	30	0.04	0.1	0.4	0.7
Ta ppm	30	-0.07	0.2	0.5	0.8
Th ppm	30	-0.09	0.05	0.1	0.2
Ti %	30	-0.00006	0.00009	0.0003	0.0005
U ppm	30	0.3	0.5	1	2
V ppm	30	-0.01	0.008	0.02	0.04
Y ppm	30	-0.005	0.001	0.003	0.005
Yb ppm	30	0.01	0.004	0.01	0.02
Zn ppm	30	-0.05	0.03	0.1	0.2

Bibliography

- Canadian Certified Reference Materials Project, 1995, Certificate of Analysis, TILL-1, TILL-2, TILL-3, and TILL-4.
- Crock, J.G., Lichte, F.E., and Briggs, P.H., 1983, Determination of elements in National Bureau of Standards geological reference materials SRM 278 obsidian and SRM 688 basalt by inductively coupled plasma-atomic emission spectroscopy: *Geostandards Newsletter*, v. 7, no. 2, p. 335-340.
- Govindaraju, K., ed., 1994, 1994 Compilation of working values and sample description of 383 geostandards: *Geostandards Newsletter*, v. 18, Special Issue, 158 p.
- Lichte, F.E., Golightly, D., and Lamothe, P.J., 1987, Inductively coupled plasma-atomic emission spectrometry, *in* Baedecker, P., ed., *Methods for Geochemical Analysis: U.S. Geological Survey Bulletin 1770*, p. B1-B10.
- National Bureau of Standards (now National Institute of Standards and Technology), 1977, Certificate of analysis: U.S. Department of Commerce, Washington, D.C.
- National Institute of Standards and Technology, 1992, Certificate of analysis: U.S. Department of Commerce, Gaithersburg, MD.
- Perkin-Elmer Corp., Optima 3000 Software, version 2.0, 1999.
- Thompson, M. and Walsh, J.N., 1983, A handbook of inductively coupled plasma spectrometry, p. 16-36.