

9.0 PROCEDURE

9.1 General:

Before operating the instrument for the first time, it is required that the analyst should familiarize him- or herself with the operating controls and other aspects of the analysis. Standard Operating Procedure (SOP) and the PE Plasma Operation Manual should be studied. Proficiency will be mandated of all analysts performing analyses.

9.2 Sample Handling:

9.2.1 Collection:

All samples received and handled by this section will meet the requirements of the Policy and Procedure Manual of the Water Quality Lab pertaining to handling and preservation criteria.

Samples will be logged by the Sample Custodian and assigned a unique log number. Pertinent information, i.e. location, date and time of collection, collector's name, sample type, date received by this section, person's name receiving the sample, preservation used, and any special remarks concerning the sample will be entered into the computer.

A batch number will be assigned by this section. All batches will be numbered with reference to protocol and year, for example pretreatment batches will be assigned PRE010102, representing the date prepped Jan 1, 2002.

9.2.2 Preservation and Preparation:

All samples will be prepared following Spectroscopy's Standard Operating Procedures Sample Preparation Manual. Samples will be preserved with 2.5 ml HNO₃ prior to admitting to lab. It will be the responsibility of sample custodian to check for pH. See QA Manual Section 1.0.

9.2.3 Storage:

Samples that may have legal ramifications due to not being within regulatory compliance limits are to be stored in a secure area but for no longer than their allowed holding times. Confer with the supervisor before discarding any samples.

9.2.4 Holding time:

All samples will be stored in the spectroscopy refrigerator until samples are prepped. Analysis will be performed within the time specified in the protocol, not to be held longer than 3 months for metal analyses.

9.3 Functions of Controls and Display of ICP Instrument

9.3.1 Checking Status:

The system informs you of the current status of the instrument in the following ways:

- ! **Instrument Diagnostics windows:** This window displays information about the spectrometer environment, spectrometer status, detector status, and plasma generator status. To open the Instrument Diagnostics window: In the **System** menu, click on **Diagnostics**. Click on the **Update** button in this window to update the information.
- ! **System Monitor status icons:** The icons in the System Monitor window give a graphical representation of the current condition of the instrument. To open the System Monitor window: In the **System** menu, click on **System Monitor**.
- ! **Error Messages:** Error messages are displayed in the Message Log. To view the message log: In the **System** menu, click on **Message Log**.
- ! **Indicator Lights:** The indicator lights on the front of the instrument (Plasma On, System Power, and System Ready) indicate general operating status.

9.3.2 System Monitor Status Icons:

To open the System Monitor windows:

- < In the **System** menu, click on **System Monitor**. See Figure 9.1.

Icon Color

Meaning

Green	The system component is powered on.
Yellow	For the plasma torch icon, this indicates that the ignition process has started. For other icons, yellow indicates a caution status. A problem has occurred within this part of the system. Read the Message Log for more details. This

problem requires attention since it may eventually cause the system to shut down.

Red A problem has occurred within this part of the system. Read the Message Log for more details. This problem needs to be fixed before you can operate the instrument.

Gray This part of the system is functioning properly.

9.3.3 Information Display Windows:

- ! **Spectra Display Window:** View spectra and associated information, as each spectrum is generated during an analysis.

- ! **The Results Display Window:** View a list of the analytical results as they are generated from an analysis.

- ! **Calibration Display Window:** View the calibration curves
 - < during Method Development to ascertain an appropriate curve fitting algorithm
 - < during routine analysis to provide a visual check of the calibrations.

- ! **Examine Spectra Window:** View spectra that have been saved to a data set so that you can:
 - < modify method parameters directly from the Examine window
 - < manipulate the appearance of the spectra to aid in method development and set up spectral displays for reports.

- ! **Analytical Sequence Window:** View the samples listed according to the analytical sequence.

- ! **Autosampler Loading List:** View the samples listed according to the autosampler location number.

9.4 System Initialization and Setup:

9.4.1 Interlocks:

Interlocks are designed to ensure operator safety and protect the instrument from damage. The main system interlocks are described below.

The following interlocks must be satisfied in order to ignite the plasma. If any of these interlocks is interrupted while the plasma is on, the plasma will automatically be shut down.

- ! The front and side doors on the sample compartment must be closed before you can ignite the plasma.
- ! Argon pressures for the torch must be correct.
- ! Cooling water must be flowing to the RF coil.

The following interlocks are for the detectors. If any of these interlocks is interrupted while the system is operating, the detectors will be shut down.

- ! The argon purge for the detectors and nitrogen purge for the spectrometer optics, monitored by flow switches, must be functioning properly.
- ! The temperature of the optics housing must be 36-40C.
- ! The temperature of the detector must be less than -30C.

9.4.2 Initialization Steps:

The Main Instrument switch is used to turn on the spectrometer (and is normally left on). When the switch is first turned on, a system clock checks to see how long the instrument has been off. Based on this time, the system determines what startup routines must be carried out before the instrument can be operated and displays the wait time on the computer (13 to 73 minutes). The software sends instructions to the instrument for the required initialization sequence. If the software is not on initially, the instrument will wait up to 10 minutes for instructions; otherwise it carries out its startup routine independently. The RF generator, which is switched on by the RF generator standby switch, takes about one minute to warm up.

Once the plasma has been ignited, you should wait one hour for the system to stabilize before running samples.

The spectrometer initialization sequence is as follows:

1. As soon as the Main Instrument switch is turned on, the spectrometer begins a high argon and high nitrogen purge. If the instrument has been off for more than thirty minutes, the instrument begins a A cold start \cong , followed by a A hot start \cong . If the instrument had been off for less than thirty minutes (and was previously in the ready state), the instrument begins a A hot start \cong .

2. The **A cold start** initialization sequence lasts 60 minutes. During this sequence, the instrument:
 - a. runs a high nitrogen purge of the optical system for 60 minutes.
 - b. checks to make sure that the thermal plenum is at the required temperature (36-40C) and checks functioning of slit motor and transfer optics motors.
 - c. checks the temperature inside the optics system.

3. The **A hot start** initialization sequence lasts 13 minutes. During this sequence, the instrument:
 - a. performs a 10-minute high argon and high nitrogen purge and checks the slit, and transfer optics motors.
 - b. turns on the cooling water, the thermoelectric coolers for the detectors, and then turns on the power to the detectors and switches to a low argon and nitrogen purge.
 - c. begins monitoring temperature of the detectors.
 - d. checks for proper functioning of the shutter, checks for dark current, turns on the mercury lamp and checks the mercury peak.

4. The spectrometer sends a **>system ready** message to the computer. Provided the software and the RF generator are on, the green **System Ready** indicator on the front of the instrument lights up.

9.4.3 Hardware Setup:

The following hardware settings are controlled by the software.

RF Power: Power levels can be adjusted in 5-watt increments. **Recommended setting 1450 watts.**

Plasma and auxiliary argon flow rates: Flow rates can be automated during the analysis with specific flow rates for each element if desired. Plasma argon is adjustable in 1 L/min increments. Auxiliary argon is adjustable in 0.1 L/min increments. **Recommended setting 0.5 L/min Auxiliary Flow.**
Recommended setting 15 L/min Plasma Flow.

Nebulizer argon flow rate: The flow rate is automatically controlled using a mass flow controller in .001 L/min increments. **Recommended setting 0.75 L/min.**

Pump rate: By specifying the desired flow rate (adjustable in 0.1 mL/min increments) and the tubing diameter, the software calculates the pump speed. A fast pump speed can be set for the read delay or rinse cycle. **Recommended setting 1.5 mL/min.**

Equilibration Time: **Recommended setting 15 sec.**

Nitrogen purge rate: A high or low rate can be specified.

Mercury wavelength recalibration: An automatic recalibration can be specified to occur at defined intervals (from every 3 minutes to 999 hrs).

Dark current measurement: A shutter under software-control can be closed to block light from reaching the detector in order to measure dark current.

The software also has a built-in optimization function for plasma viewing. The area of the plasma viewed by the optical system can be adjusted horizontally (in 0.4 mm steps) and vertically (in 1 mm steps) using a moveable transfer optic. A software function lets you automatically adjust the horizontal alignment.

9.5 Setting Up the Instrument:

9.5.1 Checking the System:

Proper setup and routine maintenance are required to keep your instrument in proper working condition so that you will get acceptable performance.

! **Exhaust Vent**

- < Check that your vent system is switched on, is working properly, and is not blocked.

! **Argon Supply**

- < Be sure that an adequate supply of argon is available and is connected to the system.
- < Check that the cylinder valve is open and the regulator for the argon outlet pressure is set within the following range: 70-120 psig.
- < Check for leaks at the gas connections on the back of the instrument.

! **Nitrogen purge gas supply**

- < Be sure that an adequate supply of nitrogen is available and connected to the system.
- < Check that the cylinder valve is open and the regulator for the nitrogen outlet pressure is set within this range: 30-120 psig.
- < Check for leaks at the gas connections on the back of the instrument.

! **Shear gas supply**

- < For compressed air, check that the shut-off valve on the air dryer filter is open and the pressure for the air compressor is set within this range: 75-100 psig.
- < Check for leaks at the gas connections on the back of the instrument.

! **Cooling water supply**

- < Check that the cooling water supply is connected to the back of the instrument.
- < Using a chiller, check the electrical connections. Check that the External Pressure Reducer has been set to 45 +/-5 psig. Be sure the chiller is turned on and has been filled. Periodically check the liquid level.
- < Check for any leaks in the plumbing connections.

! **Sample compartment**

- < Check that the spill tray is in place underneath the Quick-Change Torch module.
- < Inspect the purge window that fits into the ceramic purge extension. The window should not be cloudy or dirty.

! **Sample introduction system**

- < Inspect the torch, glassware, and aerosol injector tube. The glassware should be clean. Small amounts of deposits are acceptable in most cases. Also check for signs of melting.
- < Check that the RF coil is clean. This prevents arcing across the coils.

Leakage of air into any part of the torch, nebulizer or spray chamber may cause ignition problems. Therefore, be sure to regularly check the following four items:

Be sure that all gas fittings to the torch are finger-tight.

Check that the nebulizer/end cap is tightly secured to the spray chamber.

Check that the sample capillary tubing is attached to the nebulizer sample inlet. The tubing should be clean and in good condition.

Check that the drain fitting is secured on the spray chamber drain. A loose-fitting drain can cause pressure leaks and consequent plasma instability.

! **Peristaltic pump and tubing**

- < Check that the pump tubing is in good condition, has no flat spots, and is correctly installed around the pump head.

- < Replace tubing daily.
- < Check that the pump rollers are clean and can move freely.
- ! **Drains**
- < Check that the spray chamber drain is properly set up on the pump so that waste is pumped out of the spray chamber. Replace the drain tubing if it has deteriorated.
- < Check the spray chamber drain and the drain underneath the sample shelf lead to the drain bottle. Empty the drain bottle if necessary. Dispose of waste properly.

- ! **Autosampler**
- < Install the tray in the autosampler and place the wash beaker in the rear left of the tray (location).
- < Check that the sampling probe is installed at the correct height and the probe capillary is attached to the pump tubing for the sample.

- ! **Printer**
- < Check that the printer has an adequate supply of paper to print your results.

9.5.2 Starting the Instrument and Accessories:

It is recommended that you leave the Main Instrument switch on even when the instrument is not in use. If the Main Instrument switch is turned off, the system will need to reinitialize when power is switched back on.

The RF Generator Standby switch may also be left on when the instrument is not in use. It is only necessary to turn off this switch if you are replacing the RF coil.

To conserve the argon purge gas for the detectors, you can use the Auto Shutdown/Startup feature. Automatic Shutdown is an alternative to turning off the power to the main instrument, which also then shuts off the purge gas.

A. Starting up the instrument:

1. If you have not already done so, turn on the gases and the cooling water chiller.
2. If the Main Instrument switch has been turned off, turn it on.
3. Close the doors to the sample compartment. Be sure that the front door is fully closed and secured.

Be sure to close and secure the doors to the sample compartment before turning on or off the RF Generator Standby switch.

4. Turn the RF Generator Standby switch on if it is not already on. With both the Main Instrument and RF Generator switch on, the blue System Power light goes on.

5. Turn on the autosampler if it is not already on.

B. Starting the Computer and Software

1. Switch on the computer, monitor, and printer.

2. **For Window 95:** click on the **Start** button. Select the **ICP WinLab** program.

3. Verify instrument configuration:

< Check that the autosampler type and tray are correct.

< Check that the peristaltic pump and tubing type are correct.

< Check that the correct plasma viewing position is set. Change if required, using the **Change to Axial** or **Change to Radial** commands in the **System** menu for the desired position.

C. Setting Up and Starting the Peristaltic Pump

! **Setting up the peristaltic pump:**

1. Check that the sample tubing and the drain tubing leading from the spray chamber are properly set up on the pump.

2. If the pump tubing is new, gently stretch it. Position the clips on the tubing in front of the tubing stops.

3. Replace the tubing clamps for each channel and swing the cam levers over to apply tension to the clamps.

4. Start the pump before igniting the plasma to purge air out of the sample capillary tubing.

! **Starting the pump:**

1. Place the sample capillary tubing in a container of a 1% nitric acid solution.

2. Display the Plasma Control window by clicking on the Plasma icon or by clicking on **Plasma Control** in the Tools menu. Click on the **Pump** button to turn on the pump.

3. Adjust the flow rate to 1mL/min. Before adjusting the rate, the Override Method box must be checked.

4. If necessary, adjust the tension on the pump tubing, one channel at a time. For the sample tubing, gradually tighten the adjustment screw until the liquid flows smoothly without bubbles. The drain tubing should have a segmented flow of liquid leading to the drain bottle. Bubbles in the drain tubing are normal.

D. Igniting the Plasma

The process requires approximately 90 seconds. Argon gas flows through the torch and spray chamber, purging the sample introduction system of air. With the argon continuing to flow, power is applied to the RF coil. After the RF power has ramped to approximately 1300 watts, a high voltage spark is injected into the argon flow causing the argon to ionize. The free electrons that are created then interact with the applied RF field to cause further argon ionization and for a plasma.

When you ignite the plasma, be sure to observe it closely through the viewing window. If the plasma is unstable, immediately click the Plasma switch to OFF in the Plasma Control window or press the red Emergency Plasma Off button above the sample compartment.

1. Observe the following precautions when igniting the plasma:
2. If you have not yet started the pump, place the sample capillary tubing in a container of 1% HNO₃. The pump automatically shuts off at the beginning of the ignition sequence and is restarted at the end of the ignition sequence.
3. Click the **Plasma** switch to **On** in the Plasma Control window to turn on the plasma. A successful ignition is indicated by a colored plasma icon on the switch in the Plasma Control window. The yellow **Plasma On** light on the front of the instrument will be on.
4. Immediately examine the plasma through the viewing window.

Note: During the ignition process, messages appear in the lower left corner of the screen to inform you of the status. The first time the plasma is ignited, a message reads **AInitial Purge** and the system counts down 75 seconds. This message is followed by **ASetting nebulizer flow**. After a **AFinal Purge** of 15 seconds, the plasma is ignited. This is followed by a **ASettling Time** of 15

seconds to give the plasma time to stabilize. At the end of this process, a message reads APlasma has ignited.

E. Examining the Plasma:

After igniting the plasma, you must examine to determine whether it is stable or unstable. A stable plasma will be situated just above the inner quartz tube in the torch and will have a bright discharge of the shape as seen Figure 9.2. An unstable plasma is indicated by a doughnut-shaped glow or discharge ring that forms around the top of the torch.

! To proceed if the plasma is stable

1. Wait minimum of 45 minutes before running samples. This ensures accuracy in your results, since it allows the temperature of the sample introduction system to fully stabilize.
2. Proceed to Optimizing and Verifying Performance.

! To proceed if the plasma is unstable

1. Click the **Plasma** switch to **Off** in the Plasma Control window or press the red **Emergency Plasma Off** button located above the sample compartment.

An unstable plasma is usually caused by air leaking into the system. To correct this problem, do the following:

1. Check that the gas fittings on the torch are finger-tight.
2. If the fittings are properly tightened, then check the argon gas connections at the back of the instrument.
3. Repeat the ignition procedure. Click the **Plasma** switch to **On**.

! Correcting Unsuccessful Ignition

If the plasma fails to ignite during the ignition sequence, the plasma switch will go to the off position, indicating that ignition has failed but may be attempted again. Check the following:

! Sampling technique

- < Be sure you are using the correct sampling technique settings for your application. Click on the Sampling button in the toolbar to display a dialog where you can select the sampling technique with the proper plasma conditions.

! Ignitor

- < Open the doors to the sample compartment. Check that the ignitor cable is plugged in.

! **Air Leaks**

- < An unstable plasma is usually caused by air leaking into any part of the torch, nebulizer, or spray chamber. To correct check the following:
 1. Check that the gas fittings on the torch are finger-tight. Do not over tighten them.
 2. If the gas fittings on the torch are properly tightened, check the argon gas connections at the back of the instrument.

! **Drain:**

- < Check that the fitting on the spray chamber drain is secure. Be sure that the pump is properly draining the spray chamber and that the drain liquid is not backing up into the spray chamber.

! **Nebulizer end cap:**

- < Check that it is tightly secured to the spray chamber.

! **Sample capillary and tubing:**

- < Check that one end is attached to the nebulizer and that the other end is immersed in solution.

9.6 Optimizing and Verifying Performance:

9.6.1 Performing the Background Equivalent Concentration (BEC) Test:

In performing the BEC test, you determine what concentration of analyte is equivalent to the plasma background at the analyte wavelength. This provides an indication of the sensitivity of the instrument. To do the test you first perform a Hg realign from the mnBEC work space. Second calibration by running a blank followed by a standard. The calibration establishes the relation between emission intensity and concentration. A calibration also takes into account the dark current noise. Once the calibration is complete, you close the shutter, thereby blocking the light from the plasma source from reaching the detector. At this point, you take a reading (at Zero intensity). This extrapolates the calibration line backward until it intercepts the concentration axis. The resulting concentration value is negative. By changing the sign to positive, you obtain the BEC value. See Figure 9.3.

1. Check that the method you created for the performance tests is open. If not, in the **File** menu, click on **Open Method**.

2. Click on the **Manual** icon to open the Manual Analysis Control window,.
3. In the Manual Analysis Control window, select **Print Log** to disable, click on the box to remove **Save Log**.
4. Aspirate the blank and click on **Analyze Blank**.
5. Aspirate the 1ppm manganese solution and click on **Analyze Standard**. This completes the calibration.
6. Aspirate the rinse solution.
7. In the **Tools** menu, click on **Spectrometer Control**.
8. In the Spectrometer Control window, close the shutter by clicking on the **Shutter** button. Check that the status of the Shutter changes to AClosed≡ in the status display.
9. Click on the Manual Analysis Control window to bring it to the front.
10. While continuing to aspirate the rinse solution, click on **Analyze Sample**. This will give you results in concentration units.
11. Click on the Spectrometer Control window to bring it to the front. Click on the **Shutter** button to change the shutter from Closed to the Auto position, which returns it to instrument control. Check that the status of the Shutter changes to AAuto≡ in the status display.
12. By closing the **Manual & Spectrometer Control** windows the results are sent to the printer.

Check the results. By taking the negative value and making it a positive value, this gives you the BEC. It should be less than or equal to 0.04 mg/L. Attempt a second run, if it is not within required value refer to Troubleshooting section 11.

9.6.2 Performing the Precision Test:

The coefficient of variation (CV) test expresses the short-term precision for several measurements for a strong emission line. This is often used as an indicator of noise associated with sample introduction.

1. Check that the method you created for the performance tests is open.
2. If you have just run the BEC test, the Manual Analysis Control window should be displayed. Click on the **Manual** icon.
3. In the Method Editor, check on the **Spectrometer, Read Time, Replicates** page and specify 10 replicates if not already on file.
4. In the Manual Analysis Control window, disable **Override Method** to use the Read Delay specified in the method. To disable it, click on the box to remove the check mark.
5. Aspirate the 1 ppm manganese. Select autosampler position 160, specific location for 1 ppm manganese.

6. In the Manual Analysis Control window, click on **Analyze Sample**.
7. Close **Manual** window to send results to printer.
Check the results. The RSD should be less than 1.0%. Also attempt a second run. If value is not acceptable refer to the Troubleshooting section 11.

9.6.3 Sodium bullet test:

This test allows you to visualize the sample flow in the plasma, so that you can check that the sample introduction system is working correctly.

1. Aspirate a 1000 mg/L solution of sodium.
2. Examine the plasma through the viewing window in the sample compartment door.

A yellow-orange bullet should be visible in the center of the discharge, and should extend from the base of the discharge to about 2-3 mm above the top of the RF coil. If the bullet height is unsatisfactory, adjust the nebulizer argon flow in the Plasma Control windows.

If no bullet appears or the bullet is faint:

- # Check that your sample does contain sodium at the required concentration.
- # Check that sample is being pumped to the nebulizer.
- # Check that the drain is being pumped properly.

If the above checks fail, turn off the plasma, then check the following:

- # Check that the nebulizer end camp is connected tightly.
- # Check the nebulizer spray pattern.
- # Check that the injector is not clogged.

9.6.4 Spectral Interference Check (SIC) Solutions:

When interelement corrections are applied, SIC solutions are needed containing concentrations of the interfering elements at level that will provide an adequate test of the correction factors. Use the following PE interference check standards:

- 1.) Interferents Alone---PEN930-0226
Method-dilute PEN930-0226 by a Factor of 10

Ca	Mg	Al	Fe
----	----	----	----

2.) Combination of interferences and analytes

Method-dilute 10ml of Interferant solution PEN930-0226 + 1ml of analyte solution PEN930-0227 to 100ml with 5% HNO₃.

Be	Ba
Mn	
V	

3.) Alternate Interferant A

Method-10ml PEN930-0226 + 10ml PEN930-0228 and dilute to 100ml with 5% HNO₃.

V	Mn
---	----

4.) Solution of interference and analytes

Method-10ml PEN930-0226 + 10ml PEN930-0228 + 1ml PEN930-0227 + 1ml PEN930-0229 dilute to 100ml

B	Na	Al	Mo
Ca	Si		
Fe	Mg		

Note: Record all data in Performance Logbook, calculate all %Recoveries.

9.7 Sample Analysis:

9.7.1 Setting up the Autosample Tray:

Before arranging samples in the autosampler using the autosampler loading list and/or analytical sequence; you must have entered all information for your samples and solutions in your method and in the sample information file.

X To load samples and solutions in the autosampler using the **Automated Analysis Control** window.

1. In the **Automated Analysis Control** window, click on the **Analyze** tab.

2. In the **Analysis** menu, click on **Analytical Sequence** to see a list of the locations of your samples and solutions. Print sequence on file.
 3. Load samples and solutions in an empty tray according to preselected sequence.
 4. Place each sample in autotray position corresponding to sequence on file.
- X To load samples and solutions in the autosampler when the **Manual Analysis** window is open.
1. In the **Analysis** menu, click on **Autosampler Loading List** to see a list of the locations of your samples and solutions.
 2. Load samples and solutions in an empty tray.
 3. In the **Analysis** menu, click on **AS Controls Load Tray**.
 4. Place each sample in autotray position corresponding to sequence on file.
 5. Click on **Go to Wash Loc.** In the Manual Analysis Control window.

9.7.2 To perform automated analyses:

- < See the following steps:

Step	What it does
Examining the Analytical Sequence	Before analyzing samples, check the sequence and make any changes.
Calibrating and Analyzing Samples	Solutions and samples are analyzed. A calibration is performed according to specifications.
Stopping and Restarting	During the analysis, you may need to stop operations and restart the analysis.
Analyzing Additional Samples	If you have completed an analysis and need to analyze additional samples.

9.7.3 Examining the Analytical Sequence:

Examine the analytical sequence before analyzing samples.

1. In the **Analysis** menu, click on **Analytical Sequence**. Verify that your sequence is correct. If it is correct, skip steps 2 and 3 below.
2. To make changes to the sequence for unknown samples, QC samples, and matrix check samples, go to the Sample Information Editor and make the changes in the Sample Information file.
3. To make changes to the sequence for calibration solutions, go to the Method

Editor Calib.

4. To analyze a selected group of samples that are a subset of samples in the sample information file, type the locations on the Automated Analysis Control window.

9.7.4 Calibrating and Analyzing Samples:

According to your calibration choices, select the appropriate steps below. All controls (**Analyze All**, **Calibrate**, and **Analyze Samples**) are located on the Automated Analysis Control Analyze page. You select the calibration parameters in the Method Editor Calib pages and view the calibration curves in the Calibration Display window or view calibration data in the Results Display window.

A. Calibrating

Select a calibration option and analyze mode

To do this:

To generate a new calibration and continue with samples (overrides Initial Calib page of the Method Editor)

To generate a new calibration, examine it, and continue with samples

To recalibrate

To analyze samples when a calibration is automatically recalled with the method (as directed on the Initial Calib page of the Method Editor)

To manually recall a calibration and then analyze samples

Use this:

Click on **Analyze All**.

Click on **Calibrate**. Examine the calibration and recalibrate if desired. When ready click on **Analyze Samples**.

Click on **Calibrate**

Click on **Analyze Samples** to analyze samples using this calibration.

In the **Analysis** menu, click on **Recall Calibration...** Select the results data set that contains the desired calibration. Then, click on **Analyze Samples**.

B. Stopping and Restarting

When you stop an analysis, the autosampler probe does not automatically go back to the wash, but remains where you stopped it.

Select stop and restart options from the following:

To do this:

To stop an analysis

Use these steps:

1. In the Automated Analysis Control window Analyze page, click on the button that you used to start the analysis, press **F8 (Cancel)**, or select **Cancel Analysis...** from the Analysis menu.
2. In the Stopping and Analytical Sequence dialog, select the option that indicates how you want to stop.

To change Set Up Option after an analysis has begun.

1. Stop the analysis as described above.
2. Click on the tab for the Set Up page
3. Change set up options. You can select whether or not to save data, print a log of the results, or perform an auto wavelength realign. You can also program an automatic shutdown.

To restart analysis from a selected solution

1. In the Automated Analysis Control window, click on the button that you originally used to start the analysis.
2. In the Continuing an Analytical Sequence dialog, select where you want the analysis to continue.

To restart an analysis from beginning

In the Automated Analysis Control window, from the click on **Reset** and then on a button (**Analyze All**, **Calibrate**, or **Analyze Samples**) the starts an analysis.

C. Analyzing Additional Samples:

If you have completed an analysis using an autosampler and need to analyze some additional samples, use the following procedure.

For automated analyses using the Automated Analysis Control Window:

1. Open the method you were previously using in the Method Editor window (if it is not already open). Click on the **Calib** tab in the Method Editor. Remove the calibration blank and all calibration standards listed on the Ids and Location page. To remove them, you must click on the row label to select the entire row. Then, in the **Edit** menu, click on **Clear**.
2. Click on the **Initial Calib, Special Options** tab (also in the Calib section of the method). Select the option for **Use calibration curves from previous method**.
3. Save the method under a different name. This will keep the original method intact. In the **File** menu, click on **Save As... Method**. Type a name for the new method.
4. In the Automated Analysis Control window, remove the methods that are currently listed. Enter the name of the new method you just created.
5. Select the autosampler locations for the additional samples that you want to analyze.
6. Click on the **Analyze** tab in the Automated Analysis Control window.
7. Click on **Analyze All**.

Note: The system will use the calibration from the previous method and analyze the additional samples.

9.8 Expression of Data:

9.8.1 Generating Reports:

If you saved data in a Results data set (as requested in an analysis control window), you can use the WinLab Reporter to generate reports.

A. To start the WinLab Reporter:

- < In the **File** menu, click on **Utilities... Reported**.

B. To start the WinLab Library Manager:

Archive result data sets on a routine basis to avoid filling up space on the hard disk and potentially slowing down instrument operation. Library maintenance procedures include archiving, deleting, and packing data sets. Exit software and enter thru Optima 3000 file, click on **Libman** icon.

9.8.2 Laboratory Information Management System (LIMS):

All laboratory data is entered via the **SQLLIMS** system as discussed in the Policy and Procedure Manual section 4.0. Laboratory analysts are responsible for the data and result entry of the data they produce each day.

When available the result entry using **ARE** will be used for instrument generated data. When this system is not available manual data entry will be used.

9.8.3 Reporting MDL's- See Attachment for all current metals MDLS's:

See Attachment 1.1 section 1 for all current MDL=s.

9.8.4 Result Approval and Verification:

Result approval is performed at the TASK level only. First stamp will be the analyst who performed the test and this analyst will enter data. A second analyst other than the analyst entering the data will verify and approve all data. Verification will be done 7 days from entry of data. Data out of QA/QC requirements supervisor must approve data prior to entry.

9.8.5 Changes in Data :

Changes in data will be performed at the TASK level only. All changes in LIMS must be performed by supervisor level only.

9.9 Shut Down:

There are two ways to extinguish the plasma. You can extinguish the plasma manually or you can extinguish it automatically using Automatic Shutdown.

A. To automatically extinguish the plasma and put the spectrometer on standby:

If you have selected Automatic Shutdown and indicated how you want the system to shut down, the system will automatically shut down and flush the sample introduction system. The system next extinguishes the plasma and puts the spectrometer into standby if that option was selected. The software remains on.

B. To manually extinguish the plasma:

1. Flush the sample introduction system for five minutes with the plasma on. During this five minutes, if you analyzed aqueous solutions, flush with 1% HNO₃.

2. Extinguish the plasma by clicking on the **Plasma Off** switch in the Plasma Control Window.

C. To manually put the spectrometer on standby:

1. In the **System** menu, click on **Auto Startup/Shutdown**.

2. In the dialog, select **Shutdown enabled...** and **immediately on OK**. Be sure **Put Spectrometer into Standby** is selected.

3. If the plasma has already been extinguished, be sure that the **Wash** option is not selected and that **Turn off Plasma and Pump** is not selected.